# Vitassay qPCR

## H. pylori ClariRes

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96 <sup>TM</sup> Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

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| Termocicladores con bloque de bajo perfil    | Termocicladores con bloque de alto perfil   |
|--|---|
| <b>Agilent Technologies</b>                  | <b>Abbott</b>                               |
| AriaMx Real-Time PCR System                  | Abbott m2000 RealTime System                |
| <b>Applied Biosystems</b>                    | <b>Applied Biosystems</b>                   |
| 7500 Fast Real-Time PCR System               | 7300 Real-Time PCR System                   |
| 7500 Fast Dx Real-Time PCR System            | 7500 Real-Time PCR System                   |
| QuantStudio™ 12K Flex 96-well Fast           | 7900 HT Real-Time PCR System                |
| QuantStudio™ 6 Flex 96-well Fast             | ABI PRISM 7000                              |
| QuantStudio™ 7 Flex 96-well Fast             | ABI PRISM 7700                              |
| QuantStudio™ 5 Real-Time PCR System          | QuantStudio™ 12K Flex 96-well               |
| QuantStudio™ 3 Real-Time PCR System          | QuantStudio™ 6 Flex 96-well                 |
| StepOne Plus™ Real-Time PCR System           | QuantStudio™ 7 Flex 96-well                 |
| StepOne™ Real-Time PCR System                | QuantStudio™ 5 Real-Time PCR                |
| ViiA™ 7 Fast Real-Time PCR System            | QuantStudio™ 3 Real-Time PCR System         |
| <b>BIONEER</b>                               | ViiA™ 7 Real-Time PCR                       |
| Exicycler™ 96                                | <b>Analytik Jena Biometra</b>               |
| <b>Bio-Rad</b>                               | TOptical                                    |
| CFX96™ Real-Time PCR Detection System        | qTOWER 2.0                                  |
| Mini Opticon™ Real-Time PCR Detection System | <b>BIONEER</b>                              |
| <b>Cepheid</b>                               | Exicycler™ 96                               |
| SmartCycler®                                 | <b>Bio-Rad</b>                              |
| <b>Qiagen</b>                                | CFX96™ Deep Well Real-Time PCR Detection    |
| Rotor-Gene® Q                                | iCycler iQ™ Real-Time PCR                   |
| <b>Roche</b>                                 | iCycler iQ™ 5 Real-Time PCR                 |
| LightCycler® 480 Real-Time PCR System        | MyiQ™ Real-Time PCR Detection System        |
| LightCycler® 96 Real-Time PCR System         | MyiQ™ 2Real-Time PCR Detection System       |
| Cobas z480 Analyzer                          | <b>Cepheid</b>                              |
|  | SmartCycler®                                |
|  | <b>DNA-Technology</b>                       |
|  | DTlite Real-Time PCR System*                |
|  | DTprime Real-time Detection Thermal Cycler* |
|  | <b>Eppendorf</b>                            |
|  | Mastercycler™ ep <i>realplex</i>            |
|  | <b>Qiagen</b>                               |
|  | Rotor-Gene® Q                               |
|  | <b>Stratagene / Agilent Technologies</b>    |
|  | Mx3000P™ Real Time PCR System               |
|  | Mx3005P™ Real Time PCR System               |

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## Intended use

Vitassay qPCR *H. pylori* ClariRes allows the detection of *Helicobacter pylori* and clarithromycin (CLR) resistance in *H. pylori* by real-time PCR in human stool samples and gastric tissue biopsies. The product is intended for use in the diagnosis of *Helicobacter pylori* and clarithromycin (CLR) resistance in *H. pylori* infections alongside clinical data of the patient and other laboratory tests outcomes. The assay for detecting CLR resistance in *H. pylori* is based on detection of point mutations in the 23S rRNA gene.

## References

Vitassay qPCR *H. pylori* ClariRes 4x8-well strip, low profile 7041034

Vitassay qPCR *H. pylori* ClariRes 4x8-well strip, high profile 7042034

## Materials/reagents provided

| Reference             | Reagent/Material                                  | Colour | Amount         |
|-----------------------|---|--------|----------------|
| 7041S034/<br>7042S034 | <i>H. pylori</i> ClariRes strips low/high profile | -      | 4x8-well strip |
| 7C034                 | <i>H. pylori</i> ClariRes Positive Control        | red    | 1 vial         |
| 7001A                 | PCR grade water                                   | white  | 1 vial x 1 mL  |
| 7002B                 | Resuspension Buffer                               | green  | 1 vial x 1 mL  |
| 7003N                 | Negative control                                  | yellow | 1 vial x 1 mL  |
| 7004O                 | Optical caps                                      | -      | 4x8 cap strip  |

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)

- Filter tips
- Powder-free disposal gloves

## Summary

*Helicobacter* species have a wide host range and can be subdivided into two major lineages, the gastric species (*H. pylori*, *H. felis*, *H. mustelae*, *H. acinonychis*, *H. heilmannii*) and the enterohepatic (nongastric) species (*H. hepaticus*). Among them, *H. pylori* and *H. heilmannii* infect both humans and result in gastritis and even in ulcer disease.

*Helicobacter pylori* is a common bacterium infecting about half of world's population, with higher prevalence in developing countries, where *H. pylori* could infect up to 80% of the population, than in developed ones. *H. pylori* is associated with the development of gastrointestinal disorders as chronic gastritis, peptic ulcer, and gastric adenocarcinoma. *H. pylori* is also involved in the development of other extra-gastric disorders such as mucosa-associated lymphoid tissue lymphoma (MALT), idiopathic thrombocytopenic purpura, vitamin B12 deficiency, and iron deficiency. Eradication of *H. pylori* could help in the management of these *H. pylori*-associated disorders. For the last two decades, the recommended treatment for *H. pylori* eradication is the standard triple therapy, using a proton pump inhibitor or ranitidine bismuth citrate, combined with clarithromycin and amoxicillin or metronidazole.

Due to the increase in the prevalence of *H. pylori* resistance to antibiotics, triple therapy with clarithromycin is no longer the best treatment for *H. pylori*, especially in some areas where the local resistance to this antibiotic is higher than 20%.

Recently, a polymerase chain reaction has been developed as an alternative tool for the detection of clarithromycin resistance. This technique accurately assesses mutations in the peptidyl transferase region encoded in domain V of *H. pylori* 23S ribosomal RNA gene conferring clarithromycin resistance.

## Principle of the test

Vitassay qPCR *H. pylori* ClariRes test is based on the real-time amplification of a conserved region of the *ureB* gene for the identification of *H. pylori* and 23S rRNA gene for the identification of Clarithromycin resistance and/or Clarithromycin wild-type sequence in the 23S rRNA. After DNA isolation, the presence of *Helicobacter Pylori* and/or Clarithromycin resistance and/or Clarithromycin wild-type sequence in the 23S rRNA is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *H. pylori* ClariRes test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an

internal control allows the detection of a possible reaction inhibition. The amplification of the Point mutations in the 23S rRNA gene of *H. pylori* (A2142G and A2143G), which confer resistance to Clarithromycin DNA target sequence is detected through the FAM channel, *H. pylori* DNA target in ROX channel, Clarithromycin wild-type sequence in the 23S rRNA gene DNA target in HEX, VIC or JOE channel whereas the internal control (IC) in Cy5 channel (depending on the equipment used select the proper detection channel, see Attached II).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## Procedures

### Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® RSC Blood DNA Kit, Maxwell® 16 instrument (Promega).

ZP02004 MagPurix Tissue DNA Extraction Kit and ZP02006 MagPurix Bacterial DNA Extraction Kit (biopsies specimens), using the MagPurix 12A instrument (Zinexts Life Science Corp)

ZP02011 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit A y ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B (stool samples), using the MagPurix 12A instrument (Zinexts Life Science Corp).

Invisorb® Spin Universal Kit (Stratec)\*.

NucleoMag® Pathogen (Macherey Nagel).

NucleoSpin RNA Virus (Macherey Nagel).

MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.).

NX-48 Stool DNA Kit using the Nextractor® NX-48 system, (Genolution).

\*In order to improve the yield and quality of bacterial DNA from stool or biopsies specimens, we recommend the sample pretreatment with Lysozyme at 37°C as described on “Instructions for use” for bacterial DNA isolation. In addition, usage of small elution volumes may raise the DNA/RNA concentration.

## **Positive control preparation**

Reconstitute the lyophilized *H. pylori* ClariRes Positive Control (red tube) with 100 µL of PCR grade water (transparent tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

## **Reaction setup**

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.

- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

### Programme your thermocycler

Set your thermocycler following the conditions below:

| Step                                      | Temperature | Time   | Cycles |
|---|-------------|--------|--------|
| Initial denaturation                      | 95°C        | 2 min  | 1      |
| Denaturation                              | 95°C        | 10 sec | 45     |
| Annealing/Extension<br>(Data collection*) | 63°C        | 50 sec |        |

Set the fluorescence data collection during the extension step (\*) through the FAM (Clarithromycin resistance), ROX (*H. pylori*), HEX, JOE or VIC channels (Clarithromycin wild-type sequence in the 23S rRNA) and Cy5 channel (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none (Attached II).

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

#### Positive control

The positive controls used in each run, must show an amplification curve for FAM (Clarithromycin resistance), ROX (*H. pylori*), HEX, JOE or VIC channels (Clarithromycin wild-type sequence in the 23S rRNA), which validates the reaction.

#### Negative control

The negative controls included in each run, must show the absence of signal in ROX, HEX and FAM which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

| Clarithromycin resistance (FAM) | <i>H. pylori</i> (ROX) | Clarithromycin wild-type sequence (HEX) | Internal control (Cy5) | Negative Control | Positive Control | Interpretation   |
|---------------------------------|------------------------|---|------------------------|------------------|------------------|--|
| +                               | +                      | +                                       | +/-                    | -                | +                | <i>H. pylori</i> positive, Clarithromycin resistance positive and Clarithromycin wild-type sequences in the 23S rRNA positive <sup>a</sup> |
| +                               | +                      | -                                       | +/-                    | -                | +                | <i>H. pylori</i> positive and Clarithromycin resistance positive; Clarithromycin wild-type sequence in the 23S rRNA negative <sup>b</sup>  |
| -                               | +                      | +                                       | +/-                    | -                | +                | <i>H. pylori</i> positive and Clarithromycin wild-type sequence in the 23S rRNA positive, Clarithromycin resistance negative <sup>c</sup>  |
| +                               | -                      | +                                       | +/-                    | -                | +                | <i>H. pylori</i> negative; Clarithromycin resistance positive and Clarithromycin wild-type sequences in the 23S rRNA positive              |
| -                               | -                      | -                                       | +                      | -                | +                | <i>H. pylori</i> negative, Clarithromycin resistance negative and Clarithromycin wild-type sequences in the 23S rRNA negative              |
| -                               | -                      | +                                       | +/-                    | -                | +                | Clarithromycin wild-type sequence in the 23S rRNA positive   |
| +                               | -                      | -                                       | +/-                    | -                | +                | Clarithromycin resistance positive   |
| -                               | +                      | -                                       | +/-                    | -                | +                | Experiment fail  |

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

<sup>a</sup> Presence of Clarithromycin resistant strains and wild-type simultaneously in the same clinical sample. Similar fluorescent levels in HEX and FAM channels should be seen.

On the contrary, if differences are noted between the fluorescent signal from both channels look at the possibility of interpretation b or c.

<sup>b</sup> Small level of fluorescence can be observed in HEX channels.

<sup>c</sup> Small level of fluorescence can be observed in FAM channels.

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

206 human gastric tissue biopsies from symptomatic patients were tested by Real Time PCR using Vitassay qPCR *H. pylori* ClariRes and RIDA@GENE *Helicobacter pylori* (R-biopharm) for *H. pylori*. *H. pylori* was detected in 99 samples by both kits. In addition to this, by *H. pylori* + *Clarithromycin resistance* Real Time PCR was able to detect 3 more positives samples (The low amount of DNA template of these samples is below the detection limit of the method used).

102 human gastric tissue biopsies from symptomatic patients were tested by Real Time PCR using Vitassay qPCR *H. pylori* ClariRes and RIDA@GENE *Helicobacter pylori* (R-biopharm) for Clarithromycin resistance. Clarithromycin resistance was detected in 23 samples by both kits. Among these, there are 2 samples that are sensitive strains to Clarithromycin according to Vitassay qPCR *H. pylori* ClariRes and resistant strains to Clarithromycin according to R-biopharm. They are sensitive strains to Clarithromycin according to sequencing, confirming our results.

102 human gastric tissue biopsies from symptomatic patients were tested by Real Time PCR using Vitassay qPCR *H. pylori* ClariRes and RIDA@GENE *Helicobacter pylori* (R-biopharm) for Clarithromycin wild-type sequences in the 23S rDNA. Clarithromycin wild-type sequences in the 23S rDNA was detected in 77 samples by both kits. Among these, there are 2 samples that are sensitive strains to Clarithromycin according to Vitassay qPCR *H. pylori* ClariRes and resistant strains to Clarithromycin according to R-

biopharm. They are sensitive strains to Clarithromycin according to sequencing, confirming our results.

Vitassay qPCR *H. pylori* ClariRes shows a high sensitivity and specificity to detect *H. pylori* and clarithromycin resistance in *H. pylori*.

### **Analytical sensitivity**

The analytical sensitivity was determined by analysis of 10-fold dilution series of *H. pylori*, clarithromycin resistance in *H. pylori* and Clarithromycin wild-type sequences in the 23S rDNA templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  DNA copies per reaction.

### **Analytical specificity**

The analytical specificity for Vitassay qPCR *H. pylori* ClariRes was tested within the panel of different microorganisms.

No cross-reactivity was seen between any of the species, except *H. heilmanni* whose presence can also cause various chronic gastrointestinal diseases in humans:

| Cross-reactivity assay                            |  |                                       |
|---|--|---------------------------------------|
| <i>Helicobacter felis</i> *                       | <i>Campylobacter upsaliensis</i>                       | <i>Candida albicans</i>               |
| <i>Helicobacter hepaticus</i> *                   | <i>Campylobacter hyointestinalis</i> *                 | <i>Arcobacter butzleri</i>            |
| <i>Helicobacter cinaedi</i> *                     | <i>Proteus vulgaris</i>                                | <i>Pseudomonas aeruginosa</i>         |
| <i>Helicobacter heilmannii</i> *                  | <i>Citrobacter freundii</i>                            | <i>Enterococcus faecalis</i>          |
| <i>Shigella flexneri</i>                          | <i>Staphylococcus aureus</i> subsp. <i>aureus</i>      | <i>Yersinia enterocolitica</i> O:3*   |
| <i>Shigella dysenteriae</i>                       | <i>Serratia liquefaciens</i>                           | <i>Yersinia enterocolitica</i> O:9    |
| <i>Salmonella typhi</i>                           | <i>Vibrio parahaemolyticus</i>                         | <i>Bacteroides fragilis</i>           |
| <i>Salmonella paratyphi A</i>                     | <i>Clostridium difficile</i>                           | Adenovirus serotypes 40/41            |
| <i>Salmonella paratyphi B</i>                     | <i>Clostridium difficile</i> 027*                      | Adenovirus serotypes 1-5              |
| <i>Salmonella typhimurium</i>                     | <i>Clostridium perfringens</i>                         | Adenovirus serotype 8/15/31           |
| <i>Salmonella bongori</i>                         | <i>Enterotoxigenic Escherichia coli</i>                | Rotavirus A                           |
| <i>Salmonella enteritidis</i>                     | <i>Enterohemorrhagic Escherichia coli</i>              | Norovirus Genotypes I and II          |
| <i>Salmonella enterica</i> subsp. <i>enterica</i> | <i>Enteroinvasive Escherichia coli</i>                 | Astrovirus Genotype I-VIII            |
| <i>Salmonella pullorum</i>                        | <i>Enteropathogenic Escherichia coli</i>               | Sapovirus                             |
| <i>Salmonella gallinarum</i>                      | <i>Klebsiella oxytoca</i>                              | <i>Entamoeba histolytica</i>          |
| <i>Campylobacter lari</i> *                       | <i>Listeria monocytogenes</i>                          | <i>Entamoeba dispar</i>               |
| <i>Campylobacter fetus</i> *                      | <i>Aeromonas caviae</i>                                | <i>Cryptosporidium parvum/hominis</i> |
| <i>Campylobacter coli</i> *                       | <i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i> * | <i>Giardia intestinalis</i> *         |
| <i>Campylobacter jejuni</i> subsp. <i>jejuni</i>  | <i>Blastocystis hominis</i>                            | <i>Dientamoeba fragilis</i>           |

\*These microorganisms showed amplification signal in Clarithromycin wild-type sequences in the 23S rRNA channel.

The specificity of the clarithromycin resistance assay was confirmed by testing a panel consisting of different antimicrobial resistant organism. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

| Cross-reactivity assay   |   |
|--|---|
| Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) N315   | <i>Citrobacter braakii</i> with VIM-1 gene  |
| Methicillin-resistant <i>Staphylococcus aureus</i> ST398   | <i>Citrobacter freundii</i> -complex isolate with KPC-3 and VIM-4 genes   |
| mecC Methicillin-resistant <i>Staphylococcus aureus</i>  | <i>H. pylori</i> Clarithromycin resistant (23S rRNA A2146G)   |
| cMRSA isolate (oxa <sup>R</sup> , PVL-positive, spa:t 310)   | <i>H. pylori</i> Clarithromycin resistant (23S rRNA A2147G)   |
| <i>Serratia marcescens</i> isolate with OXA-48 gene  | VanA-type <i>Enterococcus avium</i>   |
| <i>Klebsiella pneumonia</i> isolate with SHV-1 (non-ESBL), KPC-3, and OXA-48 genes                           | VanA- type <i>Enterococcus faecium</i> LMG16165 strain  |
| <i>Klebsiella pneumonia</i> isolate with TEM-1 (non-ESBL), SHV-1 (non-ESBL), CTX-M-2 (ESBL), and KPC-2 genes | VanA- type <i>Enterococcus faecium</i> IOWA 1   |
| <i>Escherichia coli</i> with OXA-244 gene  | VanB- type <i>Enterococcus faecium</i> IOWA 2   |
| <i>Escherichia coli</i> with TEM-1 (non-ESBL) and IMP-1 genes  | VanA-type <i>Enterococcus faecalis</i>  |
| <i>Enterobacter cloacae</i> with SHV-12 (ESBL), CTX-M-9 (ESBL) and OXA-48 genes                              | VanB-type <i>E. faecalis</i> (Andrewes and Horder) Schleifer and Kilpper-Balz   |
| <i>Enterobacter cloacae</i> with TEM-1 (non ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL) and NDM-1 genes            | VanC and VanB- types <i>Enterococcus gallinarum</i> ENT20120142   |
| <i>Enterobacter cloacae</i> -complex with a NDM-7 gene   | VanC type- <i>Enterococcus gallinarum</i> (Bridge and Sneath 1982) Collins, Jones, Farrow, Kilpper-Balz and Schleifer 1984 VP |

## Analytical reactivity

Analytical reactivity of Vitassay qPCR *H. pylori* ClariRes for *Helicobacter pylori* was evaluated against *Helicobacter pylori* strain J99 and Sydney, showing positive results.

Analytical reactivity of Vitassay qPCR *H. pylori* ClariRes for Clarithromycin resistance was evaluated against *H. pylori* strains wild-type and harbour point mutations A2142G or A2143G, showing negative and positives results, respectively.

## Compatibles real-time PCR equipment

Vitassay qPCR *H. pylori* ClariRes has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96 <sup>TM</sup> Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of *H. pylori*, Clarithromycin resistance and Clarithromycin wild-type sequences in the 23S rRNA infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with human faecal samples and human gastric tissue biopsies. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *H. pylori*, Clarithromycin resistance in *H. pylori* and Clarithromycin wild-type

sequences in the 23S rDNA, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

| Low profile Block Thermocyclers              | High profile Block Thermocyclers           |
|--|--|
| <b>Agilent Technologies</b>                  | <b>Abbott</b>                              |
| AriaMx Real-Time PCR System                  | Abbott m2000 RealTime System               |
| <b>Applied Biosystems</b>                    | <b>Applied Biosystems</b>                  |
| 7500 Fast Real-Time PCR System               | 7300 Real-Time PCR System                  |
| 7500 Fast Dx Real-Time PCR System            | 7500 Real-Time PCR System                  |
| QuantStudio™ 12K Flex 96-well Fast           | 7900 HT Real-Time PCR System               |
| QuantStudio™ 6 Flex 96-well Fast             | ABI PRISM 7000                             |
| QuantStudio™ 7 Flex 96-well Fast             | ABI PRISM 7700                             |
| QuantStudio™ 5 Real-Time PCR System          | QuantStudio™ 12K Flex 96-well              |
| QuantStudio™ 3 Real-Time PCR System          | QuantStudio™ 6 Flex 96-well                |
| StepOne Plus™ Real-Time PCR System           | QuantStudio™ 7 Flex 96-well                |
| StepOne™ Real-Time PCR System                | QuantStudio™ 5 Real-Time PCR               |
| ViiA™ 7 Fast Real-Time PCR System            | QuantStudio™ 3 Real-Time PCR System        |
| <b>BIONEER</b>                               | ViiA™ 7 Real-Time PCR                      |
| Exicycler™ 96                                | <b>Analytik Jena Biometra</b>              |
| <b>Bio-Rad</b>                               | TOptical                                   |
| CFX96™ Real-Time PCR Detection System        | qTOWER 2.0                                 |
| Mini Opticon™ Real-Time PCR Detection System | <b>BIONEER</b>                             |
| <b>Cepheid</b>                               | Exicycler™ 96                              |
| SmartCycler®                                 | <b>Bio-Rad</b>                             |
| <b>Qiagen</b>                                | CFX96™ Deep Well Real-Time PCR Detection   |
| Rotor-Gene® Q                                | iCycler iQ™ Real-Time PCR                  |
| <b>Roche</b>                                 | iCycler iQ™ 5 Real-Time PCR                |
| LightCycler®480 Real-Time PCR System         | MyiQ™ Real-Time PCR Detection System       |
| LightCycler®96 Real-Time PCR System          | MyiQ™ 2Real-Time PCR Detection System      |
| Cobas z480 Analyzer                          | <b>Cepheid</b>                             |
|  | SmartCycler®                               |
|  | <b>DNA-Technology</b>                      |
|  | DTlite Real-Time PCR System                |
|  | DTprime Real-time Detection Thermal Cycler |
|  | <b>Eppendorf</b>                           |
|  | Mastercycler™ ep <i>realplex</i>           |
|  | <b>Qiagen</b>                              |
|  | Rotor-Gene® Q                              |
|  | <b>Stratagene / Agilent Technologies</b>   |
|  | Mx3000P™ Real Time PCR System              |
|  | Mx3005P™ Real Time PCR System              |

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

| REAL-TIME PCR THERMOCYCLER           | Vitassay CHANNEL | DETECTION CHANNEL | OBSERVATIONS  |
|--------------------------------------|------------------|-------------------|---|
| <b>Bio-Rad CFX96™</b>                | FAM              | FAM               |   |
|                                      | HEX              | HEX               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |
| <b>ABI 7500 Applied Biosystems</b>   | FAM              | FAM               | Passive reference option ROX is none  |
|                                      | HEX              | VIC               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |
| <b>Roche Lightcycler®480II</b>       | FAM              | 465/510           | Colour Compensation required  |
|                                      | HEX              | 533/580           |   |
|                                      | ROX              | 533/610           |   |
|                                      | Cy5              | 618/660           |   |
| <b>Smartcycler® Cepheid</b>          | FAM              | Channel 1         |   |
|                                      | HEX              | Channel 2         |   |
|                                      | ROX              | Channel 3         |   |
|                                      | Cy5              | Channel 4         |   |
| <b>Abbott m2000rt</b>                | FAM              | FAM               |   |
|                                      | HEX              | VIC               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |
| <b>Mx3000P™ Mx 3005P™ Stratagene</b> | FAM              | FAM               | Passive reference option ROX is none  |
|                                      | HEX              | VIC               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |
| <b>AriaMx Agilent</b>                | FAM              | FAM               |   |
|                                      | HEX              | HEX               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |
| <b>Rotor-Gene®Q Qiagen</b>           | FAM              | Green             | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition". |
|                                      | HEX              | Yellow            |   |
|                                      | ROX              | Orange            |   |
|                                      | Cy5              | Red               |   |
| <b>Exicycler™ 96 BIONEER</b>         | FAM              | FAM               |   |
|                                      | HEX              | HEX               |   |
|                                      | ROX              | ROX               |   |
|                                      | Cy5              | Cy5               |   |

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

| Thermocycler   | Vitassay channel | Exposure values |
|--|------------------|-----------------|
| <b>DTlite Real-Time PCR System (DNA-Technology)</b>                | FAM              | 500             |
|  | HEX              | 500             |
|  | ROX              | 500             |
|  | Cy5              | 500             |
| <b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b> | FAM              | 500             |
|  | HEX              | 1000            |
|  | ROX              | 1000            |
|  | Cy5              | 1000            |

## Bibliography/Bibliografía

1. Novel real-time PCR assay for detection of *Helicobacter pylori* infection and simultaneous clarithromycin susceptibility testing of stool and biopsy specimens. Schabereiter-Gurtner C, Hirschl AM, Dragosics B, Hufnagl P, Puz S, Kovách Z, Rotter M, Makrithathis A. J Clin Microbiol, 2004; 42(10): 4512-4518.
2. Urease activity and urea gene sequencing of coccoid forms of *H. pylori* induced by different factors. Can F, Karahan C, Dolapci I, Demirbilek M, Tekeli A, Arslan H. Curr Microbiol 2008; 56(2): 150-155.
3. Diagnosis of *Helicobacter pylori*: what should be the gold standard? Patel SK, Pratap CB, Jain AK, Gulati AK, Nath G. World J Gastroenterol 2014; 20(36): 12847-12859.
4. Diagnosis of *Helicobacter pylori* using the rapid urease test. Uotani T, Graham DY. Ann Transl Med 2015; 3(1): 9.
5. *Helicobacter pylori* clarithromycin resistance detected by Etest and TaqMan real-time polymerase chain reaction: a comparative study. Monno R, Giorgio F, Carmine P, Soleo L, Cinquepalmi V, Ierardi E. APMIS 2012; 120(9):712-7.
6. *Helicobacter pylori* treatment: antibiotics or probiotics. Goderska K, Agudo Pena S, Alarcon T. Appl Microbiol Biotechnol. 2017; Oct 26. doi: 10.1007/s00253-017-8535-7.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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