Vitassay qPCR

Rhinovirus & Enterovirus

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Rhinovirus & Enterovirus allows the detection and differentiation of Rhinovirus and Enterovirus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Rhinovirus and ENterovirus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Rhinovirus & Enterovirus 4x8 -well strip, low profile 7041041

Vitassay qPCR Rhinovirus & Enterovirus 4x8-well strip, high profile 7042041

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S041/ 7042S041	Rhinovirus & Enterovirus strips low/high profile	-	4x8-well strip
7C041	Rhinovirus & Enterovirus Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

Summary

Rhinoviruses, like Enteroviruses, are non-SARS viruses (SARS, acute and severe respiratory syndrome) related as pathogens for the respiratory system that cause upper respiratory tract infections or more commonly known as the common cold. In the case of Rhinoviruses, they are found in more than 35% of human Coronavirus (HCoV) infections, while enteroviruses are less abundant (around 7%). The frequency of suffering from this infection increases in immunosuppressed patients.

The symptoms of the infection, after the incubation period of 2 to 5 days, is a fever of 1 to 5 days. In the most severe cases it can cause shortness of breath, gastrointestinal symptoms, nausea, vomiting, fatigue and myalgia and have a longer duration. Transmission is through droplets in a cough or sneeze or being exposed to contaminated surfaces.

Since there are so many pathogens found in a SARS chart related to HCoV, a systematic characterization technique with excellent sensitivity is necessary. It is there where real-time PCR takes on its importance and vital role in the development of this function.

Principle of the test

Vitassay qPCR Rhinovirus & Enterovirus test is based on the real-time amplification of conserved fragments of the *5'UTR* gene for Rhinovirus and Enterovirus. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increment in the fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Rhinovirus & Enterovirus test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of Rhinovirus RNA target sequence is detected through the FAM channel, internal control in Cy5 channel and Enterovirus is detected through the HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Total Nucleic Acid Isolation (TNAI) kit (Roche).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized Rhinovirus & Enterovirus Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Rhinovirus), Cy5 (Internal Control) and HEX, JOE or VIC channels (Enterovirus). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3029P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for Rhinovirus (FAM), Enterovirus (HEX, JOE o VIC) and Internal Control (Cy5), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, HEX, JOE o VIC and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Rhinovirus (FAM)	Enterovirus (HEX)	Internal Control (Cy5)	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	Rinovirus y Enterovirus Positive
+	-	+/-	-	+	Rinovirus Positive, Enterovirus Negative
-	+	+/-	-	+	Rinovirus Negative, Enterovirus Positive
-	-	+	-	+	Rinovirus y Enterovirus Negative
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 187 throat swabs from symptomatic patients were analyzed by real-time PCR using Vitassay qPCR Rhinovirus & Enterovirus and the results were compared with those obtained with CLART® PneumoVir DNA array (Genomica).

Rhinovirus was detected in 87 positive samples by both kits. Apart from these samples, Vitassay qPCR Rhinovirus & Enterovirus detected 13 more positive samples and CLART® PneumoVir DNA array (Genomica) detected 2 more positive samples. These 15 samples (The low amount of template RNA detected in this sample is below the detection limit of the method used) were evaluated by an additional commercial molecular detection method (FastTrack Respiratory Pathogens 33 or FastTrack Respiratory Pathogens 21 (FastTrack Diagnosis) of which resulted in 10 of 15 positive samples.

Enterovirus was detected in 78 positive samples by both kits. Apart from these samples, Vitassay qPCR Rhinovirus & Enterovirus detected 2 more positive samples and CLART® PneumoVir DNA array (Genomica) detected 2 more positive samples. These 4 samples (The low amount of template RNA detected in this sample is below the detection limit of the method used) were evaluated by an additional commercial molecular detection method (FastTrack Respiratory Pathogens 33 or FastTrack Respiratory Pathogens 21 (FastTrack Diagnosis) of which resulted in 3 of 4 positive samples.

These results indicate the high sensitivity and specificity to detect Rhinovirus and Enterovirus using the molecular diagnostic test Vitassay qPCR Rhinovirus & Enterovirus.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Rhinovirus and Enterovirus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Rhinovirus and Enterovirus was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Cross-reactivity testing		
<i>Bordetella bronchiseptica</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza B/Florida/04/06 Virus
<i>Bordetella holmesii</i>	Methicillin-resistant <i>Staphylococcus aureus</i>	Influenza B/Phuket/3073/2013 Virus
<i>Bordetella parapertussis</i>	<i>Moraxella catarrhalis</i>	Human 1, 2, 3 and 4 parainfluenza Virus
<i>Bordetella pertussis</i>	<i>Neisseria meningitidis</i> serogrupo A	Respiratory Syncytial Virus A y B
<i>Cryptococcus gatti</i> Z156	Influenza A/New Caledonia/20/99(H1N1) Virus	Human A and B metapneumovirus Virus
<i>Chlamydomphila pneumoniae</i>	Influenza A/California/7/2009(H1N1)pdm09 Virus	Human Coronavirus 229E y NL63
<i>E. coli</i> 0. 1285;O18:H7:K1	Influenza A/Michigian/45/2015 (H1N1)pdm09 Virus	MERS Coronavirus
<i>Haemophilus influenzae</i> MinnA	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 Virus	Human Adenovirus 2 cepa Adenoid 6
<i>Legionella bozemanii</i>	Influenza A/Perth/16/2009(H3N2) Virus	Adenovirus humano 5
<i>Legionella micdadei</i>	Influenza A/Thüringen/5/17 (H3N2) Virus	Bocavirus
<i>Legionella dumoffii</i>	Influenza A/Switzerland/9715293/2013 (H3N2) Virus	Citomegalovirus cepa AD-169
<i>Legionella pneumophila</i>	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) Virus	Parechovirus Tipo 3
<i>Legionella longbeachae</i>	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) Virus	HSV-1 MacIntyre strain
<i>Listeria monocytogenes</i> Serotipo 1/2b	Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) Virus	HSV-2 MS
<i>Mycoplasma pneumoniae</i>	Influenza A/Anhui/1/2013 (H7N9) Virus	HHV6 Z29 strain
<i>Streptococcus pneumoniae</i> Z022	Influenza B/Brisbane/60/2008-like Virus	Varicella-Zoster Virus Ellen
<i>Streptococcus agalactiae</i> Z019		

Analytical reactivity

Analytical reactivity of Vitassay qPCR Rhinovirus & Enterovirus for Rhinovirus was evaluated against strains: HRV type A 8, HRV type A 16, HRV type A 30, HRV type A 49, HRV type A 90, type B 5, HRV type B 72 and HRV type C, showing positive result.

Analytical reactivity of Vitassay qPCR Rhinovirus & Enterovirus for Enterovirus was evaluated against strains: Enterovirus D68 and A71, Echovirus Type 7, 11, 30 and 68, Coxsackievirus types A9, A21, A24, B3, and B4, showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Rhinovirus & Enterovirus has been validated on the following equipments:

- 7500 Fast Real-Time PCR System (Applied Biosystems)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)
- Cobas z480 (Roche)

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis Rhinovirus and Enterovirus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with samples in transport medium, plasma, throat swabs, CSF and stool samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.

- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Duc H. Do, Stella Laus, Amy Leber, Mario J. Marcon, Jeanne A. Jordan, Judith M. Martinand Robert M. Wadowsky. A One-Step, Real-Time PCR Assay for Rapid Detection of Rhinovirus. *Journal of Molecular Diagnostics* 2010;102-108.
2. Riikka Osterback, Tuire Tevaluoto, Tiina Ylinen, Ville Peltola, Petri Susi, Timo Hyypiä Matti Warisa. Simultaneous Detection and Differentiation of Human Rhino- and Enteroviruses in Clinical Specimens by Real-Time PCR with Locked Nucleic Acid Probes. *Journal of Clinical Microbiology* 2013; 3960-3967.
3. Julien Dupouey, Laetitia Ninove, Vanessa Ferrier, Odile Py, Céline Gazin, Laurence Thirion-Perrier and Xavier de Lamballerie. Molecular detection of human rhinoviruses in respiratory samples:a comparison of Taqman probe-, SYBR green I- and BOXTO-based real-time PCR assays. *Virology Journal* 2014, 11:31.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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