Vitassay qPCR

SARS-CoV-2

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Intended use

Vitassay qPCR SARS-CoV-2 allows the qualitative detection of 2019 Novel Coronavirus (SARS-CoV-2) by real-time RT-PCR in respiratory samples. The product is intended for use in the diagnosis of SARS-CoV-2 infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR SARS-CoV-2 4 x 8-well strip, low profile	7041046
Vitassay qPCR SARS-CoV-2 4 x 8-well strip, high profile	7042046
Vitassay qPCR SARS-CoV-2 8 x 8-well strip, low profile	7081046
Vitassay qPCR SARS-CoV-2 8 x 8-well strip, high profile	7082046
Vitassay qPCR SARS-CoV-2 96-well plate, low profile	7091046
Vitassay qPCR SARS-CoV-2 96-well plate, high profile	7092046

Materials/reagents provided

Reagents provided in references 7041046, 7042046, 7081046 and 7082046:

Reference	Reagent/Material	Colour	Amount
7041S046/ 7042S046	SARS-CoV-2 strips	-	4/8 x8-well strip
7C046	SARS-CoV-2 Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4/8 x 8 cap strip

Reagents provided in references 7091046 and 7092046:

Reference	Reagent/Material	Colour	Amount
7091P046/ 7092P046	SARS-CoV-2 Plate	-	1 plate
7C046	SARS-CoV-2 Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

In December 2019 in Wuhan (Hubei, China) certain people who lived or worked near a seafood market presented pneumonia of unknown cause. The sequencing of the respiratory samples determined the presence of a new type of coronavirus, 2019 new coronavirus (2019-nCoV or SARS-CoV-2).

Coronaviruses are non-segmented positive stranded RNA viruses that belong to the *Coronaviridae* family. Six coronavirus species that cause human diseases are known: 229E, OC43, NL63, HKU1, SARS-CoV, and, MERS-CoV. The first four cause symptoms of a common cold, while the last two are zoonotic and cause more severe complications.

Transmission of this virus is occurring from person to person, even during the asymptomatic incubation period. 2019-nCoV causes severe respiratory disease such as that produced by SARS-CoV. Pneumonia is the main associated disease, although some patients have developed pulmonary edema, severe pneumonia, acute respiratory distress syndrome or multi-organ failure and death. It also causes less common symptoms such as headache, sore throat, diarrhea and vomiting. The symptoms of 2019-nCoV can appear in just two days or up to 14 days after exposure, the most common cough, fever, dyspnea and myalgia.

The samples recommended by the WHO for the identification of 2019-nCoV are those from the lower respiratory tract (sputum, bronchoalveolar lavage or endotracheal aspirates) and/or upper respiratory tract specimens (nasopharyngeal and oropharyngeal

swabs). If this is not possible, samples of the upper respiratory tract, blood, urine and feces will be taken.

The 2019-nCoV diagnosis was made by detecting the conventional causes of early pneumonia and by massive sequencing or real-time RT-PCR methods.

Principle of the test

Vitassay qPCR SARS-CoV-2 is based on the real-time amplification of specific conserved fragments of the *ORF1ab* and *N* genes encoded by the SARS-CoV-2. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of SARS-CoV-2 is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR SARS-CoV-2 is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence *ORF1ab* gene is detected through the FAM channel, the amplification of the target sequence *N* gene is detected through the ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Close protective pouches of reagents promptly with the zip seal after each use
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.

- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (ROCHE).

Invisorb® Spin Universal Kit (Stratec)

MagDEA Dx SV kit, using magLEAD® 12gC (Precision System Science Co.).

MagCore® Viral Nucleic Acid Extraction kit, using MagCore® HF16 automated Nucleic Acid Extractor System.

EZ1 Virus Mini Kit, using EZ1 instrument (Qiagen).

mSample Preparation Systems RNA, using the Abbott m2000 RealTime System (Abbott Molecular).

Positive control preparation

Reconstitute the lyophilized SARS-CoV-2 Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*ORF1ab* gene), ROX (*N* gene) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for SARS-CoV-2 in FAM and ROX channel, which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for SARS-CoV-2 in FAM and ROX channel, which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly and to review all the parameters and the sigmoid shape of the obtained curve. It is recommended to repeat the test.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and / or the results obtained in other diagnostic tests.

The result interpretation is summarized in the following table:

<i>ORF1ab</i> gene (FAM)	<i>N</i> gene (ROX)	Internal Control HEX	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	SARS-CoV-2 Positive
-	-	+	-	+	SARS-CoV-2 Negative
+	-	+/-	-	+	SARS-CoV-2 Positive*
-	+	+/-	-	+	SARS-CoV-2 Presumptive Positive* ^φ
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) Positive: Amplification signal

(-) Negative: No amplification signal

* Amplification of a single gene may occur due to some samples have a RNA copy number close to or below the detection limit, a difference in the amplification performance of the target regions, presence of mutations in the sequences of the target genes, or other factors.

^φ In the case of obtaining a positive result only for the target gene *N*, the result will be interpreted as SARS-CoV-2 Presumptive Positive. In the case that is clinically convenient, it is recommended to confirm the result by a reference laboratory to discard other unknown zoonotic coronaviruses.

To consider a result as positive, the Ct value obtained after the analysis of the sample must be less than 38.

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 100 respiratory samples (nasopharyngeal swab) from symptomatic patients with suspicion of COVID-19, were tested by Real Time PCR using Vitassay qPCR SARS-CoV-2 and the results were compared with those obtained using a molecular detection method by the Spanish National Reference Center (ISCIII). SARS-CoV-2 was detected in the 2 positive samples by Vitassay qPCR SARS-CoV-2.

Another clinical evaluation of Vitassay qPCR SARS-CoV-2 was conducted with 80 anonymized respiratory samples from the biobank of the National Center for Microbiology (NCM)(Spain), were tested by Real Time PCR using Vitassay qPCR SARS-CoV-2 and the results were compared with those obtained using a molecular detection method by the Spanish National Reference Center (ISCIII). SARS-CoV-2 was detected in the 41 positive samples by Vitassay qPCR SARS-CoV-2. A discordant sample was obtained since it was positive only for the *N* gene and it was considered presumptive positive.

These results indicate the high sensitivity and specificity to detect SARS-CoV-2 using the molecular diagnostic test Vitassay qPCR SARS-CoV-2.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of SARS-CoV-2 template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction for *ORF1ab* and *N* genes.

Analytical specificity

The analytical specificity for SARS-CoV-2 was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

SARS-CoV-2		
<i>Bordetella pertussis</i>	<i>Mycoplasma pneumoniae</i>	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus
<i>Bordetella parapertussis</i>	<i>Mycobacterium tuberculosis</i> not rifampin resistant	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella holmesii</i>	<i>Pneumocytis jirovecii</i>	Influenza B/Brisbane/60/2008 virus
<i>Bordetella bronchiseptica</i>	<i>Streptococcus pneumoniae</i> Z022	Influenza B/Florida/04/06 virus
<i>Chlamydia caviae</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza B/Phuket/3073/2013 virus
<i>Chlamydia psittaci</i> genotypes A and C	Influenza A/New Caledonia/20/99(H1N1) virus	Human Adenovirus Types 1-5, 8, 15, 31, 40 and 41
<i>Chlamydophila pneumoniae</i> CM-1	Influenza A/California/7/2009(H1N1) virus	Human Bocavirus
<i>Haemophilus influenzae</i> MinnA	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Human coronavirus 229E, OC43, NL63 and HKU1
<i>Legionella bozemanii</i>	Influenza A/Singapore/GP1908/2015 IVR-180 (H1N1)pdm09 virus	MERS Coronavirus
<i>Legionella dumoffii</i>	Influenza A/Victoria/210/2009 (H3N2) virus	SARS Coronavirus Strain Frankfurt 1
<i>Legionella longbeache</i>	Influenza A/Thüringen/5/2017 (H3N2) virus	Human metapneumovirus A and B
<i>Legionella micdadei</i>	Influenza A/Switzerland/9715293/2013 (H3N2) virus	Human parainfluenza 1, 2, 3 and 4 viruses
<i>Legionella pneumophila</i>	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	Respiratory Syncytial virus (RSV) A and B
<i>Moraxella catarrhalis</i>	Influenza A/ South Australia/55/2014, IVR-175 (H3N2) virus	Human rhinovirus type C

Analytical reactivity

The reactivity of Vitassay qPCR SARS-CoV-2 was confirmed by the real-time amplification using RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1 as template.

Compatibles real-time PCR equipment

Vitassay qPCR SARS-CoV-2 has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}

- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)¹
- SmartCycler® (Cepheid)¹

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of SARS-CoV-2 infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with nasopharyngeal swab and oropharyngeal swab. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- Vitassay qPCR SARS-CoV-2, designed for the specific detection of SARS-CoV-2, mostly showed no significant combined homologies with the human genome, other coronaviruses, or human microflora that would predict potential false positive RT-qPCR results. Even though, an alignment performed with the oligonucleotide primer and probe sequences of *N* gene showed sequence homology with a few bat and pangolin coronavirus genomes.
- Detection of SARS-CoV-2 may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including RNA isolation); c) RNA degradation during specimen shipping, storage, and/or preparation; d) potential mutations of the target regions of the SARS-CoV-2 genome covered by this test which may result in RNA being undetectable when pathogen load below the limit of detection for the assay; e) the presence of retrotranscription and/or Real Time amplification inhibitors or other types of

interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs related to COVID-19 was not performed); f) failure to follow the manufacturer's instructions and procedures.

- A negative result for *ORF1ab* gene and the only amplification of the *N* gene may occur in specimens with low pathogen load, due to subtle differences in clinical sensitivity. In case of an ambiguous result, it is recommended to seek confirmatory testing by a reference laboratory.
- Detection of viral RNA (*ORF1ab* and/or *N* gene sequences) may not indicate the presence of viable and/or infectious virus or that SARS-CoV-2 is the causative agent for clinical symptoms.
- Negative results do not preclude infection with the SARS-CoV-2 virus and should not be the sole basis of a patient treatment/management decision. Specimens types for the identification of SARS-CoV-2 and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the virus.
- If the clinical data of the patient, laboratory tests and epidemiological studies suggest possible infection with SARS-CoV-2, and other respiratory illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Zhu, N., Zhang, D., Wang, W., Li, X., Yang, B., Song, J., ... & Niu, P. (2020). A novel coronavirus from patients with pneumonia in China, 2019. *New England Journal of Medicine*.
2. Chen, N., Zhou, M., Dong, X., Qu, J., Gong, F., Han, Y., ... & Xia, J. A. (2020). Epidemiological and Clinical Characteristics of 99 Cases of 2019-Novel Coronavirus (2019-nCoV) Pneumonia in Wuhan, China. *The Lancet*.
3. Huang, C., Wang, Y., Li, X., Ren, L., Zhao, J., Hu, Y., ... & Cheng, Z. (2020). Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *The Lancet*.
4. Lu, R., Zhao, X., Li, J., Niu, P., Yang, B., Wu, H., ... & Tan, W. (2020). Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *The Lancet*.
5. Rothe, C., Schunk, M., Sothmann, P., Bretzel, G., Froeschl, G., Wallrauch, C., ... & Hoelscher, M. (2020). Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*.
6. Na Zhu, Ph.D., Dingyu Zhang, M.D., Wenling Wang *et al.* (2019) A novel coronavirus from patients with pneumonia in China. *New England Journal of Medicine*.
7. Corman *et al* (2020). Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *European communicable disease bulletin* 2020;25(3).
8. Daniel K W Chu, Yang Pan, Samuel M S Cheng, Kenrie P Y Hui, Pavithra Krishnan, Yingzhi Liu, Daisy Y M Ng, Carrie K C Wan, Peng Yang, Quanyi Wang, Malik Peiris, Leo L M Poon, Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clinical Chemistry*.
9. Ding-feng Lv, Qi-ming Ying, Yue-song Weng, Chi-bin Shen, Jin-guo Chu, Jing-ping Kong, Ding-he Sun, Xiang Gao, Xing-bei Weng and Xue-qin Chen. Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clinica Chimica Acta*.
10. Centers of Disease Control and Prevention (CDC). 2019 Novel Coronavirus, Symptoms. Available from <https://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html> Accessed February 2020.
11. World Health Organization. Global Surveillance for human infection with novel coronavirus (2019-nCoV). Interim guidance v3. 31 January 2020. Available from [https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-\(2019-ncov\)](https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-(2019-ncov)) Accessed February 2020.
12. World Health Organization. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases. Interim guidance. 14 January 2020. Available from <https://www.who.int/publications-detail/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed February 2020.

13. World Health Organization. MERS situation update. December 2019. Available from <http://applications.emro.who.int/docs/EMCSR246E.pdf?ua=1&ua=1> Accessed February 2020.
14. World Health Organization. Novel Coronavirus(2019-nCoV). Situation Report – 14. Data as reported by 3 February 2020. Available from https://www.who.int/docs/default-source/coronaviruse/situation-reports/20200203-sitrep-14-ncov.pdf?sfvrsn=f7347413_2 Accessed February 2020.
15. World Health Organization. Novel Coronavirus (2019-nCoV) technical guidance: Laboratory testing for 2019-nCoV in humans. Available from <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/laboratory-guidance> Accessed February 2020.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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