# Vitassay qPCR

**BK & JC**

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## Intended use

Vitassay qPCR BK & JC allows the detection and differentiation of BK and JC Virus by real-time PCR in clinical samples. The product is intended for use in the diagnosis of BK and/or JC Virus infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR BK & JC 4x8 -well strip, low profile 7041047

Vitassay qPCR BK & JC 4x8-well strip, high profile 7042047

## Materials/reagents provided

Code	Reagent/Material	Colour	Amount
7041S047/ 7042S047	BK & JC strips low/high profile	-	4x8-well strip
7C047	BK & JC Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yelow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Summary

Polyomaviruses are ubiquitous, small, non-enveloped DNA, species-specific viruses which belong to the family *Papovaviridae*. The two most commonly known human polyomaviruses BK virus (BKV) and JC virus (JCV), were firstly described in the 1970s.

Human polyomavirus viruses, infect immunocompetent individuals without specific signs or symptoms. They are able to establish subclinical and persistent infections and may reactivate from latency in case of immunosuppression.

Primary infection with both BK and JC virus typically occurs in childhood, with 50% of children having seroconversion for BK virus by 3-4 years of age, and for JC virus by 10-15 years of age. Primary infection is usually asymptomatic or associated with mild upper respiratory symptoms. By adulthood, 80% of people have been infected by either BK virus, JC virus, or both. [1]

JCV reactivation in immunocompromised patients may cause a demyelinating disease of the central nervous system (CNS), named progressive multifocal leukoencephalopathy (PML), whereas BKV reactivation may cause hemorrhagic cystitis and polyomavirus associated nephropathy (PVAN) in renal transplant patients.

BK virus can be isolated in the urine of both symptomatic and asymptomatic patients and plasma. The detection of JC virus is carried out using cerebrospinal fluid (CSF) samples and plasma.

PCR Detection method has demonstrated to be a sensitive and specific diagnosis tool for the detection of BK Virus and JC Virus.

## Principle of the test

Vitassay qPCR BK & JC test is based on the real-time amplification of a conserved region of the *VP1* gene for BK Virus and *T* gene for JC Virus. After DNA isolation, the presence of BK Virus and JC Virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR BK & JC test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the BK Virus DNA target sequence is detected through the Cy5 channel and JC Virus DNA in FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex II).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## Procedures

### Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit, utilizando Maxwell® 16 instrument (Promega).

Invisorb® Spin Universal Kit (Stratec).

### Positive control preparation

Reconstitute the lyophilized BK & JC Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

### Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the Cy5 (BK Virus), FAM (JC Virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

## Positive control

The positive controls used in each run, must show an amplification curve for Cy5 (BK Virus) and FAM (JC Virus), which validates the reaction.

## Negative control

The negative controls included in each run, must show the absence of signal in Cy5 and FAM which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

JC Virus (FAM)	BK Virus (Cy5)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	JC Virus and BK Virus Positives
-	-	+	-	+	JC Virus and BK Virus Negative
+	-	+/-	-	+	JC Virus Positive and BK Virus Negative
-	+	+/-	-	+	BK Virus Positive and JC Virus Negative
-	-	-	-	+	Invalid
+	+	+	+	-	Invalid

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

A total of 112 clinical samples (urine and plasma) from different panels of QCMD programs were analyzed using Vitassay qPCR BK & JC.

Vitassay qPCR BK & JC detected BK virus in 41 positive samples and JC virus in 38 positive samples. The results indicated were compared with the final reports of these programs.

The results show a high sensitivity and specificity to detect BK and JC Viruses using Vitassay qPCR BK & JC.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of BK and JC Viruses templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  DNA copies per reaction.

### Analytical specificity

The analytical specificity for BK and JC Viruses was tested within the panel of different microorganisms.

No cross-reactivity of BK and JC Viruses was seen between any of the species:

Cross-reactivity assay		
<i>E. coli</i> 0.1285;O18:H7:K1	<i>Acinetobacter baumannii</i>	HSV-2 MS
<i>Listeria innocua</i> Serotype 6a/strain CCUG 15531	<i>Candida albicans</i>	HHV6 strain Z29
<i>Listeria ivanovii</i> Serovar 5/strain CCUG 15528	<i>Aspergillus fumigatus</i> pte	HHV6 Type A
<i>Listeria monocytogenes</i> Serotype 1/2b	<i>Pneumocystis jirovecii</i> Type A1	HHV6 Type B
<i>Listeria monocytogenes</i> Serovar 4b/Strain CIP 59.53	<i>Pneumocystis jirovecii</i> g885652	Varicella-Zoster Virus Ellen
<i>Bacteroides fragilis</i>	<i>Staphylococcus aureus</i> subsp <i>aureus</i>	<i>Toxoplasma gondii</i> Type II
<i>Pseudomonas aeruginosa</i>	<i>Plasmodium falciparum</i> 3D7	Hepatitis A
<i>Yersinia enterocolitica</i> O:3	<i>Treponema pallidum</i>	Parvovirus
<i>Klebsiella oxytoca</i>	Epstein-Barr virus	Adenovirus 40
<i>Enterobacter cloacae</i> Serotype Cloaca A	Citomegalovirus strain AD-169	Adenovirus 41
<i>Enterobacter aerogenes</i> Serotype Cloaca B	HSV-1 strain MacIntyre	

### Analytical reactivity

Analytical reactivity of Vitassay qPCR BK & JC for BK Virus was evaluated against BK Virus Type Ib-2 and BK Virus Type IV, showing positive results.

Analytical reactivity of Vitassay qPCR BK & JC for JC Virus was evaluated against JC Virus Type 1A and JC Virus Type 2B, showing positive results.

### **Compatibles real-time PCR equipment**

Vitassay qPCR BK & JC has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

### **Limitations**

- This test provides a presumptive diagnosis of BK and JC Viruses infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with human urine, CSF and plasma samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different parasites, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers
<b>Agilent Technologies</b>
AriaMx Real-Time PCR System
<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
<b>BIONEER</b>
Exicycler™ 96
<b>Bio-Rad</b>
CFX96™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
<b>Cepheid</b>
SmartCycler®
<b>Qiagen</b>
Rotor-Gene® Q
<b>Roche</b>
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
<b>Abbott</b>
Abbott m2000 RealTime System
<b>Applied Biosystems</b>
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR
<b>Analytik Jena Biometra</b>
TOptical
qTOWER 2.0
<b>BIONEER</b>
Exicycler™ 96
<b>Bio-Rad</b>
CFX96™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™ 5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
<b>Cepheid</b>
SmartCycler®
<b>DNA-Technology</b>
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
<b>Eppendorf</b>
Mastercycler™ ep <i>realplex</i>
<b>Qiagen</b>
Rotor-Gene® Q
<b>Stratagene / Agilent Technologies</b>
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene®Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## **Bibliography/Bibliografía**

1. Michelle Pinto, Simon Dobson. BK and JC virus: A review. *Journal of Infection* (2014) 68, 52-58.
2. Malihe Fathi Haghighi, Noorossadat Seyyedi, Ali Farhadi, Farahnaz Zare, Leila Kasraian, Gholam Reza Refiei Dehbidi, Reza Ranjbaran, Abbas Behzad-Behbahani. Polyomaviruses BK and JC DNA infection in peripheral blood cells from blood donors. *BRAZ J INFECT DIS* 2019;23(1):22-26.
3. Deirdre Sawinski and Simin Goral. BK virus infection: an update on diagnosis and treatment. *Nephrol Dial Transplant* (2015) 30: 209–217.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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