Vitassay qPCR

West Nile Virus

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- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cyclers (DNA-Technology)
- DTLite Real Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR West Nile Virus allows the qualitative detection of West Nile virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of West Nile virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR West Nile Virus 4 x 8-well strip, low profile 7041004

Vitassay qPCR West Nile Virus 4 x 8-well strip, high profile 7042004

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S004/ 7042S004	West Nile Virus strips high/low profile	-	4x8-well strip
7C004	West Nile Virus Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

West Nile virus (WNV), has become one of the main cause of viral encephalitis. The virus was originally isolated in 1937 in Uganda, and later produced epidemic outbreaks in Asia, Europe, Australia and the USA. It is a member of the *flavivirus* genus and belongs to the Japanese encephalitis antigenic complex of the family *Flaviviridae*. Flaviviruses are single-strand positive-sense RNA viruses. Though there are seven known genetic lineages of WNV, lineages 1 and 2 are responsible for the major epidemics in humans and animals.

WNV is maintained in nature in a mosquito-bird-mosquito transmission cycle. Mosquitoes of the genus *Culex* are generally considered the principal vectors of WNV, in particular *Cx. Pipiens*, while birds are the reservoir hosts of the virus. Human transmission results mostly from mosquito bites. Occasional cases of WNV transmission have been reported through transplantation of infected organs, breast-feeding and blood transfusions from asymptomatic infected individuals. WNV infection is asymptomatic in around 80% of cases, while about 20% of people develop West Nile fever, showing a flu-like illness. 1% of infected individuals may develop severe West Nile disease, characterized by neuroinvasive conditions including meningitis, encephalitis, acute flaccid paralysis, and even death.

WNV neuroinvasive disease is usually diagnosed by identification of WNV-specific IgM antibodies in cerebrospinal fluid (CSF) or blood serum by enzyme-linked immunosorbent assay (ELISA). Heterologous cross-reactions induced by infection with or vaccination against other flaviviruses may occur, making it necessary to confirm ELISA results by plaque-reduction neutralization assays. As an alternative, PCR-based diagnostic tools offer the possibility of specifically identifying WNV-RNA in different tissue specimens.

Principle of the test

Vitassay qPCR West Nile Virus is based on the real-time amplification of a conserved sequence of genomic region 5'UTR encoded by the West Nile virus genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of West Nile virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR West Nile Virus is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (Qiagen).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized West Nile Virus Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (West Nile virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for West Nile virus, which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for West Nile virus, which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

West Nile Virus FAM	Internal Control HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	Virus West Nile Positive
-	+	-	+	Virus West Nile Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 20 samples from QCMD panels were evaluated by real-time PCR using: Vitassay qPCR West Nile Virus. 14 positive samples for West Nile virus were detected

using the Vitassay qPCR West Nile Virus. The results were compared with the final reports. Vitassay West Nile qPCR shows high sensitivity and specificity for detecting West Nile virus.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of West Nile virus template ranging from 10⁷ to 10¹ copies/rxn. This assay has a detection limit of ≥10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for West Nile virus was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

West Nile Virus		
Zika virus strain MR 766	Dengue 2 virus strain New Guinea C	St Louis Encephalitis virus strain 17D
Chikungunya virus strain S27 Petersfield	Dengue 3 virus strain H87	Yellow Fever virus strain 17D
Dengue 1 virus strain Hawaii	Dengue 4 virus strain H241	

Analytical reactivity

The reactivity of Vitassay qPCR West Nile was confirmed by the real-time amplification using West Nile virus strain H160/99, Heja, Ug37 and NY99 as template.

Compatibles real-time PCR equipment

Vitassay qPCR West Nile Virus has been validated on the following equipments:

- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTLite Real Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of West Nile virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with serum and urine samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by West Nile virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers
Agilent Technologies
AriaMx Real-Time PCR System
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Cepheid
SmartCycler®
Qiagen
Rotor-Gene® Q
Roche
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
Abbott
Abbott m2000 RealTime System
Applied Biosystems
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR
Analytik Jena Biometra
TOptical
qTOWER 2.0
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™ 5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
Cepheid
SmartCycler®
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Eppendorf
Mastercycler™ep <i>realplex</i>
Qiagen
Rotor-Gene® Q
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Genetic and serologic properties of West Nile virus associated with an epidemic, Yap State, Micronesia, 2007. Lanciotti RS, Kosoy OL, Laven JJ, Velez JO, Lambert AJ, Johnson AJ, Stanfield SM, Duffy MR. Emerg Infect Dis. 2008 Aug;14(8):1232-9.
2. Quantitative real-time PCR detection of West Nile virus and evaluation with field-caught mosquitoes. Faye O, Faye O, Diallo D, Diallo M, Weidmann M, Sall AA. Virol J. 2013 Oct 22;10:311.
3. Detection of West Nile virus in urine. Gourinat AC, O'Connor O, Calvez E, Goarant C, Dupont-Rouzeyrol M. Emerg Infect Dis. 2015 Jan;21(1):84-6.
4. Centers for Disease Control and Prevention. West Nile Virus (<http://www.cdc.gov/WestNile/index.html>).
5. World Health Organization. West Nile virus and potential complications (<http://www.who.int/emergencies/WestNile-virus/en/>).

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.

Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastecycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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