Vitassay qPCR

Bordetella pertussis + B. parapertussis+ B. holmesii

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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| Termocicladores con bloque de bajo perfil | Termocicladores con bloque de alto perfil |
|--|---|
| Agilent Technologies | Abbott |
| AriaMx Real-Time PCR System | Abbott m2000 RealTime System |
| Applied Biosystems | Applied Biosystems |
| 7500 Fast Real-Time PCR System | 7300 Real-Time PCR System |
| 7500 Fast Dx Real-Time PCR System | 7500 Real-Time PCR System |
| QuantStudio™ 12K Flex 96-well Fast | 7900 HT Real-Time PCR System |
| QuantStudio™ 6 Flex 96-well Fast | ABI PRISM 7000 |
| QuantStudio™ 7 Flex 96-well Fast | ABI PRISM 7700 |
| QuantStudio™ 5 Real-Time PCR System | QuantStudio™ 12K Flex 96-well |
| QuantStudio™ 3 Real-Time PCR System | QuantStudio™ 6 Flex 96-well |
| StepOne Plus™ Real-Time PCR System | QuantStudio™ 7 Flex 96-well |
| StepOne™ Real-Time PCR System | QuantStudio™ 5 Real-Time PCR |
| ViiA™ 7 Fast Real-Time PCR System | QuantStudio™ 3 Real-Time PCR System |
| BIONEER | ViiA™ 7 Real-Time PCR |
| Exicycler™ 96 | Analytik Jena Biometra |
| Bio-Rad | TOptical |
| CFX96™ Real-Time PCR Detection System | qTOWER 2.0 |
| Mini Opticon™ Real-Time PCR Detection System | BIONEER |
| Cepheid | Exicycler™ 96 |
| SmartCycler® | Bio-Rad |
| Qiagen | CFX96™ Deep Well Real-Time PCR Detection |
| Rotor-Gene® Q | iCycler iQ™ Real-Time PCR |
| Roche | iCycler iQ™ 5 Real-Time PCR |
| LightCycler® 480 Real-Time PCR System | MyiQ™ Real-Time PCR Detection System |
| LightCycler® 96 Real-Time PCR System | MyiQ™ 2Real-Time PCR Detection System |
| Cobas z480 Analyzer | Cepheid |
| | SmartCycler® |
| | DNA-Technology |
| | DTlite Real-Time PCR System* |
| | DTprime Real-time Detection Thermal Cycler* |
| | Eppendorf |
| | Mastercycler™ ep <i>realplex</i> |
| | Qiagen |
| | Rotor-Gene® Q |
| | Stratagene / Agilent Technologies |
| | Mx3000P™ Real Time PCR System |
| | Mx3005P™ Real Time PCR System |

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Intended use

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* allows the detection and differentiation of *Bordetella pertussis*, *Bordetella parapertussis* y/o *Bordetella holmesii* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Bordetella pertussis*, *Bordetella parapertussis* y/o *Bordetella holmesii* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* 4x8 -well strip, low profile 7041020

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* 4x8-well strip, high profile 7042020

Materials/reagents provided

| Reference | Reagent/Material | Colour | Amount |
|-----------------------|--|--------|------------------|
| 7041S020/ 7042S020 | <i>Bordetella pertussis</i> + <i>B. parapertussis</i> + <i>B. holmesii</i> strips low/high profile | - | 4 x 8-well strip |
| 7C020 | <i>Bordetella</i> Positive Control | red | 1 vial |
| 7001A | PCR grade water | white | 1 vial x 1 mL |
| 7002B | Resuspension Buffer | green | 1 vial x 1,8 mL |
| 7003N | Negative control | yellow | 1 vial x 1 mL |
| 7004O | Optical caps | - | 4 x 8 cap strip |

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

The genus *Bordetella* bacteria are gram negative bacilli which are isolated from the upper respiratory tract. *Bordetella pertussis* is the main causative agent of whooping cough, although *Bordetella parapertussis*, can also produce the disease, but with a milder symptomatology. Additionally, *Bordetella bronchiseptica* and *Bordetella holmesii* have been identified as causing respiratory tract diseases with clinical manifestations similar to that of the pertussis syndrome but mainly in immunocompromised individuals

The estimated incidence of this disease worldwide ranges from 30 to 50 million cases causing about 195,000 deaths a year. Most cases occur in developing countries where vaccination programs do not exist. In developed countries, despite the dramatic decrease in the incidence of this disease since the introduction of vaccination in the 1950s, a global resurgence of the disease has been observed in the last two decades.

Isolation of *B. pertussis* in culture is definitive for diagnosis, and although the low sensitivity (50%) remains the diagnostic reference method. PCR allows rapid diagnosis and enhances culture sensitivity.

Principle of the test

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* test is based on the real-time amplification of a conserved region of the insertion sequence *IS481* for the identification of *Bordetella pertussis*/*Bordetella holmesii*, *Bordetella holmesii* for the insertion sequence *hIS1001* and *Bordetella parapertussis* for the insertion sequence *pIS1001*. After DNA isolation, the presence of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Bordetella pertussis*/*Bordetella holmesii* DNA target sequence is detected through the FAM channel, *Bordetella holmesii* DNA target

in ROX channel and *Bordetella parapertussis* DNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega)

EZ1 Virus Mini Kit (Qiagen)

RIDA® Xtract (r-Biopharm)

Invisorb® Spin Universal Kit (Stratec)

Positive control preparation

Reconstitute the lyophilized Bordetella Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

| Step | Temperature | Time | Cycles |
|---|-------------|--------|--------|
| Initial denaturation | 95°C | 2 min | 1 |
| Denaturation | 95°C | 10 sec | 45 |
| Annealing/Extension (Data collection*) | 60°C | 50 sec | |

Set the fluorescence data collection during the extension step (*) through the FAM (*Bordetella pertussis/Bordetella holmesii*), ROX (*Bordetella holmesii*), Cy5 (*Bordetella parapertussis*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Bordetella pertussis/Bordetella holmesii*), ROX (*Bordetella holmesii*) and Cy5 (*Bordetella parapertussis*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

| <i>Bordetella pertussis/Bordetella holmesii</i> | <i>Bordetella holmesii</i> | <i>Bordetella parapertussis</i> | Internal control | Negative control | Positive control | Interpretation |
|---|----------------------------|---------------------------------|------------------|------------------|------------------|--|
| + | - | - | +/- | - | + | <i>Bordetella pertussis</i> Positive, <i>B. holmesii</i> and <i>B. parapertussis</i> Negatives |
| + | + | - | +/- | - | + | <i>Bordetella holmesii</i> Positive, <i>B. pertussis</i> and <i>B. parapertussis</i> Negatives |
| - | - | + | +/- | - | + | <i>Bordetella parapertussis</i> Positive, <i>B. holmesii</i> and <i>B. pertussis</i> Negatives |
| - | - | - | + | - | + | <i>B. pertussis</i> , <i>B. holmesii</i> and <i>B. parapertussis</i> Negatives |
| + | + | + | + | + | + | Experiment fail |
| - | - | - | - | - | - | Experiment fail |

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* was evaluated with the different QCMD panels from *Bordetella pertussis*. These panels consist of a total of 30 clinical specimens dissolved in saline or amies transport medium. The results were compared with those of EQA programme final report. Vitassay qPCR correctly detected all the samples (16 *Bordetella pertussis* positives, 3 *B. parapertussis* positives and 3 *B. holmesii* positives).

In addition, Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* was evaluated with 187 respiratory samples in a Multicenter Evaluation through collaboration with different entities. The results were compared with those obtained by three commercial real-time PCR kits: *Bordetella pertussis* & *parapertussis* primer / probes sets (Cepheid), RIDA@GENE *Bordetella* real-time PCR (R-biopharm) and *Bordetella pertussis*, *B. parapertussis* and *B. holmesii* (BioGx, USA).

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* and the other three commercial kits detected 50 positive samples for *Bordetella pertussis*.

For *B. parapertussis*, 3 samples were detected by Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* while the other three kits only detected 1 positive sample. The Vitassay kit detects *B. pertussis* and *B. parapertussis* in 2 clinical samples diagnosed as *B. pertussis* by the Cepheid assay. This co-infection cannot be ruled out.

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* and the other three commercial kits detected 3 positive samples for *Bordetella holmesii*.

The results show high sensitivity and specificity to detect *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* using Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity for *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

| Cross-reactivity assay | | |
|--|---|-------------------------------------|
| <i>Legionella bozemanii</i> | <i>Haemophilus influenzae</i> | Influenza B/Phuket/3073/2013 |
| <i>Legionella micdadei</i> | <i>Moraxella catarrhalis</i> | Respiratory syncytial virus (RSV) |
| <i>Legionella dumoffii</i> | Influenza A/California/7/2009(H1N1) | Human metapneumovirus A and B |
| <i>Legionella longbeachae</i> | Influenza A/Perth/16/2009(H3N2) | Human parainfluenza 1, 2, 3, 4 |
| <i>Legionella pneumophila</i> | Influenza A/New Caledonia/20/99(H1N1) virus | Human coronavirus 229E |
| <i>Mycoplasma pneumoniae</i> | Influenza A/Switzerland/9715293/2013 | Human rhinovirus |
| <i>Streptococcus pneumoniae</i> | Influenza A/Turkey/Germany R2485+86/2014 | Human Adenovirus 2 strain Adenoid 6 |
| <i>Staphylococcus aureus</i> subsp. <i>aureus</i> | Influenza B/Brisbane/60/2008 | Human Adenovirus 5 |
| Methicillin-resistant <i>Staphylococcus aureus</i> | Influenza B/Florida/04/06 | <i>Bordetella bronchiseptica</i> |

Analytical reactivity

Analytical reactivity of Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* was evaluated against specific clinical strains of each species, *Bordetella pertussis* (*Bordetella pertussis* (Bergey et al. 1923) Moreno-Lopez 1952 CECT 7974 strain), *Bordetella parapertussis* (*B. parapertussis* clinical strain isolated which was confirmed by molecular diagnosis assays (SmartCycler *Bordetella pertussis/parapertussis* assay (Cepheid) and RIDA@GENE *Bordetella* (r-Biopharm)) and *Bordetella holmesii* (*Bordetella holmesii* Weyant et al. 1995 DSM 13416 strain), showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with nasopharyngeal aspirates and pernasal swabs samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different bacteria, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

| Low profile Block Thermocyclers | High profile Block Thermocyclers |
|--|--|
| Agilent Technologies | Abbott |
| AriaMx Real-Time PCR System | Abbott m2000 RealTime System |
| Applied Biosystems | Applied Biosystems |
| 7500 Fast Real-Time PCR System | 7300 Real-Time PCR System |
| 7500 Fast Dx Real-Time PCR System | 7500 Real-Time PCR System |
| QuantStudio™ 12K Flex 96-well Fast | 7900 HT Real-Time PCR System |
| QuantStudio™ 6 Flex 96-well Fast | ABI PRISM 7000 |
| QuantStudio™ 7 Flex 96-well Fast | ABI PRISM 7700 |
| QuantStudio™ 5 Real-Time PCR System | QuantStudio™ 12K Flex 96-well |
| QuantStudio™ 3 Real-Time PCR System | QuantStudio™ 6 Flex 96-well |
| StepOne Plus™ Real-Time PCR System | QuantStudio™ 7 Flex 96-well |
| StepOne™ Real-Time PCR System | QuantStudio™ 5 Real-Time PCR |
| ViiA™ 7 Fast Real-Time PCR System | QuantStudio™ 3 Real-Time PCR System |
| BIONEER | ViiA™ 7 Real-Time PCR |
| Exicycler™ 96 | Analytik Jena Biometra |
| Bio-Rad | TOptical |
| CFX96™ Real-Time PCR Detection System | qTOWER 2.0 |
| Mini Opticon™ Real-Time PCR Detection System | BIONEER |
| Cepheid | Exicycler™ 96 |
| SmartCycler® | Bio-Rad |
| Qiagen | CFX96™ Deep Well Real-Time PCR Detection |
| Rotor-Gene® Q | iCycler iQ™ Real-Time PCR |
| Roche | iCycler iQ™ 5 Real-Time PCR |
| LightCycler®480 Real-Time PCR System | MyiQ™ Real-Time PCR Detection System |
| LightCycler®96 Real-Time PCR System | MyiQ™ 2Real-Time PCR Detection System |
| Cobas z480 Analyzer | Cepheid |
| | SmartCycler® |
| | DNA-Technology |
| | DTlite Real-Time PCR System |
| | DTprime Real-time Detection Thermal Cycler |
| | Eppendorf |
| | Mastercycler™ ep <i>realplex</i> |
| | Qiagen |
| | Rotor-Gene® Q |
| | Stratagene / Agilent Technologies |
| | Mx3000P™ Real Time PCR System |
| | Mx3005P™ Real Time PCR System |

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

| REAL-TIME PCR THERMOCYCLER | Vitassay CHANNEL | DETECTION CHANNEL | OBSERVATIONS |
|--------------------------------------|------------------|-------------------|---|
| Bio-Rad CFX96™ | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| ABI 7500 Applied Biosystems | FAM | FAM | Passive reference option ROX is none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Roche Lightcycler®480II | FAM | 465/510 | Colour Compensation required |
| | HEX | 533/580 | |
| | ROX | 533/610 | |
| | Cy5 | 618/660 | |
| Smartcycler® Cepheid | FAM | Channel 1 | |
| | HEX | Channel 2 | |
| | ROX | Channel 3 | |
| | Cy5 | Channel 4 | |
| Abbott m2000rt | FAM | FAM | |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Mx3000P™ Mx 3005P™ Stratagene | FAM | FAM | Passive reference option ROX is none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| AriaMx Agilent | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Rotor-Gene® Q Qiagen | FAM | Green | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition". |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Exicycler™ 96 BIONEER | FAM | FAM | |
| | HEX | JOE | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

| Thermocycler | Vitassay channel | Exposure values |
|--|-------------------------|------------------------|
| DTlite Real-Time PCR System (DNA-Technology) | FAM | 500 |
| | HEX | 500 |
| | ROX | 500 |
| | Cy5 | 500 |
| DTprime Real-time Detection Thermal Cycler (DNA-Technology) | FAM | 500 |
| | HEX | 1000 |
| | ROX | 1000 |
| | Cy5 | 1000 |

Bibliography/Bibliografía

1. Novel Multitarget Real-Time PCR Assay for Rapid Detection of *Bordetella* Species in Clinical Specimens Kathleen M. Tatti, Kansas N. Sparks, Kathryn O. Boney, Maria Lucia Tondella *J Clin Microbiol.* 2011 Dec; 49(12): 4059–4066.
2. Detection and Differentiation of *Bordetella* spp. by Real-Time PCR. Christoph Koidl, Michael Bozic, Anja Burmeister, Markus Hess, Egon Marth, Harald H. Kessler. *J Clin Microbiol.* 2007 Feb; 45(2): 347–350.
3. Real-Time PCR Assay Targeting IS481 of *Bordetella pertussis* and Molecular Basis for Detecting *Bordetella holmesii*. Udo Reischl, Norbert Lehn, Gary N. Sanden, Mike J. Loeffelholz *J Clin Microbiol.* 2001 May; 39(5): 1963–1966).

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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