Vitassay qPCR

Flu A + Flu B + RSV

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RSV

<i>Bordetella pertussis</i>	Enterovirus Coxsackievirus A9, A21, A24, B3 yB4	Virus Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus
<i>Bordetella bronchiseptica</i>	Virus parainfluenza humano 1, 2, 3 y 4	Virus Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) virus
<i>Bordetella holmesii</i>	Metapneumovirus humano A y B	Virus Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella parapertussis</i>	Coronavirus humano 229E, OC43 y NL63	Virus Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus
<i>Legionella bozemanii</i>	Rhinovirus humano	Virus Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2)
<i>Legionella micdadei</i>	Bocavirus	Virus Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2)
<i>Legionella dumoffii</i>	Adenovirus humano tipos 1-5, 8, 31, 40 y 41	Virus Influenza A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1
<i>Legionella longbeachae</i>	MERS Coronavirus	Virus Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4
<i>Legionella pneumophila</i>	Virus Influenza A/New Caledonia/20/99(H1N1) virus	Virus Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2)
<i>Mycoplasma pneumoniae</i>	Virus Influenza A/California/7/2009(H1N1)pdm09-like virus	Virus Influenza A/South Australia/55/2014
<i>Streptococcus pneumoniae</i>	Virus Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1)	Virus Influenza A/Uruguay/716/2007 (H3N2)(NYMC-175C)
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Virus Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Virus Influenza B/Brisbane/60/2008-like virus
Methicillin-resistant <i>Staphylococcus aureus</i>	Virus Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	Virus Influenza B/Netherlands/207/06 virus
<i>Haemophilus influenzae</i> MinnA	Virus Influenza A/Perth/16/2009(H3N2)-like virus	Virus Influenza B/Florida/04/06 virus
<i>Moraxella catarrhalis</i>	Virus Influenza A/Thüringen/5/17 (H3N2) virus	Virus Influenza B/Phuket/3073/2013 virus
<i>Chlamydia caviae</i>	Virus Influenza A/Switzerland/9715293/2013 (H3N2) virus	Virus Influenza B/Netherlands/2518/2016 (clade 1A) virus
<i>Chlamydia psittaci</i> genotipos A y C	Virus Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	Virus Influenza B/Netherlands/365/2016 (clade 3)
<i>Chlamydophila pneumoniae</i>	Virus Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)	Virus Influenza B/Colorado/6/2017
Enterovirus 68 y 71	Virus Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a)	Virus Influenza B/Maryland/15/2016
Enterovirus Echovirus tipos 11 y 30	Virus Influenza A/Hong Kong/213/2003 (H5N1) virus	

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Vitassay qPCR Flu A+ Flu B + RSV ha sido evaluado para Influenza A con las siguientes cepas: A/New Caledonia/20/99(H1N1), A/California/7/2009(H1N1)pdm09-like, A/Netherlands/1250/2016 (H1N1)pdm09 (clade 6B.1), A/Michigan/45/2015 (H1N1)pdm09, A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09, A/Perth/16/2009(H3N2)-like, A/Thüringen/5/17 (H3N2), A/Switzerland/9715293/2013 (H3N2), A/Hong Kong/4801/2014, NYMC X-263B (H3N2), A/Netherlands/398/2014 (H3N2) (clade 3C.3a), A/Netherlands/2393/2015 (H3N2) (clade 3C.2a), A/Hong Kong/213/2003 (H5N1), A/Turkey/Germany R2485+86/2014 (H5N8), A/DE-SH/Reiherente/AR8444/ 2013 (H5N8), A/Mallard/Netherlands/2/2009 (H7N7), A/Anhui/1/2013 (H7N9), A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2), A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2), A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1, A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4, A/Singapore/INFIMH-16-0019/2016 (H3N2), A/South Australia/55/2014 y A/Uruguay/716/2007 (H3N2)(NYMC-175C), obteniéndose un resultado positivo para todas ellas.

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ / CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Flu A+ Flu B + RSV allows the detection and differentiation of Flu A, Flu B and Human Respiratory Syncytial Virus (RSV) by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Flu A, Flu B and RSV infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Flu A+ Flu B +RSV 4x8 -well strip, low profile 7041027

Vitassay qPCR Flu A+ Flu B +RSV 4x8-well strip, high profile 7042027

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S027/ 7042S027	Flu A + Flu B + RSV Virus strips low/high profile	-	4 x 8-well strip
7C027	Flu A + Flu B + RSV Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

Summary

Influenza A and B viruses, which belong to the *Orthomyxoviridae* family, cause the majority of viral lower respiratory tract infections; elderly and compromised individuals are especially at risk of developing severe illness and complications. Influenza viruses are a significant cause of morbidity and mortality worldwide.

Seasonal influenza epidemics impose a heavy burden on society, with 3–5 million cases and 250-500 000 deaths worldwide every year. The resulting economic impact is large and includes both direct and indirect costs. Also, influenza A can be divided into numerous subgroups, but the most important ones are the ones that can cause illness in humans, them being, H1N1, H1N2, H2N2, H3N2 and H5N1.

Influenza can be transmitted either by maintaining direct contact with infected individuals, contaminated objects or by inhalation of contaminated aerosols. After an incubation period of one to two days, the illness has an abrupt onset, with chills, fever, myalgia, headache, and anorexia. Pneumonia is the most important complication of influenza. There may be primary viral pneumonia, secondary bacteria pneumonia, or a combination of both.

The human respiratory syncytial virus (RSV) is an enveloped, single stranded linear RNA genome virus. RSV, which belongs to the genus *Orthopneumovirus*, is divided into two major groups, A and B, based on antigenic and genomic differences.

RSV infection most commonly causes a cold-like illness. But it can also cause bronchitis, croup, and lower respiratory infections like bronchiolitis and pneumonia. Of every 100 infants and young children with RSV infection, 25 to 40 (25% to 40%) will show signs of pneumonia or bronchiolitis.

RSV is transmitted via large nasopharyngeal secretion droplets coming from infected individuals. These droplets enter via the mucus membranes of the eyes, nose and mouth following close contact, or self-inoculation after touching contaminated surfaces.

Several different types of laboratory tests are available for diagnosis of Influenza or RSV infection. Laboratory diagnosis of influenza has been accomplished using cell culture, rapid antigen testing, and, more recently, real-time PCR, concretely qPCR that targets the *M1* gene. Whereas for RSV, antigen detection tests and cultures are generally reliable when diagnosing young children, but they are less useful in older children and adults. Use of highly sensitive RT-PCR assays, such as the ones that target the *N*-gene, should be considered, particularly, when testing older children and adults, as they may have low viral loads in the obtained respiratory specimens.

Principle of the test

Vitassay qPCR Flu A+ Flu B + RSV test is based on the real-time amplification of specific conserved fragments of the *M1* gene for Flu A and B and specific conserved fragments of the N gene for RSV. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Flu A+ Flu B + RSV test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Influenza A RNA target sequence is detected through the FAM channel, Influenza B virus RNA target in ROX channel and RSV virus RNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

RIDA® Xtract (r-Biopharm).

Total Nucleic Acid Isolation (TNAI) Kit (ROCHE).

Invisorb® Spin Universal Kit (Stratec)

Positive control preparation

Reconstitute the lyophilized Flu A+ Flu B + RSV Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Influenza A), ROX (Influenza B), Cy5 (RSV) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (Influenza A), ROX (Influenza B) and Cy5 (RSV), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Flu A	Flu B	RSV	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Flu A, Flu B and RSV Positive
-	-	-	+	-	+	Flu A, Flu B and RSV Negative
+	-	-	+/-	-	+	Flu A Virus Positive, Flu B and RSV Negative
+	+	-	+/-	-	+	Flu A and Flu B Viruses Positive, and RSV Negative
+	-	+	+/-	-	+	Flu A and RSV Positive, and Flu B Virus Negative
-	+	-	+/-	-	+	Flu B Virus Positive, Flu A and RSV Negative
-	+	+	+/-	-	+	Flu B and RSV Positive, Flu A Virus Negative
-	-	+	+/-	-	+	RSV Positive, Flu A and Flu B Viruses Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

(+) Positive: Amplification signal

(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 256 throat swabs from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Flu A+ Flu B +RSV and the results were compared with those obtained using CLART® PneumoVir DNA array assay (Genomica). Influenza A virus was detected in 105 samples with Vitassay qPCR Flu A+ Flu B + RSV and in 107 samples using CLART PneumoVir DNA assay. These two samples of difference could not be detected with Vitassay qPCR Flu A+ Flu B + RSV because the amount of template RNA is below the detection limit of the method used.

Influenza B virus was detected in 51 samples with Vitassay qPCR Flu A+ Flu B + RSV and in 52 samples using CLART PneumoVir DNA assay. The amount of RNA template of this non-detected sample was below the detection limit of the method used.

Finally, both molecular diagnostic assays detected 41 positive samples for RSV.

These results indicate the high sensitivity and specificity to detect Influenza A, Influenza B and RSV viruses using the molecular diagnostic test Vitassay qPCR Flu A + Flu B + RSV.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Influenza A, Influenza B and RSV templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Influenza A and Influenza B virus was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
Influenza A		
<i>Bordetella pertussis</i>	<i>Haemophilus influenzae</i> MinnA	Bocavirus
<i>Bordetella bronchiseptica</i>	<i>Moraxella catarrhalis</i>	Human Adenovirus types 1-5, 8, 31, 40 and 41
<i>Bordetella holmesii</i>	<i>Chlamydia caviae</i>	MERS Coronavirus
<i>Bordetella parapertussis</i>	<i>Chlamydia psittaci</i> genotype A and C	Influenza B/Brisbane/60/2008-like virus
<i>Legionella bozemanii</i>	<i>Chlamydomytila pneumoniae</i>	Influenza B/Netherlands/207/06 virus
<i>Legionella micdadei</i>	Enterovirus 68 and 71	Influenza B/Florida/04/06 virus
<i>Legionella dumoffii</i>	Enterovirus Echovirus types 11 and 30	Influenza B/Phuket/3073/2013 virus
<i>Legionella longbeachae</i>	Enterovirus Coxsackievirus A9, A21, A24, B3 and B4	Influenza B/Netherlands/2518/2016 (clade 1A) virus
<i>Legionella pneumophila</i>	Human parainfluenza 1, 2, 3 and 4 viruses	Influenza B/Netherlands/365/2016 (clade 3)
<i>Mycoplasma pneumoniae</i>	Human metapneumovirus A and B	Influenza B/Colorado/6/2017
<i>Streptococcus pneumoniae</i>	Human coronavirus 229E, OC43 and NL63	Influenza B/Maryland/15/2016
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Human rhinovirus	Respiratory syncytial virus (RSV)
Methicillin-resistant <i>Staphylococcus aureus</i>		

Cross-reactivity assay

Influenza B

<i>Bordetella pertussis</i>	Enterovirus 68 and 71	Influenza A/Switzerland/9715293/2013 (H3N2) virus
<i>Bordetella bronchiseptica</i>	Enterovirus Echovirus types 11 and 30	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus
<i>Bordetella holmesii</i>	Enterovirus Coxsackievirus A9, A21, A24, B3 and B4	Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)
<i>Bordetella parapertussis</i>	Human parainfluenza 1, 2, 3 and 4 viruses	Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a)
<i>Legionella bozemanii</i>	Human metapneumovirus A and B	Influenza A/Hong Kong/213/2003 (H5N1) virus
<i>Legionella micdadei</i>	Human coronavirus 229E, OC43 and NL63	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus
<i>Legionella dumoffii</i>	Human rhinovirus	Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) virus
<i>Legionella longbeachae</i>	Bocavirus	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Legionella pneumophila</i>	Human Adenovirus types 1-5, 8, 31, 40 and 41	Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus
<i>Mycoplasma pneumoniae</i>	MERS Coronavirus	Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2)
<i>Streptococcus pneumoniae</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2)
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza A/California/7/2009(H1N1)pdm09-like virus	Influenza A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1
Methicillin-resistant <i>Staphylococcus aureus</i>	Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1)	Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4
<i>Haemophilus influenzae</i> MinnA	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2)
<i>Moraxella catarrhalis</i>	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	Influenza A/South Australia/55/2014
<i>Chlamydia caviae</i>	Influenza A/Perth/16/2009(H3N2)-like virus	Influenza A/Uruguay/716/2007 (H3N2)(NYMC-175C)
<i>Chlamydia psittaci</i> genotypes A and C	Influenza A/Thüringen/5/17 (H3N2) virus	Respiratory syncytial virus (RSV)
<i>Chlamydophila pneumoniae</i>		

Cross-reactivity assay

RSV

<i>Bordetella pertussis</i>	Enterovirus Coxsackievirus A9, A21, A24, B3 and B4	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus
<i>Bordetella bronchiseptica</i>	Human parainfluenza 1, 2, 3 and 4 viruses	Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) virus
<i>Bordetella holmesii</i>	Human metapneumovirus A and B	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella parapertussis</i>	Human coronavirus 229E, OC43 and NL63	Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus
<i>Legionella bozemanii</i>	Human rhinovirus	Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2)
<i>Legionella micdadei</i>	Bocavirus	Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2)
<i>Legionella dumoffii</i>	Human Adenovirus types 1-5, 8, 31, 40 and 41	Influenza A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1
<i>Legionella longbeachae</i>	MERS Coronavirus	Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4
<i>Legionella pneumophila</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2)
<i>Mycoplasma pneumoniae</i>	Influenza A/California/7/2009(H1N1)pdm09-like virus	Influenza A/South Australia/55/2014
<i>Streptococcus pneumoniae</i>	Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1)	Influenza A/Uruguay/716/2007 (H3N2)(NYMC-175C)
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Influenza B/Brisbane/60/2008-like virus
Methicillin-resistant <i>Staphylococcus aureus</i>	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	Influenza B/Netherlands/207/06 virus
<i>Haemophilus influenzae</i> MinnA	Influenza A/Perth/16/2009(H3N2)-like virus	Influenza B/Florida/04/06 virus
<i>Moraxella catarrhalis</i>	Influenza A/Thüringen/5/17 (H3N2) virus	Influenza B/Phuket/3073/2013 virus
<i>Chlamydia caviae</i>	Influenza A/Switzerland/9715293/2013 (H3N2) virus	Influenza B/Netherlands/2518/2016 (clade 1A) virus
<i>Chlamydia psittaci</i> genotypes A and C	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	Influenza B/Netherlands/365/2016 (clade 3)
<i>Chlamydomydia pneumoniae</i>	Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)	Influenza B/Colorado/6/2017
Enterovirus 68 and 71	Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a)	Influenza B/Maryland/15/2016
Enterovirus Echovirus types 11 and 30	Influenza A/Hong Kong/213/2003 (H5N1) virus	

Analytical reactivity

Analytical reactivity of Vitassay qPCR Flu A+ Flu B + RSV for Influenza A was evaluated against strains: A/New Caledonia/20/99(H1N1) virus, A/California/7/2009(H1N1)pdm09-like virus, A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1), A/Michigan/45/2015 (H1N1)pdm09 virus, A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus, A/Perth/16/2009(H3N2)-like virus, A/Thüringen/5/17 (H3N2) virus, A/Switzerland/9715293/2013 (H3N2) virus, A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus, A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a), A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a), A/Hong Kong/213/2003 (H5N1) virus, A/Turkey/Germany R2485+86/2014 (H5N8) virus, A/DE-SH/Reiherente/AR8444/2013 (H5N8) virus, A/Mallard/Netherlands/2/2009 (H7N7) virus, A/Anhui/1/2013 (H7N9) virus, A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2), A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2), A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1, A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4, A/Singapore/INFIMH-16-0019/2016 (H3N2), A/South Australia/55/2014 and A/Uruguay/716/2007 (H3N2)(NYMC-175C), showing positive result.

Analytical reactivity of Vitassay qPCR Flu A+ Flu B + RSV for Influenza B was evaluated against strains: B/Brisbane/60/2008-like virus, B/Netherlands/207/06 virus, B/Florida/04/06 virus, B/Phuket/3073/2013 virus, B/Netherlands/2518/2016 (clade 1A) virus, B/Netherlands/365/2016 (clade 3), B/Colorado/6/2017 and B/Maryland/15/2016, showing positive result.

Analytical reactivity of Vitassay qPCR Flu A+ Flu B + RSV for RSV was evaluated against Human Respiratory Syncytial Virus (RSV) A and B, showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Flu A+ Flu B has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Flu A, Flu B and RSV infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swab samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™/ CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Methods for molecular surveillance of influenza. Wang R, Taubenberger JK. Expert review of anti-infective therapy. 2010;8(5):517-527.
2. WHO information for molecular diagnosis of influenza virus—update. World Health Organization. Available: [http://www.who.int/influenza/gisrs_laboratory/molecular_ diagnosis/en/](http://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/). Accessed 2015 Dec 30.
3. Kuo RL, Yang SL, Liu YC, Chen LT, Mok CK, Kuo SM, Shih SR, Tsao KC. Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. J Virol Methods. 2014 Nov; 208:41-6.
4. Bawage SS, Tiwari PM, Pillai S, Dennis V, Singh SR. Recent advances in diagnosis, prevention, and treatment of human respiratory syncytial virus. Adv Virol. 2013;2013:595768.
5. Falsey AR. Respiratory syncytial virus infection in adults. Semin Respir Crit Care Med. 2007 Apr;28(2):171-81.
6. Van Woensel JB, Kimpen JL, Brand PL. Respiratory tract infections caused by respiratory syncytial virus in children. Diagnosis and treatment. Minerva Pediatr. 2001 Apr;53(2):99-106.
7. I.H. Brown, et al. Multiple genetic reassortment of avian and human influenza A viruses in European pigs, resulting in the emergence of an H1N2 virus of novel genotype. Journal of General Virology, 1998; 79: 2947-2955.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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