# Vitassay qPCR

## Japanese Encephalitis

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## Intended use

Vitassay qPCR Japanese Encephalitis allows the qualitative detection of Japanese encephalitis virus by real-time RT-PCR in blood or CSF samples. The product is intended for use in the diagnosis of Japanese encephalitis infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Japanese Encephalitis 4 x 8-well strip, low profile 7041038

Vitassay qPCR Japanese Encephalitis 4 x 8-well strip, high profile 7042038

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S038/ 7042S038	Japanese Encephalitis strips	-	4x8-well strip
7C038	Japanese Encephalitis Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Summary

Japanese encephalitis (JE) virus, a flavivirus, is closely related to West Nile and St. Louis encephalitis viruses. JE virus is transmitted to humans through the bite of infected *Culex* species mosquitoes, particularly *Culex tritaeniorhynchus*.

The virus is maintained in a cycle between mosquitoes and vertebrate hosts, primarily pigs and wading birds. Humans are incidental or dead-end hosts, because they usually do not develop high enough concentrations of JE virus in their bloodstreams to infect feeding mosquitoes.

JE virus transmission occurs primarily in rural agricultural areas, often associated with rice production and flooding irrigation. In some areas of Asia, these conditions can occur near urban centers.

Less than 1% of people infected with Japanese encephalitis (JE) virus develop clinical illness. In persons who develop symptoms, the incubation period (time from infection until illness) is typically 5-15 days. Initial symptoms often include fever, headache, and vomiting. Mental status changes, neurologic symptoms, weakness, and movement disorders might develop over the next few days. Seizures are common, especially among children.

## Principle of the test

Vitassay qPCR Japanese Encephalitis is based on the real-time amplification of specific conserved fragments of the *NS2A* gene encoded by the Japanese Encephalitis genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of Japanese Encephalitis is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Japanese Encephalitis is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and RNA extraction**

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

QIAamp Viral RNA Mini Kit, QIAcube instrument (Qiagen).

Invisorb® Spin Universal Kit (Stratec).

### **Positive control preparation**

Reconstitute the lyophilized Japanese Encephalitis Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

## Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Japanese Encephalitis) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

## Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

### Positive control

The positive controls used in each run, must show an amplification curve for Japanese Encephalitis, which validates the reaction.

## Negative control

The negative controls included in each run, must show the absence of signal for Japanese Encephalitis which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Japanese Encephalitis FAM	Internal Control HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	Japanese Encephalitis Positive
-	+	-	+	Japanese Encephalitis Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

Vitassay qPCR Japanese Encephalitis was evaluated with QCMD programs of tropical viruses (West Nile Virus and Dengue) of the years 2014, 2017 and 2018. These programs contained 6 positive samples for Japanese encephalitis and were detected correctly. The results were compared with the final reports.

These results indicate that Vitassay qPCR Japanese Encephalitis shows a high sensitivity and specificity for detecting Japanese encephalitis virus.

## Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Japanese Encephalitis template ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

## Analytical specificity

The analytical specificity for Japanese Encephalitis was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Japanese Encephalitis		
<i>Candida albicans</i>	Parechovirus Type 3	Human herpesvirus 6 (HHV6) Type B
<i>Escherichia coli</i> 0.1285;O18:H7:K1	Citomegalovirus strain AD169	Zika Virus (African strain)
<i>Haemophilus influenzae</i> MinnA	Varicella Zoster Virus Ellen	Zika Virus (Asian strain PF13/251013-18)
<i>Enterococcus durans</i>	Varicella Zoster Virus OKA	Zika Virus (French Polynesian strain 11474/16)
<i>Enterococcus faecalis</i>	Coxsackievirus A24	Chikungunya virus strain S27 Petersfield
<i>Enterococcus faecium</i> serotype 11	Coxsackievirus A9	Dengue 1 virus strain Hawaii
<i>Cryptococcus gatti</i> Z156	Coxsackievirus B3	Dengue 2 virus strain New Guinea C
<i>Listeria monocytogenes</i> serovar 1/2b	Echovirus 30	Dengue 3 virus strain H87
<i>Listeria monocytogenes</i> serovar 4b	Echovirus Type 11	Dengue 4 virus strain H241
<i>Listeria innocua</i> serovar 6a	Enterovirus 68	St Louis Encephalitis virus strain 17D
<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	Enterovirus 71	West Nile virus strain NY99
<i>Neisseria meningitidis</i> serogroup A	Epstein Barr virus	West Nile virus Heja
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Herpes simplex virus 1 (HSV1) MacIntyre strain	West Nile virus Ug37
<i>Streptococcus agalactiae</i> Z019	Herpes simplex virus 2 (HSV2) MS	Tick-Borne encephalitis virus (cepa Neudorfl)
<i>Streptococcus pneumoniae</i> Z022	Human herpesvirus 6 (HHV6) strain Z29	Yellow Fever Virus (17D strain)
<i>Plasmodium falciparum</i> (3D7 strain)	Human herpesvirus 6 (HHV6) Type A	

## Analytical reactivity

The reactivity of Vitassay qPCR Japanese Encephalitis was confirmed by the real-time amplification using Japanese Encephalitis as template.

## Compatibles real-time PCR equipment

Vitassay qPCR Japanese Encephalitis has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- StepOne™ Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTIite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of Japanese Encephalitis infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with blood samples and CSF. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Japanese Encephalitis, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	<b>Eppendorf</b>
	Mastercycler™ ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene®Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

### Bibliography/Bibliografía

1. <https://www.cdc.gov/japaneseencephalitis/>

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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