# Vitassay qPCR

## HCoV

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**Intended use**

Vitassay qPCR HCoV allows the detection and differentiation of Coronavirus 229E, NL63, OC43 and/or HKU1 subtypes by real-time RT-PCR in respiratory samples. The product is intended for use in the diagnosis of Coronavirus 229E, NL63, OC43 and/or HKU1 subtypes infections alongside clinical data of the patient and other laboratory tests outcomes.

**References**

- Vitassay qPCR HCoV 4x8 -well strip, low profile 7041040
- Vitassay qPCR HCoV 4x8-well strip, high profile 7042040

**Materials/reagents provided**

Reference	Reagent/Material	Colour	Amount
7041S040/ 7042S040	HCoV strips low/high profile	-	4 x 8-well strip
7C040	VITASSAY qPCR HCoV Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

**Transport and storage**

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

**Additional equipment and material required**

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

**Summary**

Coronaviruses are enveloped RNA viruses and are subdivided in 6 groups: HCoV-229E, HCoV-NL63, HCoV-OC43, HCoV-HKU1, severe acute respiratory syndrome coronavirus (MERS-CoV) and Middle East respiratory syndrome coronavirus (SARS-CoV). Coronaviruses were first identified in the mid-1960s as one of the main causes of the common cold. The symptoms of the infection, once past the incubation period, is fever of 1 to 5 days. In the most severe cases it can cause shortness of breath, gastrointestinal symptoms, nausea, vomiting, fatigue and myalgia and have a longer duration.

The transmission of the virus can occur between people by inhaling droplets of water from coughs or sneezes or by indirect transmission when in contact with contaminated materials such as equipment, surfaces and devices. The incubation period of the disease is 5 to 12 days.

The three strains have very similar symptoms and the low effectiveness of the microbiological methods make real-time PCR one of the most effective methods to improve the low specificity of the other methods. This technique makes it possible to detect the replicase gene in the HKU1 strain and the nucleoprotein gene in strains 229E, NL63 and OC43.

**Principle of the test**

Vitassay qPCR HCoV test is based on the real-time amplification of specific conserved fragments of a conserved region of the *nucleoprotein* gene for strains 229E, OC43, NL63 and HKU1 and the as targets, using specific primers and a fluorescent-labelled probes. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately by polymerase chain reaction. The presence of the RNA viral is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR HCoV test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected in the FAM (Coronavirus 229E), ROX (Coronavirus NL63), Cy5 (Coronavirus OC43) and HEX, VIC or JOE channels (Coronavirus HKU1). Choose channel according to the equipment used, see Annex II.

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## Procedures

### Specimen collection, processing and RNA extraction

For pre-treatment and RNA isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (Roche).

Invisorb® Spin Universal Kit (Strattec)

### Positive control preparation

Reconstitute the lyophilized HCoV Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

### Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Coronavirus 229E), ROX (Coronavirus NL63), Cy5 (Coronavirus OC43) and HEX, VIC or JOE channels (Coronavirus HKU1). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

### **Analysis and interpretation of results**

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

#### **Positive control**

The positive controls used in each run, must show an amplification curve for FAM (Coronavirus 229E), ROX (Coronavirus NL63), Cy5 (Coronavirus OC43) and/or HEX, VIC or JOE channels (Coronavirus HKU1), which validates the reaction.

#### **Negative control**

The negative controls included in each run, must show the absence of signal in FAM (Coronavirus 229E), ROX (Coronavirus NL63), Cy5 (Coronavirus OC43) and HEX, VIC or JOE channels (Coronavirus HKU1), which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Coronavirus 229E (FAM)	Coronavirus HKU1 (HEX)	Coronavirus NL63 (ROX)	Coronavirus OC43 (Cy5)	Negative Control	Positive Control	Interpretation
+	+	+	+	-	+	Coronavirus 229E, HKU1, NL63 and OC43 Positive
-	-	-	-	-	+	Coronavirus 229E, HKU1, NL63 and OC43 Negative
+	-	-	-	-	+	Coronavirus 229E Positive, Coronavirus HKU1, NL63 and OC43 Negative
-	+	-	-	-	+	Coronavirus HKU1 Positive, Coronavirus 229E, NL63 and OC43 Negative
-	-	+	-	-	+	Coronavirus NL63 Positive, Coronavirus 229E, HKU1 and OC43 Negative
-	-	-	+	-	+	Coronavirus OC43 Positive, Coronavirus 229E, HKU1 and NL63 Negative
+	+	-	-	-	+	Coronavirus 229E and HKU1 positive, Coronavirus NL63 and OC43 Negative
+	-	+	-	-	+	Coronavirus 229E and NL63 positive, Coronavirus HKU1 and OC43 Negative
+	-	-	+	-	+	Coronavirus 229E and OC43 positive, Coronavirus HKU1 and NL63 Negative
-	+	+	-	-	+	Coronavirus HKU1 and NL63 positive, Coronavirus 229E and OC43 Negative
-	+	-	+	-	+	Coronavirus HKU1 and OC43 positive, Coronavirus 229E and NL63 Negative
-	-	+	+	-	+	Coronavirus NL63 and OC43 positive, Coronavirus 229E and HKU1 Negative
+	+	+	-	-	+	Coronavirus 229E, HKU1 and NL63 Positive, Coronavirus OC43 Negative
-	+	+	+	-	+	Coronavirus HKU1, NL63 and OC43 Positive, Coronavirus 229E Negative
+	-	+	+	-	+	Coronavirus 229E, NL63 and OC43 Positive, Coronavirus HKU1 Negative
+	+	-	+	-	+	Coronavirus 229E, HKU1 and OC43 Positive, Coronavirus NL63 Negative
+	+	+	+	+	+	Experiment Fail
-	-	-	-	-	-	Experiment Fail

**(+) Positive:** Amplification signal

**(-) Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

### Quality Control

An Internal Control (IC) is included in each reaction for confirming the appropriate performance of the molecular diagnostic technique. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

A total of 45 clinical samples from the EQA panels were tested by Vitassay qPCR HCoV.

Vitassay pPCR HCoV detected Coronavirus 229E in 7 of the 45 clinical samples tested.

Vitassay qPCR HCoV detected Coronavirus HKU1 in 1 of the 45 clinical samples tested.

Vitassay qPCR HCoV detected Coronavirus NL63 in 11 of the 45 clinical samples tested.

Vitassay qPCR HCoV detected Coronavirus OC43 in 18 of the 45 clinical samples tested.

The result of Vitassay qPCR HCoV kit were concordant for all the samples. The results show a high sensitivity and specificity to detect Coronavirus 229E, NL63, OC43 and HKU1 using the molecular diagnostic test Vitassay qPCR HCoV.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of each virus templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

### Analytical specificity

The analytical specificity for Coronavirus 229E, NL63, OC43 y HKU1 was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity testing		
<i>Bordetella pertussis</i>	<i>Moraxella catarrhalis</i>	Virus Influenza A/Anhui/1/2013 (H7N9)
<i>Bordetella parapertussis</i>	<i>Haemophilus influenzae MinnA</i>	Virus Influenza B/Brisbane/60/2008-like
<i>Bordetella holmesii</i>	<i>Chlamydomphila pneumoniae</i>	Virus Influenza B/Florida/04/06
<i>Bordetella bronchiseptica</i>	<i>Chlamydia caviae</i>	Virus Influenza B/Phuket/3073/2013
<i>Legionella bozemanii</i>	<i>Chlamydia psittaci</i> genotypes A and C	Human parainfluenza 1, 2, 3 and 4 viruses
<i>Legionella micdadei</i>	Virus Influenza A/New Caledonia/20/99(H1N1)	Human metapneumovirus A and B
<i>Legionella dumoffii</i>	Virus Influenza A/California/7/2009(H1N1)	Human coronavirus 229E, OC43 and NL63
<i>Legionella pneumophila</i>	Virus Influenza A/Michigan/45/2015 (H1N1)pdm09	MERS Coronavirus
<i>Legionella longbeachea</i>	Virus Influenza A/Perth/16/2009(H3N2)	Human rhinovirus
<i>Mycoplasma pneumoniae</i>	Virus Influenza A/Thüringen/5/17 (H3N2)	Human Adenovirus
<i>Streptococcus pneumoniae</i>	Virus Influenza A/Switzerland/9715293/2013 (H3N2)	Respiratory Syncytial virus (RSV)
<i>Staphylococcus aureus subsp. aureus</i>	Virus Influenza A/Turkey/Germany R2485+86/2014 (H5N8)	Human Bocavirus
Influenza A/Singapore/GP1908/2015 virus	Influenza A/Hong Kong/4801/2014(H3N2) virus	Virus Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8)

### Analytical reactivity

The reactivity of de Vitassay qPCR HCoV was evaluated against the strains: Human coronavirus 229E, NL63, OC43 and HKU1, showing positive result in all of them.

### Compatibles real-time PCR equipment

Vitassay qPCR HCoV has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>

- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cyclers (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® (Qiagen)<sup>1</sup>
- SmartCycler® (Cepheid)<sup>1</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using Applied Biosystems with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

### **Limitations**

- This test provides a presumptive diagnosis of Coronavirus 229E, NL63, OC43 and HKU1 infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with RNA extracted from throat samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different viruses, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

**Attached I: Compatibility of the most common real-time PCR equipment**

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers
<b>Agilent Technologies</b>
AriaMx Real-Time PCR System
<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
<b>BIONEER</b>
Exicycler™ 96
<b>Bio-Rad</b>
CFX96™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
<b>Cepheid</b>
SmartCycler®
<b>Qiagen</b>
Rotor-Gene® Q
<b>Roche</b>
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
<b>Abbott</b>
Abbott m2000 RealTime System
<b>Applied Biosystems</b>
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR
<b>Analytik Jena Biometra</b>
TOptical
qTOWER 2.0
<b>BIONEER</b>
Exicycler™ 96
<b>Bio-Rad</b>
CFX96™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
<b>Cepheid</b>
SmartCycler®
<b>DNA-Technology</b>
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cyclor
<b>Eppendorf</b>
Mastercycler™ep <i>realplex</i>
<b>Qiagen</b>
Rotor-Gene® Q
<b>Stratagene / Agilent Technologies</b>
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

**Attached II: Detection channels of most common real time PCR equipment**

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

**Attached III: Optical measurement exposure setting**

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

<b>Thermocycler</b>	<b>Vitassay channel</b>	<b>Exposure values</b>
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## Bibliography/Bibliografía

1. Ryan K. Dare, Alicia M. Fry, Malinee Chittaganpitch, Pathom Sawanpanyalert, Sonja J. Olsen, and Dean D. Erdman. Human Coronavirus Infections in Rural Thailand: A Comprehensive Study Using Real-Time Reverse-Transcription Polymerase Chain Reaction Assays, *The Journal of Infectious Diseases* 2007; 196:1321–8
2. Narrative review of Middle East respiratory syndrome coronavirus (MERS-CoV) infection: updates and implications for practice, Abbas Al Mutair and Zainab Ambani. *Journal of International Medical Research*, 0(0) 1–6
3. Susan S. Chiu, Kwok Hung Chan, Ka Wing Chu, See Wai Kwan, Yi Guan, Leo Lit Man Poon and J. S. M. Peiris. Human Coronavirus NL63 Infection and Other Coronavirus Infections in Children Hospitalized with Acute Respiratory Disease in Hong Kong, China. *Clinical Infectious Diseases* 2005; 40:1721–9 2005 by the Infectious Diseases Society of America
4. Kate E. Templeton, Sitha A. Scheltinga, Willian C. J. F. M. van den Eeden, A. Willy Graffelman, Peterhans J. van den Broek and Eric C. J. Claas. Improved Diagnosis of the Etiology of Community-Acquired Pneumonia with Real-Time Polymerase Chain Reaction. *Clinical Infectious Diseases* 2005; 41:345–51

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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