# Vitassay qPCR

## Primary Atypical Pneumonia

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## Intended use

Vitassay qPCR Primary Atypical Pneumonia allows the detection and differentiation of *Chlamydophila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila* by real-time PCR in respiratory samples. This test is intended to be used for research purposes, without any medical objective is not regarded as devices for performance evaluation.

## References

Vitassay qPCR Primary Atypical Pneumonia 4 x 8-well strip, low profile	7041049R
Vitassay qPCR Primary Atypical Pneumonia 4 x 8-well strip, high profile	7042049R
Vitassay qPCR Primary Atypical Pneumonia 8 x 8-well strip, low profile	7081049R
Vitassay qPCR Primary Atypical Pneumonia 8 x 8-well strip, high profile	7082049R
Vitassay qPCR Primary Atypical Pneumonia 96-well plate, low profile	7091049R
Vitassay qPCR Primary Atypical Pneumonia 96-well plate, high profile	7092049R

## Materials/reagents provided

Reagents provided in references 7041049R, 7042049R, 7081049R and 7082049R:

Reference	Reagent/Material	Colour	Amount
7041S049R/ 7042S049R	Primary Atypical Pneumonia strips	-	4/8 x 8-well strip
7C049R	Primary Atypical Pneumonia Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4/8 x 8 cap strip

Reagents provided in references 7091049R and 7092049R:

Reference	Reagent/Material	Colour	Amount
7091P049R/ 7092P049R	Primary Atypical Pneumonia Plate	-	1 plate
7C049R	Primary Atypical Pneumonia Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

## Summary

Pneumonia is an infection of the lungs that can cause mild to severe illness in people of all ages. Community-acquired pneumonia (CAP) is an infection acquired in the community, outside of the hospital setting, which is associated with high rates of morbidity and mortality worldwide. Intracellular bacteria are common causes of CAP. Some intracellular pathogens that are well established as causes of pneumonia are *Legionella pneumophila*, *Mycoplasma pneumoniae* and *Chlamydia pneumoniae*.

The genus *Legionella* includes Gram-negative microorganisms living in natural and artificial water systems. To date, about 60 species of *Legionella* are known. *Legionella pneumophila* (*Lpn*) is the species most frequently associated with human disease and includes 16 serogroups.

The most common form of transmission of *Legionella* is inhalation of contaminated aerosols. Aerosolization of contaminated water supplies allows the bacteria to be inhaled into the human lung, where *L. pneumophila* can be phagocytosed by alveolar macrophages and replicate intracellularly. There is no direct human-to-human transmission.

Most *Legionella* infections are thought to be asymptomatic or lead to a self-limiting and mild respiratory disease (non-pneumonic form). However, some individuals, particularly the elderly and those with weakened immune systems or chronic lung diseases, are at

increased risk of developing severe pneumonia (Legionnaires' disease). Of the reported cases 75–80% are over 50 years and 60–70% are male.

The non-pneumonic form has an incubation period from a few and up to 48 hours and the main symptoms are fever, chills, headache, malaise and muscle pain (myalgia). The main symptoms of Legionnaires' disease are fever, mild cough, phlegm, shortness of breath, loss of appetite, headache, malaise and lethargy. Some patients may also have muscle pain, diarrhoea and confusion. Symptoms usually begin 2 to 10 days after being exposed to the bacteria, but it can take longer.

With over 100 different species, the genus *Mycoplasma* is a unique bacterium that lacks a cell wall and causes a wide range of symptoms and infections. *Mycoplasma pneumoniae* (MP) is the most common pathogenic species of the genus *Mycoplasma* infecting humans. MP is an atypical respiratory bacterium that mainly affects the lining of the respiratory system (throat, lungs, windpipe) in children and adults of all ages. MP causes up to 40% of community-acquired pneumonia (CAP) in children and about 18% of infections in patients requiring hospitalization.

*M. pneumoniae* spread from person to person by coughing or sneezing, which creates small respiratory droplets in the air that contain the bacteria which other people then breathe in. Once in the body, symptoms usually begin between 1 to 4 weeks.

The most common type of illness, especially in children, is tracheobronchitis, commonly known as a chest cold. Common symptoms of a chest cold include sore throat, fatigue, fever, slowly worsening cough that can last for weeks or months and headache. Experts have estimated that up to 1 in 10 people who get ill from MP get pneumonia. Common symptoms of pneumonia include cough that may produce mucus, fever, chills, shortness of breath, chest pain and fatigue.

*Chlamydophila pneumoniae* (*C. pneumoniae*), also known as *Chlamydia pneumoniae*, is a type of bacteria of the genus *Chlamydia* that can cause respiratory tract infections, such as CAP. *C. pneumoniae* accounts for 6 to 20% of CAP cases, depending on several factors such as setting of the studied population, age group examined, and diagnostic methods used. The bacteria affect the lining of the respiratory tract including the throat, windpipe, and lungs and can cause chronic infection which might contribute to chronic conditions, such as asthma, arthritis and atherosclerosis.

The most common form of transmission is direct human-to-human. People spread *C. pneumoniae* by coughing or sneezing, which creates small respiratory droplets that contain the bacteria. Other people in close contact then breathe in the bacteria. People can also get sick if they touch something with droplets from a sick person on it and then touch their mouth or nose.

People of all ages can get sick from *C. pneumoniae*. It most commonly infects people for the first time when they are school-aged children or young adults. However, older adults are at increased risk for severe disease caused by *C. pneumoniae* infection, including pneumonia and reinfection is most common in older adults.

In general, *C. pneumoniae* infection is a mild illness that most commonly causes an upper respiratory tract infection which symptoms can include a sore throat, runny or stuffy nose, fatigue, low-grade fever, hoarseness or loss of voice, slowly worsening cough and headache. *C. pneumoniae* can also cause lower respiratory tract infections like bronchitis and lung infections like pneumonia. Symptoms usually begin 3 to 4 weeks after exposure.

### **Principle of the test**

Vitassay qPCR Primary Atypical Pneumonia test is based on the real-time amplification of specific conserved fragments of the *argR* gene (*Chlamydomphila pneumoniae*), *CARD5* (*Mycoplasma pneumoniae*) and *mip* (*L. pneumophila*). It is detected by an increment in the fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' nuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Primary Atypical Pneumonia is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Legionella pneumophila* is detected through the FAM channel, the amplification of the *Chlamydomphila pneumoniae* is detected through the ROX channel, the amplification of the *Mycoplasma pneumoniae* is detected through the Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

### **Precautions**

- For Research Use Only.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.

- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and DNA extraction**

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

RIDA® Xtract (r-Biopharm).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).

ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

NX-48 Bacterial DNA Kit, using the Nexttractor® NX-48 system (Genolution).

MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.).

Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (Roche).

Invisorb® Spin Universal Kit (Stratec).

### **Positive control preparation**

Reconstitute the lyophilized Primary Atypical Pneumonia Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

### Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (*Legionella pneumophila*), HEX, JOE o VIC (Internal Control), ROX (*Chlamydomphila pneumoniae*) and Cy5 (*M. pneumoniae*). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none (Attached II).

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer’s instructions.

For a valid diagnostic test run, the following control conditions must be met:

**Positive control**

The positive controls used in each run, must show an amplification curve for FAM (*Legionella pneumophila*), HEX, JOE o VIC (Internal Control), ROX (*Chlamydomphila pneumoniae*) and Cy5 (*M. pneumoniae*), which validates the reaction.

**Negative control**

The negative controls included in each run, must show the absence of signal in FAM, ROX, HEX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table

<i>L. pneumophila</i> (FAM)	<i>C. pneumoniae</i> (ROX)	<i>M. pneumoniae</i> (Cy5)	Internal Control HEX	Negative Control	Positive Control	Interpretation
-----------------------------	----------------------------	----------------------------	----------------------	------------------	------------------	----------------

+	+	+	+/-	-	+	<i>L. pneumophila</i> , <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Positive
-	-	-	+	-	+	<i>L. pneumophila</i> , <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Negative
+	-	-	+/-	-	+	<i>L. pneumophila</i> Positive, <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Negative
+	+	-	+/-	-	+	<i>L. pneumophila</i> and <i>C. pneumoniae</i> Positive, and <i>M. pneumoniae</i> Negative
+	-	+	+/-	-	+	<i>L. pneumophila</i> and <i>M. pneumoniae</i> Positive, and <i>C. pneumoniae</i> Negative
-	+	-	+/-	-	+	<i>C. pneumoniae</i> Positive, <i>L. pneumophila</i> and <i>M. pneumoniae</i> Negative
-	+	+	+/-	-	+	<i>C. pneumoniae</i> and <i>M. pneumoniae</i> Positive, <i>L. pneumophila</i> Negative
-	-	+	+/-	-	+	<i>M. pneumoniae</i> Positive, <i>L. pneumophila</i> and <i>C. pneumoniae</i> Negative
-	-	-	-	-	-	Experiment fail
+	+	+	+	+	+	Experiment fail

**(+) Positive:** Amplification signal

**(-) Negative:** No amplification signal

If no amplification signal is detected in the negative control and / or no signal is present in the negative control, the experiment is considered unsuccessful.

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

### **Quality Control**

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

### **Preliminary results of the test characteristics**

#### **Clinical sensitivity and specificity**

A total of 90 respiratory samples were analysed by Vitassay strip qPCR Primary Atypical Pneumonia and (FTD Bacterial pneumonia CAP (Fast Track Diagnostics).

Vitassay qPCR Primary Atypical Pneumonia detected *C. pneumoniae* in 5 positive samples, *M. pneumoniae* in 4 positive samples y *L. pneumophila* in 1 positive sample.

1 discordant sample (samples determined at the limit of detection) were obtained in *M. pneumoniae*.

These results show the high sensitivity and specificity to detect *C. pneumoniae*, *M. pneumoniae* and *L. pneumophila* using the molecular diagnostic test Vitassay qPCR Primary Atypical Pneumonia.

#### **Analytical sensitivity**

The analytical sensitivity was determined by analysis of 10-fold dilution series of the template DNA of the different pathogens ( $10^7$ - $10^1$  copies / reaction). This assay has a limit of detection of  $\geq 10$  copies of DNA per reaction.

#### **Analytical specificity**

The analytical specificity for Primary Atypical Pneumonia was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

## Analytical reactivity

Primary Atypical Pneumonia		
<i>Bordetella pertussis</i>	<i>Moraxella catarrhalis</i>	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella parapertussis</i>	<i>Haemophilus influenzae</i>	Influenza B/Brisbane/60/2008-like virus
<i>Bordetella holmesii</i>	<i>Chlamydia caviae</i>	Influenza B/Florida/04/06 virus
<i>Bordetella bronchiseptica</i>	<i>Chlamydia psittaci</i> genotypes A and C	Influenza B/Phuket/3073/2013 virus
<i>Legionella bozemanii</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Human parainfluenza 1, 2, 3 and 4 viruses
<i>Legionella micdadei</i>	Influenza A/California/7/2009(H1N1) virus	Human metapneumovirus A and B
<i>Legionella dumoffii</i>	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Human coronavirus 229E
<i>Legionella longbeache</i>	Influenza A/Perth/16/2009(H3N2) virus	MERS Coronavirus
<i>Legionella hackeliae</i>	Influenza A/Thüringen/5/17 (H3N2) virus	Human rhinovirus
<i>Streptococcus pneumoniae</i>	Influenza A/Switzerland/9715293/2013 (H3N2) virus	Human Adenovirus
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus	Respiratory Syncytial virus (RSV A and RSV B)
<i>Methicillin-resistant Staphylococcus aureus</i>	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	Human Bocavirus

The reactivity of Vitassay qPCR Primary Atypical Pneumonia for *Chlamydomphila pneumoniae* was confirmed by the real time amplification using *Chlamydomphila pneumoniae*, as template.

The reactivity of Vitassay qPCR Primary Atypical Pneumonia for *Mycoplasma pneumoniae* was confirmed by the real time amplification using *Mycoplasma pneumoniae* as template.

The reactivity of Vitassay qPCR Primary Atypical Pneumonia for *L. pneumophila* was confirmed by the real time amplification using *Legionella pneumophila* sg1 (ST47) and (ST62), sg2-14, sg3, and sg6, *Legionella pneumophila* subsp. *pneumophila* Brenner et al. 1979 strain Philadelphia-1, sg 1, *Legionella pneumophila* subsp. *pneumophila* Brenner et al. strain Togus-1, sg 2, *Legionella pneumophila* subsp. *pneumophila* Brenner et al. strain Bloomingtoon-2, sg 3, *Legionella pneumophila* subsp. *pneumophila* Brenner et al. strain Concord 3, sg 8, *Legionella pneumophila* subsp. *fraseri* Brenner et al. 1989 strain

Los Angeles-1, sg 4, *Legionella pneumophila* subsp. fraseri Brenner et al. 1989 strain Dallas 1E, sg 5, as templates.

### **Compatibles real-time PCR equipment**

Vitassay qPCR Primary Atypical Pneumonia has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTIite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>

I: For Rotor-Gene® Q thermocycler the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

### **Limitations**

- The test must be used for RUO purposes (Exclusive Use in Research), without any medical objective.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positives due to cross-contamination with different pathogens, either by samples containing high concentrations of target template DNA or by carryover contamination from PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cyclers
	<b>Eppendorf</b>
	Mastercycler™ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTIite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## Bibliography/Bibliografía

1. Lanks CW, Musani AI, Hsia DW. (2019) Community-acquired Pneumonia and Hospital-acquired Pneumonia. *Med Clin North Am.* 103(3):487-501.
2. Cillóniz C, Torres A, Niederman M, et al. (2016). Community-acquired pneumonia related to intracellular pathogens. *Intensive Care Med.* 42(9):1374-1386.
3. Centers for Disease Control and Prevention. Available from <https://www.cdc.gov/pneumonia/> & <https://www.cdc.gov/pneumonia/atypical/> Accessed Jun 2020
4. World Health Organization. Legionellosis. Available from <https://www.who.int/news-room/fact-sheets/detail/legionellosis> Accessed Jun 2020
5. Liu X, Shin S. (2019). Viewing Legionella pneumophila Pathogenesis through an Immunological Lens. *J Mol Biol.* 431(21):4321-4344.
6. Montagna MT, De Giglio O, Napoli C, et al. (2018). Control and prevention measures for legionellosis in hospitals: A cross-sectional survey in Italy. *Environ Res.* 166:55-60.
7. Chaudhry R, Ghosh A, Chandolia A. (2016). Pathogenesis of Mycoplasma pneumoniae: An update. *Indian J Med Microbiol.* 34(1):7-16.
8. Lanao AE, Chakraborty RK, Pearson-Shaver AL. (2020). Mycoplasma Infections. In: StatPearls. Treasure Island (FL): StatPearls Publishing; May 29.
9. Medjo B, Atanaskovic-Markovic M, Radic S, Nikolic D, Lukac M, Djukic S. (2014). Mycoplasma pneumoniae as a causative agent of community-acquired pneumonia in children: clinical features and laboratory diagnosis. *Ital J Pediatr.* 40:104. Published Dec 18.
10. Blasi F, Tarsia P, Aliberti S. (2009). Chlamydia pneumoniae. *Clin Microbiol Infect.* 15(1):29-35.
11. Blasi F, Tarsia P, Aliberti S, Cosentini R, Allegra L. (2005). Chlamydia pneumoniae and Mycoplasma pneumoniae. *Semin Respir Crit Care Med.* 26(6):617-624.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for components and reagents/ Símbolos utilizados para componentes y reactivos

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