Vitassay qPCR

Malaria 5

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Malaria 5 allows the detection and differentiation of *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and/or *Plasmodium knowlesi* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and/or *Plasmodium knowlesi* alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Malaria 5, 8x8-well strip, low profile	7041050
Vitassay qPCR Malaria 5, 8x8-well strip, high profile	7042050

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S050A/ 7042S050A	<i>P. malariae</i> , <i>P. knowlesi</i> & <i>P. ovale</i> strips low/high profile	-	4 x 8-well strip
7041S050B/ 7042S050B	<i>P. falciparum</i> + <i>P. vivax</i> strips low/high profile	-	4 x 8-well strip
7C050	Malaria 5 Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	8 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes

- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Malaria is a life-threatening disease caused by parasites that are transmitted to people through the bites of infected female *Anopheles* mosquitoes. The mosquitoes that spread malaria are found in Africa, Central and South America, parts of the Caribbean, Asia, Eastern Europe, and the South Pacific. Malaria causes high levels of morbidity and mortality in human beings worldwide. In 2018, there were an estimated 228 million cases of malaria worldwide. The estimated number of malaria deaths stood at 405 000 in 2018. There is no malaria vaccine.

Malaria is caused by parasites belonging to the *Plasmodium* genus. Five species of *Plasmodium* are known to cause disease in humans: *P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae*, and *P. knowlesi*. Two of these species – *Plasmodium falciparum* and *Plasmodium vivax* – pose the greatest threat.

P. falciparum is lethal and responsible for severe disease pathology and the majority of deaths due to malaria, especially in sub-Saharan Africa. In 2018, *P. falciparum* accounted for 99.7% of estimated malaria cases in the WHO African Region, 50% of cases in the WHO South-East Asia Region, 71% of cases in the Eastern Mediterranean and 65% in the Western Pacific. If not treated within 24 hours, *P. falciparum* malaria can progress to severe illness, often leading to death.

P. vivax is the most widespread of all the malaria species and can cause severe, even fatal infections. *P. vivax* is the most common malaria parasite causing clinical disease outside of Africa. *P. vivax* is highly prevalent across endemic Asia and South America. *P. vivax* is the predominant parasite in the WHO Region of the Americas, representing 75% of malaria cases. *P. vivax* malaria is associated with severe morbidity and with mortality.

P. ovale was shown to be two distinct species (*P. ovale curtisi* and *P. ovale wallikeri*), which only differ by a shorter latency period in *P. ovale wallikeri* and genetic sequence differences. These two sympatric organisms are difficult to distinguish since they present with the same clinical syndrome and respond to the same therapy. The two sympatric species of *P. ovale* produce dormant liver stages that can relapse months after the initial infection and are also associated with severe and life-threatening syndromes.

P. malariae is the most benign form of malaria infection with several distinct clinical features. *P. malariae* is present in most malaria-endemic settings and associated with greatly prolonged parasite carriage, severe anemia, and renal impairment. *P. malariae* is often considered to cause a chronic malaria they may last decades.

P. knowlesi is found in a limited distribution in Malaysian/Indonesian Borneo with cases reported in other southeast Asian countries. *P. knowlesi* is not unique among the nonhuman vertebrate malarial parasites that have been transmitted to humans. The disease presents like other malarial parasites with fever/chills and headache with uncommon features like nausea/vomiting, myalgia/arthralgia, upper respiratory symptoms, and jaundice. Although rare, fatal complications of *P. knowlesi* have occurred and do so with higher frequency than seen in *P. vivax* and *P. falciparum* proportionally.

Malaria symptoms usually appear within 7 to 30 days but can take up to one year to develop. In a non-immune individual, symptoms usually appear 10–15 days after the infective mosquito bite. The first symptoms – fever, headache, and chills but no serious organ disturbance – may be mild and difficult to recognize as malaria. This can be caused by all the *Plasmodium* spp. Severe malaria, on the other hand, is the most dangerous form of the disease with numerous complications such as severe anemia and multiple organ damage, including the brain (cerebral malaria), lungs, and kidneys and is mainly produced by *P. falciparum*. Nonetheless, it can also be caused by *P. vivax* or *P. knowlesi* at a much lower frequency.

Children under 5 years of age are the most vulnerable group affected by malaria; in 2018, they accounted for about two thirds of all malaria deaths worldwide. Children with severe malaria frequently develop one or more of the following symptoms: severe anemia, respiratory distress in relation to metabolic acidosis, or cerebral malaria. In adults, multi-organ failure is also frequent. Moreover, in malaria endemic areas, people may develop partial immunity, allowing asymptomatic infections to occur.

Malaria diagnostic tools have expanded beyond the conventional microscopic examination of Giemsa-stained blood films. Contemporary and innovative techniques have emerged, mainly the rapid diagnostic tests (RDT) and other molecular diagnostic methods such as PCR, qPCR and loop-mediated isothermal amplification (LAMP). Early and accurate diagnosis of malaria is of paramount importance to ensure appropriate administration of treatment, reduce disease and prevent deaths.

Principle of the test

Vitassay qPCR Malaria 5 test is based on the real-time amplification of a conserved target region of the 18S rRNA gene for *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and/or *Plasmodium knowlesi* in whole blood samples. After DNA extraction, the presence of viruses is detected by an increase in the fluorescence observed during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' nuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the

fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Malaria 5 test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip (*P. malariae*, *P. knowlesi* & *P. ovale* strips low/high profile) contains the multiplex reaction mix for the detection of *P. malariae*, *P. knowlesi* and/or *P. ovale*. The second strip (*P. falciparum* + *P. vivax* strips low/high profile) contains the multiplex reaction mix for the detection of *P. falciparum* and/or *P. vivax*. *P. malariae* and *P. falciparum* DNA targets are amplified and detected in FAM channel, *P. knowlesi* and *P. vivax* DNA targets are amplified and detected in ROX channel, the internal control (IC) DNA targets is amplified and detected in HEX, VIC or JOE channel (depending on the equipment used) and *P. ovale* DNA targets are amplified and detected in Cy5 channel.

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used. For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

- Invisorb® Spin Universal Kit (Stratec).
- MagDEA Dx SV kit, utilizando magLEAD® 6gC instrument (Precision System Science Co.).
- EZ1 DNA Blood Kit (QIAGEN), utilizando BioRobot EZ1 DSP Workstation (QIAGEN).

Positive control preparation

Reconstitute the lyophilized Malaria 5 Positive Control (red tube) in the 200 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*P. malariae* and *P. falciparum*), HEX, JOE or VIC (Internal Control (IC)), ROX (*P. knowlesi* and *P. vivax*), and Cy5 channels (*P. ovale*). If you use the Applied Biosystems 7500 Fast Real-Time PCR System or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*P. malariae* and *P. falciparum*), HEX, JOE or VIC (Internal Control), ROX (*P. knowlesi* and *P. vivax*) and Cy5 (*P. ovale*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX, HEX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

Interpretation of results for *P. malariae*, *P. knowlesi* & *P. ovale* strips:

The results of the first strip containing the multiplex reaction mixture for the detection of *P. malariae*, *P. knowlesi* and *P. ovale* are interpreted as follows:

<i>P. malariae</i>	<i>P. knowlesi</i>	<i>P. ovale</i>	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>P. malariae</i> , <i>P. knowlesi</i> and <i>P. ovale</i> Positives
-	-	-	+	-	+	<i>P. malariae</i> , <i>P. knowlesi</i> and <i>P. ovale</i> Negatives
+	-	-	+/-	-	+	<i>P. malariae</i> Positive, <i>P. knowlesi</i> and <i>P. ovale</i> Negatives
+	+	-	+/-	-	+	<i>P. malariae</i> and <i>P. knowlesi</i> Positives, and <i>P. ovale</i> Negative
+	-	+	+/-	-	+	<i>P. malariae</i> and <i>P. ovale</i> Positives, and <i>P. knowlesi</i> Negative
-	+	-	+/-	-	+	<i>P. knowlesi</i> Positive, <i>P. malariae</i> and <i>P. ovale</i> Negatives
-	+	+	+/-	-	+	<i>P. knowlesi</i> and <i>P. ovale</i> Positives, <i>P. malariae</i> Negative
-	-	+	+/-	-	+	<i>P. ovale</i> Positive, <i>P. malariae</i> and <i>P. knowlesi</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

Interpretation of results for *P. falciparum* + *P. vivax* strips:

The results of the second strip containing the multiplex reaction mixture for the detection of *P. falciparum* and *P. vivax* are interpreted as follows:

<i>P. falciparum</i>	<i>P. vivax</i>	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	<i>P. falciparum</i> and <i>P. vivax</i> Positives
-	-	+	-	+	<i>P. falciparum</i> and <i>P. vivax</i> Negatives
+	-	+/-	-	+	<i>P. falciparum</i> Positive and <i>P. vivax</i> Negative
-	+	+/-	-	+	<i>P. vivax</i> Positive and <i>P. falciparum</i> Negative
-	-	-	-	+	Experiment fail
+	+	+	+	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If amplification signal is detected in the negative control and / or no signal is present in the positive control, the experiment is considered unsuccessful.

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

To analyze sensitivity and clinical specificity, a comparative clinical trial was carried out and 226 blood samples were evaluated. This comparative study was carried out with the Vitassay qPCR Malaria 5 and FTD Malaria differentiation (Fast Track Diagnostics) molecular diagnostic tests.

Vitassay qPCR Malaria 5 detected *P. malariae* in 2 samples as well as the FTD Malaria differentiation test (Fast Track Diagnostics). Both tests detected *P. ovale* in 2 positive samples. Vitassay qPCR Malaria 5 detected *P. vivax* in 5 samples as well as the FTD Malaria differentiation test (Fast Track Diagnostics). For *P. falciparum* a total of 55 samples were positive for both tests.

Sensitivity and specificity for Vitassay qPCR Malaria 5 compared with FTD Malaria differentiation (Fast Track Diagnostics) are shown in next table:

Microorganism	SE (%)	SP (%)
<i>P. malariae</i>	>99	>99
<i>P. ovale</i>	>99	>99
<i>P. vivax</i>	>99	>99
<i>P. falciparum</i>	>99	>99

These results show the high sensitivity and specificity to detect *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and/or *Plasmodium knowlesi* using the molecular diagnostic test Vitassay qPCR Malaria 5.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of the template DNA of the different pathogens (10^7 - 10^1 copies / reaction). This assay has a limit of detection of ≥ 10 copies of DNA per reaction.

Analytical specificity

The analytical specificity for Malaria pathogens was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Cross-reactivity testing		
<i>Anaplasma marginale</i>	Chikungunya virus strain S27 Petersfield	<i>Theileria annulata</i>
<i>Aspergillus fumigatus</i>	Chikungunya virus Martinique isolate	Tick Borne Encephalitis virus (TBEV) strain Neudorfl
<i>Bartonella henselae</i> strain Houston-1	Chikungunya virus strain F24	<i>Toxoplasma gondii</i> (Type II)
<i>Borrelia hermsii</i>	Dengue 1 virus strain Hawaii	<i>Treponema phagedenis</i>
<i>Borrelia lusitanae</i>	Dengue 2 virus strain New Guinea C	<i>Trypanosoma cruzi</i>
<i>Borrelia valaisiana</i>	Dengue 3 virus strain H87	Usutu Virus
<i>Borrelia azfeli</i> strain P-Ko/1984	Dengue 4 virus strain H241	West Nile Virus strain NY99
<i>Borrelia bavariensis</i>	Japanese encephalitis	West Nile Virus Heja
<i>Borrelia bisetti</i>	Japanese Encephalitis Nakayama	West Nile Virus Ug37
<i>Borrelia burgdorferi</i> sensu stricto strain IRS	<i>Leptospira</i>	Yellow Fever Virus strain 17D
<i>Borrelia burgdorferi</i> sensu stricto strain B31	<i>Pneumocystis jirovecii</i>	Yellow Fever Virus French Neurotropic
<i>Borrelia garinii</i>	Chikungunya virus WHO IS (R91064)	Zika Virus strain 11474/16 (French Polynesia)
<i>Borrelia japonica</i>	<i>Rickettsia conorii</i> . strain Moroccan	Zika Virus strain 11468/16(French Polynesia)
<i>Borrelia miyamotoi</i>	Rift Valley Fever virus AR21229	Zika Virus (African)
<i>Borrelia spielmanii</i>	Rift Valley Fever virus MP12	Zika Virus strain PF13/251013-18 (Asian)
<i>Coxiella burnetii</i> strain Nine Mile Q	St Louis Encephalitis virus	Zika Virus strain FB-GWUH-2016

Analytical reactivity

Vitassay qPCR Malaria 5 has been evaluated with the following strains *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and *Plasmodium knowlesi*, obtaining a positive result for all of them.

Compatibles real-time PCR equipment

Vitassay qPCR Malaria 5 has been validated on the following equipments:

- Cobas z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^I
- SmartCycler® (Cepheid) ^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with whole blood samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Centers for Disease Control and Prevention. Available from: <https://wwwnc.cdc.gov/travel/diseases/malaria> Accessed Jul 2020
2. World Health Organization. Malaria. <https://www.who.int/news-room/fact-sheets/detail/malaria> Accessed Jul 2020
3. Milner DA, Jr. (2018). Malaria Pathogenesis. Cold Spring Harbor perspectives in medicine. 8(1), a025569.
4. Garrido-Cardenas JA, González-Cerón L, Manzano-Agugliaro F, Mesa-Valle C. (2019). Plasmodium genomics: an approach for learning about and ending human malaria. Parasitol Res. 118(1):1-27.
5. Amir A, Cheong FW, De Silva JR, Lau YL. (2018). Diagnostic tools in childhood malaria. Parasit Vectors. 11(1):53.
6. Lover AA, Baird JK, Gosling R, Price RN. (2018). Malaria Elimination: Time to Target All Species. Am J Trop Med Hyg. 99(1):17-23.
7. Meibalan E, Marti M. (2017). Biology of Malaria Transmission. Cold Spring Harb Perspect Med. 7(3):a025452.
8. Menkin-Smith L, Winders WT. (2019). Malaria (Plasmodium Vivax). In: StatPearls. Treasure Island (FL): StatPearls Publishing; May 4, 2019.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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