Vitassay qPCR

MAYV

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Intended use

Vitassay qPCR MAYV allows the qualitative detection of Mayaro virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Mayaro virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR MAYV 4 x 8-well strip, low profile	7041053
Vitassay qPCR MAYV 4 x 8-well strip, high profile	7042053

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S053/ 7042S053	MAYV strips high/low profile	-	4x8-well strip
7C053	MAYV Positive Control	Red	1 vial
7001A	PCR grade water	White	1 vial x 1 mL
7002B	Resuspension Buffer	Green	1 vial x 1,8 mL
7003N	Negative control	Yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Mayaro virus (MAYV) is an enveloped single-stranded RNA virus which belongs to the genus *Alphavirus*, the family *Togaviridae*, and was discovered incidentally in blood samples of rural workers in Trinidad and Tobago in 1954. The presence of MAYV has been mainly reported in Latin America, with Brazil having the highest number of reported cases of Mayaro fever. MAYV is also likely to circulate in the Caribbean.

The genome of MAYV is a single-stranded RNA with positive polarity and a length of approximately 11.5 kilobases (kb). Its genomic structure forms two open reading frames, one coding for structural proteins and the other coding for non-structural proteins. The latter are expressed as a polyprotein that is cleaved both during and after translation into four non-structural proteins (nsP1, nsP2, nsP3, and nsP4), whereas the structural genes generate the six structural proteins (C, E1, E2, E3, 6K, and transframe). The viral particles are icosahedral, comprising a closed nucleocapsid in a compact envelope. Phylogenetic studies have classified MAYV strains into three major genotypes, genotype D (causing most cases, widely dispersed in South Americans), L (limited in Brazil) and N (newly described in Peru).

Human infections with MAYV result in an acute febrile illness that is difficult to differentiate from dengue or chikungunya. The acute phase is characterized by a short, transient viremia (3–4 days), followed by an incubation period of 7–12 days, in which the systemic symptoms become evident. In most cases, MAYV disease is nonspecific, mild, and self-limited. The main symptoms include chills, fever, gastrointestinal manifestations, eye pain, myalgia, and arthralgia. Acute symptoms generally last three to five days in most patients with Mayaro disease. Arthralgia can last for months to years. Arthralgia and myalgia account for 50–89% and 75% of infected patients, respectively. Severe complications can occur due to MAYV infection, among which are chronic polyarthritis, myocarditis, hemorrhagic and neurological manifestations and even death.

Mayaro virus is primarily transmitted by tree dwelling *Haemagogus* species mosquitoes but MAYV has been isolated from other genera of mosquitoes, including *Culex*, *Mansonia*, *Aedes*, *Psorophora*, and *Sabethes*. MAYV transmission occurs when female mosquitoes mainly *Haemagogus janthinomys*, acquire blood from an infected host with high concentrations of MAYV in its blood. Like other alphaviruses, MAYV can infect, replicate, and disseminate in both vertebrate and invertebrate hosts. After a bite from an infected vector, MAYV spreads via the blood vessels into the body of the susceptible host. The virus replicates in white blood cells and spreads to bones, muscles, and joints via the main sites of replication (spleen).

Accurate diagnosis of MAYV infection in human communities is rendered difficult by the existence of coinfections and due to the high prevalence of asymptomatic cases. The diagnosis is often confirmed by serological testing, but the use of serological methods

can be affected by cross-reactivity and the window period. Virus is detectable by culture or reverse-transcription PCR (RT-PCR) within the first 3-5 days of illness.

Principle of the test

Vitassay qPCR MAYV is based on the real-time amplification of specific conserved fragment of the *nsp1* gene encoded by the Mayaro virus genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of Mayaro virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity that uses two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to extend the primers, the probe is hydrolyzed by its 5'-3' exonuclease activity causing the spatial separation of the fluorophore and the quencher. The resulting increase in fluorescent signal is proportional to the amount of amplified product in the sample and is detected by a real-time PCR equipment.

Vitassay qPCR MAYV is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do

not return samples, equipment and reagents to the area in which the previous step was performed.

- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national sanitary waste legislation. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use.

Positive control preparation

Reconstitute the lyophilized MAYV Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative (yellow tube) and positive (red tube) controls and add them to the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Mayaro virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for Mayaro virus (FAM), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for Mayaro virus (FAM), which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Mayaro Virus	Internal Control	Negative Control	Positive Control	Interpretation
+	+/-	-	+	Mayaro Virus Positive
-	+	-	+	Mayaro Virus Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction must be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Mayaro virus template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Mayaro virus was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Virus Mayaro		
Chikungunya Virus strain S27 Petersfield	Yellow fever virus strain 17D	St Louis Encephalitis virus strain 17D
Dengue 2 Virus strain New Guinea C	Zika Virus strain MR 766	Dengue 1 Virus strain Hawaii
West Nile Virus strain H160/99	Dengue 3 Virus strain H87	Dengue 4 Virus strain H241
West Nile Virus Heja	West Nile Virus Ug37	

Analytical reactivity

The reactivity of Vitassay qPCR MAYV was confirmed by the real-time amplification using Mayaro virus as template, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR MAYV has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Mayaro virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Mayaro virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Sun J, Wu D. Mayaro virus, a regional or global threat? *Travel Med Infect Dis.* 2019 Jul 25:101462.
2. Diagne CT, Bengue M, Choumet V, Hamel R, Pompon J, Missé D. Mayaro Virus Pathogenesis and Transmission Mechanisms. *Pathogens.* 2020 Sep 8;9(9):738.
3. Waggoner JJ, Rojas A, Mohamed-Hadley A, de Guillén YA, Pinsky BA. Real-time RT-PCR for Mayaro virus detection in plasma and urine. *J Clin Virol.* 2018 Jan; 98:1-4.
4. Acosta-Ampudia Y, Monsalve DM, Rodríguez Y, Pacheco Y, Anaya JM, Ramírez-Santana C. Mayaro: an emerging viral threat? *Emerg Microbes Infect.* 2018 Sep 26;7(1):163.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.

Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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