Vitassay qPCR

Carbapenemase-generating Enterobacteriaceae

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Intended use

Vitassay qPCR Carbapenemase-generating Enterobacteriaceae allows the qualitative detection and differentiation of the main carbapenemase-encoding genes (*NDM*, *VIM*, *OXA*, *KPC* and/or *IMP*) by real-time PCR directly from rectal swabs and from bacterial isolates of clinical samples. The product is intended for use in the diagnosis of carbapenem-resistant *Enterobacteriaceae* infections alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR Carbapenemase-generating Enterobacteriaceae 8x8-well strip, low profile	7041057
Vitassay qPCR Carbapenemase-generating Enterobacteriaceae 8x8-well strip, high profile	7042057
Vitassay qPCR Carbapenemase-generating Enterobacteriaceae 16x8-well strip, low profile	7081057
Vitassay qPCR Carbapenemase-generating Enterobacteriaceae 16x8-well strip, high profile	7082057

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S057A/ 7042S057A	Carbapenemase-generating Enterobacteriaceae I strips low/high profile	-	4/8 x 8-well strip
7041S057B/ 7042S057B	Carbapenemase-generating Enterobacteriaceae II strips low/high profile	-	4/8 x 8-well strip
7C057	Carbapenemase-generating Enterobacteriaceae Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	8/16 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

The emergence and spread of drug-resistant pathogens that have acquired new resistance mechanisms, leading to antimicrobial resistance, continues to threaten our ability to treat common infections. Carbapenems are a potent class of β -lactam antibiotics that are often used as “last-line agents” or “antibiotics of last resort” when infected patients become severely ill or are suspected of harboring resistant bacteria. However, the emergence of carbapenem-resistant Enterobacteriaceae (CRE) has become a major public health concern. CRE are an emerging cause of Health care-associated infections which have a significant impact on morbidity, mortality and quality of life and represent an economic burden at the societal level. Since the identification of carbapenem-resistant Enterobacteriaceae in the 1990s, CRE has spread worldwide.

The current CRE definition by CDC is Enterobacteriaceae that test resistant to at least one of the carbapenem antibiotics (ertapenem, meropenem, doripenem, or imipenem) or produce a carbapenemase (an enzyme that can make them resistant to carbapenem antibiotics). Some strains are innately resistant to carbapenems, while others contain mobile genetic elements (plasmids, transposons) that result in the production of carbapenemase enzymes (carbapenemases), which break down most beta-lactam antibiotics, including carbapenems.

Enterobacteriaceae comprise a large family of gram-negative bacteria that commonly colonize humans and sometimes cause disease. They are facultative anaerobes, frequently residing in the intestinal tract, as part of the normal gut microbiota. Enterobacteriaceae are a common cause of infections in both community and healthcare settings. Many different types of Enterobacteriaceae can develop resistance, including *Klebsiella pneumoniae* and *Escherichia coli* (*E. coli*). The widespread use of carbapenems for empiric and directed treatment of severe infections has resulted in the emergence of carbapenem-hydrolyzing β -lactamases, also known as carbapenemases, as the most well-recognized mechanism of resistance to carbapenems. Many of the genes coding for carbapenemases are located in plasmids which can be transferred between bacteria and hence this type of resistance is more likely to spread. The

carbapenemases, according to Ambler's classification, belong to different classes of β -lactamases: class A, serine β -lactamases (*Klebsiella pneumoniae* carbapenemases (KPC)); class B, also named metallo- β -lactamase, (Verona integron-encoded metallo- β -lactamase (VIM), New Delhi metallo- β -lactamase (NDM), and metallo- β -lactamase active on imipenem (IMP)); class D, oxacillinases (OXAs; which is represented in Enterobacteria by OXA-48-like enzymes).

Although resistance patterns may be geographically distinct, the threat remains that any new resistance mechanism may spread rapidly to other areas. KPC serine carbapenemases are present more frequently in the Americas, the Mediterranean countries, and China, whereas NDM metallo- β -lactamases are more prevalent in the Indian subcontinent and Eastern Europe. OXA-48 was first described among clinical isolates of *K. pneumoniae* from Turkey and has subsequently disseminated in Europe and the Mediterranean region, but it is found less frequently in the Americas. IMP metallo- β -lactamase was first identified in Japan and have remained continuing contributors to carbapenem resistance in other regions of Southeast Asia. Although they have not spread extensively throughout the rest of the world, IMP carbapenemases are being reported more frequently in Middle Eastern countries. VIM metallo-carbapenemases were first identified in Italy and France and initially spread rapidly throughout southern Europe, but today, VIM metallo- β -lactamases are found globally, mainly in *K. pneumoniae* and *E. cloacae* complex.

CRE can cause infections in almost any body part, including bloodstream infections, ventilator-associated pneumonia, intra-abdominal abscesses, and infections that involve the urinary tract. Carbapenemase-producing Enterobacteriaceae carriers are often asymptomatic. CRE are usually spread person to person through contact with infected or colonized people, particularly contact with wounds or stool. In healthcare settings, this contact can occur via the hands of healthcare personnel, or through contaminated medical equipment. Treatment with long-term β -lactam antibiotics, fluoroquinolones, mechanical ventilation, parenteral nutrition, permanent urinary catheters, organ transplantation and ICU admission pose independent risk factors. The risk is also increased in those patients with diabetes and immunocompromised. Additionally, sink drains and toilets are increasingly recognized as an environmental reservoir and CRE transmission source.

Clinical laboratories can perform phenotypic tests for carbapenemase production (Carba NP, Modified Carbapenem Inactivation Method (mCIM), Modified Hodge Test (MHT)) or molecular assays for the presence of a carbapenemase gene such as polymerase chain reaction (PCR).

NOTE: We will refer to carbapenemase-encoding genes as follows: $bla_{NDM} = NDM$, $bla_{VIM} = VIM$, $bla_{KPC} = KPC$, $bla_{OXA-48 \text{ like}} = OXA$ and $bla_{IMP} = IMP$ (IMP-1 cluster + IMP-8 cluster)*

*IMP-1 cluster includes: IMP-1, IMP-4, IMP-3, IMP-6, IMP-10, IMP-25, IMP-30, IMP-34, IMP-40, IMP-42, IMP-52, IMP-55. IMP-8 cluster includes: IMP-2, IMP-8, IMP-16, IMP-19, IMP-20, IMP-22, IMP-24.

Principle of the test

Vitassay qPCR Carbapenemase-generating Enterobacteriaceae is based on the real-time amplification of a conserved region of the *NDM* gene, *VIM* gene, *OXA-48* and *OXA-48-like* genes, *KPC* gene and *IMP* gene. After DNA extraction, the presence of carbapenemases is detected by an increase in the fluorescence observed during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its 5'-3' exonuclease activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Carbapenemase-generating Enterobacteriaceae is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip (Carbapenemase-generating Enterobacteriaceae I strips low/high profile) contains the multiplex reaction mix for the detection of *NDM* and/or *VIM* carbapenemase-encoding genes and contains the internal control (IC). The second strip (Carbapenemase-generating Enterobacteriaceae II strips low/high profile) contains the multiplex reaction mix for the detection of *OXA*, *KPC* and/or *IMP* carbapenemase-encoding genes. *NDM* and *OXA* are amplified and detected in FAM channel, *VIM* and *KPC* are amplified and detected in ROX channel, the internal control (IC) is amplified and detected in HEX, VIC or JOE channel (depending on the equipment used) and *IMP* is amplified and detected in Cy5 channel.

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.

- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Make sure to use a well for determining *NDM* and *VIM* genes and another well for determining *OXA*, *KPC* and *IMP* genes. Be careful not to mix them throughout the process.

Procedures

Specimen collection, processing, and DNA extraction

Patient samples must be collected, transport and storage according to appropriate laboratory guidelines. For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use. For rectal swabs and bacterial isolates from clinical samples, the assay has been validated with the following extraction kits:

GXT NA Extraction kit in GenoXtract® system (Hain).

MagDEA Dx SV Kit, using the magLEAD® 6gC or 12gC instrument (Precision System Science Co.).

Invisorb® Spin Universal Kit (Invitex).

Positive control preparation

Reconstitute the lyophilized Carbapenemase-generating Enterobacteriaceae Positive Control (red tube) with 200 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative (yellow tube) or positive (red tube) controls and add them to the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*NDM* and *OXA* genes), ROX (*VIM* and *KPC* genes), Cy5 (*IMP* gene) and HEX, JOE, or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied

Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The results' analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve ($Ct \leq 40$) in FAM, ROX, Cy5 and HEX, VIC, or JOE channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence ($Ct \geq 40$ or no signal) in FAM, ROX, and Cy5 channels, which validates the reaction.

The Internal Control (IC) should show an amplification signal ($Ct \leq 40$) in Positive and Negative Controls wells.

The experiment seems to be failed if there is amplification signal in negative control or signal' absence in the positive control for any channel. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Carbapenemase-generating Enterobacteriaceae I			Carbapenemase-generating Enterobacteriaceae II			Interpretation
NDM (FAM)	VIM (ROX)	Internal Control (HEX)	OXA (FAM)	KPC (ROX)	IMP (Cy5)	
+	-	+/-*	-	-	-	NDM-expressing Enterobacteriaceae
-	+	+/-*	-	-	-	VIM-expressing Enterobacteriaceae
-	-	+/-*	+	-	-	OXA-expressing Enterobacteriaceae
-	-	+/-*	-	+	-	KPC-expressing Enterobacteriaceae
-	-	+/-*	-	-	+	IMP-expressing Enterobacteriaceae
-	-	+#	-	-	-	Carbapenemase-generating Enterobacteriaceae not detected [#]
-	-	- [#]	-	-	-	Invalid [#]

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct ≥40 or no signal)

* Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

[#] In the case of negative Carbapenemase-generating Enterobacteriaceae target genes detection, IC must show an amplification signal with Ct less than 35. If there is a signal' absence or Ct value ≥ 35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

Table shows only the most representative results that can be expected with the Vitassay qPCR Carbapenemase- generating Enterobacteriaceae kit. One patient sample can simultaneously contain several carbapenemase-encoding genes.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The clinical sensitivity and specificity of Vitassay qPCR Carbapenemase-generating Enterobacteriaceae was evaluated using directly rectal swabs and from bacterial isolates from clinical samples.

In the two studies carried out in the Reference Laboratory of Catalonia, the extraction was carried out with the GXT NA extraction kit using the GenoXtract® system (Hain), following the instructions for use and the thermal cycler used was the Applied BioSystem 7500 RT PCR System.

In both studies the results obtained with Vitassay qPCR Carbapenemase-generating Enterobacteriaceae were compared with phenotypic characterization followed by the LightMix® modular carbapenemase assay. The results obtained were the following:

Reference Laboratory of Catalonia (bacterial isolates from clinical samples)								
Target	TP	TN	FP	FN	SE	SP	PPV	NPV
<u>Global:</u> Carbapenemase-generating <i>Enterobacteriaceae</i>	80	17	0	0	1 (0.95-1)	1 (0.80-1)	1 (0.95-1)	1 (0.80-1)
OXA gene	33	64	0	0	1 (0.89-1)	1 (0.94-1)	1 (0.89-1)	1 (0.94-1)
VIM gene	23	74	0	0	1 (0.85-1)	1 (0.95-1)	1 (0.85-1)	1 (0.95-1)
NDM gene	19	78	0	0	1 (0.82-1)	1 (0.95-1)	1 (0.82-1)	1 (0.95-1)
KPC gene	2	95	0	0	1 (0.16-1)	1 (0.96-1)	1 (0.16-1)	1 (0.96-1)
IMP gene	7	90	0	0	1 (0.59-1)	1 (0.96-1)	1 (0.59-1)	1 (0.96-1)

TP = True Positive, TN = True Negative, FP = False Positive, FN = False Negative, SE = Sensibility, SP = Specificity, PPV = Positive Predictive Value, NPV = Negative Predictive Value.

Reference Laboratory of Catalonia (rectal swabs)								
Target	TP	TN	FP	FN	SE	SP	PPV	NPV
<u>Global:</u> Carbapenemase- generating <i>Enterobacteriaceae</i>	34	16	0	0	1 (0.91-1)	1 (0.79-1)	1 (0.91-1)	1 (0.79-1)
<i>OXA</i> gene	13	37	0	0	1 (0.75-1)	1 (0.90-1)	1 (0.75-1)	1 (0.90-1)
<i>VIM</i> gene	11	39	0	0	1 (0.71-1)	1 (0.91-1)	1 (0.71-1)	1 (0.91-1)
<i>NDM</i> gene	10	40	0	0	1 (0.69-1)	1 (0.91-1)	1 (0.69-1)	1 (0.91-1)
<i>KPC</i> gene	5	45	0	0	1 (0.47-1)	1 (0.92-1)	1 (0.47-1)	1 (0.92-1)
<i>IMP</i> gene	0	50	0	0	n.a*	1 (0.91-1)	n.a*	1 (0.91-1)

* Due to all analyzed samples were negative for *IMP* gene, the analytical sensitivity of the test could not be performed.

These results show high agreement to detect Carbapenemase-producing *Enterobacteriaceae* using the molecular diagnostic test Vitassay qPCR Carbapenemase-generating *Enterobacteriaceae*.

Analytical sensitivity

The analytical sensitivity was determined from serial dilutions (1:10) of standards of the *NDM*, *VIM*, *OXA*, *KPC* and *IMP* genes (10^7 - 10^1 copies/reaction). This assay has a detection limit of 10 DNA copies per reaction for the *NDM*, *KPC* and *IMP* genes, 50 DNA copies per reaction for *OXA* and 100 DNA copies per reaction for *VIM*.

Analytical specificity

The analytical specificity for the Carbapenemase-producing *Enterobacteriaceae* detection was confirmed by testing a panel composed of the following microorganisms, where no cross-reactivity was seen between any of the species except with the target pathogens detected by each assay:

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Analytical reactivity

The reactivity of Vitassay qPCR Carbapenemase-generating Enterobacteriaceae was evaluated against TEM-1 (non-ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL) and NDM-1 producing *Enterobacter cloacae* isolate; NDM-7 producing *Enterobacter cloacae-complex* isolate; VIM-1 producing *Citrobacter braakii* isolate; KPC-3 and VIM-4 producing *Citrobacter freundii-complex* isolate; SHV-12 (ESBL), CTX-M-9 (ESBL) and OXA-48 producing *Enterobacter cloacae* isolate; OXA-244 producing *Escherichia coli* isolate; SHV-1 (non-ESBL), KPC-3, and OXA-48 producing *Klebsiella pneumoniae* isolate; OXA-48 producing *Serratia marcescens* isolate; TEM-1 (non-ESBL), SHV-1 (non-ESBL), CTX-M-2 (ESBL), and KPC-2 producing *Klebsiella pneumoniae* isolate and TEM-1 (non-ESBL) and IMP-1 producing *Escherichia coli* isolate, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR Carbapenemase-generating Enterobacteriaceae has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- 7500 RT PCR System (Applied Biosystems)
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^I
- LightCycler 480 II (Roche)

I: For Rotor-Gene® Q thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a preliminary diagnosis of carbapenem-resistant Enterobacteriaceae infection. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with rectal swabs and bacterial isolates of clinical samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; DNA from clinical specimens must be extracted properly.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- There is a possibility of false positive results due to cross-contamination by Carbapenemase-generating Enterobacteriaceae Positive Control, samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA isolation); c) DNA degradation during specimen shipping, storage, and/or preparation; d) pathogen load in the specimen below the limit of detection for the assay; e) the presence of Real Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent drug-resistant infections or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures.
- DNA detection may not indicate the presence of viable and/or infectious pathogens or that those pathogens are the causative agents for clinical symptoms.
- Negative results do not preclude Carbapenemase-producing Enterobacteriaceae infection and should not be the sole basis of a patient treatment/management decision. Optimum specimen types for the Carbapenemase-producing Enterobacteriaceae identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple

specimens from the same patient at different time points, which may increase the probability of detecting the pathogens.

- If the patient clinical data, laboratory tests and epidemiological studies suggest possible Carbapenemase-producing Enterobacteriaceae infection, and other drug-resistant related illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.
- A bioinformatics tool (NCBI BLAST) analysis of the primer alignments showed that this kit detects the following members for each family of the target carbapenemase genes:

Family	Family Members
NDM	1, 2, 3, 4, 5, 6, 7, 9, 11, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29
VIM	1, 2, 3, 4, 6, 8, 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 23, 24, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 39, 40, 41, 42, 43, 44, 45, 46, 48, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 62, 63, 64, 65, 66, 67, 68, 70
OXA	48, 162, 163, 181, 199, 204, 232, 244, 245, 247, 252, 370, 405, 416, 438, 439, 484, 505, 514, 515, 517, 519, 538, 546, 547, 566, 567, 788, 793, 833, 894
KPC	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 51, 52, 54, 56
IMP	1, 4, 5, 6, 8, 9, 10, 11, 12, 13, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 28, 29, 30, 33, 34, 37, 38, 39, 40, 41, 42, 44, 45, 49, 52, 53, 55, 56, 58, 59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Bio Molecular Systems	Exicycler™ 96
Mic Real Time PCR Cyclor	Bio-Rad
Cepheid	CFX96™ Deep Well Real-Time PCR Detection
SmartCycler®	iCycler iQ™ Real-Time PCR
Precision System Science Co., Ltd.	iCycler iQ™5 Real-Time PCR
geneLEAD VIII System	MyiQ™ Real-Time PCR Detection System
Qiagen	MyiQ™ 2Real-Time PCR Detection System
Rotor-Gene® Q	Bio Molecular Systems
Roche	Mic Real Time PCR Cyclor
LightCycler @480 Real-Time PCR System	Cepheid
LightCycler @96 Real-Time PCR System	SmartCycler®
Cobas z480 Analyzer	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cyclor*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Precision System Science Co., Ltd.
	geneLEAD VIII System
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Color Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquaring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliografía


1. World Health Organization. Antimicrobial resistance. Available from: <https://www.who.int/news-room/fact-sheets/detail/antimicrobial-resistance> Accessed May 2021.
2. Chen HY, Jean SS, Lee YL, Lu MC, Ko WC, Liu PY, Hsueh PR. (2021) Carbapenem-Resistant Enterobacterales in Long-Term Care Facilities: A Global and Narrative Review. *Front Cell Infect Microbiol.* Apr 23;11:601968.
3. Hughes LD, Aljawadi A, Pillai A. (2019) An overview of carbapenemase producing enterobacteriaceae (CPE) in trauma and orthopaedics. *J Orthop.* Jul 2;16(6):455-458.
4. Ramsamy Y, Misana KP, Allam M, Amoako DG, Abia ALK, Ismail A, Singh R, Kisten T, Han KS, Muckart DJJ, Hardcastle T, Suleman M, Essack SY. (2020) Genomic Analysis of Carbapenemase-Producing Extensively Drug-Resistant *Klebsiella pneumoniae* Isolates Reveals the Horizontal Spread of p18-43_01 Plasmid Encoding *bla*_{NDM-1} in South Africa. *Microorganisms.* Jan 17;8(1):137.
5. World Health Organization. Guidelines for the prevention and control of carbapenem-resistant Enterobacteriaceae, *Acinetobacter baumannii* and *Pseudomonas aeruginosa* in health care facilities. Available from: <https://www.who.int/infection-prevention/publications/guidelines-cre/en/> Accessed May21.
6. Centers of Disease Control and Prevention (CDC). CRE Technical Information. Available from: <https://www.cdc.gov/hai/organisms/cre/technical-info.html> Accessed May21
7. Centers of Disease Control and Prevention (CDC). Clinicians: Information about CRE. Available from: <https://www.cdc.gov/hai/organisms/cre/cre-clinicians.html> Accessed May21
8. Centers of Disease Control and Prevention (CDC). Patients: Information about CRE. Available from: <https://www.cdc.gov/hai/organisms/cre/cre-patients.html> Accessed May21
9. Fontana C, Angeletti S, Mirandola W, Cella E, Alessia L, Zehender G, Favaro M, Leoni D, Rose DD, Gherardi G, Florio L, Salemi M, Andreoni M, Sarmati L, Ciccozzi M. (2020) Whole genome sequencing of carbapenem-resistant *Klebsiella pneumoniae*: evolutionary analysis for outbreak investigation. *Future Microbiol.* Feb;15:203-212.
10. Logan LK, Weinstein RA. (2017) The Epidemiology of Carbapenem-Resistant Enterobacteriaceae: The Impact and Evolution of a Global Menace. *J Infect Dis.* Feb 15;215(suppl_1):S28-S36.
11. Bush K, Bradford PA. (2020) Epidemiology of β -Lactamase-Producing Pathogens. *Clin Microbiol Rev.* Feb 26;33(2):e00047-19.

12. Centers of Disease Control and Prevention (CDC). Carbapenemase Producing Carbapenem-Resistant Enterobacteriaceae (CP-CRE) 2018 Case Definition. Available from: <https://wwwn.cdc.gov/nndss/conditions/carbapenemase-producing-carbapenem-resistant-enterobacteriaceae/case-definition/2018/>
Accessed May21

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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