# Vitassay qPCR

## Sexually Transmitted Infections

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\*The Notified Body 0318 is only involved in the identification of *Chlamydia Trachomatis* included in 7041S035A/7042S035A.

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## Referencias

Vitassay qPCR Sexually Transmitted Infections 8x8-well strip, low profile	7041035
Vitassay qPCR Sexually Transmitted Infections 8x8-well strip, high profile	7042035
Vitassay qPCR Sexually Transmitted Infections 16x8-well strip, low profile	7081035
Vitassay qPCR Sexually Transmitted Infections 16x8-well strip, high profile	7082035

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## Intended use

Vitassay qPCR Sexually Transmitted Infections allows the detection and differentiation of *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and/or *Mycoplasma hominis* by real-time PCR in urogenital, endocervical and urine samples from patients with signs and symptoms of sexually transmitted diseases (STDs). The product is intended for use in the diagnosis of, Sexually Transmitted Infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Sexually Transmitted Infections 8x8-well strip, low profile	7041035
Vitassay qPCR Sexually Transmitted Infections 8x8-well strip,high profile	7042035
Vitassay qPCR Sexually Transmitted Infections 16x8-well strip, low profile	7081035
Vitassay qPCR Sexually Transmitted Infections 16x8-well strip,high profile	7082035

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S035A/ 7042S035A	<i>Chlamydia trachomatis</i> <i>Mycoplasma genitalium</i> & <i>Neisseria gonorrhoeae</i> strip low/high profile	-	4/8 tiras de 8 pocillos
7041S035B/ 7042S035B	<i>Trichomonas vaginalis</i> , <i>Ureaplasma urealyticum</i> , <i>Ureaplasma parvum</i> & <i>Mycoplasma hominis</i> strip low/high profile	-	4/8 tiras de 8 pocillos
7C035	Sexually Transmitted Infections Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	8/16 x 8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.

- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

### **Additional equipment and material required**

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

### **Summary**

Sexually transmitted diseases (STDs) are very common and are transmitted from one person to another through intimate physical contact and sexual activity that includes vaginal, oral, and anal sex. STDs can mostly be prevented by not having sex. STDs do not always cause symptoms, so it is possible to have an infection and not know it.

Chlamydia, a sexually transmitted infection is caused by the bacterium *Chlamydia trachomatis*. The majority of genital chlamydial infections in both males and females are asymptomatic. When symptoms do occur, lower urogenital tract infection can manifest as cervicitis in females and urethritis in males and females. Whether symptomatic or asymptomatic, untreated chlamydia can ascend to the upper genital tract. In males, this can cause epididymitis, which is not thought to be an important cause of long-term sequelae. However, in females, upper tract infection can result in pelvic inflammatory disease (PID), a spectrum of clinical disorders involving infection and inflammation of the uterus, fallopian tubes, ovaries, or adjacent peritoneum.

*Mycoplasma genitalium* is a sexually transmitted pathogen that causes no gonococcal urethritis, and recent evidence indicates that it increases the risk for cervicitis, pelvic inflammatory disease, preterm delivery, and spontaneous abortion.

*Neisseria gonorrhoeae* is a sexually transmitted pathogen that causes infections at the anatomic site of exposure (e.g., urethra, cervix, pharynx, and rectum). Urethral infections often cause discharge and pain with urination; cervical, pharyngeal, and rectal infections are frequently asymptomatic. In women, gonorrhoea is a major cause of pelvic inflammatory disease, which can lead to severe reproductive health complications (e.g., infertility, ectopic pregnancy, and chronic pelvic pain). Infrequently, disseminated gonococcal infection can result in localized septic arthritis, endocarditis, and meningitis.

Gonococcal infections can increase the risk for sexual transmission of human immunodeficiency virus (HIV).

*Trichomonas vaginalis* infection is one of the most common sexually transmitted diseases (STDs) in the world. In women it can be found in the vagina and in the urethra, while in men it can be found in the urethra, the prostate, and the epididymis. *Trichomonas vaginalis* infection has been associated with vaginitis, cervicitis and urethritis, premature rupture of membranes and premature delivery in pregnant women. *Trichomonas vaginalis* infection has also been associated with an increased risk of HIV acquisition and transmission in women.

*Ureaplasma parvum* and *Ureaplasma urealyticum* are part of the genital flora of men and women and are present in almost 70% of sexually active population. These bacteria cause inflammation and lead to chorioamnionitis, preterm deliveries, and premature rupture of membranes.

*Mycoplasma hominis* colonizes the lower urogenital tract and is associated with urogenital infections, particularly bacterial vaginosis, and non-gonococcal urethritis. In adults, bacteremia, septic arthritis, osteitis, endocarditis, mediastinitis, brain abscesses and respiratory infections have been described. It is also involved in extra genital infections, such as postpartum or post-abortion fever, in post-cesarean wound infections or after a hysterectomy. In neonates, it can cause meningitis, brain abscesses and eye infections.

### **Principle of the test**

Vitassay qPCR Sexually Transmitted Infections test is based on the real-time amplification of specific conserved fragments of the *T. vaginalis*-specific 2-kb repeated sequence gene (*Trichomonas vaginalis*), *ureasa* gene (*Ureaplasma urealyticum* and *Ureaplasma parvum*), *yidC* gene (*Mycoplasma hominis*), *porA* and *Opa* genes (*Neisseria gonorrhoeae*), a region within ORF2 of the chlamydial plasmid (*Chlamydia trachomatis*) and *MgPa adhesin* gene (*Mycoplasma genitalium*). It is detected by an increment in the fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' nuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Sexually Transmitted Infections test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip

(*C. trachomatis*, *M. genitalium* & *N. gonorrhoeae* ref #7041S035/7042S035) contains the multiplex reaction mix for the detection of *Chlamydia trachomatis*, *Mycoplasma genitalium* and *Neisseria gonorrhoeae* as well as the internal control (IC). The second strip (*T. vaginalis*, *U. urealyticum*, *U. parvum* & *M. hominis* ref #7041S035B/7042S035B) contains the multiplex reaction mix for the detection of *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*. *Chlamydia trachomatis* and *Trichomonas vaginalis* DNA targets are amplified and detected in FAM channel, *Neisseria gonorrhoeae* and *Ureaplasma parvum* DNA targets are amplified and detected in ROX channel, *Mycoplasma genitalium* and *Ureaplasma urealyticum* DNA targets are amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2) and *Mycoplasma hominis* and the internal control (IC) DNA targets are amplified and detected in Cy5 channel.

## Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Make sure to use a well for determining *Neisseria gonorrhoeae*, *Chlamydia trachomatis* and *Mycoplasma genitalium* and another well for determining *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*. Be careful not to mix them throughout the process.

- Vitassay qPCR Sexually Transmitted Infections has only been validated with the equipment mentioned in Section Compatibles real-time PCR equipment of these Instructions for Use. The rest of the compatible Real Time PCR instrument indicated in Annex 1 is based on bibliographic data.

## **Procedures**

### **Specimen collection, processing, and DNA extraction**

Clinical specimens (urogenital specimens, urine samples, rectal and endocervical samples) should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend using fresh samples.

For longer storage, the samples should be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenize sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).

ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

US121 NX-48 Urine/Swab DNA Kit using the Nextractor® NX-48 (Genolution Inc.).

MagDEA® Dx SV, using magLEAD® 6gC (Precision System Science)).

NucleoSpin® RNA Virus (MACHEREY-NAGEL).

NucleoMag® Pathogen (MACHEREY-NAGEL).

STARMag 96x4 Universal Cartridge Kit (Seegene, Korea)

Microlab STAR Let automatic extraction system (Hamilton, Switzerland)

Invisorb® Spin Universal Kit (Invitex).

## Positive control preparation

Reconstitute the lyophilized Sexually Transmitted Infections Positive Control (red tube) with 200 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

## Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (*Chlamydia trachomatis* y *Trichomonas vaginalis*), HEX, JOE o VIC (*Mycoplasma genitalium* and *Ureaplasma urealyticum*), ROX (*Neisseria gonorrhoeae* and *Ureaplasma parvum*) and Cy5 (Internal Control (IC) and *Mycoplasma hominis*). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none (Attached II).

## Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

### Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Chlamydia trachomatis* and *Trichomonas vaginalis*), HEX, JOE or VIC (*Mycoplasma genitalium* and *Ureaplasma urealyticum*), ROX (*Neisseria gonorrhoeae* and *Ureaplasma parvum*) and Cy5 (Internal Control and *Mycoplasma hominis*), which validates the reaction.

### Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX, HEX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The results of the first strip containing the multiplex reaction mixture for the detection of *Chlamydia trachomatis*, *Mycoplasma genitalium* and *Neisseria gonorrhoeae* are interpreted as follows:

- *Chlamydia trachomatis*. A sample is considered positive for this pathogen if it shows a signal in the FAM channel and a Ct less than 40, and the internal control shows or not a signal.
- *Mycoplasma genitalium*. A sample is considered positive for this pathogen if it shows a signal in the HEX, JOE or VIC channel and a Ct less than 40, and the internal control shows or not a signal.
- *Neisseria gonorrhoeae*. A sample is considered positive for this pathogen if it shows a signal in the ROX channel and a Ct less than 40, and the internal control shows or not a signal.
- If no amplification curve is detected above the threshold value and the internal control is present, the sample is considered negative.

The results of the second strip containing the multiplex reaction mixture for the detection of *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis* are interpreted as follows:

- *Trichomonas vaginalis*. A sample is considered positive for this pathogen if it shows a signal in the FAM channel and a Ct less than 40.
- *Ureaplasma urealyticum*. A sample is considered positive for this pathogen if it shows a signal in the HEX, JOE or VIC channel and a Ct less than 40.
- *Ureaplasma parvum*. A sample is considered positive for this pathogen if it shows a signal in the ROX channel and a Ct less than 40.
- *Mycoplasma hominis*. A sample is considered positive for this pathogen if it shows a signal in the Cy5 channel and a Ct less than 40.
- If no amplification curve is detected above the threshold value, the sample is considered negative.

If amplification signal is detected in the negative control and / or no signal is present in the positive control, the experiment is considered unsuccessful.

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction must be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

A total of 948 samples (first-void urines, rectal swabs, urethral swabs, endocervical swabs/exudate, vaginal swabs and pharyngeal swabs) were analysed by Vitassay strip qPCR Sexually Transmitted Infections and Allplex™ STI Essential Assay (Seegene).

The first Vitassay strip qPCR Sexually Transmitted Infections (*Chlamydia trachomatis*, *Mycoplasma genitalium* and *Neisseria gonorrhoeae*) detected *Chlamydia trachomatis* in 68 positive samples, *Mycoplasma genitalium* in 28 positive sample and *Neisseria gonorrhoeae* in 19 positive samples.

The second Vitassay strip qPCR Sexually Transmitted Infections (*Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*) detected *Trichomonas vaginalis* in 20 positive samples, *Ureaplasma urealyticum* in 166 positive samples, *Ureaplasma parvum* in 402 positive samples and *Mycoplasma hominis* in 183 positive samples.

46 discordant samples were obtained (samples close to the detection limit). Of the 46 discordant samples, 35 were confirmed by another molecular detection kit FTD Urethritis plus (Fast Track Diagnostics). The remaining discrepant samples cannot be re-analysed due to lack of sample.

Sensitivity, specificity, PPV and NPV values for Vitassay qPCR Sexually Transmitted Infections compared with Allplex™ STI Essential Assay (Seegene) are shown in next table:

Microorganism	SE (%)	SP (%)	PPV (%)	NPV (%)
<i>C. trachomatis</i>	97.14	99.6	95.77	99.77
<i>M. genitalium</i>	93.33	99.89	96.55	99.78
<i>N. gonorrhoeae</i>	95	99.89	95	99.79
<i>T. vaginalis</i>	100	99.68	86.96	100
<i>U. urealyticum</i>	95.95	99.23	96.51	99.1
<i>U. parvum</i>	98.77	97.6	96.87	99.06
<i>M. hominis</i>	100	99.74	98.92	100

These results show the high sensitivity and specificity to detect *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis* using the molecular diagnostic test Vitassay qPCR Sexually Transmitted Infections.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of the template DNA of the different pathogens ( $10^7$ - $10^1$  copies / reaction). This assay has a limit of detection of  $\geq 10$  copies of DNA per reaction.

### Analytical specificity

The analytical specificity for Sexually Transmitted Infections was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Sexually Transmitted Infections		
<i>Candida parapsilosis</i>	<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i>
<i>Candida tropicalis</i>	<i>Escherichia Coli</i>	<i>Serratia marcescens</i>
<i>Candida glabrata</i>	<i>Gardnerella vaginalis</i>	<i>Staphylococcus aureus</i>
<i>Candida krusei</i>	<i>Haemophilus influenzae</i> MinnA	<i>Stenotrophomonas maltophilia</i>
<i>Candida dubliniensis</i>	<i>Haemophilus ducreyi</i> class 1	<i>Streptococcus agalactiae</i>
<i>Candida albicans</i>	<i>Klebsiella oxytoca</i>	<i>Streptococcus pneumoniae</i>
<i>Aspergillus fumigatus</i>	<i>Klebsiella pneumoniae</i>	<i>Treponema pallidum</i>
<i>Bacteroides fragilis</i>	<i>Listeria Ivanovii</i>	Cytomegalovirus
<i>Acinetobacter baumannii</i>	<i>Listeria monocytogenes</i>	Herpes simplex virus 1
<i>Enterobacter cloacae</i>	<i>Listeria Innocua</i>	Herpes simplex virus 2
<i>Enterobacter aerogenes</i>	<i>Neisseria meningitidis</i>	Human Papillomavirus 16
<i>Enterococcus faecalis</i>	<i>Proteus mirabilis</i>	Human Papillomavirus 18

### Analytical reactivity

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Chlamydia trachomatis* was confirmed by the real time amplification using *Chlamydia trachomatis* (SW and LGV), Swedish strain, and Genovar F, and Serovars D, E, I, K and J, as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Mycoplasma genitalium* was confirmed by the real time amplification using *Mycoplasma genitalium* as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Neisseria gonorrhoeae* was confirmed by the real time amplification using *Neisseria gonorrhoeae* strains St 49226 and Lvl Ng PorA, as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Trichomonas vaginalis* was confirmed by the real time amplification using *Trichomonas vaginalis* as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Ureaplasma urealyticum* was confirmed by the real time amplification using *Ureaplasma urealyticum* as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Ureaplasma parvum* was confirmed by the real time amplification using *Ureaplasma parvum* as template.

The reactivity of Vitassay qPCR Sexually Transmitted Infections for *Mycoplasma hominis* was confirmed by the real time amplification using *Mycoplasma hominis* as template.

## Compatibles real-time PCR equipment

Vitassay qPCR Sexually Transmitted Infections has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of Sexually Transmitted Infections. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with urogenital specimens and endocervical specimens for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*, rectal specimens (rectal-vaginal specimens and rectal samples) for *Neisseria gonorrhoeae*, *Ureaplasma urealyticum* and *Mycoplasma hominis*; and urine samples for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positives due to cross-contamination with different pathogens, either by samples containing high concentrations of target template DNA or by carryover contamination from PCR products from previous reactions

- ON 0318 only intervenes in the conformity assessment of the test for *Chlamydia trachomatis*. The scope of the CE certification covers the detection of *Chlamydia trachomatis* from urine samples and urogenital and endocervical swabs. The detection of *Chlamydia trachomatis* in pharyngeal, semen, serum and rectal, urethral and vaginal swabs are out of the scope of certification by ON 0318. The rest of the pathogens have self-certification for CE marking.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cyclor
	<b>Eppendorf</b>
	Mastercycler™ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTIite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## Bibliography/Bibliografía

1. G. Kusdian et al. The biology of *Trichomonas vaginalis* in the light of urogenital tract infection. *Molecular and Biochemical Parasitology Journal* (2015).
2. F. Kazumasa et al. Rapid and simple detection of *Ureaplasma* species from vaginal swab samples using a loop-mediated isothermal amplification method. *American Journal of Reproductive Immunology* 2017; e12771.
3. CDC Grand Rounds: Chlamydia prevention: challenges and strategies for reducing disease burden and sequelae. *Centers for Disease Control and Prevention (CDC)*.
4. D. Férandon et al. Development of a real-time PCR targeting the *yidC* gene for the detection of *Mycoplasma hominis* and comparison with quantitative culture. *Clinical Microbiology and Infection* 2011; 17(2): 155-159.
5. JB. Slifirski et al. *Mycoplasma genitalium* Infection in Adults Reporting Sexual Contact with Infected Partners, Australia, 2008–2016. *Emerging Infectious Diseases*. 2017;23(11):1826-1833.
6. RD. Kirkcaldy et al. *Neisseria gonorrhoeae* Antimicrobial Susceptibility Surveillance. The Gonococcal Isolate Surveillance Project, 27 Sites, United States, 2014. *MMWR Surveill Summ* 2016;65(No. SS-7):1–19.
7. <http://www.who.int/mediacentre/factsheets/fs110/es/>
8. <https://www.cdc.gov/std/>

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.




Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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