Vitassay qPCR

Pneumocystis jirovecii

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR *Pneumocystis jirovecii* allows the qualitative and quantitative detection of *Pneumocystis jirovecii* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Pneumocystis jirovecii* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR *Pneumocystis jirovecii* 4 x 8-well strip, low profile 7041026

Vitassay qPCR *Pneumocystis jirovecii* 4 x 8-well strip, high profile 7042026

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S026/ 7042S026	<i>Pneumocystis jirovecii</i> strips low/high profile	-	4x8-well strip
7C026	<i>Pneumocystis jirovecii</i> Positive Control	red	1 vial
7Q026	<i>Pneumocystis jirovecii</i> Quantitative standard	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips

- Powder-free disposal gloves

Summary

Pneumocystis pneumonia (PCP) is a serious illness caused by the fungus *Pneumocystis jirovecii*. PCP is one of the most frequent and severe opportunistic infections in people with weakened immune systems, particularly people with HIV/AIDS. Although people with HIV/AIDS are less likely to get PCP today than in recent years, PCP is still a significant public health problem.

The following groups are at risk for *P. jirovecii* pneumonia: persons with HIV infection with low CD4 count, persons with primary immune deficiencies, persons receiving long-term immunosuppressive regimens, persons with hematologic and nonhematologic malignancies, persons with severe malnutrition.

The presentation of PCP in HIV-positive patients is well-known and consists of a triad of dyspnea, fever, and cough, whereas the presentation of PCP in HIV-negative patients is atypical and consists of a sudden outbreak, O₂ desaturation, and a rapid lethal outcome without therapy.

The laboratory diagnosis of Pneumocystis pneumonia (PCP), still relies on tinctorial and/or immunofluorescent staining of bronchoalveolar lavage (BAL) fluid samples. Nucleic acid amplification tests can overcome the difficulties of microscopic examination. Real-time quantitative PCR (qPCR) reduces the risk of false-positive results because of the closed-tube nature of the amplification process, and the resulting data are semi-quantitatives.

Principle of the test

Vitassay qPCR *Pneumocystis jirovecii* is based on the real-time amplification of specific conserved fragments of the *mt large-subunit* rRNA gene of *Pneumocystis jirovecii* genome. After DNA extraction, *P. jirovecii* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *Pneumocystis jirovecii* is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used). In addition, quantification of *P. jirovecii* DNA can be achieved by generating a standard curve using *P. jirovecii* quantitative standard provide with the kit.

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

EZ1 Mini Kit, EZ1 instrument (Qiagen)

Invisorb® Spin Universal Kit (Stratec).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega)

Positive control and Quantitative Standard preparation

Reconstitute the lyophilized *Pneumocystis jirovecii* Positive Control and the Quantitative Standard (red tubes) with 100 µl of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly.

For *Pneumocystis jirovecii* Positive Control, after first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Standard Curve Preparation

To perform a quantitative assay, a standard curve must be prepared by serial dilution of *Pneumocystis jirovecii* Quantitative Standard (which contains approximately 2×10^7 copies/ μl *). Prepare the standard curve dilution series as follows:

- Pipette 90 μl of PCR grade water into 6 microcentrifuge tubes.
- Add 10 μl of *Pneumocystis jirovecii* Quantitative Standard to the first tube to get a standard with about 2×10^6 copies/ml. Mix by vortex and centrifuge briefly.
- Add 10 μl of standard with $\sim 2 \times 10^6$ copies/ μl to the second tube to get a standard with approximately 2×10^5 copies/ μl . Mix by vortex and centrifuge briefly.
- Repeat the previous step sequentially to complete the dilution series for standards with about 2×10^4 , 2×10^3 , 2×10^2 and 2×10^1 copies/ μl .

* For more details see “Certificate of Analysis”, where lot number and DNA copy number of *Pneumocystis jirovecii* Quantitative Standard are detailed.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 μl of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 μl of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*Pneumocystis jirovecii*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (Attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Pneumocystis jirovecii*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for FAM (*Pneumocystis jirovecii*), which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

Qualitative analysis

The result interpretation is summarized in the following table:

<i>Pneumocystis jirovecii</i> FAM	Internal Control HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	<i>Pneumocystis jirovecii</i> Positive
-	+	-	+	<i>Pneumocystis jirovecii</i> Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quantitative analysis

Using standard curve dilution series results a standard curve for quantitative analysis can be generated:

$$Ct = m \log(Q) + b \quad \text{Where: } Ct = \text{Threshold Cycle}$$

$$m = \text{Slope}$$

$$Q = \text{Concentration}$$

$$b = \text{Intercept}$$

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{(Ct-b)/m}$$

The concentration of the samples is displayed in copies/μl and refers to the concentration in the Eluted DNA not to the original clinical sample. To determine the target concentration of the original sample, take into account the dilutions of the extraction procedure and PCR Set-up.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The clinical performance of Vitassay qPCR *Pneumocystis jirovecii* was evaluated using 14 bronchoalveolar lavage from patients with suspicious of *Pneumocystis* pneumonia. We obtained a positive result for 11 samples and a negative result for 3 of them with Vitassay qPCR *Pneumocystis jirovecii*. The results were compared with other real-time PCR commercial kit (FTD *Pneumocystis jirovecii*, fast-track DIAGNOSTICS). The concordance between both tests was 100%.

In addition, Vitassay qPCR *Pneumocystis jirovecii* was evaluated with the EQA programme: QCMD 2017 *Pneumocystis jirovecii* pneumonia DNA (PCPDNA17). This programme consisted on 10 samples with different concentrations of the studied pathogen or were negative to it. Vitassay qPCR *Pneumocystis jirovecii* detected correctly all samples from the studied panel.

These results show the high sensitivity and specificity to detect *Pneumocystis jirovecii* using the molecular diagnostic kit Vitassay qPCR *Pneumocystis jirovecii*.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Pneumocystis jirovecii* template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity for *Pneumocystis jirovecii* was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay	
<i>Staphylococcus aureus</i>	Influenza A/California/7/2009(H1N1) virus
<i>Bordetella pertussis</i>	Influenza A/Perth/16/2009(H3N2) virus
<i>Fluoribacter bozemanæ</i>	Influenza A/New Caledonia/20/99(H1N1) virus
<i>Legionella micdadei</i>	Influenza A/Switzerland/9715293/2013 virus
<i>Legionella dumoffii</i>	Influenza A/Turkey/Germany R2485+86/2014 virus
<i>Legionella longbeachæ</i>	Influenza B/Brisbane/60/2008-like virus
<i>Legionella pneumophila</i>	Influenza B/Florida/04/06 virus
<i>Mycoplasma pneumoniae</i>	Influenza B/Phuket/3073/2013 virus
Methicillin-resistant <i>Staphylococcus aureus</i>	Human Adenovirus 5
<i>Haemophilus influenzae</i>	Human Adenovirus 2 strain Adenoid 6
<i>Moraxella catarrhalis</i>	Human coronavirus 229E
<i>Streptococcus pneumoniae</i>	Human parainfluenza 1, 2, 3 and 4 viruses
Respiratory syncytial virus (RSV)	Human metapneumovirus A and B
Human rhinovirus	

Analytical reactivity

The reactivity of Vitassay qPCR *Pneumocystis jirovecii* was confirmed by the real-time amplification using *Pneumocystis jirovecii* Type 1A, *Pneumocystis jirovecii* g885652 and *Pneumocystis jirovecii* j888023 as template.

Compatible real-time PCR equipment

Vitassay qPCR *Pneumocystis jirovecii* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Pneumocystis jirovecii* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with bronchoalveolar lavage (BAL) specimens. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from respiratory clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Pneumocystis jirovecii*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. *Pneumocystis jirovecii* pneumonia in patients with and without human immunodeficiency virus infection. Su YS, Lu JJ, Perng CL, Chang FY. *J Microbiol Immunol Infect.* 2008 Dec. 41(6):478-82.
2. *Pneumocystis* pneumonia in patients with HIV infection: clinical manifestations, laboratory findings, and radiological features. Fujii T, Nakamura T, Iwamoto A. *J Infect Chemother.* 2007 Feb. 13(1):1-7.
3. Detection of *Pneumocystis jirovecii* by two staining methods and two quantitative PCR assays. Rohner P, Jacomo V, Studer R, Schrenzel J, Graf JD. *Infection.* 2009 Jun;37(3):261-5.
4. Cryptococcal vs *Pneumocystis (carinii) jiroveci* pneumonia (PCP): Clinical & Microbiology differential diagnostic considerations. Cunha BA, Schoch PE, Barbari N. *Infect Dis Pract.* 2006. 30:514-517.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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