# Vitassay qPCR

## E. coli pathotypes

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## Intended use

Vitassay qPCR E. coli pathotypes allows the specific identification and differentiation of Enterohemorrhagic *Escherichia coli* (EHEC), Shiga toxin-producing *Escherichia coli* (STEC), Enteropathogenic *Escherichia coli* (EPEC), Enterotoxigenic *Escherichia coli* (ETEC) and Enteroinvasive *Escherichia coli* (EIEC)/Shigella in human stool samples from patients with suspected gastrointestinal infection. This product is intended to aid in the diagnosis of *Escherichia coli* infections, alongside the patient's clinical data and other laboratory tests outcomes.

## References

Vitassay qPCR E. coli pathotypes 8x8-well strip, low profile	7041064
Vitassay qPCR E. coli pathotypes 8x8-well strip, high profile	7042064

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S064A/ 7042S064A	EHEC, EPEC & EIEC strips low/high profile	-	4 x 8-well strip
7041S064B/ 7042S064B	ETEC + EIEC strips low/high profile	-	4 x 8-well strip
7C064	<i>E. coli</i> pathotypes Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	8 x 8-cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Summary

*Escherichia coli* (*E. coli*) is a facultative anaerobic Gram-negative bacteria member of the bacterial family of *Enterobacteriaceae*. It is the most prevalent commensal inhabitant in human gastrointestinal tract and warm-blooded animals, as well as one of the most important pathogens. The bacterium is able to grow under both aerobic and anaerobic condition, making it an important host organism in biotechnology. The niche of commensal *E. coli* is the mucous layer of the mammalian colon.

This commensal *E. coli* strains cause disease when the host is immunocompromised or where the normal gastrointestinal barriers are breached. It is a major cause of morbidity and mortality in infants and young children. Three general clinical syndromes can result from infection with one of these pathotypes: enteric/diarrheal disease, urinary tract infections (UTIs) and sepsis/meningitis.

The major categories of diarrheagenic *E. coli* strains include enteropathogenic *E. coli* (EPEC), enterotoxigenic *E. coli* (ETEC), enteroinvasive *E. coli* (EIEC) and enterohaemorrhagic *E. coli* (EHEC).

- EPEC was the first pathotype of *E. coli* to be described. These strains cause diarrhea primarily in children, particularly under conditions of poor hygiene, as well as in animals. Clinically, EPEC illness is characterized by acute diarrhea, fever, malaise, and vomiting.
- ETEC are the most common pathogens causing travelers' diarrhea with mild to severe watery diarrhea in humans of all ages. It is also the most common causes of diarrheal illness in military deployed to endemic areas. Some of the symptoms of infection include vomiting, none or mild dehydration, watery stool, and fever.
- EIEC is frequently cause of watery diarrhea and occasionally dysentery in both children and adults. EIEC strains are closely related to *Shigella* spp. EIEC has similar characteristics to *Shigella* (and sometimes similar serotypes) with a milder virulence and a higher infectious dose required. The symptoms are very similar to shigellosis, such as severe bloody, mucous diarrhea, fever, and abdominal cramps.

- EHEC is a typically food-borne pathogen causing hemorrhagic colitis or hemorrhagic uremic syndrome (HUS). EHEC strains produce Shiga-like toxins (named Shiga toxin product *E. coli*, STEC) similar to those produced by *Shigella dysenteriae* making them the most virulent diarrhoeagenic *E. coli* known to date. HUS symptoms are acute renal failure, thrombocytopenia, and microangiopathic hemolytic anemia, which occurs in 2-15% of cases. Mortality from HUS is high (3-17% and a significant number suffer a range of permanent disabilities, including chronic renal insufficiency, hypertension, and neurological deficits.

Nowadays, antimicrobial resistance is a major and increasing global healthcare problem. Most of bacteria have evolved and transmitted antimicrobial resistance to other species. The excessive consumption of antibiotics and the inappropriate use are among factors which further accelerated this phenomenon. In Europe, antimicrobial resistance in Gram-negative bacteria is on the rise, particularly in *E. coli*, making difficult the treatment of several serious infections.

### **Principle of the test**

Vitassay qPCR *E. coli* pathotypes is based on the real-time amplification of a conserved region of *stx1*, *stx2*, *eae*, *lt*, *st1a*, *st1b* e *ipaH* genes for *E. coli* EHEC, *E. coli* EPEC, *E. coli* EIEC/Shigella and *E. coli* ETEC. After nucleic acid isolation, the pathogens presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR *E. coli* pathotypes is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of *E. coli* EHEC, EPEC and EIEC/Shigella. After the amplification reaction, *stx1* and *stx2* genes are detected in FAM channel, the *ipaH* gene is detected in ROX channel, *eae* gene is detected in Cy5 channel and the Internal Control (IC) is detected in HEX, VIC, or JOE channels (depending on the equipment used). The second strip contains the multiplex reaction mix for the detection of *E. coli* ETEC and EIEC/Shigella. After the amplification reaction, *lt* gene is detected in FAM channel, *ipaH* gene is detected in ROX channel, *st1a* and *st1b* genes are detected in Cy5

channel and the Internal Control (IC) is detected in HEX, VIC, or JOE channels (depending on the equipment used).

## **Precautions**

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Make sure to use a well for determining *E. coli* EHEC, STEC, EPEC, and EIEC/Shigella and another well for determining *E. coli* ETEC and EIEC/Shigella. Be careful not to mix them throughout the process.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing, and DNA extraction**

The Vitassay qPCR *E. coli* pathotypes has been tested in human stool samples. Other sample types must be validated by the user.

Overall, samples should be collected and labelled appropriately in clean containers and be processed as soon as possible to guarantee the test quality. It is recommended to use

fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples.

To preserve samples for a long time, we recommend freeze at  $-20^{\circ}\text{C}$ . The samples must be fully thawed to room temperature and homogenized thoroughly before they can be used in the test. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

For nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

- Invisorb® Spin Universal Kit (Invitex Molecular).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).

### **Positive control preparation**

Reconstitute the lyophilized *E. coli* pathotypes Positive Control (red tube) with 200  $\mu\text{L}$  of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at  $-20^{\circ}\text{C}$ .

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

### **Reaction setup**

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15  $\mu\text{L}$  of Resuspension buffer (green tube) and add them into each well.
- Pipette 5  $\mu\text{L}$  of DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (*stx1*, *stx2* and *It* genes), ROX (*ipaH* gene), HEX, JOE or VIC (Internal Control), and Cy5 (*ae* and *st1a* and *st1b* genes). If you use the Applied Biosystems 7500 Fast Real-Time PCR System, and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none (See Attached II).

### Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

For a valid diagnostic test run, the following control conditions must be met:

#### Positive control

The positive control used in each run must show an amplification curve ( $Ct \leq 40$ ) in FAM, ROX, and Cy5 channels, which validates the reaction.

#### Negative control

The negative control included in each run must show signal' absence ( $Ct > 40$  or no signal) in FAM, ROX, and Cy5, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal ( $Ct \leq 40$ ) in positive control and negative control wells.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following tables:

**Results' interpretation for EHEC, EPEC & EIEC 8-well strips:**

EHEC, EPEC & EIEC				Interpretation
<i>stx1/stx2</i> genes (FAM)	<i>eae</i> gene (Cy5)	<i>ipaH</i> gene (ROX)	Internal Control (HEX)	
+	+	-	+/- <sup>1</sup>	EHEC detected
+	-	-	+/- <sup>1</sup>	STEC (EHEC) detected
-	+	-	+/- <sup>1</sup>	EPEC detected
-	-	+	+/- <sup>1</sup>	EIEC/ <i>Shigella</i> detected <sup>3</sup>
+	-	+	+/- <sup>1</sup>	<i>Shigella dysenteriae</i> Type 1 detected <sup>3</sup>
-	-	-	+ <sup>2</sup>	Targets not Detected <sup>2</sup>
-	-	-	- <sup>2</sup>	Test failure <sup>2</sup>

**(+) Positive:** Amplification signal (Ct ≤40)

**(-) Negative:** No amplification signal (Ct >40 or no signal)

**Results' interpretation for ETEC + EIEC 8-well strips:**

ETEC + EIEC				Interpretation
<i>ft</i> gene (FAM)	<i>st1a/st1b</i> genes (Cy5)	<i>ipaH</i> gene (ROX)	Internal Control (HEX)	
+	+	-	+/- <sup>1</sup>	ETEC detected
+	-	-	+/- <sup>1</sup>	ETEC detected
-	+	-	+/- <sup>1</sup>	ETEC detected
-	-	+	+/- <sup>1</sup>	EIEC/ <i>Shigella</i> detected <sup>3</sup>
-	-	-	+ <sup>2</sup>	Targets not Detected <sup>2</sup>
-	-	-	- <sup>2</sup>	Test failure <sup>2</sup>

**(+) Positive:** Amplification signal (Ct ≤40)

**(-) Negative:** No amplification signal (Ct >40 or no signal)

<sup>1</sup> Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

<sup>2</sup> In the case of negative *E. coli* target genes detection, IC must show an amplification signal with Ct ≤35. If there is a signal' absence or Ct value >35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

<sup>3</sup> If a patient sample presents an amplification curve for *ipaH* gene (*Shigella*) in ROX channel in EHEC, EPEC & EIEC and/or ETEC + EIEC reaction mixes, it will be considered positive for *Shigella*/EIEC. In addition, if the *Shigella* positive sample (ROX channel, *ipaH* gene) is also FAM positive in EHEC, EPEC & EIEC reaction mix (*stx1/stx2* genes), *Shigella* will be subtyping as *Shigella dysenteriae* type 1.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to review all the parameters, and the sigmoid shape of the curve. It is also recommended to repeat the test, preferably in duplicate.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

## Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

In one study performed at company facilities, 25 fecal samples from symptomatic patients with clinical suspicion of *E. coli* infection were analyzed with the EHEC, EPEC & EIEC reaction mix. The obtained results were confirmed with that obtained by the commercial RIDA®GENE EHEC/EPEC Real-Time multiplex PCR kit (R-Biopharm) and are shown below:

- 6 positive samples for STEC
- 1 positive sample for *Shigella dysenteriae* Type 1
- 2 positive samples for EHEC
- 4 positive samples for EPEC

In another study also performed at company facilities, 87 fecal samples from patients were analyzed with the ETEC + EIEC reaction mix. The obtained results were confirmed

with that obtained by the commercial RIDA®GENE ETEC/EIEC Real-Time multiplex PCR kit (R-Biopharm) and are shown below:

- 2 positive samples for ETEC
- 5 positive samples for EIEC/Shigella

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10<sup>7</sup> to 10<sup>1</sup> copies/reaction. This assay has a detection limit of ≥10 DNA copies per reaction for *stx1*, *stx2*, *eae*, *ipaH*, *lt*, *st1a* and *st1b* genes.

### Analytical specificity

The analytical specificity for the *Escherichia coli* pathogen variants detection was tested within a panel consisting of the following different microorganisms, where no cross-reactivity was observed between almost any of the species except the targeted pathogens of each assay:

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter coli</i>	-	<i>Candida albicans</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-	<i>Arcobacter butzleri</i>	-
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter</i> <i>upsaliensis</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Helicobacter heilmannii</i>	-	<i>Proteus vulgaris</i>	-	<i>Enterococcus faecalis</i>	-
<i>Shigella flexneri</i>	-/+	<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3	-
<i>Shigella dysenteriae</i>	-/+	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Salmonella typhi</i>	-	<i>Serratia liquefaciens</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella paratyphi A</i>	-	<i>Listeria monocytogenes</i>	-	Adenovirus serotypes 40/41	-
<i>Salmonella paratyphi B</i>	-	<i>Clostridium difficile</i>	-	Rotavirus A	-
<i>Salmonella typhimurium</i>	-	<i>Clostridium perfringens</i>	-	Norovirus Genotypes I and II	-
<i>Salmonella bongori</i>	-	<i>Escherichia coli</i> <i>Enterotoxigénica</i>	-/+	Astrovirus Genotypes I-VIII	-
<i>Salmonella enteritidis</i>	-	<i>Blastocystis hominis</i>	-	Sapovirus	-
<i>Salmonella enterica</i> subsp. <i>enterica</i>	-	<i>Escherichia coli</i> <i>Enterohemorrágica</i>	-/+	<i>Entamoeba histolytica</i>	-
<i>Salmonella pullorum</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Entamoeba dispar</i>	-
<i>Salmonella gallinarum</i>	-	<i>Vibrio parahaemolyticus</i>	-	<i>Cryptosporidium</i> <i>parvum/hominis</i>	-
<i>Campylobacter lari</i>	-	<i>Aeromonas caviae</i>	-	<i>Giardia intestinalis</i>	-
<i>Campylobacter fetus</i>	-	<i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i>	-	<i>Dientamoeba fragilis</i>	-

## Analytical reactivity

The Vitassay qPCR *E. coli* pathotypes reactivity was evaluated against *E. coli* strains of each pathotype and different *Shigella* species which contain the following virulence factors: *stx1* and/or *stx2* and *eae* genes (EHEC), *stx1* and/or *stx2* genes (STEC), *eae* gene (EPEC), *lt* and/or *st1a* and/or *st1b* genes (ETEC) and *ipaH* gene (EIEC, *S. dysenteriae*, *S. flexneri*, *S. sonnei* and *S. boydii*), showing positive results.

## Compatibles real-time PCR equipment

Vitassay qPCR *E. coli* pathotypes has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>1</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)

<sup>1</sup>: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- All obtained results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay has been validated with DNA extracted from human fecal samples. The use of other samples has not been established.
- The test quality depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by *E. coli* pathotypes Positive Control, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 <sup>(1)</sup>
AriaDx Real-Time PCR System	<b>Applied Biosystems</b>
<b>Applied Biosystems</b>	7300 Real-Time PCR System <sup>(3) (1)</sup>
7500 Fast Real-Time PCR System <sup>(1) (2)</sup>	7500 Real-Time PCR System <sup>(1)</sup>
7500 Fast Dx Real-Time PCR System <sup>(1) (2)</sup>	7900 Real-Time PCR System <sup>(3)</sup>
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 <sup>(3)</sup>
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 <sup>(3)</sup>
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System <sup>(4)</sup>	QuantStudio™ 3 Real-Time PCR System <sup>(3)</sup>
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
<b>Azure Biosystems</b>	<b>Agilent Technologies</b>
Azure Cielo 3 <sup>(5)</sup>	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
<b>BIONEER</b>	<b>Analytik Jena</b>
Exicycler™ 96	qTOWER
<b>Bio-Rad</b>	<b>BIONEER</b>
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System <sup>(5)</sup>	<b>BIOER</b>
<b>Roche</b>	QuantGene 9600
LightCycler®480 Real-Time PCR System <sup>(1) (6)</sup>	<b>Bio-Rad</b>
LightCycler®96 Real-Time PCR System <sup>(1)</sup>	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer <sup>(1) (6)</sup>	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System <sup>(5)</sup>
	My iQ™ 2 Real-Time PCR Detection System <sup>(5)</sup>
<b>Special Formats <sup>(7)</sup></b>	<b>DNA-Technology</b>
<b>Bio Molecular Systems</b>	DTlite Real-Time PCR System <sup>(8)</sup>
Mic Real Time PCR Cyclers	DTprime Real-time Detection Thermal Cyclers <sup>(8)</sup>
<b>Cepheid</b>	<b>Eppendorf</b>
SmartCycler®	Mastercycler™ ep <i>realplex</i>
<b>Precision System Science Co., Ltd.</b>	<b>Qiagen</b>
geneLEAD VIII System	QIAquant 96
<b>Qiagen</b>	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

## Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Applied Biosystems ABI 7500</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Roche Cobas z 480</b>	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
<b>Cepheid Smartcycler®</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Agilent Technologies Mx3000P™ Mx 3005P™</b>	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Bio Molecular Systems Mic Real Time PCR Cycler</b>	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Agilent Technologies AriaMx</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Qiagen Rotor-Gene® Q</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>BIONEER Exicycler™ 96</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

\*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

## Bibliography/Bibliografía

- Allocati N, Masulli M, Alexeyev MF and Di Ilio, Carminie. (2013). Escherichia coli in Europe: An Overview. International journal of environmental research and public health. 10(12), 6235-6254.
- Belotserkovsky I and Sansonetti PJ. (2018). Shigella and Enteroinvasive Escherichia Coli. Current topics in microbiology and immunology. 416:1-16.
- Chen HD and Frankel G. (2005). Enteropathogenic Escherichia coli: unravelling pathogenesis. FEMS microbiology reviews. 29(1):83-98.
- Fleckenstein JM and Kuhlmann FM. (2019). Enterotoxigenic Escherichia coli Infections. Current infectious disease reports. 21(3):9.
- Gomes T A, Elias WP, Scaletsky IC, Guth BE, Rodrigues JF, Piazza RM, Ferreira LC, and Martinez MB. (2016). Diarrheagenic Escherichia coli. Brazilian journal of microbiology. Brazilian Society for Microbiology. 47(Suppl 1):3–30.
- Kaper JB, Nataro JP and Mobley HL. (2004). Pathogenic Escherichia coli. Nature reviews. Microbiology. 2(2):123-140.
- Qadri F, Svennerholm AM, Faruque AS and Sack RB. (2005). Enterotoxigenic Escherichia coli in developing countries: epidemiology, microbiology, clinical features, treatment, and prevention. Clinical microbiology reviews. 18(3):465-483.
- Welinder-Olsson C and Kaijser B. (2005). Enterohemorrhagic Escherichia coli (EHEC). Scandinavian Journal of Infectious Diseases. 37(6-7):405-416.

## Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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