Vitassay qPCR

BK

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Intended use

Vitassay qPCR BK allows the identification and quantitative detection of DNA from BK virus by Real-Time PCR in human clinical samples. This product is intended to aid in the diagnosis of BK virus infections in transplant patients, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes. The test can be used to assess the need and monitoring the treatment of the infection in these patients.

References

Vitassay qPCR BK 4x8-well strip, low profile	7041072
Vitassay qPCR BK 4x8-well strip, high profile	7042072

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S072/ 7042S072	BK strips low/high profile	-	4 x 8-well strip
7Q0721- 7Q0725	BK Quantitative standards	red	5 vials (1 vial per concentration)
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended BK Quantitative Standards are stable at room temperature during the time needed to dispense 96-wells, but they should be stored at 4 °C for up to 5 months, or at -20°C or lower for longer period. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves
- 1st WHO International Standard for BK Virus DNA (NIBSC code 14/212), or another reference material (optional, when using a different workflow and in case the final user needs to calculate its own conversion factor).

Summary

BK polyomavirus is a circular, double-stranded DNA virus from the *Polyomaviridae* family and *Betapolyomavirus* genus. BK and JC virus are the more representatives of polyomavirus. BK polyomavirus can be divided into six subtypes or genotypes. Genotype I is the most frequent worldwide (80%), followed by Genotype IV (15%). More than 80% of adults are seropositive for BK polyomavirus, and infection usually is acquired during childhood. This is because BK polyomavirus remains latent in genitourinary system, mostly in urothelial and renal tubular epithelial cells.

The mode of transmission is thought to be by person to person contact and via contaminated surfaces, food, and water. Primary infection is usually asymptomatic, and if there are some symptoms, they are mild such as fever, upper respiratory symptoms, and transient cystitis. A progressive multifocal leukoencephalopathy, meningitis and encephalitis, retinitis, pneumonitis, prostate cancer, HIV-associated salivary gland disease, and renal carcinoma are pathologies less common associated with BK virus. In the case of immunocompromised individuals, BK virus can be associated with nephropathy, ureteral stenosis, late-onset hemorrhagic cystitis and even bladder cancer. The renal spectrum begins with viruria and ends with graft failure. BK polyomavirus could also cause hemorrhagic cystitis in patients with allogeneic stem cell transplantation.

Regardless extensive research, the risk factors for BK polyomavirus nephropathy remain controversial. Male gender and extremes of age, among other things, have been suggested as possible risk factors. The widely held belief is that BK polyomavirus infects people who are immunocompromised. In bone marrow transplant recipients, the BK polyomavirus is prevalent, also BK polyomavirus affects other non-kidney solid organ transplants such as heart, liver, and lung.

As BK virus infections develops asymptotically, regular monitoring of the viral load through quantitative PCR is the most effective technique for monitoring viruria and viraemia, as well as early diagnosis and treatment of BK nephropathy. Furthermore, measuring BKV mRNA allows for the detection of ongoing viral transcription.

Principle of the test

Vitassay qPCR BK is based on real-time amplification of a conserved region of the *vp1* gene. After DNA extraction, the viruses presence is detected by an increase in fluorescence observed during the reaction, following hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect the accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to extend the primers, the probe is hydrolysed by its 5'-3' exonuclease activity resulting in spatial separation of the fluorophore and the quencher. The resulting increase in fluorescent signal is proportional to the amount of amplified product in the sample and is detected by real-time PCR equipment.

Vitassay qPCR BK is a ready-to-use test that contains in each well all the necessary reagents in a stabilised format to perform real-time PCR. In addition, an endogenous internal control allows the extraction process control and the detection of a possible inhibition reaction. The assay uses a human housekeeping gene as an endogenous internal control (IC) (*RNase P* gene present in human DNA) that is expected to be present in all nucleated human cells. After the amplification reaction, the BK virus is detected in the FAM channel and the Endogenous Internal control (IC) is detected on HEX, VIC, or JOE channels (depending on the equipment used).

In addition, this kit contains 5 lyophilized vials with different standardized concentrations of BK virus specific DNA (BK Quantitative Standards) for the quantification of BK virus DNA in clinical samples. These quantitative standards are calibrated against the 1st WHO International Standard for BK virus DNA (NIBSC code: 14/212), ensuring the quality of the results and the comparison with results from other laboratories.

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.

- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Do not mix BK Quantitative Standards from different kits and/or lots. Each quantitative standard set has been validated for each lot of master mix with which it is supplied.
- Inaccurate volumes due to imprecise pipetting or malpractices may result in an erroneous quantification. Use calibrated micropipettes and proper tips and follow the test procedure of this instructions of use.
- Standard curve should be included in the same run of samples to be quantify and should be generated by running each BK Quantitative Standard in triplicate.
- In case of using **1st WHO International Standard for BK Virus DNA** (NIBSC code 14/212), note that **this standard stability can be affected by some practices**. It is recommended to reconstitute the standard and prepared aliquots for single use only by maintaining the cold chain since freeze-thawing cycles can affect the reference material integrity.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR BK has been tested on human EDTA-plasma samples. Other sample types should be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 2 hours. For long term transport (more than 2 hours), we recommend shipping at -20°C or lower. The samples can be stored frozen at -20°C or lower (at -80°C ideally) for conservation. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guideline:

- Miller JM, et al. (2018). A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2018 Update by the Infectious Diseases Society of America and the American Society for Microbiology. Clin Infect Dis. Aug 31;67(6):e1-e94.

DNA extraction

For the sample's pretreatment follow the instructions for use of the extraction kit used. Vitassay qPCR BK has been validated with EDTA-plasma samples using the following reference workflow (A), which is the one applied in the analytical and clinical validation, and the one that serves as the basis for the interpretation of results:

MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.) and CFX96™ Real-Time PCR Detection System (Bio-Rad).

Reference workflow for Conversion Factor (CF) calculation

Before BK virus quantification from clinical samples, an extraction-amplification procedure (workflow) for which the CF is to be calculated must be decided in advance and maintained for further analysis. Due to BK virus concentration must be express in International Units (IU), and the quantitative results are obtained in copies per mL after interpolation in the linear calibration curve, a conversion factor (CF) is required to be calculated to relate copies/mL values with IU/mL data. **The conversion factor (CF) depends on the tested sample type and the used workflow.**

The conversion factor value for the reference workflow was calculated using the 1st WHO International Standard for BK Virus DNA (NIBSC code: 14/212) and established at **0.420 IU/copies**.

IU/mL = copies/mL x 0.420

For instance, a viral load determined at 10,000 copies/mL in a sample would correspond to 4,200 IU/mL.

Other workflows were validated using Vitassay qPCR BK, even if they are not considered reference workflows:

Workflow	Extraction	Amplification	CF (only applicable with EDTA-plasma samples)
B	MagDEA Dx SV kit - magLEAD® 12gC instrument (Precision System Science Co.)	NEOS-96 Real-Time PCR System (Linear)	CF = 0.50 IU/copies
C	MagDEA Dx SV kit - magLEAD® 12gC instrument (Precision System Science Co.)	7500 Fast Real-Time PCR System (Applied Biosystems)	CF = 0.54 IU/copies
D	Maxwell® RSC 16 Viral TNA/ Maxwell® RSC (Promega)	CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)	CF = 0.71 IU/copies
E	QIAamp® MiniElute® Virus Spin Kit (QIAGEN)	CFX96™ Real-Time PCR Detection System (Bio-Rad)	CF = 0.925 IU/copies

Each laboratory must establish its own conversion factor. If another workflow is set, then the CF value must be calculated by the user.

Standard curve preparation

This kit contains 5 lyophilized vials (7Q0721-7Q0725) with different standardized concentrations of BK virus specific DNA (BK Quantitative Standards) for the quantification of BK virus viral DNA in clinical samples. BK Quantitative Standards can be used to generate a standard curve with which to determine the amount of BK virus viral DNA in a sample in IU/mL.

These Quantitative Standards were calibrated against the 1st WHO International Standard for BK Virus DNA (NIBSC code: 14/212), thus ensuring the results quality and comparison with those from other laboratories. The lyophilized BK Quantitative Standards panel of

Vitassay qPCR BK are ready to use after reconstitution with PCR grade water, and the BK Quantitative Standards have the following concentrations:

BK Quantitative Standard	Concentration (copies/mL)
7Q0721	10^7
7Q0722	10^6
7Q0723	10^5
7Q0724	10^4
7Q0725	10^3

To prepare the standard curve to perform a quantitative assay, reconstitute each lyophilized BK Quantitative Standard (red tube) with 100 μ L of PCR grade water (white tube), changing the micropipette tip between vials. To ensure a complete resuspension, vortex the tube thoroughly (20 sec).

After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C .

BK Quantitative Standards contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples. It is recommended that the BK Quantitative Standards be rehydrated and vortexed at the time of dispensing the samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that at least one replicate of each BK Quantitative Standard and one replicate of negative control must be included in each run for each assay. *The standard curve should be generated by adding each BK Quantitative Standard in triplicate and should be included on the same run as the samples to be quantified.
- Peel off protective aluminium seal from the strips.
- Pipette 15 μ L of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 μ L of DNA sample, negative control (yellow tube) or each reconstituted BK Quantitative Standard (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (BK virus), and HEX, JOE, or VIC (Endogenous Internal Control). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, and the extraction procedure check the endogenous internal control (IC) signal emission.

Use the BK Quantitative Standards amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

BK Quantitative Standards

The BK Quantitative Standards included in each run should show an amplification curve in FAM channel, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence in FAM, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the BK Quantitative Standards for any channel. The assay should be repeated.

The endogenous internal control (IC) should show no signal in the BK Quantitative Standards and negative control wells. In Negative Control, it is possible that non-specific amplifications in late Ct values (Ct >35) may appear due to human manipulation of the process, since human *RNAse P* gene is a human housekeeping gene that is expected to be present in all nucleated human cells.

From the BK Quantitative Standards dilution series, a standard curve for quantitative analysis can be generated using the formula:

$$Ct = m \log (Q) + b$$

Where: Ct = Threshold Cycle

m = Slope

Q = Concentration

b = Intercept

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{(Ct-b)/m}$$

The obtained quantification results are obtained in the concentration units of the calibration set.

NOTE: Standard curve should be included in the same run of samples to be quantify and should be generated by running each of BK Quantitative Standard in triplicate.

The quantification results are valid if the generated standard curve reaches the following control parameter values:

Efficiency: 80-120%.

R ² : ≥ 0,98

The concentration of the "Positive Sample" refers to the concentration of the eluted DNA after extraction, not of the original clinical sample. **To determine the concentration of the original sample, take into account the dilutions corresponding to the extraction and PCR preparation.**

Once the controls results have been validated, in case the workflow used was the reference workflow with EDTA-plasma samples extracted with MagDEA Dx SV Kit/magLEAD® 12 gC instrument (Precision System Science Co.) and amplified with

CFX96™ Real-Time PCR Detection System (Bio-Rad), use the following table to analyze the samples results:

BK Virus (FAM)	Endogenous Internal Control (HEX)	Interpretation	
< 102,31 IU/mL	+/- ¹	Valid	Non-quantifiable BK Virus positive sample
102,31 a 420 IU/mL	+/- ¹	Valid	Quantifiable BK Virus positive sample by extrapolation ²
420 a 4,2 x 10 ⁶ IU/mL	+/- ¹	Valid	Quantifiable BK Virus positive sample in the quantification range
> 4,2 x 10 ⁶ IU/mL	+/- ¹	Valid	Quantifiable BK Virus positive sample by extrapolation ³
-	+ ⁴	Valid	Target not detected⁴
-	- ⁴	Invalid	Test Failed⁴

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (No signal)

¹Sometimes, the Endogenous Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Endogenous Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

²BK virus target gene positive below the quantification range let quantification of BK virus viral load by extrapolation of the standard curve, but results could not be accurate.

³BK virus target gene positive above the quantification range let quantification of BK virus viral load by extrapolation of the standard curve, but results could not be accurate. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100.

⁴In the case of negative BK virus target genes detection, endogenous IC must show an amplification signal with Ct ≤35. The Ct value could be very variable due to the Endogenous Internal Control is a human *housekeeping* gene that should be present in all human nucleated cells in the original sample. If there is a signal' absence or Ct value > 35 of the Endogenous Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Endogenous Internal Control (IC) is included in each reaction. Moreover, a negative control and BK Quantitative Standards must be included in each assay to interpret and quantify the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR BK kit's clinical performance was assessed using plasma leftovers collected in EDTA-containing tubes from patients with compatible symptoms of BK virus infection. After the extraction method and the analysis with the routine method, which is a commercial CE-IVD real-time PCR assay, the BK virus DNA was quantitated using the quantification standards providing in the kit and results was given by the corresponding software. Of the total collected samples, 93 samples could be quantifiable, and they were later analysed with Vitassay qPCR BK and correlation between the results obtained with the reference methods and Vitassay assay were calculated using the Spearman's rank correlation statistical analysis, showing a strong coefficient of correlation ($r_s=0.987$; $p<0.001$).

A qualitative comparative analysis, using 161 specimen remnants, was also carried out, indicating a sensitivity value of 0,99 (0,96-0,99) and a specificity value of 1(0,96-1). In conclusion, Vitassay qPCR BK shows a great consistency to detect and quantify the BK virus presence in clinical samples.

Linear range

The Vitassay qPCR BK kit's linearity was determined using specific and synthetic DNA belonging to BK virus, and further was confirmed with the use of reference material in plasma samples.

In the first place, a series of ten-fold dilutions containing a known concentration (ranging from 10^7 to 10^1 copies/ μ L) of specific BK virus DNA (7Q072) was amplified, and Ct values in all dilution tested were obtained using CFX96™ Real-Time PCR Detection System (Bio-Rad) (FAM channel), and the following data were obtained:

$$E= 83.3 \% ; R^2=0.999 ; \text{Slope}=-3.799 \text{ and } -\text{int}= 40.966$$

In the second place, linear range was further confirmed using negative plasma samples spiked with a known concentration of BK virus reference material (1st WHO IS for BK Virus

DNA (NIBSC code 14/212) using CFX96™ Real-Time PCR Detection System (Bio-Rad) (FAM channel).

Vitassay qPCR BK obtained amplifications in all dilutions tested. End-point fluorescence level was maintained, and no crosstalk was observed. The assays with 1st WHO IS for BK Virus DNA showed statistical significance:

$$Y = 3.6868x + 46.672; R^2 = 0.9961$$

Analytical sensitivity

For the Vitassay qPCR BK LoD determination, a dilution series of the 1st WHO International Standard for BK Virus DNA (NIBSC code 14/212) and a tentative LoD concentration of NATrol™ BK Virus External Run Controls (Zeptometrix, Ref: NATBK-ERCL) were set up.

The extraction of 71 negative plasma samples per concentration level of 1st WHO IS for BK Virus DNA reference material and 106 negative plasma samples replicates of the NATrol™ BK Virus External Run Controls reference material was performed with the MagDEA Dx SV kit / magLEAD® 12gC instrument (Precision System Science Co.) and their amplification was carried out in the CFX96™ Real-Time PCR Detection System (Bio-Rad).

Vitassay qPCR BK kit showed a detection limit of 243.6 copies/mL, equivalent to 102.31 IU/mL for BK virus (positive rate ≥ 95%) using the 1st WHO IS BK Virus DNA (NIBSC code 14/212) and Probit Model to analyze the results.

Samples with 1 st WHO IS BK Virus DNA (NIBSC code 14/212)				
BK virus concentration (copies/mL)	BK virus concentration (IU/ml)	Nº of replicates	Nº of positive replicates	Hit rate (%)
1000	420	71	71	100
800	336	71	71	100
600	252	71	70	98.6
400	168	71	69	97.2
200	84	71	68	95.8
50	21	71	37	52.1

Vitassay qPCR BK showed a detection limit of 100 copies/mL, equivalent to 36.4 IU/mL for BK virus (positive rate ≥ 95%) using NATrol™ BK Virus External Run Controls reference material. **IU** = International Units.

Samples with NATtrol™ BK Virus External Run Controls (Zeptomatrix, Ref: NATBK-ERCL)				
BK virus concentration (copies/mL)	BK virus concentration (IU/ml)	Nº of replicates	Nº of positive replicates	Hit rate (%)
100	36.4	106	106	100

Lower Limit of Quantification

The lower limit of quantification (LLoQ) is the lowest concentration value that BK virus can be quantified and at which the total error meets the acceptance criterium of being below 40%. The LLoQ for BK virus was determined to be 1000 copies/mL, equivalent to 420 IU/mL.

Analytical specificity

The analytical specificity for BK virus detection was tested within the panel of following microorganisms, where no cross-reactivity was observed with any of the species.

Cross-reactivity testing					
Parvovirus B19 genotypes 1, 2 and 3a	-	<i>Klebsiella oxytoca</i>	-	<i>Plasmodium falciparum</i> strain 3D7	-
<i>Acinetobacter baumannii</i>	-	Herpes simplex virus 1 (HSV-1) strain MacIntyre	-	<i>Pneumocystis jirovecii</i> type A1 and g885652	-
<i>Aspergillus fumigatus</i>	-	Herpes simplex virus 2 (HSV-2) strain MS	-	<i>Pseudomonas aeruginosa</i>	-
<i>Bacteroides fragilis</i>	-	Human herpesvirus 6 (HHV6) strain Z29 and types A and B	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
JC virus types 1A and 2B	-	Hepatitis A Virus	-	<i>Toxoplasma gondii</i> type II	-
<i>Candida albicans</i>	-	Cytomegalovirus strain AD-169	-	<i>Treponema pallidum</i>	-
<i>Enterobacter aerogenes</i> serotype <i>Cloacae B</i>	-	<i>Listeria monocytogenes</i> serotype 1/2b and serovar 4b	-	Varicella-Zoster Virus strain Ellen	-
<i>Enterobacter cloacae</i> serotype <i>Cloacae A</i>	-	<i>Listeria innocua</i> serotype 6a	-	<i>Yersinia enterocolitica</i> O:3	-
Epstein-Barr virus	-	<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	-	Adenovirus types 1-5, 8, 15, 31, 40 and 41	-
<i>Escherichia coli</i> strain 0.1285;O18:H7:K1	-				

Analytical reactivity

Vitassay qPCR BK kit's reactivity for BK virus was tested against DNA from BK virus, subtypes 1b-1 and 1b-2, and subtype IV as templates, showing positive results.

Limitations

- All obtained results should be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from EDTA-plasma samples. The use of other samples has not been established.
- The correct test performance depends on the sample's quality; DNA should be properly extracted from clinical samples.
- Lower levels of target below the quantification limit might be detected, but results may not be accurately quantifiable.
- Low copy numbers of the target DNA template may be detected below the detection limit, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination, either by BK Quantitative Standards during reconstitution, by samples containing high concentrations of target DNA template or by PCR products from previous reactions.
- Detection can be affected by several factors and their combinations that can lead to false negative results, including: a) inadequate specimen sampling, shipping, storage and handling; b) procedural errors (including DNA extraction); c) DNA degradation during specimen shipping, storage and/or preparation; d) mutations or polymorphisms in primer or probe binding regions may affect the detection of new or unknown BK virus genes mutations; e) pathogen load is below the detection limit for the assay; f) presence of real-time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapies, antibiotics, chemotherapeutics or immunosuppressive drugs used to prevent infection or during treatment of infection was not performed); g) failure to follow the manufacturer's instructions and procedures.
- Detection can be affected by several factors and their combinations that can lead to inaccurate viral DNA values (loss of viral load), which include: a) inadequate specimen sampling, shipping, storage and handling; b) procedural errors (including DNA extraction); c) incorrect maintenance of commonly used equipment, especially micropipettes and extraction/amplification systems; d) use of uncalibrated micropipettes and/or inadequate tips; e) failure to follow the manufacturer's instructions and procedures.
- Detection of viral DNA may not indicate the presence of viable and/or infectious virus or that these viruses are the causative agents for clinical symptoms.
- Negative results do not preclude BK virus infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same

patient at different time points, which may increase the probability of detecting the virus.

- If the patient clinical data, laboratory tests and epidemiological studies suggest possible BK virus infection, and other immunodepression-related illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.
- The conversion factor (CF) depends on the tested sample type and the used workflow (extraction and amplification). For this reason, each sample type-workflow requires the calculation of its own CF.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyler	DTprime Real-time Detection Thermal Cycler ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	


Bibliography/Bibliografía

- Blackard J.T., Davies S.M., & Laskin B.L. (2020). BK polyomavirus diversity – Why viral variation matters. *Reviews in medical virology*. 30(4):e2102.
- Furmaga, J., Kowalczyk, M., Zapolski, T., Furmaga, O., Krakowski, L., Rudzki, G., Jaroszyński, A., & Jakubczak, A. (2021). BK Polyomavirus-Biology, Genomic Variation and Diagnosis. *Viruses*, 13(8), 1502.
- Kotla, S. K., Kadambi, P. V., Hendricks, A. R., & Rojas, R. (2021). BK polyomavirus-pathogen, paradigm and puzzle. *Nephrology, dialysis, transplantation : official publication of the European Dialysis and Transplant Association - European Renal Association*, 36(4), 587–593.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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