Vitassay qPCR

Aspergillosis

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- Cobas z480 Analyzer (Roche)
- LightCycler® 480 Instrument II (Roche)
- QuantGene 9600 (Bioer Technology)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- CFX96™ Real-Time PCR Detection System (Bio-Rad).
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- NEOS-96 Real-Time PCR System (Linear)

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Intended use

Vitassay qPCR Aspergillosis allows the qualitative detection and identification of DNA from *Aspergillus fumigatus*, *Aspergillus flavus* and *Aspergillus terreus* by Real-Time PCR in BAL, BAS and sputum samples from patients suspected of aspergillosis by their healthcare professional. This product is intended to aid in the aspergillosis diagnosis, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes.

References

Vitassay qPCR Aspergillosis 4x8-well strip, low profile	7041076
Vitassay qPCR Aspergillosis 4x8-well strip, high profile	7042076

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S076/ 7042S076	Aspergillosis strips low/high profile	-	4 x 8-well strip
7C076	Aspergillosis Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes

- Vortex
- Micropipettes (1-20 μ L, 20-200 μ L)
- Filter tips
- Powder-free disposal gloves

Summary

Aspergillus is a worldwide distributed specie, whose members are saprotrophic fungi which are ubiquitous in the environment, and can be found in the soil, water and building materials. It is believed that regional differences in precipitation, humidity, and temperature contribute to the environmental prevalence of these species. For example, in Barcelona, *A. flavus* was the most common species of airborne aspergilli, whereas in Madrid, it was *A. fumigatus*, according to two studies conducted in Spain.

Aspergillosis causes a wide range of clinical conditions, and the way in which they manifest depends on both the site of infection and the host's immune response. Although up to 10% of patients may experience dissemination, the main locations which are most usually affected are the upper airway, bronchi, lung parenchyma, and contiguous structures (eg. pleural or lymph nodes). Allergic bronchopulmonary aspergillosis (ABPA) and allergic fungal rhinosinusitis are noninvasive types, whereas chronic pulmonary aspergillosis and invasive pulmonary aspergillosis (IPA) are more invasive disorders. Some nonspecific symptoms of aspergillosis are dyspnea, persistent productive cough, and wheezing, low-grade fever, weight loss, night sweats, chest pain, mucus, and hemoptysis.

A hypersensitivity reaction to *Aspergillus* spp. in the airway could lead to an allergic bronchopulmonary aspergillosis, which is primarily described in asthma and cystic fibrosis patients, although case reports have been published in other lung illnesses, as well.

Chronic pulmonary aspergillosis is often found in immunocompetent patients with underlying lung conditions. *Aspergillus* nodules and aspergillomas are two examples of manifestations.

The most serious form of pulmonary aspergillosis is invasive pulmonary aspergillosis, which is common in immunocompromised hosts. The severity of the illness is what distinguishes this etiology from others because the symptoms are tachypnea, tachycardia, and hypoxia, causing acute respiratory failure over a few days.

A.fumigatus is the most common *Aspergillus* spp. that causes invasive aspergillosis and allergic disease, although other species, such as *Aspergillus flavus* and *Aspergillus terreus*, can also cause these diseases.

One of the most common fungi with airborne conidia is *A. fumigatus*. The conidia released into the atmosphere have an alveolar-reaching diameter. In immunocompromised hosts frequently causes severe and fatal invasive infections in developed countries.

The second cause of human invasive aspergillosis is *A. flavus*, which can produce potent mycotoxins. It is also transmitted by airborne conidia, which easily disperse by air movements and possibly by insects.

Nearly 4% of all cases of invasive aspergillosis are caused by *A. terreus*, an increasingly well-known opportunistic fungi which is usually associated with dissemination and poor outcome. A growing number of reports have linked its infections to both acute and chronic aspergillosis.

Aspergillus spp. usually grow on routine media. The gold standard for their diagnosis is via histopathological examination and culture of a surgical lung biopsy. Due to the fact that fungal culture is slow and insensitive, and histopathology is not organism specific, a good alternative for *Aspergillus* spp. diagnostic is the use of molecular diagnostic, which is a faster and more accurate method.

Principle of the test

Vitassay qPCR Aspergillosis is based on the real-time amplification of a conserved region of *ITS1* of *A. flavus* and *A. terreus*, and *23S rRNA* of *A. fumigatus*. After DNA isolation, the pathogens presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Aspergillosis is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. After the amplification reaction, *Aspergillus flavus* is detected in FAM channel, *Aspergillus fumigatus* is detected in HEX, VIC, or JOE channels (depending on the equipment used), *Aspergillus terreus* is detected in ROX channel and the Internal Control (IC) is detected in Cy5 channel.

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR Aspergillosis has been tested in respiratory samples (bronchoalveolar lavages (BAL), bronchoalveolar aspirates (BAS), and sputum). Other sample types should be validated by the user.

Overall, respiratory samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the

test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 2 hours or at 2°C to 8°C for up to 5 days. For long term transport (more than 5 days), we recommend shipping at -20°C or lower. The samples can be stored at 2°C to 8°C for up to 5 days or frozen at -20°C or lower (at -80°C ideally) for a long time. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- Miller JM, et al. (2018). A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2018 Update by the Infectious Diseases Society of America and the American Society for Microbiology. Clin Infect Dis. Aug 31;67(6):e1-e94.

DNA extraction

For nucleic acid isolation from respiratory samples, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, empleando el equipo Maxwell® RSC 16 instrument (Promega)

MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.)

*The main problem with the use of PCR for the diagnosis of aspergillosis is DNA extraction. This is because they have a complex and robust cell wall that is resistant to lysis, therefore, DNA extraction is impaired by the methods most commonly used with bacteria. Different protocols for the extraction of fungal DNA have been described in the literature (Griffiths et al., 2006). A pre-treatment of the sample consisting of a 10-minute incubation at 56 °C in a mixture of Lysis Buffer with Proteinase K has been shown to be very useful for the detection of *Aspergillus*.

Positive control preparation

Reconstitute the lyophilized Aspergillosis Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*M. tuberculosis* complex), and HEX, JOE, or VIC (Internal Control). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve (Ct ≤40) in FAM, HEX and ROX channel, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence (Ct >40 or no signal) in FAM, HEX and ROX channel, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in positive control and negative control wells.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Aspergillosis				Interpretation
<i>Aspergillus flavus</i> (FAM)	<i>Aspergillus fumigatus</i> (HEX)	<i>Aspergillus terreus</i> (ROX)	Internal Control (Cy5)	
+	-	-	+/- ¹	<i>A. flavus</i> DNA Detected
-	+	-	+/- ¹	<i>A. fumigatus</i> DNA Detected
-	-	+	+/- ¹	<i>A. terreus</i> DNA Detected
+	+	-	+/- ¹	<i>A. flavus</i> and <i>A. fumigatus</i> DNA Detected
+	-	+	+/- ¹	<i>A. flavus</i> and <i>A. terreus</i> DNA Detected
-	+	+	+/- ¹	<i>A. fumigatus</i> and <i>A. terreus</i> DNA Detected
+	+	+	+/- ¹	All targets' DNA Detected
-	-	-	+ ²	Targets not detected ²
-	-	-	- ²	Failed test

(+) **Positive:** Amplification signal (Ct ≤40)
 (-) **Negative:** No amplification signal (Ct >40 or no signal)

¹ Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

² In the case of negative pathogens target genes detection, IC must show an amplification signal with Ct ≤ 35. If there is a signal' absence or Ct value >35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR Aspergillosis was tested with clinical samples collected from male and female adult patients with clinical manifestations of a low respiratory tract infection. A total of six hundred and eighty remnants of sputum, bronchoalveolar lavages and aspirates, previously diagnosed as positive or negative for *Aspergillus*, among other microorganisms, were analysed. The initial diagnosis using the laboratory gold standard (mass spectrometry, culture, microscopy, and molecular methods) for diverse pathogens, showed the presence of 43 positive samples for different species of *Aspergillus*, along with gram-negative and positive bacteria and other fungi and yeast.

Comparative analysis between the reference method and Vitassay qPCR, and after assessing the discrepancies, showed that Vitassay was able to detect 3/680 positive samples for *A. flavus* (true positive). No false negatives or false positives were obtained for this target. Hence, specificity and sensitivity values for *A. flavus* were of 1(0.54-1) and 1 (0.99-1), respectively.

Regarding *A. fumigatus*, 41/680 samples were identified as positive for both reference method and Vitassay qPCR. No discrepancies were found, obtaining sensitivity and specificity values for this target of 1(0.91-1) and 1 (0.99-1), respectively.

Finally, for *A. terreus*, 9 samples were positive in Vitassay qPCR and the reference method. After analysing the discrepant results, one false negative sample was found for this target, leading to a sensitivity of 0.9 (0.55-0.99) and specificity of 1 (0.99-1), for this target.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of 4×10^{-3} CFU per DNA reaction for *Aspergillus flavus*, 4×10^{-2} CFU per reaction for *Aspergillus fumigatus*, and 8×10^{-2} CFU per reaction for *Aspergillus terreus* (positive rate $\geq 95\%$), on serum samples, and 8×10^{-2} CFU per reaction for *Aspergillus flavus*, 0.1 CFU per reaction for *Aspergillus fumigatus* and 4×10^{-2} CFU per reaction for *Aspergillus terreus* (positive rate of $\geq 95\%$), on BAL samples.

Analytical specificity

The analytical specificity for the different pathogen detection, was tested within the panel of following microorganisms, where no cross-reactivity was observed with any of the species:

Cross-reactivity testing					
Adenovirus humano tipos 1-5, 8, 15, 31, 40 y 41	-	<i>Haemophilus influenzae</i> <i>MinnA</i>	-	<i>Legionella bozemanii</i> serovar 1	-
<i>Aspergillus niger</i>	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Legionella dumoffii</i>	-
Bocavirus	-	Influenza A/California/7/2009(H1N1) pdm09 virus	-	<i>Legionella longbeachae</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella micdadei</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Singapore/GP1908/201 5, IVR-180 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i> serogrupo 1	-
<i>Bordetella parapertussis</i>	-	Influenza A/Victoria/210/2009 (H3N2) virus	-	Metapneumovirus humano A1, A2 y B	-
<i>Bordetella pertussis</i>	-	Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	<i>Moraxella catarrhalis</i>	-
<i>Chlamydia caviae</i>	-	Influenza A/Switzerland/9715293/20 13 (H3N2) virus	-	<i>Mycoplasma pneumoniae</i>	-
<i>Chlamydia psittaci</i> genotipos A y C	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> no resistente a rifampicina	-
<i>Chlamydomydia pneumoniae</i> CM-1	-	Influenza A/South Australia/55/2014, IVR- 175 (H3N2) virus	-	Virus parainfluenza humano 1, 2, 3 y 4	-
Coronavirus humano 229E, OC43, NL63	-	Influenza A/DE- SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Pneumocystis jirovecii</i> Tipo A1 y g885652	-
MERS-CoV	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	Rinovirus humano tipo C	-
Enterovirus 68 y 71	-	Influenza B/Brisbane/60/2008 virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Enterovirus Echovirus 30	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> Z022	-
Enterovirus Coxsackievirus A24, A9 y B3	-	Influenza B/Phuket/3073/2013 virus	-	Virus Respiratorio Sincitial tipos A y B	-

Analytical reactivity

The Vitassay qPCR Aspergillosis kit's reactivity for *Aspergillus* species was evaluated against extracted DNA from the following strains, showing positive results:

Aspergillus flavus Link 1809 cepa *MUCL 19006*

Aspergillus fumigatus Fresenius 1863 cepa DSM 819

Aspergillus terreus var. *terreus* Thom 1918 cepa DSM 1958

Compatibles real-time PCR equipment

Vitassay qPCR Aspergillosis has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- LightCycler® 480 Instrument II (Roche)
- QuantGene 9600 (Bioer Technology)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- CFX96™ Real-Time PCR Detection System (Bio-Rad).
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- DTIite Real-Time PCR System (DNA-Technology)
- NEOS-96 Real-Time PCR System (Linear)

¹: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from respiratory samples (bronchoalveolar lavages, bronchoalveolar aspirates and sputum). The use of other samples has not been established.
- The correct test performance depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by Aspergillosis Positive Control during its reconstitution, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA isolation); c) DNA degradation during specimen shipping, storage, and/or preparation; d)

potential mutations of the target regions of the *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* covered by this test which may result in DNA being undetectable e) pathogen load below the limit of detection for the assay; f) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); g) failure to follow the manufacturer's instructions and procedures.

- Detection of DNA may not indicate the presence of viable and/or infectious bacteria or that these are the causative agents for clinical symptoms.
- Negative results do not preclude *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *Aspergillus flavus*, *Aspergillus terreus* and/or *Aspergillus fumigatus* infection, and other respiratory illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽⁴⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽⁴⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽⁴⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ^{(3) (4)}	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER ⁽⁶⁾
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyclers	DTprime Real-time Detection Thermal Cyclers ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Qiagen	Qiagen
Rotor-Gene® Q	QIAquant 96

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.


Bibliography/Bibliografía

- Cadena, J., Thompson, G. R., 3rd, & Patterson, T. F. (2021). Aspergillosis: Epidemiology, Diagnosis, and Treatment. *Infectious disease clinics of North America*, 35(2), 415–434.
- El-Baba, F., Gao, Y., & Soubani, A. O. (2020). Pulmonary Aspergillosis: What the Generalist Needs to Know. *The American journal of medicine*, 133(6), 668–674.
- Griffiths, L. J., Anyim, M., Doffman, S. R., Wilks, M., Millar, M. R., & Agrawal, S. G. (2006). Comparison of DNA extraction methods for *Aspergillus fumigatus* using real-time PCR. *Journal of medical microbiology*, 55(Pt 9), 1187–1191.
- Hedayati, M. T., Pasqualotto, A. C., Warn, P. A., Bowyer, P., & Denning, D. W. (2007). *Aspergillus flavus*: human pathogen, allergen and mycotoxin producer. *Microbiology (Reading, England)*, 153(Pt 6), 1677–1692.
- Lass-Flörl, C., Dietl, A. M., Kontoyiannis, D. P., & Brock, M. (2021). *Aspergillus terreus* Species Complex. *Clinical microbiology reviews*, 34(4), e0031120.
- Latgé J. P. (1999). *Aspergillus fumigatus* and aspergillosis. *Clinical microbiology reviews*, 12(2), 310–350.
- Powers-Fletcher, M. V., & Hanson, K. E. (2016). Molecular Diagnostic Testing for *Aspergillus*. *Journal of clinical microbiology*, 54(11), 2655–2660.
- van de Veerdonk, F. L., Gresnigt, M. S., Romani, L., Netea, M. G., & Latgé, J. P. (2017). *Aspergillus fumigatus* morphology and dynamic host interactions. *Nature reviews. Microbiology*, 15(11), 661–674.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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