Vitassay qPCR

MTBC/NTM

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Intended use

Vitassay qPCR MTBC/NTM allows the detection and differentiation of the *M. tuberculosis* complex (MTBC), the genus mycobacterium, and the specie *M. tuberculosis* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of the *M. tuberculosis* complex (MTBC), the genus mycobacterium, and the specie *M. tuberculosis* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR MTBC/NTM 4x8-well strip, low profile 7041048

Vitassay qPCR MTBC/NTM 4x8-well strip, high profile 7042048

Materials/reagents provided

Reference	Reagent/Material	Colour	Quantity
7041S048/ 7042S048	MTBC/NTM strips low/high profile	-	4 x 8-well strip
7C048	MTBC/NTM Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Mycobacterium tuberculosis infections constitute a major public health problem due to extensive antibiotic resistance in some strains. Mycobacteria are subdivided into three groups: The *Mycobacterium tuberculosis* complex that produces tuberculosis, the non-tuberculous mycobacteria called NTM or MOTT (Mycobacteria Other Than Tuberculosis) that can cause a broad spectrum of diseases and *Mycobacterium leprae* that produces leprosy. The differentiation of *Mycobacterium tuberculosis* complex from NTM is of primary importance for infection control and choice of antimicrobial therapy.

Tuberculosis (TB) is an infectious chronic respiratory disease caused by *Mycobacterium tuberculosis* complex (MTBC). The bacteria usually attack the lungs (pulmonary TB), but TB bacteria can also attack any part of the body such as the kidney, spine, and brain. TB disease in the lungs or throat can be infectious. This means that the bacteria can be spread to other people. TB in other parts of the body, such as the kidney or spine, is usually not infectious.

TB bacteria are spread through the air from one person to another. When a person with TB disease of the lungs or throat coughs, sneezes, spits, speaks, or sings, they propel the TB germs into the air. People nearby may breathe in the TB germs and become infected. Once in the body, the bacteria can settle in the lungs and begin to grow. From there, they can move through the blood to other parts of the body, such as the kidney, spine, and brain. People with TB disease are most likely to spread it to people they spend time with every day. This includes family members, friends, and coworkers or schoolmates.

Not everyone infected with TB bacteria becomes sick. TB bacteria can live in the body without making you sick. This is called latent TB infection. People with latent TB infection have no symptoms, do not feel sick and cannot spread TB bacteria to others. In most people who breathe in TB bacteria and become infected, the body can fight the

bacteria to stop them from growing. As a result, two TB-related conditions exist: latent TB infection (LTBI) and TB disease.

Tuberculosis mostly affects adults in their most productive years. However, all age groups are at risk. According to the WHO report, almost one third of the worldwide population has been infected by *M. tuberculosis*. However, only 10% of infected individuals develop an active disease state with the appearance of clinical symptoms, suggesting that the immune system can control the infection and prevent active disease in most of the population. People with active TB can infect 5–15 other people through close contact over the course of a year. About one-quarter of the world's population has latent TB.

For people, whose immune systems are weak, especially those with HIV infection, the risk of developing TB disease is much higher than for people with normal immune systems. Without proper treatment, 45% of HIV-negative people with TB on average and nearly all HIV-positive people with TB will die. The risk of active TB is also greater in persons suffering from other conditions that impair the immune system like undernutrition.

Symptoms of TB disease depend on where in the body the TB bacteria are growing. Common symptoms of TB disease include weakness or fatigue, chills, lack of appetite, weight loss, fever and night sweats. TB disease in the lungs may cause symptoms such as chest pain, prolonged cough and coughing up blood or sputum. Symptoms of TB disease in other parts of the body depend on the area affected.

Early, rapid and accurate identification of *M. tuberculosis* and the determination of drug susceptibility is essential for the treatment and management of this disease.

Principle of the test

Vitassay qPCR MTBC/NTM test is based on the real-time amplification of a conserved region of the 16S rRNA for the identification of the genus mycobacterium, the amplification of the insertion sequences IS6110 and IS1081 for the identification of *M. tuberculosis* complex and the amplification of a fragment of the TbD1 region for the identification of *M. tuberculosis* specie. After DNA isolation, the presence of the *M. tuberculosis* complex (MTBC), the genus mycobacterium, and the specie *M. tuberculosis* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is

proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR MTBC/NTM test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of two regions of the Insertion sequences IS6110 and IS1081 are detected through the FAM channel, a fragment of the 16S rRNA is detected through the ROX channel and a fragment of TbD1 region is detected through the Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex II).

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.

- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR MTBC/NTM has been tested on mycobacterial culture, clinical isolates of mycobacteria and decontaminated sputum (with positive or negative acid-fast bacilli (AFB) smear). Other sample types must be validated by the user.

Overall, respiratory samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 2 hours or at 4°C for up to 24 hours. For long term transport (more than 24 hours), we recommend shipping at -20°C or lower. The samples can be stored at 4°C for up to 24 hours or frozen at -20°C or lower (at -80°C ideally) for long-term conservation. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The respiratory specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- CDC (Specimen collection guidelines. Website: <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>)
- IDSA (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94)
- García-Lechuz Moya JM, et al. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Procedimientos en

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For sputum samples, the recommended amount of sputum is approximately 200µL that can be pipetted.

DNA extraction

For nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

MagDEA Dx SV kit, using magLEAD® 12gC (Precision System Science Co.).

To improve the performance and quality of bacterial DNA extracted from sputum samples, pretreatment of the sample with N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) is recommended. After decontamination, the cell pellets must be resuspended in a maximum final volume of 1 ml of phosphate buffer. Furthermore, the use of small elution volumes (50-100 µL) can increase the DNA concentration.

Positive control preparation

Reconstitute the lyophilized MTBC/NTM Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative control (yellow tube) or reconstituted positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Insertion sequences IS6110 and IS1081 - *M. tuberculosis* complex), ROX (16S rRNA - genus mycobacterium), Cy5 (TbD1 region - *M. tuberculosis* specie) and HEX, JOE, or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run, must show an amplification curve ($Ct \leq 40$) for FAM, ROX and Cy5, which validates the reaction.

Negative control

The negative control included in each run, must show the signal absence ($Ct > 40$ or no signal) in FAM, ROX and Cy5 which validates the reaction.

The Internal Control (IC) should show an amplification signal ($Ct \leq 40$) in positive and negative controls wells.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

Note that to consider a positive result, the fluorescence level must be at least 20% of the total fluorescence value.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Sequences IS6110 and IS 1081 (FAM)	16S rRNA (ROX)	TbD1 Region (Cy5)	Internal Control (HEX)	Interpretation	
+	+	+	+/- ¹	Valid	<p><i>M. tuberculosis</i> specie Positive (If Ct FAM < Ct ROX) ³</p> <p>or NTM Positive (If Ct ROX < Ct FAM and Cy5) ³</p>
+	- #	+	+/- ¹	Valid	<p><i>M. tuberculosis</i> specie Positive (If Ct FAM > 35) ³</p> <p>If Ct FAM < 35, repeat the assay depending on the available material:</p> <ul style="list-style-type: none"> a) obtain a new specimen, re-extract, and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample. <p>After retesting one time, if the result is maintained, the sample should be considered positive for <i>M. tuberculosis</i> specie</p>
+	- #	- #	+/- ¹	Valid	<p>MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct FAM > 35) ³</p> <p>If Ct FAM < 35, repeat the assay depending on the available material:</p> <ul style="list-style-type: none"> a) obtain a new specimen, re-extract, and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat qPCR with the same isolated DNA sample.

					After retesting one time, if the result is maintained, the sample should be considered positive for MTBC strains different from <i>M. tuberculosis</i>
+	+	-	+/- ¹	Valid	MTBC strains different from <i>M. tuberculosis</i> Positive (If Ct FAM = Ct ROX and/or if Ct FAM < Ct ROX) ³ or NTM Positives (If Ct FAM > Ct ROX) ³
-	+ [†]	-	+/- ¹	Valid	NTM Positive (In case Cy5 amplifications appear with a Ct > 37, the sample should be considered as NTM Positive)
-	-	-	+ ²	Valid	Targets not detected ²
-	-	-	- ²	Invalid	Test failed ²

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct >40 or no signal)

NTM = non-tuberculous mycobacteria

MTBC= *M. tuberculosis* complex

†: Due to environmental NTM, amplification curves with late Ct may appear in the ROX channel.

#: Due to the amplified targets, the FAM signal is stronger than the ROX or ROX and Cy5 signals. Preferential amplification curve can be observed in this channel and absence in ROX channel or ROX and Cy5. This result also indicates the presence of a different MTBC strain from the *M. tuberculosis* or *M. tuberculosis* species.

¹ Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

² In the case of negative genus *Mycobacterium*, the *M. tuberculosis* complex (MTBC) and the specie *M. tuberculosis* target genes detection, IC must show an amplification signal with Ct ≤35. If there is a signal' absence or Ct value > 35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR MTBC/NTM was evaluated with DNA extracted from 52 clinical strains (12 NTM strains, 13 MTBC strains, 25 *M. tuberculosis* strains of different lineages and sublineages, and 2 samples with mixed strains). Vitassay qPCR MTBC/NTM detected 14 positive NTM samples, 15 positive MTBC samples and 24 positive *M. tuberculosis* samples. The results obtained with Vitassay qPCR MTBC/NTM were compared with those obtained with the reference methods in the diagnosis of active tuberculosis (culture) and are shown below:

Target	TP	TN	FP	FN	SE	SP
NTM	14	38	0	0	1 (0.76-1)	1 (0.9-1)
MTBC	15	37	0	0	1(0.78-1)	1(0.9-1)
<i>M. Tuberculosis</i>	24	27	0	1	0.96 (0.79-0.99)	1 (0.87-1)

In addition, Vitassay qPCR MTBC/NTM was also evaluated with 130 decontaminated sputum samples (49 from patients with bacterial respiratory infection, 30 from patients with respiratory NTM infection and 51 from patients with tuberculous infection). Vitassay qPCR MTBC/NTM detected 22 NTM positive samples, 52 MTBC positive samples and 47 *M. tuberculosis* positive samples. These results were compared with those obtained

with the reference method consisting of characterization with culture and identification of NTM species with INNO-LiPA® MYCOBACTERIA v2 (Fujirebio) in patients with NTM respiratory infection, and in patients diagnosed with TB, MTBC was identified with SD BIOLINE TB Ag MPT64 Rapid test (Abbott). The obtained results are shown below:

Target	TP	TN	FP	FN	SE	SP
NTM	22	91	9	8	0.73 (0.54-0.87)	0.91(0.83-0.95)
MTBC	52	77	0	1	0.98 (0.89-1)	1(0.95-1)
<i>M. Tuberculosis</i>	47	78	1	4	0.92(0.81-0.97)	0.98 (0.93-1)

These results show the high sensitivity and specificity of the molecular diagnostic kit Vitassay qPCR MTBC/NTM.

Analytical sensitivity

The analytical sensitivity was determined from serial dilutions (1:10) of the template DNA of the different targets (10^7 - 10^1 copies/reaction). This assay has a detection limit of 10 DNA copies per reaction for the 16S rRNA region, the *IS6110* and *IS1081* insertion sequences and the TbD1 region (positivity rate $\geq 95\%$).

Analytical specificity

The analytical specificity for MTBC/NTM detection was confirmed by testing a panel composed of different microorganisms. No cross-reactivity was observed between any of the species:

Cross-reactivity assay		
<i>Bordetella pertussis</i>	<i>Mycoplasma pneumoniae</i>	Influenza A/Hong Kong/4801/2014, NYMC X-263B
<i>Bordetella parapertussis</i>	Human Adenovirus 1-5, 8, 15, 31, 40 and 41	Influenza A/Michigan/45/2015 (H1N1) pdm09 virus
<i>Bordetella holmesii</i>	Human Bocavirus	Influenza A/New Caledonia/20/99(H1N1) virus
<i>Bordetella bronchiseptica</i>	Human Coronavirus 229E, OC43 and NL63	Influenza A/Victoria/210/2009(H3N2) virus
<i>Chlamydia caviae</i>	Human Metapneumovirus A and B	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1) pdm09 virus
<i>Chlamydophila pneumoniae</i> CM-1	Human parainfluenza Virus 1, 2, 3 and 4	Influenza A/South Australia/55/2014, IVR-175)
<i>Chlamydia psittaci</i> genotypes A and C	Human Rhinovirus	Influenza A/Switzerland/9715293/2013 (H3N2) virus
<i>Haemophilus influenzae</i> MinnA	MERS Coronavirus	Influenza A/Thüringen/5/17 (H3N2) virus (clade 3C2a.1)
<i>Legionella bozemanii</i>	<i>Streptococcus pneumoniae</i> Z022	Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus
<i>Legionella micdadei</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza B/Brisbane/60/2008 virus
<i>Legionella dumoffii</i>	Respiratory Syncytial Virus (RSV) type A and B	Influenza B/Florida/04/06 virus
<i>Legionella pneumophila</i>	Influenza A/Anhui/1/2013 (H7N9) virus	Influenza B/Phuket/3073/2013 virus
<i>Legionella longbeachae</i>	Influenza A/California/7/2009(H1N1) pdm09 virus	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus
<i>Moraxella catarrhalis</i>		

Analytical reactivity

Vitassay qPCR MTBC/NTM kit's analytical reactivity for the specie *M. tuberculosis* was evaluated against DNA extracted from reference strains (*Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331 and *Mycobacterium tuberculosis* X004439) and from clinical strains containing *M. tuberculosis* lineages 1, 2, 3, and lineage 4 (T sublineage, LAM sublineage, and X sublineage), showing positive results.

Vitassay qPCR MTBC/NTM kit's analytical reactivity for the *Mycobacterium* genus was evaluated against DNA extracted from reference strains (*Mycobacterium bovis* AF

2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium africanum*, *Mycobacterium microti*, *Mycobacterium abscessus*, *Mycobacterium avium* subsp. *avium*, *Mycobacterium chelonae*, *Mycobacterium fortuitum*, *Mycobacterium gordonae*, *Mycobacterium kansasii*, *Mycobacterium marinum*, *Mycobacterium peregrinum*, *Mycobacterium smegmatis*, *Mycobacterium scrofulaceum*, *Mycobacterium xenopi*), from clinical strains (*Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), *Mycobacterium africanum* lineage 5 and 6, *Mycobacterium avium* complex (MAC), *Mycobacterium avium* subsp. *paratuberculosis*, *Mycobacterium avium* subsp. *silvaticum*, *Mycobacterium chelonae* complex (group III- *M. abscessus*), *Mycobacterium intracellulare*) and from clinical samples (*Mycobacterium celatum*, *Mycobacterium chimaera-intracellulare*), showing positive results.

Vitassay qPCR MTBC/NTM kit's analytical reactivity for MTBC strains was evaluated against DNA extracted from reference strains (*Mycobacterium bovis* AF 2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium africanum* and *Mycobacterium microti*), as well as from clinical strains (*Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), and *Mycobacterium africanum* lineage 5 and 6), showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR MTBC/NTM has been validated on the following equipment:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)

¹: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from mycobacterial culture, clinical isolates of mycobacteria and decontaminated sputum (with positive or negative acid-fast bacilli (AFB) smear). The use of other samples has not been established.

- The correct test performance depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination, either by the MTBC/NTM Positive Control during its reconstitution, by samples containing high concentrations of target DNA or by PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA extraction); c) DNA degradation during specimen shipping, storage, and/or preparation; d) pathogen load below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown variants of *M. tuberculosis* complex (MTBC), the genus *Mycobacterium*, and the specie *M. tuberculosis*.
- DNA detection may not indicate the presence of viable and/or infectious bacteria or that these bacteria are the causative agents for clinical symptoms.
- The negative results do not preclude *M. tuberculosis* complex (MTBC), the genus mycobacterium, and the specie *M. tuberculosis* infection and should not be the sole basis of a patient treatment/management decision. Optimum specimen types for these pathogen's identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *M. tuberculosis* complex (MTBC), the genus *Mycobacterium*, and the specie *M. tuberculosis* infection, and other respiratory diseases have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type, and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Excycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Excycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™ 5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cycler	DTprime Real-time Detection Thermal Cycler ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.



Bibliography/Bibliografía

1. Wenping Gong, Yan Liang, and Xueqiong Wu. (2018). The current status, challenges, and future developments of new tuberculosis vaccines. *Hum Vaccin Immunother.* 14(7): 1697–1716..
2. World Health Organization. Tuberculosis. https://www.who.int/health-topics/tuberculosis#tab=tab_2 <https://www.who.int/news-room/fact-sheets/detail/tuberculosis>.
3. Centers for Disease Control and Prevention. <https://www.cdc.gov/tb/>
4. Emilia Cercenado y Rafael Cantón. (2005). *Procedimientos en Microbiología Clínica. Recomendaciones de la Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica.* 9a. Micobacterias.
5. M M Rahman, M R Rahim, A Khaled, T A Nasir, F Nasrin, M A Hasan. (2017). Molecular Detection and Differentiation of Mycobacterium Tuberculosis Complex and Non-tuberculous Mycobacterium in the Clinical Specimens by Real Time PCR. *Jul;26(3):614-620.*
6. Röltgen K, Pluschke G, Spencer JS, Brennan PJ, Avanzi C. (2020). The immunology of other mycobacteria: *M. ulcerans*, *M. leprae*. *Semin Immunopathol.* Feb 25;1-21.
7. María Cebriá-Mendoza, Rafael Sanjuán, Pilar Domingo-Calap. (2019). Directed Evolution of a Mycobacteriophage. *Antibiotics (Basel).* Apr 25;8(2):46.
8. Diana Machado, Isabel Couto, Miguel Viveiros. (2019). Advances in the Molecular Diagnosis of Tuberculosis: From Probes to Genomes. *Infect Genet Evol.* Aug; 72:93-112.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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