VA Vitassay

Real-Time PCR Kits

Vitassay qPCR *Pneumocystis jirovecii*

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Programa EQA	Variación del patógeno
INSTAND EQA programme "Fungal Genome Detection - <i>Pneumocystis jirovecii</i> " panel (November 2016, June 2017, November 2017, June 2018, May 2020)	<i>Pneumocystis jirovecii</i>
QCMD <i>Pneumocystis jirovecii</i> pneumonia (PCP) DNA EQA Programme (2016, 2017, 2018, 2019, 2020, 2022, 2023)	
QCMD <i>Pneumocystis jirovecii</i> pneumonia (PCP) DNA EQA Programme (2016, 2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> Type A1
QCMD <i>Pneumocystis jirovecii</i> pneumonia (PCP) DNA EQA Programme (2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> g885652
QCMD <i>Pneumocystis jirovecii</i> pneumonia (PCP) DNA EQA Programme (2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> j888023

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- Cobas Z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- LightCycler® 480 Real-Time PCR System (Roche Diagnostics)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR instrument (DNA-Technology)
- DTPrime Real Time PCR instrument (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I
- CFX96 Opus 96™ Real-Time PCR Detection System (Bio-Rad)

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Intended use

Vitassay qPCR *Pneumocystis jirovecii* allows the qualitative and quantitative detection of *Pneumocystis jirovecii* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Pneumocystis jirovecii* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR <i>Pneumocystis jirovecii</i> 4x8-well strip, low profile	7041026
Vitassay qPCR <i>Pneumocystis jirovecii</i> 4x8-well strip, high profile	7042026

Provided Reagents

In references 7041026 and 7042026:

Reference	Reagent/material	Description	Concentration Range	Color	Amount
7041S026/ 7042S026	Pneumocystis <i>jirovecii</i> strips low/high profile	Lyoprotectors and Stabilizers	± 6 g/100 mL *	-	4 x 8-well strips
		dNTPs	± 1 mM *		
		Primers and Probes	0,2-1 nMol/ μ L *		
		Enzymes	10-100 U/rxn *		
7C026	Pneumocystis <i>jirovecii</i> Positive Control	Non-infectious synthetic lyophilized DNA	$1,9 \times 10^4$ copies/ μ L *	red	1 vial
7Q026	<i>Pneumocystis jirovecii</i> Quantitative standard	Quantitative Standard	2×10^7 copies/ μ L *	red	1 vial
7001A	PCR grade water	RNAsa/DNAsa free water	1 g/mL	white	1 vial x 1 mL
7002B	Resuspension buffer	Saline Solution	± 13 mM	green	1 vial x 1,8 mL
		Buffer (TRIS, pH)	± 67 mM		
7003N	Negative control	RNAsa/DNAsa free water	1 g/mL	yellow	1 vial x 1 mL
7004O	Optical caps	Optical caps for sealing strips during thermal cycling	-	-	4 x 8-cap strips

* For component in stabilized format, the concentration range means after rehydration.

Material and equipment required, not provided

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or \leq -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 μ L, 20-200 μ L)
- Filter tips
- Powder-free disposal gloves
- Blank sample (See DNA extraction)

Transport and storage

- Transport and storage of the kits (before use) can be performed at 2-40°C until the expiry date stated on the label. Once in use, see storage conditions for each component. Avoid shaking during transport to prevent liquid leakage.
- The resuspended positive control and Quantitative Standard should be stored at -20°C until the expiry date stated on the kit label. The stability of the positive control and Quantitative Standard has been validated after 6 freeze-thaw cycles. However, it is recommended to distribute in aliquots to avoid repeated freeze-thaw cycles.
- For the quantitative standard, freeze-thaw cycles are only applicable for the reconstituted vial and its aliquots. For proper use of the product, dilutions of the standard curve must be discarded after each use.
- For proper use of the kit, it is recommended that unneeded wells are not rehydrated. Unused wells can be easily cut before removing the protective foil and stored in the foil bags provided with the silica gel inside and closed with the zip closure to avoid contact with light or moisture and stored at 2-40°C for the shelf life indicated on the kit label. Once the wells have been rehydrated, nucleic acid addition and qPCR should be performed immediately.
- Once opened, the Resuspension buffer, PCR grade water and Negative Control components should be stored at 2-40°C for the shelf life indicated on the kit label.
- Keep all reagents in the darkness.

Summary

Pneumocystis jirovecii (PJ) is an opportunistic fungal pathogen that can cause severe pneumonia in immunocompromised hosts. Risk factors for *Pneumocystis jirovecii* pneumonia (PJP) include HIV, organ transplant, malignancy, certain inflammatory or rheumatologic conditions, and associated therapies and conditions that result in cell-mediated immune deficiency (Weyant et al. 2021).

The signs and symptoms of this PJP are non-specific and include progressive dyspnea, low-grade fever, non-productive cough, hypoxemia, as well as diffuse dry rales on examination (Weyant et al. 2021).

The cyst form of *Pneumocystis* is spread through the air from person-to-person. In the host the cysts are received in the lungs, and in most cases their presence is asymptomatic. Cases of PJP can be due to reactivation of latent infection or de-novo person-person transmission (especially in outbreaks). This high rate of colonization and asymptomatic transmission allows for the immunocompetent to act as unknowing reservoirs for the pathogen (Weyant et al. 2021).

The prevalence of PJP has significantly decreased in people living with HIV thanks to combined antiretroviral therapy (cART) and the standard practice of prescribing prophylaxis. Even with treatment, the mortality in PJP is high, and it is higher in non-HIV patients (20–50%) than in people living with HIV (10–30%) (Weyant et al. 2021). Non-HIV patients may present with a rapidly progressing disease likely due to an intense immune response to PJP in the lungs than in HIV patients (Ibrahim et al. 2023). Although most patients will improve with appropriate treatment, some develop progressive respiratory failure. The mortality among patients requiring intensive care admission or mechanical ventilation is reported to be as high as 60% (Truong et al. 2023).

The differential diagnosis of PJP includes acute respiratory distress syndrome, tuberculosis, cytomegalovirus (CMV), *Legionella* pneumonia, mycoplasma pneumonia, other viral or bacterial pneumonias, pulmonary embolism, lymphocytic interstitial pneumonia, and coronavirus disease 2019 infection (Ibrahim et al. 2023). As the clinical signs of PJP are nonspecific, definitive diagnosis requires direct detection of the organism in lower respiratory secretions or tissue (Weyant et al. 2021).

Workup of PJP typically involves a combination of pulmonary imaging, biochemical tests, and appropriate respiratory samples for analysis (Weyant et al. 2021). Traditional diagnostic testing has relied on staining and direct visualization of the life-forms in bronchoalveolar lavage fluid. This method has proven insensitive, and invasive procedures may be needed to obtain adequate samples. Molecular methods of detection such as polymerase chain reaction (PCR), loop-mediated isothermal amplification (LAMP), and antibody-antigen assays have been developed in an effort to solve these problems. These techniques are very sensitive and have the potential to detect *Pneumocystis* life-forms in non-invasive samples such as sputum, oral washes, nasopharyngeal aspirates, and serum (Bateman et al. 2020).

Test principle

Vitassay qPCR *Pneumocystis jirovecii* is based on the real-time amplification of a conserved region of the mitochondrial *large-subunit* ribosomal RNA (*mtLSU* rRNA) gene of *Pneumocystis jirovecii* genome. After DNA extraction, *P. jirovecii* is detected by an increase in observed fluorescence during the reaction, following hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR

equipment. Additionally, quantification of *P. jirovecii* DNA can be performed by creating a standard curve from the Pneumocystis jirovecii Quantitative Standard included in the kit.

Vitassay qPCR Pneumocystis jirovecii is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition and checks whether the amplification reaction is working correctly. The detection channels of the target sequences are described as follows:

Target	Detection channel
<i>mtLSU</i> rRNA gene (<i>P. jirovecii</i>)	FAM
Internal Control (IC)	HEX, VIC or JOE (depending on the equipment used)

Precautions

- For professional *in vitro* diagnostic use.
- Do not use the kit after the expiration date.
- Do not use the kit without having read and understood the information on procedures, precautions, and limitations provided in the instructions for use.
- Do not mix components from different kits and/or lots.
- Do not use if the package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Use the strips immediately after removing the protective foil to protect the reaction mixture from sunlight. It is not recommended to leave the strips open without resuspending. Once the required wells have been rehydrated, the qPCR assay should be performed immediately.
- It is recommended to cut the required/necessary wells before removing the protective seal from the entire strip. Store the rest inside the bag with the desiccant.
- Close protective pouches of reagents promptly with the zip seal after each use (do not remove the desiccant from the pouches). Remove excess air from the bags before sealing.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Rehydrate the vials immediately after opening the Positive Control and Quantitative Standard envelope and use them in the PCR reaction or store them properly. Leaving the vial open without resuspending may cause deterioration of the reagent.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area. Once the test is completed wash your hands thoroughly.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended. A new tip should be used for each sample. It is necessary to change gloves before handling the reagents.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- The product Vitassay qPCR Pneumocystis jirovecii has only been validated with the equipment mentioned in Section Compatibles real-time PCR equipment of this Instructions for Use.
- You must ensure that the conditions are set up in the thermal cycler according to the instructions in section "Programme your thermocycler".
- Refer to Safety data sheets, available on request.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport the conditions validated by the user must be followed. Vitassay qPCR *Pneumocystis jirovecii* has been tested in respiratory samples (bronchoalveolar lavage fluid, bronchoalveolar aspirates, and sputum). Other sample types must be validated by the user.

Overall, clinical samples should be collected properly in clean containers with or without transport media (depending on sample type), labelled correctly and be processed as soon as possible to guarantee the test quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 2h or at 4°C for up to 7 days and for long term transport (more than 7 days), it is recommended shipping at -20°C or lower. The samples can be stored at 4°C for up to 7 days or frozen at -20°C or lower (ideally at -80°C) for conservation during a long time. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

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DNA extraction

For pre-treatment and nucleic acid isolation from clinical samples, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use. The assay has been validated with the following extraction kits:

- EZ1 Virus Mini Kit, using the EZ1 instrument (Qiagen)
- Invisorb® Spin Universal Kit (Invitex).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega), using the Maxwell® 16 instrument (Promega).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

For the use of DNA extraction kits other than those mentioned above, validation must be performed by the end user, following the manufacturer's instructions for use.

It should be noted that the positive and negative control should not be extracted. To monitor and control the extraction process and rule out possible contamination, a sample blank (not supplied in the kit) can be extracted, which should consist of the same matrix as the clinical samples (real or simulated matrix), and previously characterized as negative for all targets.

Positive control and quantitative standard preparation

Reconstitute the lyophilized *Pneumocystis jirovecii* Positive Control (red tube) with 400 µL of PCR grade water (white tube) and the lyophilized *Pneumocystis jirovecii* Quantitative Standard (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tubes thoroughly.

Once the *Pneumocystis jirovecii* Positive Control and the *Pneumocystis jirovecii* Quantitative Standard have been resuspended, dispense them into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C. For the quantitative standard, freeze-thaw cycles are only applicable for the reconstituted vial and its aliquots (dilutions of the standard curve should be discarded after each use).

Both components contain a high-copy number template and entail a very significant contamination risk. Therefore, we recommend opening and manipulating them in a separate area of the laboratory away from the other components and samples.

Standard Curve Preparation (quantitative assay)

This kit can be used to quantify the DNA copy number of *Pneumocystis jirovecii*. To perform a quantitative assay, a standard curve must be prepared by serial dilution of *Pneumocystis jirovecii* Quantitative Standard (red tube, which contains approximately 2×10^7 copies/µL*). The contents of this vial will be the standard with the highest concentration from which the rest of the serial dilutions will be prepared as indicated below:

- Pipette 90 µL of PCR grade water into 6 microcentrifuge tubes.
- Add 10 µL of *Pneumocystis jirovecii* Quantitative Standard to the first tube to get a standard with about 2×10^6 copies/µL. Mix by vortex and centrifuge briefly.
- Add 10 µL of standard with 2×10^6 copies/µL to the second tube to get a standard with approximately 2×10^5 copies/µL. Mix by vortex and centrifuge briefly.
- Repeat the previous step sequentially to complete the dilution series for standards with about 2×10^4 , 2×10^3 , 2×10^2 and 2×10^1 copies/µL.

* For further details refer to the “Certificate of Analysis”, for the lot number and DNA copy number of *Pneumocystis jirovecii* Quantitative Standard.

For qualitative assays, the preparation of dilutions is not necessary and it is recommended to use *Pneumocystis jirovecii* Positive Control to minimize the risk of cross-contamination.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
It is recommended not to rehydrate wells that are not needed. Unused wells can be easily cut before removing the protective foil and stored as indicated in the Transport and Storage Section.
- Peel off the protective aluminum seal from the strips to be used.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA extracted from each sample, of negative control (yellow tube) and of reconstituted positive control (red tube) for qualitative assays or dilutions for the standard curve (for quantitative assays), and add them into the corresponding wells.
- Cover the wells with the optical caps provided. It is recommended to spin down briefly.
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the annealing/extension step (*) through the FAM (*Pneumocystis jirovecii*) and HEX, JOE, or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel.

Analysis and interpretation of results

The analysis of the results is performed by the software itself of the used real-time PCR system following manufacturer’s instructions. To verify the correct operation of the amplification mix, check the internal control signal emission.

For each assay, it is recommended to set the threshold values for each channel (target) independently by the end-user. To do so, select the well corresponding to the positive control for each channel and set the threshold value within the exponential phase of the fluorescence curve and above any background signal (below the baseline). In case a threshold is set automatically by the equipment used, check, and verify that it is adjusted to the positive control or adjust it manually. Once the threshold value is set, the test samples can be interpreted.

The chemical nature of the different fluorophores can affect the threshold value of the channels, just as different signal intensities can influence the threshold value for different instruments.

In addition, it is recommended to include a blank (confirmed negative sample from the same matrix as the samples tested) for the detected targets in order to establish the baseline of the assay.

For a valid diagnostic test run, before analysing the result of the clinical samples, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve (Ct ≤40) in FAM channel (*Pneumocystis jirovecii*), which validates the reaction.

Negative control

The negative control included in each run must show the absence of signal (Ct >40 or no signal) in FAM (*Pneumocystis jirovecii*), which validates the reaction.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in the positive control and negative control wells.

The experiment is invalid if there is no signal in the positive control or there is amplification signal in the negative control, except for the HEX, JOE or VIC channel which should show amplification signal (corresponding to Internal Control). In either case, the assay must be repeated.

Qualitative analysis

Once the results of the controls have been validated, for the interpretation of the sample results **select only the channels where the targets are detected**. The result interpretation is summarized in the following table:

Negative Control	Positive Control	<i>Pneumocystis jirovecii</i> FAM	Internal Control HEX	Interpretation
-	+	+	+/- ¹	<i>Pneumocystis jirovecii</i> DNA detected
-	+	-	+ ²	<i>Pneumocystis jirovecii</i> DNA Not detected ²
-	+	-	- ²	Invalid ²

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct >40 or no signal)

¹ Sometimes, in the case of clinical samples, detection of the internal control is not necessary as the presence of a high initial copy number of the target nucleic acid may cause preferential amplification of the target nucleic acid. The internal control (IC) either shows or does not show an amplification signal (Ct ≤40 or no signal).

² In case the *Pneumocystis jirovecii* target genes are negative, the Internal Control should show an amplification signal with Ct ≤ 35. In case of no signal or Ct value > 35 of the internal control, the result is considered "invalid" and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1: 100, or to re-extract and repeat the assay to check for possible failure of the extraction procedure and/or inhibition problems.

If the result obtained is unclear or doubtful, it is necessary to check that all steps have been performed correctly, verify the correct performance of each qPCR step, check all parameters, the sigmoid shape of the curve and the fluorescence intensity. If the situation is not resolved, it is recommended to repeat the test preferably in duplicate depending on the available material: (a) obtain a new specimen, re-extract and re-test (ideal), (b) re-extract and re-test another aliquot of the same sample, or (c) repeat the qPCR with the same isolated DNA sample.

The test results should be evaluated by a healthcare professional, in the context of the medical history, clinical symptoms and other diagnostic tests.

Quantitative analysis

Using the serial dilutions of *Pneumocystis jirovecii* Quantitative Standard a standard curve can be created according to the formula:

$$\text{Ct} = m \log (\text{Q}) + b \quad \text{Where:} \quad \begin{array}{l} \text{Ct} = \text{Threshold Cycle} \\ m = \text{Slope} \\ \text{Q} = \text{Concentration} \\ b = \text{Intercept} \end{array}$$

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{-(Ct-b)/m}$$

The concentration of the positive samples is displayed in copies/μL and refers to the concentration in the Eluted DNA not to the original clinical sample. To determine the target concentration of the original sample, take into account the dilutions of the extraction procedure and PCR Set-up.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

In a first assessment, the clinical performance of Vitassay qPCR *Pneumocystis jirovecii* was evaluated using 14 bronchoalveolar lavage and sputum from patients with suspicious of *Pneumocystis pneumonia*. We obtained a positive result for 11 samples and a negative result for 3 of them. The obtained results with Vitassay qPCR *Pneumocystis jirovecii* were compared with other real-time PCR commercial kit (FTD *Pneumocystis jirovecii*, Fast-track Diagnostics). An automated nucleic acid extraction system (EZ1, Qiagen) and the CFX96™ Real-Time PCR Detection System (Bio-Rad) were used for this evaluation. The concordance between both tests was 100% showing a Specificity of 1 (0.29-1) and a Sensitivity of 1 (0.71-1).

A second evaluation was carried out on a total of 712 samples of bronchoalveolar lavages and aspirates and sputum with suspected *Pneumocystis jirovecii* infection, using the MagDEA Dx SV extraction kit with the magLEAD® 12gC (Precision System Science Co.) and the LightCycler® 480 Real-Time PCR System (Roche Diagnostics). The results obtained with Vitassay qPCR *Pneumocystis jirovecii* were 10 positive and 702 negative and were compared with the RealStar® *Pneumocystis jirovecii* PCR kit (Altona Diagnostics), showing a Specificity of 1 (0.99 -1) and a Sensitivity of 1 (0.69-1).

These results show the high sensitivity and specificity to detect *Pneumocystis jirovecii* using the molecular diagnostic kit Vitassay qPCR *Pneumocystis jirovecii*.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of the template DNA of *Pneumocystis jirovecii* (10^7 - 10^1 copies / reaction) in CFX96™ Real-Time PCR Detection System (Bio-Rad) equipment. This assay has a limit of detection of 10 copies of DNA per reaction, with a positive rate of ≥90%.

Veracity

The trueness of Vitassay qPCR *Pneumocystis jirovecii* was evaluated against reference material from different EQA programmes:

EQA Program	Pathogen variety
INSTAND EQA programme "Fungal Genome Detection - <i>Pneumocystis jirovecii</i> " panel (November 2016, June 2017, November 2017, June 2018, May 2020)	<i>Pneumocystis jirovecii</i>
QCMD <i>Pneumocystis jirovecii pneumonia</i> (PCP) DNA EQA Programme (2016, 2017, 2018, 2019, 2020, 2022, 2023)	
QCMD <i>Pneumocystis jirovecii pneumonia</i> (PCP) DNA EQA Programme (2016, 2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> Type A1
QCMD <i>Pneumocystis jirovecii pneumonia</i> (PCP) DNA EQA Programme (2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> g885652
QCMD <i>Pneumocystis jirovecii pneumonia</i> (PCP) DNA EQA Programme (2017, 2018, 2019, 2020)	<i>Pneumocystis jirovecii</i> j888023

Reference strains were analysed with the CFX96™ Real-Time PCR Detection System (Bio-Rad) and AriaMx Real-Time PCR System (Agilent Technologies).

Precision

The accuracy of Vitassay qPCR *Pneumocystis jirovecii* was evaluated by several assays: intra-assay (repeatability), inter-assay (reproducibility), inter-batch assay (intermediate precision) and inter-thermocycler assay (inter-laboratory). The Invisorb® Spin Universal Kit (Invitex) was used for sample extraction.

The repeatability (intra-assay) was assessed by testing replicates of all samples in the same run using Vitassay qPCR *Pneumocystis jirovecii* on the AriaMx Real-Time PCR System (Agilent Technologies). The obtained results were as follows:

Sample	Target	Vitassay channel	\bar{x} (Ct)	σ	CV%
Positive	<i>Pneumocystis jirovecii</i>	FAM	24.63	0.08	0.31
Low Positive	<i>Pneumocystis jirovecii</i>	FAM	28.83	0.11	0.36
LoD sample	<i>Pneumocystis jirovecii</i>	FAM	34.33	1.04	3.03
Negative sample	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	24.96	0.09	0.37
Positive Control	<i>Pneumocystis jirovecii</i>	FAM	25.45	0.11	0.45
Negative Control	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	24.94	0.11	0.44

Ct = threshold cycle, \bar{x} = arithmetic mean of Ct, σ = standard deviation, CV % = coefficient of variation, **Neg** = negative, **NA** = not applicable.

The reproducibility (inter-assay) was assessed by testing the different samples on three different days by three different operators using Vitassay qPCR *Pneumocystis jirovecii* with the AriaMx Real-Time PCR System (Agilent Technologies). The obtained results were as follows:

Sample	Target	Vitassay channel	\bar{x} (Ct)	σ	CV%
Positive	<i>Pneumocystis jirovecii</i>	FAM	32.32	0.56	1.75
Low Positive	<i>Pneumocystis jirovecii</i>	FAM	33.82	0.64	1.89
LoD sample	<i>Pneumocystis jirovecii</i>	FAM	35.84	0.58	1.61
Negative sample	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	25.15	0.32	1.27
Positive Control	<i>Pneumocystis jirovecii</i>	FAM	25.61	0.15	0.59
Negative Control	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	25.23	0.18	0.70

Ct = threshold cycle, \bar{x} = arithmetic mean of Ct, σ = standard deviation, CV % = coefficient of variation, **Neg** = negative, **NA** = not applicable.

The intermediate precision (inter-batch) was assessed by testing three replicates of different samples using three batches of Vitassay qPCR *Pneumocystis jirovecii* on the AriaMx Real-Time PCR System (Agilent Technologies). The obtained results were as follows:

Sample	Target	Vitassay channel	\bar{x} (Ct)	σ	CV%
Positive	<i>Pneumocystis jirovecii</i>	FAM	32.62	0.31	0.96
Low Positive	<i>Pneumocystis jirovecii</i>	FAM	34.33	0.78	2.26
LoD sample	<i>Pneumocystis jirovecii</i>	FAM	35.76	0.10	0.29
Negative sample	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	25.11	0.12	0.46
Positive Control	<i>Pneumocystis jirovecii</i>	FAM	25.49	0.14	0.57
Negative Control	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	25.19	0.29	1.15

Ct = threshold cycle, \bar{x} = arithmetic mean of Ct, σ = standard deviation, CV % = coefficient of variation, **Neg** = negative, **NA** = not applicable.

The inter-laboratory values (inter-thermocycler) were assessed by analysing three replicates of the same sample used for intra-assay, inter-assay, and inter-batch assays, using Vitassay qPCR *Pneumocystis jirovecii*. The DNA extracted from the samples, the positive control and the negative control were analysed together in the same run on the CFX96™ Real-Time PCR Detection System (Bio-Rad), AriaMx Real-time PCR System (Agilent Technologies) and DTprime Real-time Detection Thermal Cycler (DNA-Technology). The obtained results were as follows:

Sample	Target	Vitassay channel	\bar{x} (Ct)	σ	CV%
Positive	<i>Pneumocystis jirovecii</i>	FAM	32.78	1.33	4.07
Low Positive	<i>Pneumocystis jirovecii</i>	FAM	35.22	2.98	8.47
LoD sample	<i>Pneumocystis jirovecii</i>	FAM	35.82	0.04	0.12
Negative sample	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	23.28	2.37	10.20
Positive Control	<i>Pneumocystis jirovecii</i>	FAM	25.97	0.95	3.66
Negative Control	<i>Pneumocystis jirovecii</i>	FAM	Neg	NA	NA
	Internal Control	HEX	23.22	2.43	10.47

Ct = threshold cycle, \bar{x} = arithmetic mean of Ct, σ = standard deviation, CV % = coefficient of variation, Neg = negative, NA = not applicable.

Analytical specificity

The analytical specificity (Cross-Reactivity) for *Pneumocystis jirovecii* detection was confirmed by experimentally testing a panel composed of the following microorganisms, and no cross-reactions were observed between any of the species.

Cross-Reactivity		
Adenovirus Type 1	Coronavirus OC43	<i>Legionella longbeache</i>
Adenovirus Type 2, Strain Adenoid 6	Coronavirus SARS-CoV-22019nCoV/USAWA1/2020	<i>Legionella micdadei</i>
Adenovirus Type 3, strain GB	Coronavirus 229E	<i>Legionella pneumophila</i> Serogroup 1
Adenovirus Type 4, strain RI-67	Enterovirus Coxsackievirus A9	Human Metapneumovirus A1
Adenovirus Type 5	Enterovirus Coxsackievirus A24	Human Metapneumovirus A2
Adenovirus Type 8	Enterovirus Coxsackievirus B3	Human Metapneumovirus B2
Adenovirus Type 15, strain CH. 38	Enterovirus Echovirus 30	<i>Moraxella catharralis</i>
Adenovirus Type 31, strain 1315/63 (NIAID V-231-001-014)	Enterovirus Enterovirus 68	<i>Mycoplasma pneumoniae</i>
Adenovirus Type 40, strain Dugan	Enterovirus Enterovirus 71	<i>Mycobacterium tuberculosis</i> Mycobacterium tuberculosis not rifampin resistant
Adenovirus Type 41, strain Tak (73-3544)	<i>Haemophilus influenzae</i> MinnA	Parainfluenza 1
Rhinovirus type C	Influenza A (H3N2) virus. Strain A/Switzerland/9715293/2013 (vaccine strain)	Parainfluenza 2
<i>Bordetella bronchiseptica</i>	Influenza A (H3N2). Strain A/Thüringen/5/2017 (Clade 3C2a.1). Origin: infected MDCK-cells (lysate)	Parainfluenza 3
<i>Bordetella holmesii</i>	Influenza A (H5N8) virus. Strain A/DE-SH/Reiherente/AR8444/ 2016	Parainfluenza 4
<i>Bordetella parapertussis</i>	Influenza A (H7N9). Strain A/Anhui/1/2013.	Bocavirus
<i>Bordetella pertussis</i>	Influenza A (H1N1 pdm09). Strain A/Michigan/45/2015 (vaccine strain)	RSV A
<i>Chlamydia caviae</i>	Influenza A (H1N1) pdm09- virus. Strain A/California/7/2009/ (vaccine strain).	RSV B

<i>Chlamydia psittaci</i> Genotype A	Influenza A, B and H1N1 A/California/7/2009 (H1N1) ; A/Victoria/210/2009 (H3N2) (Similar strain to A/Perth/16/2009 (H3N2)); B/Brisbane/60/2008	<i>Streptococcus pneumoniae</i> Z022
<i>Chlamydia psittaci</i> Genotype C	Influenza A, B and H1N1 A/California/7/2009 (H1N1)PDM09; A/South Australia/55/2014, IVR-175 (Similar strain to Switzerland/9715293/2013 (H3N2)); B/Phuket/3073/2013	<i>Streptococcus pyogenes</i>
<i>Chlamydophila pneumoniae</i> Strain CM-1	Influenza A, B and H1N1 Strain similar to A/Michigan/45/2015(H1N1) pdm09 (A/Singapore/GP1908/2015, IVR-180); Strain similar to A/Hong Kong/4801/2014(H3N2) (A/Hong Kong/4801/2014, NYMC X-263B);Strain similar to B/Brisbane/60/2008 (B/Brisbane/60/2008).	<i>Streptococcus salivarius</i>
Coronavirus HKU	Influenza A/B Influenza A/New Caledonia/20/99 (H1N1), Influenza B/Florida/04/06	<i>Staphylococcus aureus</i> subsp. aureus
Coronavirus MERS	<i>Legionella bozemanii</i> Serovar 1	<i>Staphylococcus epidermidis</i>
Coronavirus NL63	<i>Legionella dumoffii</i>	

In addition, analytical specificity was analysed using public nucleotide sequence databases and search and/or alignment tools such as:

- BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>)
- Mafft Alignment (<http://mafft.cbrc.jp/alignment/server/>)
- Primer-BLAST (<http://www.ncbi.nlm.nih.gov/tools/primer-blast/>)

The bioinformatics analysis showed that the selected primers and probes should not cause false positives in the detection of *Pneumocystis jirovecii* in the presence of other organisms.

Finally, regarding PCR interferences and inhibitors, it should be noted that no endogenous or exogenous substances interfere with the Vitassay qPCR *Pneumocystis jirovecii* kit.

Analytical reactivity

The analytical reactivity (inclusivity) of Vitassay qPCR *Pneumocystis jirovecii* was experimentally confirmed by real-time amplification using DNA from the following strains from EQA programmes as template, with positive results:

Pneumocystis jirovecii Type 1A,

Pneumocystis jirovecii g885652

Pneumocystis jirovecii j888023

In addition, analytical reactivity was analysed using public nucleotide sequence databases and search and/or alignment tools such as:

- BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>)
- Mafft Alignment (<http://mafft.cbrc.jp/alignment/server/>)
- Primer-BLAST (<http://www.ncbi.nlm.nih.gov/tools/primer-blast/>)

The bioinformatics analysis showed that the selected primers and probes matched *Pneumocystis jirovecii* DNA sequences.

Compatible real-time PCR equipment

Vitassay qPCR *Pneumocystis jirovecii* has been validated on the following equipment:

- Cobas Z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ⁱⁱ

- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- LightCycler® 480 Real-Time PCR System (Roche Diagnostics)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR instrument (DNA-Technology)
- DTPRime Real Time PCR instrument (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I
- CFX96 Opus 96™ Real-Time PCR Detection System (Bio-Rad)

^I For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

^{II} When using the Applied Biosystems 7500 Fast with strips, it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

For more information regarding thermocycler compatibility based on bibliographic data and the most common detection channels contact the manufacturer (info@vitassay.com).

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Limitations

- This test provides a presumptive diagnosis of *Pneumocystis jirovecii* infection. All results obtained should be interpreted by a specialist together with the available clinical information and laboratory findings.
- This assay has been tested on bronchoalveolar lavage, bronchoalveolar aspirate and sputum samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper nucleic acids from respiratory clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Possible crosstalk might be observed in empty channels of the thermocycler if there is no target to be detected, so it is required to select only the channels where these are amplified when interpretation of results is performed.
- Occasionally, target template nucleic acids with a copy number below the detection limit can be detected, but the results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Pneumocystis jirovecii*, either by the positive control or the Quantitative standard during reconstitution, by samples containing high concentrations of target nucleic acids or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including nucleic acids extraction); c) nucleic acids degradation during specimen shipping, storage, and/or preparation; d) microorganism load is below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or

immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect the detection of new or unknown strains.

- The detection of nucleic acids does not imply the presence of viable and/or infectious microorganisms or that these microorganisms are the causative agents of clinical symptoms.
- Negative results do not preclude *Pneumocystis jirovecii* infection and should not be the sole basis for a patient management/treatment decision. The optimal specimen types for *Pneumocystis jirovecii* identification and/or the most appropriate stage of infection for collection have not yet been established. Consider collecting multiple samples from the same patient at different times, which may increase the likelihood of detecting the microorganism(s).
- If the patient's clinical data, and epidemiological studies suggest possible infection by these microorganisms, and other respiratory diseases have been discarded with laboratory tests, a false negative result might not be dismissed, and additional tests should be performed.
- The fluorescence values may vary due to multiple factors such as PCR equipment (even if it is the same model), extraction system, sample type, and pretreatment, among others.
- It is recommended to include a sample blank during the extraction and qPCR process to rule out possible contamination and to avoid missing probable false positive results.

Bibliography/Bibliografía


1. Truong J, Ashurst JV. (2023) *Pneumocystis jirovecii* Pneumonia. Jan 21. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2023 Jan–.
2. Weyant RB, Kabbani D, Doucette K, Lau C, Cervera C. (2021) *Pneumocystis jirovecii*: a review with a focus on prevention and treatment. *Expert Opin Pharmacother*. Aug;22(12):1579-1592. Epub 2021 Apr 19.
3. Ibrahim A, Chattaraj A, Iqbal Q, Anjum A, Rehman MEU, Aijaz Z, Nasir F, Ansar S, Zangeneh TT, Iftikhar A. (2023) *Pneumocystis jirovecii* Pneumonia: A Review of Management in Human Immunodeficiency Virus (HIV) and Non-HIV Immunocompromised Patients. *Avicenna J Med*. Mar 24;13(1):23-34.
4. Bateman M, Oladele R, Kolls JK. (2020) Diagnosing *Pneumocystis jirovecii* pneumonia: A review of current methods and novel approaches. *Med Mycol*. Nov 10;58(8):1015-1028.

Trademarks/Marcas registradas

All trademarks that may appear in this package insert are the property of their respective owners.

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