

For professional *in vitro* diagnostic use only.

INTENDED USE

Vitassay Clostridium difficile GDH+Toxin A+B is a rapid, immunochromatographic assay for the simultaneous qualitative detection of *Clostridium difficile* glutamate dehydrogenase (GDH), Toxin A and Toxin B in human stool samples.

Simple, non-invasive and highly sensitivity immunoassay to make a presumptive diagnosis of *Clostridium difficile* infection.

INTRODUCTION

Clostridium (Clostridioides) difficile, is a gram-positive anaerobic enterotoxigenic bacillus identified as one of the leading causes of infectious diarrhoea related to antibiotic treatment (10-15%). It is a major public health concern, as it responsible of broad outbreaks in hospitals and nursing homes. *C. difficile* is also associated with high mortality rates in the elderly (Marra et al., 2020). The main risk factors for *C. difficile* infections are age (>65 years-old), hospitalisation, being immunocompromised and/or previous infection (Abt et al., 2016).

C. difficile can be divided into two major strains: nontoxigenic and toxigenic, but only the later one causes disease in humans. These toxigenic strains produce toxin B (TcdB) alone, or in combination with toxin A (TcdA) and the pathogenicity depends on the expression of these two toxins (Burke et al., 2014). The toxins are produced and secreted in the large intestine under limit nutrient availability and undergo endocytosis by the intestinal epithelial cells. The activation of these toxins in the cytosol, produces cell necrosis and the subsequent loss of intestinal membrane integrity. This process leads to host exposure to intestinal microorganisms and the subsequent activation of the inflammatory response (CDC).

Different toxigenic strains have been reported to cause important outbreaks around the globe. In 2003, the new strain NAP1/ribotype 027 was a source of *C. difficile* epidemics in the US and Canada. It has also been described across Europe and Chile, and it is associated with higher morbidity and recurrent rates. It is of a particular interest since it has a mutation in an inhibitory gene, leading to an increase toxin A and B production (McDonald et al., 2005). In other parts of the world, as Asia, it is rare to find 027 ribotype, however 017 or 002 accounts for a large number of cases (Burke et al., 2014).

Treatment is based on antibiotic administration, including and vancomycin and fidaxomicin and metronidazole/vancomycin. However, antibiotics are not always affective and there is a 25% risk of having recurrent episodes. Innovative treatments, like FMT (faecal microbiota transplant), antibody-based therapy and

microbiota-based drugs have risen. There are two vaccines under evaluation (Guery et al., 2019).

Signs of infection might appear after taking antibiotics or being hospitalised/during hospitalisation. Common symptoms are severe diarrhoea, fever, stomach tenderness, loss of appetite and nausea (CDC). Measurements to avoid spreading must be taking, such as good hygiene and avoiding sharing towels/linen.

Systematic testing for *C. difficile* is recommended in case a diarrhoea happens in a healthcare environment, or where the common causatives of intestinal infections have been discarded.

PRINCIPLE

Vitassay Clostridium difficile GDH+Toxin A+B is a qualitative immunochromatographic assay to make a presumptive diagnosis of *Clostridium difficile* infection.

Strip A: The test line zone of the nitrocellulose membrane is pre-coated with monoclonal antibodies against GDH.

Strip B: The test line zone of the nitrocellulose membrane is pre-coated with monoclonal antibodies against Toxin A of *Clostridium difficile*.

Strip C: The test line zone of the nitrocellulose membrane is pre-coated with monoclonal antibodies against Toxin B of *Clostridium difficile*.

During the process, the sample reacts with the antibodies against GDH (strip A), Toxin A (strip B) and/or Toxin B (strip C), forming conjugates. The mixture moves upward on the membrane by capillary action. If the sample is positive, the antigens of the diluted sample react with the conjugate complex and a red line will be visible. Although the sample is positive or negative, the mixture continues to move across the membranes and the control line always appears.

The presence of a green line (in the control zone (C)) indicates that sufficient volume is added; proper flow is obtained and serves as an internal control for the reagents.

PRECAUTIONS

- For professional *in vitro* use only.
- Do not use after expiry date.
- Read the instructions for use carefully before using the test.
- Do not use the kit if the label sealing the outer carton is torn or if the bags are open or damaged on arrival.
- Do not use the tests if the desiccant material is missing or broken inside the aluminium pouch.
- Specimens should be considered potentially hazardous and should be handled in the same manner as an infectious agent, following local/national regulations. A new test should be used for each sample to avoid contamination errors.

VITASSAY

Clostridium difficile GDH+Toxin A+B

Rapid test for the qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B in human stool samples.

IUE-7715024 Ed03 June 2023



- Material exposed to the samples should also be considered potentially hazardous and should be handled in the same manner as an infectious agent, following local/national regulations.
- Do not reuse. This is a single-use device.
- Used material should be disposed of in an appropriate biohazard container after testing.
- The reagents contain a preservative called sodium azide (<0.1%). Avoid any contact with skin or mucous membranes.
- In accordance with current regulations **Vitassay Clostridium difficile GDH+Toxin A+B** does not contain substances and/or mixtures that are hazardous or present in a concentration above the limits for their declaration. The safety data sheet is available on request (not included).
- All reagents included in the kit are approved for use with **Vitassay Clostridium difficile GDH+Toxin A+B** only. Do not mix or use the components with other batches of Vitassay. Do not use with reagents from other kits or commercial assays.
- Follow Good Laboratory Practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as lab coat, surgical or appropriate mask or face shield, disposable gloves and eye protection. Take the necessary precautions during sample collection, transport, storage, handling and disposal. Each sample must be correctly and unequivocally identified to ensure proper traceability of samples.
- In case of spillage, clean thoroughly with a suitable disinfectant.
- Do not eat, drink or smoke in the workplace.
- The presence of yellow lines in the result window (control line area and test line area), before using the test, is completely normal and does not imply a failure in the functionality of the test.
- The visual interpretation of the results is done by coloured lines, the interpretation of the results should be done by a professional user without problems of visualisation and colour interpretation.
- A certificate of analysis can be provided on request (not included).

STORAGE AND STABILITY

The storage temperature of the kits should be 2-30°C.

Do not freeze.

Under these conditions, they can be used until the expiry date indicated on the kit label.

All kit components are for single use only and must remain in their primary packaging until use. The test must remain in the sealed pouch until use.

MATERIALS

MATERIAL PROVIDED	MATERIAL REQUIRED BUT NOT PROVIDED
10X Vitassay Clostridium difficile GDH+Toxin A+B <ul style="list-style-type: none"> • Instructions for use. • 10 X vials with diluent for sample dilution. 	<ul style="list-style-type: none"> • Specimen collection container. • Disposable gloves. • Timer. • Spatula.

SPECIMEN COLLECTION

Collect sufficient quantity of faeces: 1-2g or mL for liquid samples. Stool specimens should be collected in clean and dry containers.

If processed the day after collection, samples can be stored at room temperature. However, we recommend testing immediately after sample collection.

Samples can be stored for up to 7 days at 2-8°C, and also frozen at -20°C for up to one year. Samples shall be brought to room temperature before testing.

Ensure that only the necessary quantity is thawed, as freezing freeze-thaw cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

SPECIMEN PREPARATION

- 1 Remove the cap of the vial with diluent for sample dilution (figure 1) and use the spoon to collect sufficient sample quantity. For solid stool, insert the spoon in 4 different areas of the stool sample (figure 2), remove any excess sample with a spatula (figure 3), and place the spoon cap back into the vial for sample dilution (figure 4). For liquid stool, take a spoonful of the sample (figure 3) and transfer it into the vial for sample dilution.
2. Close the vial for sample dilution tightly and shake it to dilute and mix the sample with the diluent (figure 4).



Vial for sample dilution



Insert the spoon in 4 different areas of the stool.



Remove excess sample with a spatula. Liquid samples: full spoon.




Put the sample into the vial, close the cap and shake.

PROCEDURE


Allow the test, stool sample, controls and diluent to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Shake the vial with the sample vigorously to obtain a good sample dilution.
2. Remove the **Vitassay Clostridium difficile GDH+Toxin A+B** from its sealed bag just before using it (figure 5).
3. Take the vial for sample dilution containing the diluted sample (figure 6), place it inside the multiplex tube (figure 7). Screw the cap of the multiplex tube tightly (figure 8). The bottom of the vial for sample dilution will break and the diluent+sample solution reaches the sample zone of the strips (figure 9).
4. Leave the multiplex tube vertically on a flat surface and read the results at **10 minutes**. Do not read the test results later than 10 minutes.



If the test does not run due to solid particles (the sample is not homogenized), migration process can stop on one or more strips. In this case, tap the end of the multiplex tube on hard surface to allow migration to start again.



Vitassay Clostridium difficile GDH+Toxin A+B



Vial with the diluted sample inside.



Introduce the vial with the diluted sample into the multiplex.

Close the cap and the bottom of the diluent vial will break.

Reaction takes place. **Read results at 10 minutes.**

INTERPRETATION OF THE RESULTS

Strip A: GDH, Strip B: Toxin A and Strip C: Toxin B

Diagram	Result	Interpretation
	NEGATIVE Only one green line in the control zone (C) in the three strips (A,B and C).	There is no GDH, Toxin A and Toxin B presence. A negative result should be cross-checked with other laboratory techniques to confirm the results.
	POSITIVE In addition to the green line (control line C), a red line appears in each strip, test line (T).	There is GDH, Toxin A and Toxin B presence. Possible <i>C. difficile</i> infection.
	NEGATIVE Strip C (Toxin B) → green line.	There is GDH and Toxin A presence. Possible <i>C. difficile</i> infection.
	POSITIVE Strip A (GDH) → green/red lines. Strip B (Toxin A) → green/red lines.	
	NEGATIVE Strip B (Toxin A) → green line.	There is GDH, Toxin B presence. Possible <i>C. difficile</i> infection.
	POSITIVE Strip A (GDH) → green/red lines. Strip C (Toxin B) → green/red lines.	

	<p>NEGATIVE</p> <p>Strip A (GDH) → green line.</p>	<p>There is Toxin A and Toxin B presence. Possible <i>C. difficile</i> infection.</p>
	<p>POSITIVE</p> <p>Strip B (Toxin A) → green/red lines.</p> <p>Strip C (Toxin B) → green/red lines.</p>	
	<p>NEGATIVE</p> <p>Strip B (Toxin A) → green line.</p> <p>Strip C (Toxin B) → green line.</p>	<p>There is GDH presence. Possible <i>C. difficile</i> infection.</p>
	<p>POSITIVE</p> <p>Strip A (GDH) → green/red lines.</p>	
	<p>NEGATIVE</p> <p>Strip A (GDH) → green line.</p> <p>Strip C (Toxin B) → green line.</p>	<p>There is Toxin A presence. Possible <i>C. difficile</i> infection.</p> <p>If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin A and negative for GDH, the sample should be considered positive for Toxin A.</p>
	<p>POSITIVE</p> <p>Strip B (Toxin A) → green/red lines.</p>	
	<p>NEGATIVE</p> <p>Strip A (GDH) → green line.</p> <p>Strip B (Toxin A) → green line.</p>	<p>There is Toxin B presence. Possible <i>C. difficile</i> infection.</p> <p>If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin B and negative for GDH, the sample should be considered positive for Toxin B.</p>
	<p>POSITIVE</p> <p>Strip C (Toxin B) → green/red lines.</p>	
Any other results	Invalid results either A, B or C, we recommend repeating the assay using the same sample with another test.	

A very low percentage of specimens might result negative for GDH but positive for toxins.

Notes: The intensity of the red coloured test line in the result line region (T) will vary depending on the concentration of antigens in the specimen.

Positive results detailed in the above table should be followed up with additional confirmatory diagnostic procedures.

Invalid results: Total absence of any control coloured lines (green) indicates an invalid result, regardless of the appearance or not of the test lines (red). Wrong procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the problem persists, discontinue using the kit and contact your local distributor.

If the results are read later than 10 minutes after sample addition, they shall be considered invalid.

QUALITY CONTROL

Internal procedural controls are included in **Vitassay GDH+Toxin A+B**. Green lines appearing in the results window are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

- **Vitassay Clostridium difficile GDH+Toxin A+B** test must be carried out within 2 hours of opening the sealed bag.
- An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary depending on the concentration of antigens.
- The use of other samples different from human fecal samples has not been established.
- The quality of **Vitassay Clostridium difficile GDH+Toxin A+B** depends on the quality of the sample; Proper fecal specimens must be obtained.
- Positive results determine the presence of GDH, Toxin A and/or Toxin B of *Clostridium difficile* in fecal samples. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- Negative results should not be considered as conclusive; it is possible that the concentration of antigen is lower than the detection limit value. If symptoms or situation still persist a *Clostridium difficile* determination should be carried out on a sample from an enrichment culture.
- **Bloody stool samples** and/or **mucinous** stool samples can cause non-specific reactions in the test. Such positive samples

should be followed up with other diagnostic techniques to confirm the result.

- This test is a screening tool. The final diagnosis must be confirmed by a physician after a thorough evaluation of the clinical history and additional confirmatory tests.

EXPECTED VALUES

Clostridium difficile is associated with 95-100% of cases of pseudomembranous colitis, 60-75% of cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases. The other causes of antibiotic-associated diarrhoea are largely unknown.

According to ECDC, in 2016–2017, 72.0% of the CDI (*C. difficile* infection) cases with case-based data were above 64 years old and the majority (56.4%) were female. More than half (n=3 446/5 863; 58.8%) of the CDI cases had had contact with healthcare in the three months before the current healthcare admission, of which the vast majority (n=2 804/3 446; 81.4%) had had contact with a hospital. Fewer were reported to have had contact with a long-term care facility (LTCF) (n=431; 12.5%). There were an estimated 189 526 healthcare-associated (HA) CDI cases (cumulative 95% confidence interval (95%CI): 105 154–340 978) in EU/EEA countries/administrations annually in 2016–2017. In 2016–2017, the crude incidence density of CDI was 3.48 cases per 10 000 patient-days. This was higher in tertiary care hospitals (3.87 cases per 10 000 patient-days) than in secondary or primary care hospitals (3.46 and 2.28 cases per 10 000 patient-days, respectively).

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity and Specificity

Two retrospective evaluations were carried out with stool samples, comparing the results obtained by **Vitassay Clostridium difficile GDH+Toxin A+B** against results obtained with commercial assays based on the immunochromatographic technique, considered as reference methods.

In an initial study, a total of 126 samples from patients with *C. difficile* infection were compared. The results are shown below:

GDH

		Reference method (GDH)		
		Positive	Negative	Total
Vitassay Clostridium difficile GDH+Toxin A+B	Positive	51	0	51
	Negative	0	75	75
	Total	51	75	126
GDH				

Table 1. Results of Vitassay Clostridium difficile GDH+ Toxin A + B (GDH) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile GDH + Toxin A + B (GDH) Reference method (GDH)		
	Value	IC 95%
Sensitivity	100%	93.0-100%
Specificity	100%	95.2-100%
PPV	100%	93.0-100%
NVP	100%	95.2-100%

Table 2. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (GDH) kit compared to a commercial CE-IVD kit.

Toxin A

		Reference method (Toxin A+B)		
		Positive	Negative	Total
Vitassay Clostridium difficile GDH+Toxin A+B Toxin A	Positive	28	0	28
	Negative	1	97	98
	Total	29	97	126

Table 3. Results of Vitassay Clostridium difficile GDH + Toxin A + B (Toxin A) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile GDH + Toxin A + B (Toxin A) vs Reference method (Toxin A+B)		
	Value	IC 95%
Sensitivity	96.6%	92.2-99.9%
Specificity	100%	96.2-100%
PPV	100%	87.7-100%
NVP	99.0%	94.4-100%

Table 4. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (Toxin A) kit compared to a commercial CE-IVD kit.

Toxin B

		Reference method (Toxin A+B)		
		Positive	Negative	Total
Vitassay Clostridium difficile GDH+Toxin A+B Toxin B	Positive	36	1	37
	Negative	0	89	89
	Total	36	90	126

Table 5. Results of Vitassay Clostridium difficile GDH + Toxin A + B (Toxin B) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile GDH + Toxin A + B (Toxin B) vs Reference method (Toxin A+B)		
	Value	IC 95%
Sensitivity	96.6%	92.2-99.9%
Specificity	100%	96.2-100%
PPV	100%	87.7-100%
NVP	99.0%	94.4-100%

	95% IC	IC 95%
Sensitivity	100%	90.3-100%
Specificity	98.9%	94.0-100%
PPV	97.3%	85.8-99.9%
NVP	100%	95.9-100%

Table 6. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (Toxin B) kit compared to a commercial CE-IVD kit.

In a second, more recent evaluation, 250 frozen clinical specimen (stool) samples from patients with suspected *C. difficile* infection were evaluated. This test evaluated the performance of the Vitassay kit against another commercial CE-IVD kit and discrepancies were resolved using a certified nucleic acid amplification method.

GDH

		Reference method (GDH)		
		Positive	Negative	Total
Vitassay Clostridium difficile GDH+Toxin A+B GDH	Positive	86	1	87
	Negative	3	160	163
	Total	89	161	250

Table 7. Results of Vitassay Clostridium difficile GDH + Toxin A + B (GDH) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile GDH + Toxin A + B (GDH) Reference method (GDH)		
	Value	IC 95%
Sensitivity	96.6%	90.5-99.3%
Specificity	99.4%	96.6-100%
PPV	98.9%	93.0-100%
NVP	98.2%	95.2-100%

Table 8. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (GDH) kit compared to a commercial CE-IVD kit.

Toxin A

		Reference method (Toxin A+B)		
		Positive	Negative	Total
Vitassay Clostridium difficile Toxin A+B Toxin A	Positive	51	0	51
	Negative	1	198	199
	Total	52	198	250

Table 9. Results of Vitassay Clostridium difficile GDH + Toxin A + B (Toxin A) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile Toxin A + B (Toxin A) vs Reference method (Toxin A+B)		
	Value	IC 95%
Sensitivity	98.1%	89.7-100%
Specificity	100%	98.2-100%
PPV	100%	93.0-100%
NVP	99.5%	97.2-100%

	Value	IC 95%
Sensitivity	98.1%	89.7-100%
Specificity	100%	98.2-100%
PPV	100%	93.0-100%
NVP	99.5%	97.2-100%

Table 10. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (Toxin A) kit compared to a commercial CE-IVD kit

Toxin B

		Reference (Toxin A+B)		
		Positive	Negative	Total
Vitassay Clostridium difficile Toxin A+B Toxin B	Positive	40	0	40
	Negative	4	206	210
	Total	44	206	250

Table 11. Results of Vitassay Clostridium difficile GDH + Toxin A + B (Toxin B) compared to a commercial CE-IVD kit.

Vitassay Clostridium difficile Toxin A + B (Toxin B) vs Reference method (Toxin A+B)		
	Value	IC 95%
Sensitivity	90.9%	78.3-97.5%
Specificity	100%	98.2-100%
PPV	100%	91.2-100%
NVP	98.1%	95.2-99.9%

Table 12. Sensitivity, specificity, positive predictive values and negative predictive values of the Vitassay Clostridium difficile GDH +Toxin A + B (Toxin B) kit compared to a commercial CE-IVD kit.

In brief, **Vitassay Clostridium difficile GDH + Toxin A + B** has been tested on a total of 376 clinical faecal samples. The **pooled** analysis of both studies shows an overall clinical **sensitivity** for GDH of 0.97 (0.93-0.99), toxin A of 0.99 (0.97-0.99) and for toxin B of 0.95 (0.87-0.98). On the other hand, the overall clinical **specificity** calculated for GDH is 0.99 (0.97-1), toxin A is 1 (0.98-1) and for toxin B 0.99 (0.98-1).

Altogether, the results show that **Vitassay Clostridium difficile GDH+Toxin A + B** has a high sensitivity and specificity for detecting the enzyme glutamate dehydrogenase (GDH), toxins A and B of *Clostridium difficol*.

Analytical sensitivity

The Vitassay **Clostridium difficile GDH +Toxin A+B** limit of detection (LoD) is 0.39 ng/mL for *C. difficile* Glutamate dehydrogenase (GDH) antigen, 2 ng/mL for *Clostridium difficile* Toxin A and 3.12 ng/mL for Toxin B. Note that a general, reference standard has not been approved; LoD was determined using recombinant GDH (glutamate dehydrogenase) protein, and purified *C. difficile* Toxins A and B as internal standards.

Cross reactivity

No cross reactivity was detected against other gastrointestinal pathogens, other organism, substances and/or faecal markers that are occasionally present in faeces.

GDH

Adenovirus	<i>Escherichia coli</i> O111	<i>Peptostreptococcus anaerobius</i>
Astrovirus	<i>Escherichia coli</i> O26	Rotavirus
Bovine Haemoglobin	<i>Escherichia coli</i> O157	RSV
<i>Campylobacter coli</i>	<i>Clostridium difficile</i> Toxin A	<i>Salmonella paratyphi</i> A
<i>Campylobacter jejuni</i>	<i>Giardia</i>	<i>Salmonella typhi</i>
<i>Clostridium bifermentans</i>	<i>Helicobacter pylori</i>	<i>Salmonella typhimurium</i>
<i>Clostridium butyricum</i>	Human Calprotectin	<i>Shigella boydii</i>
<i>Clostridium haemolyticum</i>	Human Haemoglobin	<i>Shigella dysenteriae</i>
<i>Clostridium novyi</i>	Human Lactoferrin	<i>Shigella flexneri</i>
<i>Clostridium tetani</i>	Human Transferrin	<i>Shigella sonnei</i>
<i>Clostridium difficile</i> Toxin B	Influenza A virus	<i>Staphylococcus aureus</i>
<i>Clostridium perfringens</i>	Influenza B virus	<i>Streptococcus pneumoniae</i>
<i>Clostridium septicum</i>	<i>Salmonella enteritidis</i>	<i>Streptococcus pyogenes</i>
<i>Cryptosporidium parvum</i>	Norovirus GII	<i>Yersinia Enterocolitica</i> O:3
<i>Entamoeba histolytica</i>	Porcine haemoglobin	<i>Yersinia Enterocolitica</i> O:9

It has been observed that *C. sporogenes* (CECT 485) y *C. bolutinium* (CECT 551) can be detected with **Vitassay Clostridium difficile GDH+Toxin A + B**. Intercomparison studies have shown that *C. sordellii* (ATCC 9714) yielded a positive signal in this Vitassay kit, however the origin of this positive result could not be established. Hence, it was not considered a genuine cross-reactivity.

Toxin A+B

Adenovirus	<i>Giardia</i>	<i>Salmonella enteritidis</i>
Astrovirus	<i>Helicobacter pylori</i>	<i>Streptococcus pyogenes</i>
Human calprotectin	Bovine haemoglobin	<i>Staphylococcus aureus</i>
<i>Campylobacter coli</i>	Porcine haemoglobin	Human transferrin
<i>Campylobacter jejuni</i>	Human haemoglobin	<i>Salmonella paratyphi</i> A
<i>Clostridium difficile</i> antigen GDH	Influenza A virus	<i>Salmonella typhi</i>
<i>Clostridium difficile</i> Toxin A (strip: Toxin B)	Influenza B virus	<i>Salmonella typhimurium</i>
<i>Clostridium difficile</i> Toxin B (strip: Toxin A)	Human lactoferrin	<i>Shigella boydii</i>
<i>Clostridium perfringens</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
Coronavirus	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Cryptosporidium parvum</i>	Norovirus GI	<i>Shigella sonnei</i>

<i>Entamoeba histolytica</i>	Norovirus GII	<i>Streptococcus pneumoniae</i>
<i>Escherichia coli</i> O111	Rotavirus	<i>Yersinia enterocolitica</i> O:3
<i>Escherichia coli</i> O26	RSV	<i>Yersinia enterocolitica</i> O:9
<i>Escherichia coli</i> O157		

Interferences

Possible interferences were assessed using the following **exogenous compounds and drugs** at the given concentrations (no interferences were found):

Acetylcysteine (Fluimucil)	3.0 mg/mL	No interference
Acetyl Salicylic (Adiro)	0.3 mg/mL	No interference
Afrin (Oxymetazoline)	5.0-10-2 mg/mL	No interference
Almagato (Almax)	3.0 mg/mL	No interference
Amantadine	0.3 mg/mL	No interference
Ambroxol hydrochloride (Mucosan)	0.3 mg/mL	No interference
Amoxicillin	3.0 mg/mL	No interference
Ampicillin	3.0 mg/mL	No interference
Benzocaine (Angileptol)	0.3 mg/mL	No interference
Biotine	100.0 µg/mL	No interference
Carbocisteine (Iniston mucolítico)	3.0 mg/mL	No interference
Ciprofloxacin	0.3 mg/mL	No interference
Cloperastine (Flutox)	0.3 mg/mL	No interference
Codeine (Toseina)	0.2 mg/mL	No interference
CVS Nasal Drops (Phenylephrine)	10.0 mg/mL	No interference
CVS Nasal Spray (Cromolyn)	4 mg/mL	No interference
Dexchlorpheniramine (Polaramine)	0.3 mg/mL	No interference
Dexketoprofen trometamol (Enantyum)	0.3 mg/mL	No interference
Ebastine (Ebaste)	3.0 mg/mL	No interference
Fluticasone Propionate	5.0-10-2 mg/mL	No interference
Fosfamicin (Monurol)	3.0-10-3 mg/mL	No interference
Heparin (Hibor)	350.0 IU/mL	No interference
Hydroxyzine dihydrochloride	0.3 mg/mL	No interference
Homeopathy	1/10 dilution	No interference
Ibuprofen (Espidifen)	3.0 mg/mL	No interference
Levofloxacin	3.0 mg/mL	No interference
Loperamide hydrochloride (Fortasec)	0.15 mg/mL	No interference
Loratadine	0.3 mg/mL	No interference
Lorazepam	3.0-10-3 mg/mL	No interference
Lysine Carbocysteinat (Pectox)	3.0-10-2 mg/mL	No interference
Macrogol 3350 (Movicol)	3.0 mg/mL	No interference
Mercaptopurine	0.3 mg/mL	No interference
Metamizole (Nolotil)	5.0 mg/mL	No interference
Metronidazole	3.0 mg/mL	No interference
Mupirocin	2.5-10-2 mg/mL	No interference
Naso GEL	0.9 mg/mL	No interference
Omeprazole	2.0-10-3 mg/mL	No interference

Oseltamivir	0.3 mg/mL	No interference
Paracetamol (Dolocati)	5.0 mg/mL	No interference
Phenoxymethylpenicillin potassium	3.0 mg/mL	No interference
Phenulpropanolamine	0.3 mg/mL	No interference
Prednisone	0.3 mg/mL	No interference
Ribavirin	3.0 mg/mL	No interference
Rifampicin (Rifaldin)	0.3 mg/mL	No interference
Sore Throat Phenol spray	0.5 mg/mL	No interference
Tobramycin	0.3 mg/mL	No interference
ZICAM	1.0 mg/mL	No interference

In addition, the following human **endogenous substances** were also analysed. As before, no interferences were found:

Human blood	5.0 %	No interference
Human calprotectin	5.0 µg/mL	No interference
Human haemoglobin	5.0 µg/mL	No interference
Human lactoferrin	5.0 µg/mL	No interference
Human transferrin	0.5 µg/mL	No interference
Mucine	0.5 %	No interference

Repeatability and Reproducibility

Repeatability and reproducibility studies performed with different positive and negative samples yielded no differences between the evaluations.

Hook effect

No inhibitory hook effect was demonstrated at concentrations higher than 10^5 times the limit of detection ($10^5 \times \text{LoD}$) for GDH, nor 10^2 times the limit of detection ($10^2 \times \text{LoD}$) for toxin A, nor at concentrations 10^4 times the limit of detection ($10^4 \times \text{LoD}$) for toxin B in the **Vitassay Clostridium difficile GDH+Toxin A + B** kit.











REFERENCES

- Abt MC, McKenney PT, Pamer EG. Clostridium difficile colitis: pathogenesis and host defence. Nat Rev Microbiol. 2016 Oct;14(10):609-20. doi: 10.1038/nrmicro.2016.108. Epub 2016 Aug 30. PMID: 27573580; PMCID: PMC5109054.
- Burke KE, Lamont JT. Clostridium difficile infection: a worldwide disease. Gut Liver. 2014 Jan;8(1):1-6. doi: 10.5009/gnl.2014.8.1.1. Epub 2014 Jan 13. PMID: 24516694; PMCID: PMC3916678.
- CDC. "What Is C. Diff?" Centers for Disease Control and Prevention, Centers for Disease Control and Prevention, 20 July 2021, <https://www.cdc.gov/cdiff/what-is.html>. Accessed 05 January 2021
- Guery B, Galperine T, Barbut F. Clostridioides difficile: diagnosis and treatments. BMJ. 2019 Aug 20;366:14609. doi: 10.1136/bmj.14609. PMID: 31431428.

5. Marra AR, Perencevich EN, Nelson RE, Samore M, Khader K, Chiang HY, Chorazy ML, Herwaldt LA, Diekema DJ, Kuxhausen MF, Blevins A, Ward MA, McDanel JS, Nair R, Balkenende E, Schweizer ML. Incidence and Outcomes Associated With Clostridium difficile Infections: A Systematic Review and Meta-analysis. JAMA Netw Open. 2020 Jan 3;3(1):e1917597. doi: 10.1001/jamanetworkopen.2019.17597. PMID: 31913488; PMCID: PMC6991241.
6. McDonald LC, Killgore GE, Thompson A, Owens RC Jr, Kazakova SV, Sambol SP, Johnson S, Gerding DN. An epidemic, toxin gene-variant strain of Clostridium difficile. N Engl J Med. 2005 Dec 8;353(23):2433-41. doi: 10.1056/NEJMoa051590. Epub 2005 Dec 1. PMID: 16322603.
7. European Centre for Disease Prevention and Control. Clostridioides (Clostridium) difficile infections. Annual epidemiological report for 2016–2017. Stockholm: ECDC; 2022.



SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	<i>in vitro</i> diagnostic device		Keep dry
	Consult instructions for use		Temperature limitation
	Use by		Manufacturer
	Batch code		Contains sufficient for <n> test
DIL	Sample diluent		Catalogue number
	CE Marking		

