VA Vitassay

Real-Time PCR Kits

Vitassay qPCR Japanese Encephalitis

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Intended use

Vitassay qPCR Japanese Encephalitis allows the qualitative detection of Japanese encephalitis virus by real-time RT-PCR in blood or CSF samples. The product is intended for use in the diagnosis of Japanese encephalitis infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Japanese Encephalitis 4x8-well strip, low profile	7041038
Vitassay qPCR Japanese Encephalitis 4x8-well strip, high profile	7042038

Provided Reagents

In references 7041038 and 7042038:

Reference	Reagent/Material	Description	Concentration Range	Colour	Amount
7041S038/ 7042S038	Japanese Encephalitis strips	Lyoprotectors and Stabilizers	±6 g/100 mL *	-	4 x 8-well strips
		dNTPs	±1 mM *		
		Primers and Probes	0,2-1 nMol/μL *		
		Enzymes	10-100 U/rxn *		
7C038	Japanese Encephalitis Positive Control	Non-infectious synthetic lyophilized cDNA	1,9x10 ⁴ copies/μL*	red	1 vial
7001A	PCR grade water	RNAsa/DNAsa free water	1 g/mL	white	1 vial x 1 mL
7002B	Resuspension buffer	Saline Solution	±13 mM	green	1 vial x 1,8 mL
		Buffer (TRIS, pH)	±67 mM		
7003N	Negative control	RNAsa/DNAsa free water	1 g/mL	yellow	1 vial x 1 mL
7004O	Optical caps	Optical caps for sealing strips during thermal cycling	-	-	4 x 8-cap strips

* For component in stabilized format, the concentration range means after rehydration.

Material and equipment required, not provided

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- RNA extraction kit
- Real-time PCR instrument (thermocycler)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 μL, 20-200 μL, 100-1000 μL)
- Filter tips
- Powder-free disposal gloves
- Blank sample (See RNA extraction)
- Control material for reverse transcription step (such as reference RNA or viral RNA particles).

Transport and storage

- Transport and storage of the kits (before use) can be performed at 2-40°C until the expiry date stated on the label. Once in use, see storage conditions for each component. Avoid shaking during transport to prevent liquid leakage.
- The resuspended positive control should be stored at -20°C until the expiry date stated on the kit label. The stability of the positive control has been validated after 6 freeze-thaw cycles. However, it is recommended to distribute in aliquots to avoid repeated freeze-thaw cycles.

- For proper use of the kit, it is recommended that unneeded wells are not rehydrated. Unused wells can be easily cut before removing the protective foil and stored in the foil bags provided with the silica gel inside and closed with the zip closure to avoid contact with light or moisture and stored at 2-40°C for the shelf life indicated on the kit label. Once the wells have been rehydrated, nucleic acid addition and RT-qPCR should be performed immediately.
- Once opened, the Resuspension buffer, PCR grade water and Negative Control components should be stored at 2-40°C for the shelf life indicated on the kit label.
- Keep all reagents in the darkness.

Summary

Japanese encephalitis (JE) virus is an enveloped RNA virus of the genus *Flavivirus*, family Flaviviridae (ECDC Factsheet Japanese Encephalitis, Mar2022). JE virus exists as one phenotype and five distinguishable genotypes (Roberts A. et al, 2020). JE virus is closely related to West Nile and St. Louis encephalitis viruses and is the leading cause of vaccine-preventable encephalitis in Asia and the western Pacific (CDC Japanese Encephalitis Virus, Dec2023). JE virus transmission occurs primarily in rural agricultural areas, often associated with rice production and flooding irrigation. In some areas of Asia, these conditions can occur near urban centers (CDC Japanese Encephalitis Virus Transmission, Dec2022).

Japanese encephalitis virus is a mosquito-borne virus spread by the bite of carrier *Culex* mosquitoes, mainly *Culex tritaeniorhynchus* and *Culex vishnui* mosquitoes (Sharma KB. et al., 2021). The virus is maintained in a zoonotic cycle between mosquitoes (the vector) and vertebrate hosts, primarily pigs and wading birds (also referred to as amplifying hosts or natural reservoirs). Humans are incidental or dead-end hosts, because they usually do not develop high enough concentrations of JE virus in their bloodstreams to infect feeding mosquitoes (CDC Japanese Encephalitis Virus Transmission, Dec2022) (Roberts A. et al, 2020).

Most people infected with JE do not have symptoms or have only mild symptoms. However, a small percentage of infected people develop inflammation of the brain (encephalitis), with symptoms including headache, fever, disorientation, seizures, weakness, and coma. About 1 in 4 cases are fatal (CDC Japanese Encephalitis Virus, Dec2023). On average, 1 in 250 infected people develops a severe neuroinvasive illness. The incubation period is usually 5 to 15 days (ECDC Factsheet Japanese Encephalitis, Mar2022).

According to JEV statistics given by WHO in 2019, 24 countries in South-East Asia and Western Pacific regions have endemic JEV transmission, which exposes 3 billion people to the risk of infection (Roberts A. et al, 2020). A literature review estimates nearly 68 000 clinical cases of JE globally each year, with approximately 13 600 to 20 400 deaths (WHO Japanese Encephalitis, May 2019). The case-fatality rate among those with encephalitis can be as high as 30%, more so in children. Permanent neurologic or psychiatric sequelae can occur in 30-50% of cases (Roberts A. et al, 2020).

The detection of JEV infection is problematic due to the short duration of viraemia and the asymptomatic nature of infection (Roberts A. et al, 2020). The diagnosis of Japanese encephalitis virus infection relies on the detection of specific IgM antibodies that are present in the cerebrospinal fluid and serum of patients and confirmation by neutralization assays. Diagnosis of JEV could be performed also by viral direct detection by RT-PCR performed on blood or cerebrospinal fluid from patients in the early phase of the disease (ECDC Factsheet Japanese Encephalitis, Mar2022).

Test principle

Vitassay qPCR Japanese Encephalitis is based on the real-time amplification of specific conserved fragments of the NS5 gene encoded by the Japanese Encephalitis genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of Japanese Encephalitis is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Japanese Encephalitis is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The detection channels of the target sequences are described as follows:

Target	Detection Channel
NS5 gene (Japanese encephalitis virus)	FAM
Internal Control (IC)	HEX, VIC, or JOE (depending on the equipment used)

Precautions

- For professional in vitro diagnostic use.
- Do not use the kit after the expiration date.
- Do not use the kit without having read and understood the information on procedures, precautions, and limitations provided in the instructions for use.
- Do not mix components from different kits and/or lots.
- Do not use if the package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Use the strips immediately after removing the protective foil to protect the reaction mixture from sunlight. It is not recommended to leave the strips open without resuspending. Once the required wells have been rehydrated, the qPCR assay should be performed immediately.
- It is recommended to cut the required/necessary wells before removing the protective seal from the entire strip. Store the rest inside the bag with the desiccant.
- Close protective pouches of reagents promptly with the zip seal after each use (do not remove the desiccant from the pouches). Remove excess air from the bags before sealing.
- Do not use the product near solvents to avoid deterioration of the label.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Rehydrate the vial immediately after opening the Positive Control envelope and use it in the PCR reaction or store it properly. Leaving the vial open without resuspending may cause deterioration of the reagent.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area. Once the test is completed wash your hands thoroughly.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended. A new tip should be used for each sample. It is necessary to change gloves before handling the reagents.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- The product Vitassay qPCR Japanese Encephalitis has only been validated with the equipment mentioned in Section Compatibles real-time PCR equipment of this Instructions for Use.
- You must ensure that the conditions are set up in the thermal cycler according to the instructions in section "Programme your thermocycler".
- Refer to Safety data sheets, available on request.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport the conditions validated by the user must be followed. Vitassay qPCR Japanese encephalitis has been tested in blood and cerebrospinal fluid (CSF). Other sample types must be validated by the user.

Overall, clinical samples should be collected properly in clean containers with or without transport media (depending on sample type), labelled correctly and be processed as soon as possible to guarantee the test quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. Blood samples should be transported at 4°C for 24 hours. For longer transport (more than 24 hours), shipment at -20°C or below is recommended (Yong W., 2019). Blood samples can be stored at 4°C for up to 24 hours or frozen at -20°C or

lower (ideally at -80°C) for conservation during a long time. In the case of CSF samples, they should be stored frozen at -20°C or lower (ideally at -80°C). Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- Yong W. (eds) Biobanking. Methods in Molecular Biology, vol 1897 (2019) Humana Press, New York, NY.

RNA extraction

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use. The assay has been validated with the following extraction kits:

- QIAamp Viral RNA Mini Kit (Qiagen).
- Invisorb® Spin Universal Kit (Invitek).

For the use of RNA extraction kits other than those mentioned above, validation must be performed by the end user, following the manufacturer's instructions for use.

It should be noted that the positive and negative control should not be extracted. To monitor and control the extraction process and rule out possible contamination, a sample blank (not supplied in the kit) can be extracted, which should consist of the same matrix as the clinical samples (real or simulated matrix), and previously characterized as negative for all targets.

Positive control preparation

Reconstitute the lyophilized Japanese Encephalitis Positive Control (red tube) with the 400 μL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C .

This component contains a high-copy number template and entails a very significant contamination risk. Therefore, we recommend opening and manipulating it in a separate area of the laboratory away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
It is recommended not to rehydrate wells that are not needed. Unused wells can be easily cut before removing the protective foil and stored as indicated in the Transport and Storage Section.
- Peel off the protective aluminium seal from the strips to be used.
- Pipette 15 μL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 μL of RNA extracted from each sample, of negative control (yellow tube) and of reconstituted positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the optical caps provided. It is recommended spin down briefly.
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the annealing/extension step (*) through the FAM (Japanese Encephalitis virus) and HEX, JOE, or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel.

Analysis and interpretation of results

The analysis of the results is performed by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control signal emission.

For each assay, it is recommended to set the threshold values for each channel (target) independently by the end-user. To do so, select the well corresponding to the positive control for each channel and set the threshold value within the exponential phase of the fluorescence curve and above any background signal (below the baseline). In case a threshold is set automatically by the equipment used, check, and verify that it is adjusted to the positive control or adjust it manually. Once the threshold value is set, the test samples can be interpreted.

The chemical nature of the different fluorophores can affect the threshold value of the channels, just as different signal intensities can influence the threshold value for different instruments.

In addition, it is recommended to include a blank (confirmed negative sample from the same matrix as the samples tested) for the detected targets to establish the baseline of the assay.

For a valid diagnostic test run, before analyzing the result of the clinical samples, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve (Ct ≤40) in FAM (Japanese Encephalitis), which validates the reaction.

Negative control

The negative control included in each run must show the absence of signal (Ct >40 or no signal) in FAM (Japanese Encephalitis) which validates the reaction.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in the positive control and negative control wells.

The experiment is invalid if there is no signal in the positive control or is amplification signal in the negative control, except for the HEX, JOE or VIC channel which should show amplification signal (corresponding to Internal Control). In either case, the assay must be repeated.

Once the results of the controls have been validated, for the interpretation of the sample results **select only the channels where the targets are detected**. The result interpretation is summarized in the following table:

Negative Control	Positive Control	Japanese Encephalitis FAM	Internal Control HEX	Interpretation
-	+	+	+/- ¹	Japanese Encephalitis RNA Detected
-	+	-	+ ²	Japanese Encephalitis RNA Not detected ²
-	+	-	- ²	Invalid ²

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct >40 or no signal)

¹ Sometimes, in the case of clinical samples, detection of the internal control is not necessary as the presence of a high initial copy number of the target nucleic acid may cause preferential amplification of the target nucleic acid. The internal control (IC) either shows or does not show an amplification signal (Ct ≤40 or no signal).

² In case the Japanese Encephalitis target genes are negative, the Internal Control should show an amplification signal with Ct ≤ 35. In case of no signal or Ct value > 35 of the internal control, the result is considered "invalid" and retesting is required. It is recommended to repeat the RT-qPCR by diluting the RNA sample 1:10 and/or 1: 100, or to re-extract and repeat the assay to check for possible failure of the extraction procedure and/or inhibition problems.

If the result obtained is unclear or doubtful, it is necessary to check that all steps have been performed correctly, verify the correct performance of each RT-qPCR step, check all parameters, the sigmoid shape of the curve and the fluorescence intensity. If the situation is not resolved, it is recommended to repeat the test preferably in duplicate depending on the available material: (a) obtain a new specimen, re-extract, and re-test (ideal), (b) re-extract and re-test another aliquot of the same sample, or (c) repeat the RT-qPCR with the same isolated RNA sample.

The test results should be evaluated by a healthcare professional, in the context of the medical history, clinical symptoms and other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR Japanese Encephalitis was evaluated with QCMD programs of tropical viruses (West Nile Virus and Dengue) of the years 2014, 2017 and 2018. These programs contained 6 positive samples for Japanese encephalitis and were detected correctly. The results were compared with the final reports.

These results indicate that Vitassay qPCR Japanese Encephalitis shows high agreement to detect Japanese encephalitis virus.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Japanese Encephalitis cDNA template (from 10^7 to 10^1 copies/reaction). This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Japanese Encephalitis was tested within the panel of following microorganisms, where no cross-reactivity was observed between any of the species:

Cross-reactivity assay		
<i>Candida albicans</i>	Parechovirus Type 3	Human herpesvirus 6 (HHV6) Type B
<i>Escherichia coli</i> 0.1285;O18:H7:K1	Citomegalovirus strain AD169	Zika Virus (African strain)
<i>Haemophilus influenzae</i> MinnA	Varicella Zoster Virus Ellen	Zika Virus (Asian strain PF13/251013-18)
<i>Enterococcus durans</i>	Varicella Zoster Virus OKA	Zika Virus (French Polynesian strain 11474/16)
<i>Enterococcus faecalis</i>	Coxsackievirus A24	Chikungunya virus strain S27 Petersfield
<i>Enterococcus faecium</i> serotype 11	Coxsackievirus A9	Dengue 1 virus strain Hawaii
<i>Cryptococcus gatti</i> Z156	Coxsackievirus B3	Dengue 2 virus strain New Guinea C
<i>Listeria monocytogenes</i> serovar 1/2b	Echovirus 30	Dengue 3 virus strain H87
<i>Listeria monocytogenes</i> serovar 4b	Echovirus Type 11	Dengue 4 virus strain H241
<i>Listeria innocua</i> serovar 6a	Enterovirus 68	St Louis Encephalitis virus strain 17D
<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	Enterovirus 71	West Nile virus strain NY99
<i>Neisseria meningitidis</i> serogroup A	Epstein Barr virus	West Nile virus Heja
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Herpes simplex virus 1 (HSV1) MacIntyre strain	West Nile virus Ug37
<i>Streptococcus agalactiae</i> Z019	Herpes simplex virus 2 (HSV2) MS	Tick-Borne encephalitis virus (cepa Neudorfl)
<i>Streptococcus pneumoniae</i> Z022	Human herpesvirus 6 (HHV6) strain Z29	Yellow Fever Virus (17D strain)
<i>Plasmodium falciparum</i> (3D7 strain)	Human herpesvirus 6 (HHV6) Type A	

Analytical reactivity

The reactivity of Vitassay qPCR Japanese Encephalitis was confirmed by the real-time amplification using Japanese Encephalitis as template, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR Japanese Encephalitis has been validated on the following equipment:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I

I: For Rotor-Gene® Q thermocycler the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

For more information regarding thermocycler compatibility based on bibliographic data and the most common detection channels contact the manufacturer (info@vitassay.com).

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTLite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Limitations

- This test provides a presumptive diagnosis of Japanese Encephalitis infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with blood samples and CSF. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper nucleic acids from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Possible crosstalk phenomena might be observed in empty channels of the thermocycler (if there is no target to be detected), therefore, when interpreting the results, it is necessary to select only those channels where the targets amplify.
- Occasionally, target template nucleic acids with a copy number below the detection limit can be detected, but the results may not be reproducible.
- This test does not provide quantitative values nor indicate the number of organisms present. It is not possible to correlate the Ct values obtained in the PCR with the concentration of the sample, as these depend on the thermal cycler used and the run itself. This is a qualitative test.
- There is a possibility of false positive results due to cross-contamination by Japanese Encephalitis, either by the positive control reconstitution, by samples containing high concentrations of target nucleic acids or contamination due to PCR products from previous reactions.

- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including nucleic acids extraction); c) nucleic acids degradation during specimen shipping, storage, and/or preparation; d) microorganism load is below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect the detection of new or unknown strains.
- The detection of nucleic acids does not imply the presence of viable and/or infectious microorganisms or that these microorganisms are the causative agents of clinical symptoms.
- Negative results do not preclude Japanese encephalitis infection and should not be the sole basis for a patient management/treatment decision. The optimal specimen types for Japanese encephalitis identification and/or the most appropriate stage of infection for collection have not yet been established. Consider collecting multiple samples from the same patient at different times, which may increase the likelihood of detecting the microorganism(s).
- If the patient's clinical data, laboratory tests and epidemiological studies suggest possible infection by these microorganisms, and other diseases have been discarded, a false negative result might not be dismissed, and additional tests should be performed.
- The fluorescence values may vary due to multiple factors such as PCR equipment (even if it is the same model), extraction system, sample type, and pretreatment, among others.
- The Positive Control and Negative Control must not be extracted.
- It is recommended to include a sample blank during the extraction and RT-qPCR process to rule out possible contamination and to avoid missing probable false positive results.
- The use of an appropriate control material such as reference RNA or viral RNA particles is recommended to control the reverse transcription step. These controls are not included in the kit.
- Detection of some circulating genotypes belonging to genotype IV (GIV) of JEV may be compromised.

Bibliography/Bibliografia


1. “Centers for Disease Control and Prevention (CDC), “Japanese Encephalitis”, 07 December 2023. Website: www.cdc.gov/japaneseencephalitis/index.html.
2. Roberts A., Gandhi S. (2020). Japanese encephalitis virus: a review on emerging diagnostic techniques. *Front. Biosci. (Landmark Ed)*, 25(10), 1875–1893.
3. Sharma KB, Vrati S, Kalia M. (2021) Pathobiology of Japanese encephalitis virus infection. *Mol Aspects Med.* Oct;81:100994.
4. European Centre for Disease Prevention and Control (ECDC), “Factsheet about Japanese Encephalitis”, 16 March 2022. Website: www.ecdc.europa.eu/en/japanese-encephalitis/facts.
5. World Health Organization (WHO). Japanese Encephalitis, 09 May 2019. Website: [Japanese encephalitis \(who.int\)](http://japanese-encephalitis.who.int)
6. Centers for Disease Control and Prevention (CDC), Japanese Encephalitis Virus Transmission, 07 December 2022. Website: <https://www.cdc.gov/japaneseencephalitis/transmission/index.html>

Trademarks/Marcas registradas

All trademarks that may appear in this package insert are the property of their respective owners.

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