Vitassay

Real-Time PCR Kits

Vitassay qPCR Vaginosis

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- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^I

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Intended use

Vitassay qPCR Vaginosis allows the detection and differentiation of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Vaginosis 4x8-well strip, low profile	7041051
Vitassay qPCR Vaginosis 4x8-well strip, high profile	7042051

Reagents provided

In references 7041051 and 7042051:

Reference	Reagent/Material	Description	Colour	Amount
7041S051/ 7042S051	Vaginosis strips low/high profile	Lyoprotectors and Stabilizers	-	4 x 8-well strips
		dNTPs		
		Primers and Probes		
		Enzymes		
7C051	Vaginosis Positive Control	Non-infectious synthetic lyophilized DNA	red	1 vial
7001A	PCR grade water	RNAsa/DNAsa free water	white	1 vial x 1 mL
7002B	Resuspension buffer	Saline Solution	green	1 vial x 1,8 mL
		Buffer (TRIS, pH)		
7003N	Negative control	RNAsa/DNAsa free water	yellow	1 vial x 1 mL
7004O	Optical caps	Optical caps for sealing strips during thermal cycling	-	4 x 8-cap strips

Material and equipment required, not provided

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL, 100-1000 µL)
- Filter tips
- Powder-free disposal gloves

Transport and storage

- Transport and storage of the kits (before use) can be performed at 2-40°C until the expiry date stated on the label. Once in use, see storage conditions for each component. Avoid shaking during transport to prevent liquid leakage.
- The resuspended positive control should be stored at -20°C until the expiry date stated on the kit label. The stability of the positive control has been validated after 6 freeze-thaw cycles. However, it is recommended to distribute in aliquots to avoid repeated freeze-thaw cycles.
- For proper use of the kit, it is recommended that unneeded wells are not rehydrated. Unused wells can be easily cut before removing the protective foil and stored in the foil bags provided with the silica gel inside and closed with the zip closure to avoid contact with light or moisture and stored at 2-40°C for the shelf life indicated on the kit label. Once the wells have been rehydrated, nucleic acid addition and qPCR should be performed immediately.
- Once opened, the Resuspension buffer, PCR grade water and Negative Control components should be stored at 2-40°C for the shelf life indicated on the kit label.
- Keep all reagents in the darkness.

Summary

Vaginosis is a bacterial dysbiosis commonly affecting women of reproductive age. Described in the mid-50s by Gardner and Dukes, who initially incorrectly named the condition as *Haemophilus vaginalis* vaginitis (Gardner and Dukes, 1955). It was later seen that the causing bacteria was now from the genus *Haemophilus*, and it was renamed as *Gardnerella vaginalis*.

Nowadays the accurate aetiology of vaginosis is not known, and it is believed to be caused by an imbalance between the overgrowth of opportunistic bacteria and a decrease in the presence of *Lactobacilli*, the healthy flora. It is thought that it might be sexually transmitted, especially in women having sex with women (Bagnall and Rizzolo, 2017; Coudray and Madhivanan, 2020).

Typically, women with vaginosis present an abnormal vaginal discharge, malodor, increased vaginal pH and itching, which are also common in other vaginal infections. It has been seen that prevalence is of around 23-29% across regions, and it can be even reach more than 50% in some specific populations. (Bautista et al., 2016; Peebles et al., 2019)

Candida albicans

Candida albicans is a eukaryotic microorganism with a great adaptability due to the extraordinary capacity of undergoing a morphogenic change from a round-ovoid yeast cell to a hyphal mycelial growing organism. The yeast form (Y) is usually associated with its commensal lifestyle, whereas the hyphal (H) form is more related to pathogenicity (Cassone et al., 2014).

C. albicans is the most aetiology of Vulvovaginal candidiasis (VVC), a type of vaginosis caused by *Candida* species. Up to the 75% of women population would experience at least one episode of VVC, and the recurrent episodes occurs in around 5%. Symptoms are nonspecific and include vaginal discharge, soreness, irritation, burning, and dyspareunia. There are several risk factors to take into account, including reproductive age, pregnancy, hormone replacement therapy, antibiotic use, immunosuppression, uncontrolled diabetes, oral contraceptive pills, frequent sexual intercourse, and receptive oral sex (Brown et al., 2019).

Candida spp. not only colonise the genital tract of women; it is also found in extragenital areas, like the oral cavity and the rectum. Despite not being considered to cause classical STIs, *Candida* can be sexually transmitted.

It is necessary to identify the *Candida* species responsible for infections in all patients presenting with VVC especially those with refractory infections. Such identification will influence antimycotic agents' selection and duration of therapy. Typically, identification is performed by microscopy techniques, which is not a very specific technique. Molecular biology techniques, such as PCR, are now being used to identify *Candida* species and other sexually transmitted pathogens (Makanjuola et al., 2018).

Gardnerella vaginalis

As aforementioned, *Gardnerella vaginalis* was the first microorganism described as a causative agent of bacterial vaginosis (BV), but it can also be isolated from the vagina of healthy women. It is believed that certain strains of *G. vaginalis* are more pathogenic than others. Recent advances in BV pathogenesis research have determined the existence of 13 different genomospecies within the genus *Gardnerella* (Potter et al., 2019).

Amsel criteria is still used clinically for the diagnosis of BV. Three of four of the following (Amsel) criteria are generally considered to support the diagnosis of BV: (1) a thin, homogenous discharge, (2) pH over 4.5, (3) the presence of epithelial cells that appeared to be coated in bacteria in wet mount, and (4) the detection of a fishy odor, with or without treatment of the sample with 10% potassium hydroxide (Morrill et al., 2020). However, some of the symptoms can easily be also due to other microorganisms.

BV disease is linked with a compelling amount of obstetric and gynecological complications which includes spontaneous abortion, preterm delivery, subclinical pelvic inflammatory disease, post-Caesarean delivery wound infections, chorioamnionitis, premature rupture of membranes, postpartum endometritis, and postsurgical infections. Moreover, bacterial vaginosis is associated with an increased risk of acquisition of sexually transmitted infections including human immunodeficiency virus (HIV).

With the aim of improving the detection of *Gardnerella*, and hence, a proper diagnosis, molecular methods have been developed. Some of the qPCRs are even able to distinguish between the different bacteria species of this genus (Balashov et al., 2014) or between different others linked to BV (Singh et al., 2020)

Trichomonas vaginalis

T. vaginalis is a one-celled flagellated protozoan and human parasite causing the most prevalent non-viral sexually transmitted infection worldwide. Trichomoniasis is a widespread, global health concern and occurring at an increasing rate. For example, CDC has calculated that in the US, there were 2.6 million trichomoniasis infections in 2018 (CDC).

In females, the genital infections can display a plethora of manifestations such as including vaginitis and cervicitis. On the other hand, infections in males are generally asymptomatic, with sporadic mild urethritis or prostatitis. Importantly, during the last decades, *T. vaginalis* infection has been associated with complications such as cervical and prostate cancer, adverse pregnancy outcomes and an increased likelihood of HIV infection. Thus, a proper diagnosis is required but is a tricky matter, since the symptoms are similar to the ones caused by *N. gonorrhoeae* and *Mycoplasma genitalium*, for example (Edwards et al., 2014). An imbalance of the normal

lactobacilli vaginal flora is observed in *T. vaginalis*, as well as an increased of vaginal pH, a common feature observed in *Gardnerella* and other vaginosis-related bacteria overgrowth.

Various methods have been used to diagnose trichomoniasis, such as wet mount and culture. Nevertheless, PCR is considered to be highly sensitive and specific technology that also allows the simultaneous differentiation with other sexually transmitted microorganisms (Noh et al., 2019).

Test principle

Vitassay qPCR Vaginosis test is based on the real-time amplification of *T. vaginalis*-specific 2-kb repeated sequence for *Trichomonas vaginalis*, a conserved region of 5.8S rRNA gene for *Candida albicans* and 16S rRNA gene for *Gardnerella vaginalis*. After DNA extraction, the presence of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Vaginosis test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The detection channels of the target sequences are described as follows:

Target	Detection channel
<i>Trichomonas vaginalis</i> -specific 2-kb repeated sequence	FAM
5.8S rRNA gene (<i>Candida albicans</i>)	ROX
16S rRNA gene (<i>Gardnerella vaginalis</i>)	Cy5
Internal Control (IC)	HEX, VIC, or JOE (depending on the equipment used)

Precautions

- For professional *in vitro* diagnostic use.
- Do not use the kit after the expiration date.
- Do not use the kit without having read and understood the information on procedures, precautions, and limitations provided in the instructions for use.
- Do not mix components from different kits and/or lots.
- Do not use if the package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Use the strips immediately after removing the protective foil to protect the reaction mixture from sunlight. It is not recommended to leave the strips open without resuspending. Once the required wells have been rehydrated, the qPCR assay should be performed immediately.
- It is recommended to cut the required/necessary wells before removing the protective seal from the entire strip. Store the rest inside the bag with the desiccant.
- Close protective pouches of reagents promptly with the zip seal after each use (do not remove the desiccant from the pouches). Remove excess air from the bags before sealing.
- Do not use the product near solvents to avoid deterioration of the label/markings.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Rehydrate the vial immediately after opening the Positive Control envelope and use it in the PCR reaction or store it properly. Leaving the vial open without resuspending may cause deterioration of the reagent.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area. Once the test is completed wash your hands thoroughly.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.

- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended. A new tip should be used for each sample. It is necessary to change gloves before handling the reagents.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the test.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- The product Vitassay qPCR Vaginosis has only been validated with the equipment mentioned in Section Compatibles real-time PCR equipment of this Instructions for Use.
- You must ensure that the conditions are set up in the thermal cycler according to the instructions in section "Programme your thermocycler".
- Refer to Safety data sheets, available on request.

Procedures

Specimen collection, transport, and conservation

DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use. The assay has been validated with the following extraction kits:

- MagDEA Dx SV kit using the magLEAD® 12gC instrument (Precision System Science Co).
- ZP02014 Mag Purix Plant DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- Invisorb® Spin Universal Kit (Invitex).

Positive control preparation

Reconstitute the lyophilized Vaginosis Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains a high-copy number template and entails a very significant contamination risk. Therefore, we recommend opening and manipulating it in a separate area of the laboratory away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
It is recommended not to rehydrate wells that are not needed. Unused wells can be easily cut before removing the protective foil and stored as indicated in the Transport and Storage Section.
- Peel off protective aluminium seal from the strips to be used.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA extracted from each sample, of negative control (yellow tube) or of reconstituted positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the optical caps provided. It is recommended spin down briefly.
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the Annealing/extension step (*) through the FAM (*Trichomonas vaginalis*), ROX (*Candida albicans*), Cy5 (*Gardnerella vaginalis*) and HEX, JOE, or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel.

Analysis and interpretation of results

The analysis of the results is performed by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control signal emission.

For each assay, it is recommended to set the threshold values for each channel (target) independently by the end-user. To do so, select the well corresponding to the positive control for each channel and set the threshold value within the exponential phase of the fluorescence curve and above any background signal (below the baseline). In case a threshold is set automatically by the equipment used, check, and verify that it is adjusted to the positive control or adjust it manually. Once the threshold value is set, the test samples can be interpreted.

The chemical nature of the different fluorophores can affect the threshold value of the channels, just as different signal intensities can influence the threshold value for different instruments.

For a valid diagnostic test run, before analyzing the result of the clinical samples, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve (Ct ≤20-30) in FAM (*Trichomonas vaginalis*), ROX (*Candida albicans*), and Cy5 (*Gardnerella vaginalis*), which validates the reaction.

Negative control

The negative control included in each run must show the absence of signal (Ct >40 or no signal) in FAM (*Trichomonas vaginalis*), and ROX (*Candida albicans*), and must also show absence of signal (Ct >36 or no signal) in Cy5 (*Gardnerella vaginalis*), which validates the reaction.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in the positive control and negative control wells.

The experiment is invalid if there is no signal in the positive control or is amplification signal in the negative control, except for the HEX, JOE or VIC channel which should show amplification signal (corresponding to Internal Control). In either case, the assay must be repeated.

Once the results of the controls have been validated, for the interpretation of the sample results **select only the channels where the targets are detected**. The result interpretation is summarized in the following table:

Negative Control	Positive Control	<i>Trichomonas vaginalis</i> (FAM)	<i>Candida albicans</i> (ROX)	<i>Gardnerella vaginalis</i> (Cy5)	Internal Control (HEX)	Interpretation
-	+	+	+	+ ³	+ ¹	<i>T. vaginalis</i> , <i>C. albicans</i> and <i>G. vaginalis</i> DNA Detected
-	+	-	-	- ³	+ ²	Targets' DNA Not detected
-	+	+	-	- ³	+ ¹	<i>T. vaginalis</i> DNA Detected, <i>C. albicans</i> and <i>G. vaginalis</i> DNA Not detected
-	+	+	+	- ³	+ ¹	<i>T. vaginalis</i> and <i>C. albicans</i> DNA Detected, <i>G. vaginalis</i> DNA Not detected
-	+	+	-	+ ³	+ ¹	<i>T. vaginalis</i> and <i>G. vaginalis</i> DNA Detected, <i>C. albicans</i> DNA Not detected
-	+	-	+	- ³	+ ¹	<i>C. albicans</i> DNA Detected, <i>T. vaginalis</i> and <i>G. vaginalis</i> DNA Not detected
-	+	-	+	+ ³	+ ¹	<i>C. albicans</i> and <i>G. vaginalis</i> DNA Detected, <i>T. vaginalis</i> DNA Not detected
-	+	-	-	+ ³	+ ¹	<i>G. vaginalis</i> DNA Detected, <i>T. vaginalis</i> and <i>C. albicans</i> DNA Not detected
-	+	-	-	- ³	- ²	Invalid

Positive (+): Amplification signal (Ct ≤40), except for ^{1,2,3}.

Negative (-): No amplification signal (Ct >40 or no signal), except for ^{1,2,3}.

¹ The internal control (IC) should show an amplification signal (Ct ≤40). Sometimes, in the case of clinical samples, detection of the internal control does not occur as the presence of a high initial copy number of the target nucleic acid may cause preferential amplification of the target-specific nucleic acids. In case of one of the targets presents a high number of copies, the rest of targets remain negative, and the Internal Control does not show amplification signal, false negatives of other targets cannot be discarded. It is recommended to repeat the qPCR by diluting the DNA sample 1:10.

² In case the *T. vaginalis*, *C. albicans*, and *G. vaginalis* target genes are negative, the Internal Control should show an amplification signal with Ct ≤ 35. In case of no signal or Ct value > 35 of the internal control, the result is considered "invalid" and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1: 100, or to re-extract and repeat the assay to check for possible failure of the extraction procedure and/or inhibition problems.

³ For *G. vaginalis* a cut-off Ct of 36 has been established. Thus:

(+) Positive: Amplification signal (Ct ≤36)

(-) Negative: No amplification signal (Ct >36 or no signal)

In case of an unclear or doubtful interpretation result, it is necessary to check that all steps have been performed correctly, the extraction process used, check the instructions for use, verify the correct performance of each qPCR step, check all parameters, the sigmoid shape of the curve and the fluorescence intensity. If the situation is not resolved, it is recommended to repeat the test preferably in duplicate depending on the available material: (a) obtain a new specimen, re-extract, and re-test (ideal), (b) re-extract and re-test another aliquot of the same sample, or (c) repeat the qPCR with the same isolated DNA sample.

NOTE: A positive test result for *C. albicans* and/or *G. vaginalis* does not necessarily indicate vulvovaginal candidiasis or bacterial vaginosis (respectively), since those microorganisms can be part of the vaginal microbiota in healthy women. The test results should be evaluated by a healthcare professional, in the context of the medical history, clinical symptoms and other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR Vaginosis clinical sensitivity and specificity was tested within two evaluations.

A first study was conducted to evaluate the kit's functionality with leftovers of 232 urogenital tract samples. These samples were collected from symptomatic patients (mainly women but also men) with suspected infection or dysbiosis. After routine analysis by microscopy and broth culture technique, the leftovers were extracted with MagDEA® Dx reagents using the MagLEAD® 12gC nucleic acid extractor (PSS instruments) and analysed with the CFX96™ Touch Real-Time PCR Detection System (Bio-Rad). Discrepancies between the initial diagnosis and Vitassay qPCR were tested with a commercially available real-time PCR panel kit. The obtained results after the analysis of discrepant results were as follows:

Microorganism	OPA	TP	TN	FP	FN	SE	SP	PPV	NPV
<i>C. albicans</i>	99.13%	88	142	0	2	0.97 (0.91-0.99)	1 (0.96-1)	1 (0.94-1)	0.98 (0.94-0.99)
<i>G. vaginalis</i>	99.13%	152	78	2	0	1 (0.96-1)	0.97 (0.90-0.99)	0.98 (0.94-0.99)	1 (0.94-1)
<i>T. vaginalis</i>	99.56%	54	177	1	0	1 (0.91-1)	0.99 (0.96-0.99)	0.98 (0.89-0.99)	1 (0.97-1)

Results of Vitassay qPCR Vaginosis compared to the reference method in urogenital samples after the analysis of discrepant results.

A second study was conducted to evaluate the kit's clinical sensitivity and specificity with leftovers of 240 vaginal samples from women with suspected abnormal vaginal discharge. After routine analysis by culture, microscopy, and two molecular CE-marked commercially available kits, the leftovers were extracted with MagDEA® Dx reagents using the MagLEAD® 12gC nucleic acid extractor (PSS instruments) and analysed with the DTPrime real-time PCR instrument (DNA-Technology). Discrepancies between the initial diagnosis and Vitassay qPCR were tested with a different method from the one initially employed. The results obtained after the analysis of discrepant results were as follows:

Microorganism	Overall agreement	TP	TN	FP	FN	SE	SP	PPV	NPV
<i>C. albicans</i>	0.97 (0.95-0.99)	74	161	5	0	1 (0.95-1)	0.97 (0.93-0.99)	0.93 (0.85-0.97)	1 (0.97-1)
<i>G. vaginalis</i>	0.94 (0.90-0.96)	164	62	14	0	1 (0.97-1)	0.81 (0.71-0.90)	0.92 (0.87-0.95)	1 (0.94-1)
<i>T. vaginalis</i>	1 (0.98-1)	46	194	0	0	1 (0.92-1)	1 (0.98-1)	1 (0.90-1)	1 (0.97-1)

Results of Vitassay qPCR Vaginosis compared to the reference methods in vaginal samples after the analysis of discrepant results.

To sum up, the clinical performance of Vitassay qPCR Vaginosis was evaluated by comparing it with conventional and molecular methods performed for routine diagnosis in hospital settings. The mean values for sensitivity and specificity obtained from the clinical performance of Vitassay qPCR Vaginosis, were as follows:

Target	Conventional methods (culture +microscopy)		Molecular reference method (real time PCR)	
	SE	SP	SE	SP
<i>T. vaginalis</i>	1 (0.93-1)	0.99 (0.96-1)	1 (0.92-1)	1 (0.98-1)
<i>C. albicans</i>	0.98 (0.95-0.99)	0.97 (0.93-0.99)	1 (0.95-1)	0.97 (0.93-0.99)
<i>G. vaginalis</i>	1 (0.98-1)	0.89(0.83-0.94)	1 (0.97-1)	0.81 (0.71-0.90)

Mean values of the Vitassay qPCR Vaginosis results for *T. vaginalis*, *C. albicans* and *G. vaginalis* compared to conventional methods and molecular reference methods.

Vitassay qPCR Vaginosis represents an efficient tool for the detection of *C. albicans*, *G. vaginalis* and *T. vaginalis* DNA, as it shows a high sensitivity and specificity in genitourinary tract samples.

Analytical sensitivity

The analytical sensitivity was determined by analysis of serial dilutions of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* standards (10^7 to 10^1 copies/reaction). This assay has a detection limit of ≥ 10 DNA copies per reaction for *T. vaginalis*, *C. albicans*, and *G. vaginalis*.

Analytical specificity

The analytical specificity for *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* was confirmed by testing the panel of following microorganisms, where no cross-reactivity was observed between any of the microorganisms, except for the targeted microorganisms:

Cross-reactivity assay					
<i>Candida albicans</i>	-/+	<i>Escherichia coli</i> 0.1285; O18:H7:K1	-	Human papillomavirus genotypes HPV16 and HPV18	-
<i>Candida parapsilosis</i>	-	<i>Gardnerella vaginalis</i>	-/+	HSV-1 and HSV-2	-
<i>Candida tropicalis</i>	-	<i>Haemophilus influenzae</i>	-	<i>Serratia marcescens</i>	-
<i>Candida glabrata</i>	-	<i>Haemophilus ducrey</i>	-	<i>Staphylococcus aureus</i>	-
<i>Candida krusei</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Stenotrophomonas maltophilia</i>	-
<i>Candida dubliniensis</i>	-	<i>Klebsiella pneumoniae</i>	-	<i>Streptococcus pneumoniae</i>	-
<i>Aspergillus fumigatus</i>	-	<i>Listeria ivanovii</i>	-	<i>Streptococcus agalactiae</i>	-
<i>Bacteroides fragilis</i>	-	<i>Listeria monocytogenes</i>	-	<i>Trichomonas vaginalis</i>	-/+
<i>Acinetobacter baumannii</i>	-	<i>Listeria Innocua</i>	-	<i>Treponema pallidum</i>	-
<i>Chlamydia trachomatis</i> (LGV)	-	<i>Mycoplasma genitalium</i>	-	<i>Ureaplasma parvum</i>	-
<i>Chlamydia trachomatis</i> (SW)	-	<i>Mycoplasma hominis</i>	-	<i>Ureaplasma urealyticum</i>	-
<i>Enterobacter cloacae</i>	-	<i>Neisseria gonorrhoeae</i>	-	Hepatitis A	-
<i>Enterobacter aerogenes</i>	-	<i>Neisseria meningitidis</i>	-	Cytomegalovirus	-
<i>Enterococcus faecalis</i>	-	<i>Proteus mirabilis</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Enterococcus faecium</i>	-				

Analytical reactivity

Analytical reactivity of Vitassay qPCR Vaginosis for *C. albicans* was evaluated against *Candida albicans*, showing positive result.

Analytical reactivity of Vitassay qPCR Vaginosis for *G. vaginalis* was evaluated against *Gardnerella vaginalis*, showing positive result.

Analytical reactivity of Vitassay qPCR Vaginosis for *T. vaginalis* was evaluated against *Trichomonas vaginalis* showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Vaginosis has been validated on the following equipment:

- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)

Vitassay qPCR Vaginosis has been tested on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^I

I: For Rotor-Gene® Q thermocycler the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

For more information regarding thermocycler compatibility based on bibliographic data and the most common detection channels contact the manufacturer (info@vitassay.com).

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR instrument (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time PCR instrument (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Limitations

- The test results should be evaluated by a healthcare professional, in the context of the medical history, clinical symptoms and other diagnostic tests.
- A positive test result for *C. albicans* and/or *G. vaginalis* does not necessarily indicate vulvovaginal candidiasis or bacterial vaginosis (respectively), since *C. albicans* and *G. vaginalis* can be part of the vaginal microbiota in healthy women.
- This assay has been validated with DNA extracted from vaginal and urethral swabs. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper nucleic acids from clinical specimens must be extracted.
- Possible crosstalk phenomena might be observed in empty channels of the thermocycler (if there is no target to be detected), therefore, when interpreting the results, it is necessary to select only those channels where the targets amplify.
- Occasionally, target template nucleic acids with a copy number below the detection limit can be detected, but the results may not be reproducible.
- This test does not provide quantitative values nor indicate the number of organisms present. It is not possible to correlate the Ct values obtained in the PCR with the concentration of the sample, as these depend on the thermal cycler used and the run itself. This is a qualitative test.
- There is a possibility of false positive results due to cross-contamination by *C. albicans*, *T. vaginalis* and/or *G. vaginalis*, either by the positive control reconstitution, by samples containing high concentrations of target nucleic acids or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including nucleic acids extraction); c) nucleic acids degradation during specimen shipping, storage, and/or preparation; d) microorganism load is below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference substances (an interference study evaluating the effect of vaccines, some antiviral therapeutics, antibiotics, chemotherapeutics,

immunosuppressant drugs or antifungals used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect the detection of new or unknown strains.

- The detection of nucleic acids does not imply the presence of viable and/or infectious microorganisms or that these microorganisms are the causative agents of clinical symptoms.
- Negative results do not preclude *C. albicans*, *T. vaginalis* and/or *G. vaginalis* infection and should not be the sole basis for a patient management/treatment decision. The optimal specimen types for *C. albicans*, *T. vaginalis* and/or *G. vaginalis* identification and/or the most appropriate stage of infection for collection have not yet been established. Consider collecting multiple samples from the same patient at different times, which may increase the likelihood of detecting the microorganism(s).
- If the patient's clinical data, laboratory tests and epidemiological studies suggest possible infection by these microorganisms, and other diseases have been discarded, a false negative result might not be dismissed, and additional tests should be performed.
- The fluorescence values may vary due to multiple factors such as PCR equipment (even if it is the same model), extraction system, sample type, and pretreatment, among others.

Bibliography/Bibliografía


1. Bagnall, Paulette MPAS, PA-C; Rizzolo, Denise PA-C, PhD Bacterial vaginosis, JAAPA: December 2017 - Volume 30 - Issue 12 - p 15-21
2. Balashov SV, Mordechai E, Adelson ME, Gyax SE. Identification, quantification and subtyping of Gardnerella vaginalis in noncultured clinical vaginal samples by quantitative PCR. J Med Microbiol. 2014 Feb;63(Pt 2):162-175.
3. Bautista CT, Wurapa E, Sateren WB, Morris S, Hollingsworth B, Sanchez JL. Bacterial vaginosis: a synthesis of the literature on etiology, prevalence, risk factors, and relationship with chlamydia and gonorrhea infections. Mil Med Res. 2016 Feb 13;3:4.
4. Brown SE, Schwartz JA, Robinson CK, O'Hanlon DE, Bradford LL, He X, Mark KS, Bruno VM, Ravel J, Brotman RM. The Vaginal Microbiota and Behavioral Factors Associated With Genital Candida albicans Detection in Reproductive-Age Women. Sex Transm Dis. 2019 Nov;46(11):753-758.
5. Cassone A. (2014). Vulvovaginal Candida albicans infections: pathogenesis, immunity and vaccine prospects. BJOG : an international journal of obstetrics and gynaecology, 122(6), 785–794.
6. CDC - Trichomoniasis Statistics. Centers for Disease Control and Prevention, Centers for Disease Control and Prevention, 5 Apr. 2021, <https://www.cdc.gov/std/trichomonas/stats.htm>.
7. Coudray MS, Madhivanan P. Bacterial vaginosis-A brief synopsis of the literature. Eur J Obstet Gynecol Reprod Biol. 2020 Feb;245:143-148.
8. Edwards, T., Burke, P., Smalley, H., & Hobbs, G. (2014). Trichomonas vaginalis: Clinical relevance, pathogenicity and diagnosis. Critical Reviews in Microbiology, 42(3), 406–417.
9. Gardner HL, Dukes CD. Haemophilus vaginalis vaginitis: a newly defined specific infection previously classified non-specific vaginitis. Am J Obstet Gynecol. 1955;69(5):962–76.
10. Makanjuola O, Bongomin F, Fayemiwo SA. An Update on the Roles of Non-albicansCandida Species in Vulvovaginitis. J Fungi (Basel). 2018 Oct 31;4(4):121.
11. Morrill S, Gilbert NM, Lewis AL. Gardnerella vaginalis as a Cause of Bacterial Vaginosis: Appraisal of the Evidence From in vivo Models. Front Cell Infect Microbiol. 2020 Apr 24;10:168.
12. Noh CS, Kim SS, Park SY, Moon HS, Hong Y, Ryu JS. Comparison of Two PCR Assays for Trichomonas vaginalis. Korean J Parasitol. 2019 Feb;57(1):27-31.
13. Peebles, Kathryn; Velloza, Jennifer MPH; Balkus, Jennifer E., McClelland, R. Scott; Barnabas, Ruanne V. MBChB. High Global Burden and Costs of Bacterial Vaginosis: A Systematic Review and Meta-Analysis, Sexually Transmitted Diseases: May 2019 - Volume 46 - Issue 5 - p 304-311
14. Potter RF, Burnham CD, Dantas G. In Silico Analysis of Gardnerella Genomespecies Detected in the Setting of Bacterial Vaginosis. Clin Chem. 2019 Nov;65(11):1375-1387.
15. Singh, R., Ramsuran, V., Mitchev, N. et al. Assessing a diagnosis tool for bacterial vaginosis. Eur J Clin Microbiol Infect Dis 39, 1481–1485 (2020).

Trademarks/Marcas registradas

All trademarks that may appear in this package insert are the property of their respective owners.

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Change Control / <i>Control de cambios</i>		
Version Nº / <i>Versión Nº</i>	Changes / <i>Cambios</i>	Date / <i>Fecha</i>
00	Original versión / <i>Versión original</i>	27/08/2020
01	New format, editorial improvement and better explanation of some sections, more information regarding material, reagents, transport and conservation, bibliographic information update, addition of some precautions, recommendations, clarifications and limitations, update with follow-up evaluations, deletion of Attached I and II, addition of symbols not covered / <i>Nuevo formato, mejora editorial y mejor explicación de algunas secciones, más información sobre el material, los reactivos, el transporte y la conservación, actualización de la información bibliográfica, adición de algunas precauciones, recomendaciones, aclaraciones y limitaciones, actualización con evaluaciones de seguimiento, supresión de los Anexos I y II y adición de símbolos no contemplados.</i>	24/04/2024

Real-Time PCR Kits

VA Vitassay

identifying pathogens worldwide

Vitassay Healthcare, S.L.U

Parque Tecnológico Walqa

Ctra. N-330 Km. 566

22197 Huesca (Spain)

Ph. (+34) 974 001 193

info@vitassay.com

www.vitassay.com

F09-81 Rev.00