For professional in vitro diagnostic use only.

INTENDED USE

Vitassay Entamoeba is a rapid, immunochromatographic, one step assay for the qualitative detection of *Entamoeba histolytica* and *Entamoeba dispar* in human stool samples.

Simple, non-invasive and highly sensitivity screening assay to make a presumptive diagnosis of *Entamoeba* infection (amoebiasis).

INTRODUCTION

Entamoeba histolytica is a protozoan parasite that invades through the intestinal epithelium via a unique cell biologic process called trogocytosis, resulting in diarrhea, dysentery, and extra intestinal disease.

The genus Entamoeba contains many species, six of which are found in the human intestinal tract: *Entamoeba histolytica*, *Entamoeba dispar*, *Entamoeba moshkovskii*, *Entamoeba coli*, *Entamoeba hartmanni*, and *Entamoeba polecki*. Of these species, only *E. histolytica* is associated with pathological injuries; the others are considered to be nonpathogenic species.

Food and drink contaminated with faeces containing the cysts is a common source of infection. Most cases arise from human carriers, or cyst passers, who pass cysts in formed or semiformed stools.

Because transmission is frequently associated with contaminated food and water, young infants are not expected to develop amebiasis very often. More severe disease is associated with young age, malnutrition and immunosuppression.

The motil (trophozoite) form of *E. histolytica* inhabits the human colon where it multiplies and differentiates into cysts that are released into the environment. In turn, these cysts are responsible for transmitting the infection to another host via the fecal-oral route. The parasite invades the intestinal mucosa and causes many forms of invasive disease, including dysentery. The parasite also exhibits bloodborne spreading and causes extraintestinal lesions, mainly liver abcesses. The latter form occurs only rarely. Invasive disease occurs when virulent trophozoites disrupt the mucoepithelial barrier by crossing the mucus layer, thereby damaging intestinal cells. This damage leads to inflammation and, consequently, dysentery. Nevertheless, the majority of infections seem to be asymptomatic.

PRINCIPLE

Vitassay Entamoeba is a qualitative immunochromatographic assay for the detection of *Entamoeba histolytica* and *Entamoeba dispar* in human stool samples.

The test line zone of the nitrocellulose membrane is pre-coated with monoclonal antibodies against Entamoeba.

During the process, the sample reacts with the antibodies against Entamoeba, forming conjugates. The mixture moves upward on the membrane by capillary action. If the sample is positive, antibodies present on the membrane (test line) capture the conjugate complex and a red line will be visible. Although the sample is positive or negative, the mixture continues to move across the membranes and the green control line always appears.

The presence of this green line (in the control zone (C)) indicates that sufficient volume is added; proper flow is obtained and serves as an internal control for the reagents.

PRECAUTIONS

- For professional in vitro use only.
- Do not use after expiration date.
- Do not use the test if its pouch is damaged.
- Specimens should be considered as potentially hazardous and handle in the same manner as an infectious agent. A new test must be used for each sample to avoid contaminations errors. Single use device.
- Tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with the skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use with the **Vitassay Entamoeba**. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/35.6-86°F) on the sealed pouch.

The test is stable until the expiration date printed on the sealed pouch.

The test must remain in the sealed pouch until use. Do not freeze.

VITASSAY

Entamoeba

Rapid test for the qualitative detection of Entamoeba histolytica and Entamoeba dispar in human stool samples.

IUE-7355036 Ed01 October 2017





MATERIALS

MATERIAL PROVIDED	MATERIAL REQUIRED BUT NOT PROVIDED
 25 tests/kit Vitassay Entamoeba. Instructions for use. 25 vials with diluent for the dilution sample. 	 Specimen collection container. Disposable gloves. Timer.


SPECIMEN COLLECTION

Collect sufficient quantity of feces: 1-2g or mL for liquid samples. Stool samples should be collected in clean and dry containers.


Samples can be stored in the refrigerator $(2-8^{\circ}C/35.6-46.4^{\circ}F)$ for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C (-4°F). Samples must be brought to room temperature before testing.

SPECIMEN PREPARATION

- 1. Remove the cap of the vial with diluent for the sample dilution (figure 1).
- 2. Use the stick to collect sufficient sample quantity. For solid stool, insert the stick in 4 different areas of the stool sample, taken approx. 125mg, (figure 2), and add it into the vial with diluent for the sample dilution. For liquid stool, take 125µL of the sample using a micropipette and transfer it into the vial with diluent for the sample dilution.
- 2. Close the vial with the diluent and stool sample. Shake vigorously the vial in order to assure good sample dispersion (figure 3).



Vial for sample dilution.



ROCEDORE


Allow the test, stool sample, controls and diluent to reach room temperature $(15-30^{\circ}C/59-86^{\circ}F)$ prior to testing. Do not open pouches until the performance of the assay.

- 1. Shake the vial with the sample to obtain a good sample dilution.
- 2. Remove the **Vitassay Entamoeba** from its sealed bag just before using it.
- 3. Take the vial containing the diluted sample, cut the end of the cap (figure 4) and dispense 4 drops in the circular window marked with the letter S (figure 5).
- Read the results at **10 minutes**. Do not read the results later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window with the stick. If it does not work, dispense a drop of diluent until seeing the liquid running through the reaction zone. 4

5

W



Dispense 4 drops in the circular window marked with the letter S.

INTERPRETATION OF THE RESULTS

WA Vitassay

Entamoeba

С

T

S

		These is as Estamosta
СТ	Only one green line in the control zone (C)	There is no Entamoeba histolytica or Entamoeba dispar presence. No infection caused by Entamoeba histolytica or E.
	POSITIVE	dispar. There is presence of Entamoeba
СТ	In addition to the green line (control line C), a red line appears, (test line T)	histolytica or Entamoeba dispar presence. Infection caused by Entamoeba histolytica or Entamoeba dispar.
ANY OTHER RESULTS		Invalid result, we recommend repeating the assay using the sample with another test. Note: Wrong procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. If the symptoms or situation persist, discontinue using the test kit and contact your local distributor.



Notes: The intensity of the red colored test line in the result line zone (T) will vary depending on the concentration of antigens in the specimen.

QUALITY CONTROL

Internal procedural control is included in **Vitassay Entamoeba**. Green line appearing in the results window is an internal control (C), which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- Vitassay Entamoeba must be carried out within 2 hours of opening the sealed bag.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary depending on the concentration of antigens.
- After one week of infection, the number of parasites in feces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- The use of other samples different from human samples has not been established.
- The quality of **Vitassay Entamoeba** depends on the quality of the sample; Proper fecal specimens must be obtained.
- Positive results determine the presence of Entamoeba histolytica and/or dispar in fecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical method or by microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- Negative results should not be considered as conclusive; it is possible that the concentration of antigens is lower than the detection limit. If symptoms or situation still persist, a Entamoeba histolytica and/or dispar determination should be carried out with another technique (for example microscopy).

EXPECTED VALUES

Entamoeba histolytica infection is the third-greatest parasitic disease responsible for death in the world after malaria and schistosomiasis. It affects approximately 180 million people, of whom 40000 to 110000 die each year. Entamoeba histolytica infections are worldwide, more common in the tropics and subtropics.

Particularly exposed to this parasite are people living and traveling to the tropical and subtropical zones (Asia, Africa, India, Indonesia, Mexico, South America, South Africa). These zones climatic conditions are optimal for the Protozoan cysts, which described the external environment can survive for many days. This fact contributes to the increase in the number of infections slider Entamoeba histolytica amongst people in tropical zones.

E. histolytica differs among countries, socio-economic and sanitary conditions, and populations. It is highly endemic throughout poor and socio-economically deprived communities in the tropics and subtropics. Environmental, socio-economic, demographic, and hygiene related behavior is known to influence the transmission and distribution of intestinal parasitic infections.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using **Vitassay Entamoeba** and these results were confirmed using another commercial qPCR test (VIASURE *Entamoeba histolytica* Real Time Detection Kit and VIASURE *Entamoeba dispar* Real Time Detection Kit, CerTest).

Results were as follows:

			<i>intamoeba h</i> ime Detectio	
		Positive	Negative	Total
Vitassay Entamoeba	Positive	4	1	5
	Negative	1	107	108
	Total	5	108	113

Vitassay Entamoeba vs VIASURE <i>Entamoeba histolytica</i> Real Time Detection Kit				
	Mean Value E. histolytica 95% confidence interval			
Sensitivity	80%	28.4-99.5%		
Specificity	99%	94.9-100%		
PPV	80%	28.4-99.5%		
NPV	97%	94.9-100%		

		VIASURE <i>Entamoeba dispar</i> Real Time Detection Kit		
		Positive	Negative	Total
Vitassay Entamoeba	Positive	15	1	16
	Negative	6	107	113
Entamoeba	Total	21	108	129

Vitassay Entamoeba vs VIASURE <i>Entamoeba dispar</i> Real Time Detection Kit				
Mean Value E. dispar 95% confidence interval				
Sensitivity	71%	47.8-88.7%		
Specificity	99%	94.9-100%		
PPV	94%	69.8-99.8%		
NPV	95%	88.8-98.0%		

	VIASURE <i>Entamoeba histo</i> Real Time Detection Kit VIASURE <i>Entamoeba dispa</i> Time Detection Kit		Kit and <i>spar</i> Real	
		Positive	Negative	Total
	Positive	19	2	21
Vitassay Entamoeba	Negative	7	214	221
	Total	26	216	242

Vitassay Entamoeba vs VIASURE Entamoeba histolytica Real Time Detection Kit and VIASURE Entamoeba dispar Real Time Detection

KIL KIL				
	Mean Value E. histolytica+ E. dispar	95% confidence interval		
Sensitivity	73%	52.2-88.4%		
Specificity	99%	96.7-100%		
PPV	90%	69.6-98.8%		
NPV	97%	93.6-98.7%		

The results showed that **Vitassay Entamoeba** has a high sensitivity and specificity to detect *Entamoeba histolytica* and *Entamoeba dispar*.

Cross reactivity

No cross reactivity was detected against gastrointestinal pathogens that are occasionally present in feces:

Campylobacter jejuni	Helicobacter pylori	Shigella boydii
Campylobacter coli	Listeria monocytogenes	Shigella dysenteriae
Clostridium difficile	Salmonella enteritidis	Shigella flexneri
Escherichia coli 0157:H7	Salmonella paratyphi	Shigella sonnei
Cryptosporidium parvum	Salmonella typhi	Staphylococcus aureus
Giardia lamblia	Salmonella typhimurium	



REFERENCES

1. HANS P. VERKERKE; BLAKE HANBURY; ABDULLAH SIDDIQUE; AMIDOU SAMIE; RASHIDUL HAQUE; JOEL HERBEIN; WILLIAM A. PETRI, JR. "Multisite Clinical Evaluation of a rapid test for Entamoeba histolytica in stool". Journal of Clinical Microbiology, February 2015, Vol. 53, Number 2, pp. 493-497.

2. THIAGO DOS SANTOS COMES; MARIANA COIMBRA GARCIA; FLAVIA DE SOUZA CUNHA; HELOISA WERNECK DE MACEDO; JOSÉ MAURO PERALTA; REGINA HELENA SARAMAGO PERALTA. "Differential diagnosis of Entamoeba spp. in clinical samples using SYBR Green real-time Polymerase Chain Reaction". Hindawi Publishing Corporation – The Scientific World Journal, vol 2014, article ID 645084, 8 pages.

3. LAU, Y.L.; JAMAIAH, I.; ROHELA, M.; FONG, M.Y.; SITI, C.O.S.; SITI, F.A. "Molecular detection of Entamoeba hystolitica and Entamoeba dispar infection among wild rats in Kuala Lumpur, Malaysia. Tropical Biomedicine 31(4): 721-727 (2014).

4. MOUSTAFA ABDELAAL HEGAZI; TABAREK AHMED PATEL; BASEM SALAMA EL-DEEK. "Prevalence and characters of Entamoeba histolytica infection in Saudi infants and children admitted with diarrhea at 2 main hospitals at south Jeddah: a reemerging serious infection with unusual presentation". The Brazilian Journal of Infectious Diseases, 2013; 1 7(1):32-40.

5. DEIVIANE APARECIDA CALEGAR; BEATRIZ CORONATO NUNES; KERLA JOELINE LIMA MONTEIRO; JESSICA PEREIRA DOS SANTOS; HELENA KEIKO TOMA; TAIS FERREIRA COMES; MARLI MARIA LIMA; MARCIO NEVES BOIA; FILIPE ANIBAL CARVALHO-COSTA. "Frequency and molecular characterization of Entamoeba histolytica, Entamoeba dispar, Entamoeba moshkovskii, and Entamoeba hartmanni in the context of water scarcity in northeastern Brazil". Mem Inst. Oswaldo Cruz, Rio de Janeiro, Vol. 111(2): 114-119, February 2016.

6. PIOTR NOWAK; KATARZYNA MASTALSKA; JAKUB LOSTER. "Entamoeba Histolytica – Pathogenic Protozoan of the Large Intestine in Humans". Journal of Clinical Microbiology and Biochemical Tecnology 1(1): 010-017.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

IVD	in vitro diagnostic device	Ť	Keep dry
Ĩ	Consult instructions for use	X	Temperature limitation
Σ	Use by	444	Manufacturer
LOT	Batch code	Σ _n	Contains sufficient for <n> test</n>
DIL	Sample diluent	REF	Catalogue number

