# **VITASSAY**

## Strep. pneumoniae

Rapid test for the qualitative detection of *Streptococcus* pneumoniae in human urine samples.

IUE-7355047 Ed02 February 2018









#### For professional in vitro diagnostic use only.

## **INTENDED USE**

**Vitassay Strep. pneumoniae** is a rapid one step immunochromatographic assay for the qualitative detection of *Streptococcus pneumoniae* in human urine samples.

Simple, non-invasive and highly sensitive screening assay to make a presumptive diagnosis of *pneumoniae* caused by *Streptococcus pneumoniae* in infected humans.

#### INTRODUCTION

Streptococcus pneumoniae is a potent human pathogen. Infections leads to common diseases such as otitis media, meningitis and pneumonia, which affect several million people and is responsible for significant infant death in developing countries. Worldwide each year, there are over 14 million serious Streptococcus pneumoniae infections in children <5 years of age leading to over 800.000 deaths.

There are at least 95 capsular serotypes, but only a few causes the majority of disease. The organism produces a range of colonization and virulence factors, including the polysaccharide capsule, surface proteins and the toxin pneumolysin (PLY).

## **PRINCIPLE**

**Vitassay Strep. pneumoniae** is a qualitative immunochromatographic assay for the detection of *Streptococcus pneumoniae* in human urine samples.

The test line zone of the nitrocellulose membrane is pre-coated with monoclonal antibodies against *Streptococcus pneumoniae*.

During the process, the sample reacts with the antibodies against *Streptococcus pneumoniae*, forming conjugates. The mixture moves upward on the membrane by capillary action. If the sample is positive, antibodies present on the membrane (test line) capture the conjugate complex and a red line will be visible. Although the sample is positive or negative, the mixture continues to move across the membranes and the blue control line always appears.

The presence of this blue line (in the control zone (C)) indicates that sufficient volume is added; proper flow is obtained and serves as an internal control for the reagents.

#### **PRECAUTIONS**

- For professional in vitro use only.
- · Do not use after expiration date.
- Do not use the test if its pouch is damaged.
- Specimens should be considered as potentially hazardous and handle in the same manner as an infectious agent. A new test

- must be used for each sample to avoid contaminations errors. Single use device.
- Tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with the skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use with the Vitassay Strep. pneumoniae. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

#### STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F).

The test is stable until the expiration date printed on the sealed pouch.

The test must remain in the sealed pouch until use.

Do not freeze.

## **MATERIALS**

MATERIAL PROVIDED	MATERIAL REQUIRED BUT NOT PROVIDED
25 Tests/kits Vitassay Strep. pneumoniae.     Instructions for use.     25 Plastic pipettes.     Positive Control: Inactivated S. pneumoniae swab + pipette + testing tube.     Negative Control: S. pneumoniae negative swab + pipette + testing tube.     Reagent (sample and controls diluent)     25 Testing tubes	<ul> <li>Specimen collection container.</li> <li>Disposable gloves.</li> <li>Timer.</li> </ul>

## SPECIMEN COLLECTION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C/59-86°F) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C (35.6-46.4°F) for up to 14 days or at 10°C to -20°C (14°F to -4°F) for longer periods before testing.

When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C (35.6-46.4°F) or frozen.



## SPECIMEN PREPARATION

Allow all specimens to equilibrate to room temperature before testing.


#### **PROCEDURE**

Allow tests, urine samples and controls to reach room temperature (15-30°C/59-86°F) prior to testing.

Do not open pouches until the performance of the assay.

#### **Patient samples:**

- Use a separate testing tube or vial for each sample. Add 3 drops of urine sample. (figure 1).
- Add 1 drop of Reagent into the testing tube or vial and mix (figure 2). Homogenize the sample.
- Remove the Vitassay Strep. pneumoniae from its sealed bag just before using it.
- Use a separate pippette and device for each sample or control.
   Dispense 3 drops from the testing tube in the circular window marked with the letter S (figure 3).
- Read the results at **15 minutes**. Do not read the test results later than 15 minutes



Add 1 drop of Reagent and mix.



Dispense 3 drops in the circular window marked with the letter S.



Put the swab into the testing tube with Reagent.



## Positive and negative Swabs controls:

- Hold Reagent vertically. Add slowly 13 free falling drops of the Reagent into the testing tube (figure 1b)
- Remove the Positive Control Swab from the pouch and put the swab into the testing tube with reagent (figure 2b), mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (figure 3b)
- Remove the Vitassay Strep. pneumoniae from its sealed bag just before using.
- Use a separate pipette and device for each sample or control.
   Dispense 3 drops from the testing tube, into the circular window marked as the letter S (figure 4b).
- Read the results at **15 minutes**. Do not read the test results later than 15 minutes.

Repeat the procedure for Negative Control swab using the Reagent, same used for test and for positive control swab.

Positive and negative controls should be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.



Mix and extract as much liquid possible from the swab.











Dispense 3 drops in the circular window marked with the letter S.





testing tube.

Parque Tecnológico WALQA Ctra. N.330, Km.566 22197-Cuarte (Huesca, SPAIN)

## INTERPRETATION OF THE RESULTS

C T	NEGATIVE	There is no Streptococcus	
	Only one blue line in the control zone (C)	pneumoniae presence. No infection caused by Streptococcus pneumoniae.  Negative control result.	
	POSITIVE		
СТ	In addition to the blue line (control line C), another red line appears, test line(T)	There is presence of Streptococcus pneumoniae. Infection caused by Streptococcus pneumoniae. Positive control result.	
ANY OTHER RESULTS		Invalid result, we recommend repeating the assay using the sample with another test. <b>Note:</b> Wrong procedural techniques, deterioration of the reagents or insufficient specimen volume are mostly the main reasons for control line failure. If the symptoms or situation persist, discontinue using the test kit and contact your local distributor.	

**Notes:** The intensity of the red colored test line in the result line zone (T) will vary depending on the concentration of antigens in the specimen.

## QUALITY CONTROL

Internal procedural control is included in **Vitassay Strep. pneumoniae**. Blue line appearing in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

External Positive and Negative Controls are included in the kit. The use of positive and negative controls is recommended to assure functionality of reagents and proper performance of assay procedure.

## LIMITATIONS

- Vitassay Strep. pneumoniae must be carried out within 2 hours of opening the sealed bag.
- The use of other samples different from human urine samples has not been established.
- Positive results determine the presence of Streptococcus pneumoniae in urine samples; nevertheless, a positive result should be followed up with additional laboratory techniques to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.

- Negative results should not be considered as conclusive; it is
  possible that the concentration of antigens is lower than the
  detection limit value. If symptoms or situation still persist, a
  Streptococcus pneumoniae determination should be carried out
  from a culture or other methods.
- Excretion of *Streptococcus pneumoniae* antigen in urine may vary depending on the individual patient. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A positive result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.
- The test is compatible with acid boric used as a preservative up to a percentage of 2%.

#### **EXPECTED VALUES**

Streptococcus pneumoniae is a common cause of invasive disease and respiratory tract infections in more and less developed countries. Risk groups for diseases caused by pneumococci such as meningitis, sepsis and pneumonia; include young children, elderly people and patients with immunodeficiencies. Each year, 1 million children younger than 5 years old die from pneumonia and invasive diseases. In USA, the annual number of fatal pneumococcal infections is 40000.

## PERFORMANCE CHARACTERISTICS

#### Analytical sensitivity (detection limit)

Detection limit value of **Vitassay Strep. pneumoniae** is: 0.25 ng/mL of CWPS.

#### Clinical sensitivity and specificity

An evaluation was performed, with urine samples, comparing **Vitassay Strep. pneumoniae** and another commercial test (BinaxNOW®*Streptococcus pneumoniae* Antigen Card, Alere).

Results were as follows:

		BinaxNOW®Streptococcus pneumoniae Antigen Card (Alere)		
		Positive	Negative	Total
Vitassay Strep. pneumoniae	Positive	33	1	34
	Negative	4	73	77
	Total	37	74	111

Vitassay Strep. pneumoniae vs BinaxNOW® <i>Streptococcus pneumoniae</i> Antigen Card				
	Mean Value	95% confidence interval		
Sensitivity	89.2%	74.6-97.0%		
Specificity	98.6%	92.7-100.0%		
PPV	97.1%	84.7-99.9%		
NPV	94.8%	87.2-98.6%		

The results showed that **Vitassay Strep. pneumoniae** has a high sensitivity and specificity to detect *Streptococcus pneumoniae*.

## **Cross reactivity**

No cross reactivity was detected against other pathogens that are occasionally present in urine:

Legionella pneumophila

#### Reproducibility Study

Evaluation studies were performed to determine reproducibility of the **Vitassay Strep. pneumoniae**, including inter-day, interlaboratory, inter and intra lot, showing high reproducibility in all cases.

#### REFERENCES

- 1. WOLFMEIER, H.; RADECKE, J,; SCHOENAUER, R.; ET. AL. "Active release of pneumolysin prepores and pores by mammalian cells undergoing a Streptococcus pneumoniae attack". BIOCHIMICA ET BIOPHYSICA ACTA-GENERAL SUBJECTS 2016; 1860(11): 2498-2509.
- 2. MITCHELL, A.M.; MITCHELL, T.J. "Streptococcus pneumoniae: Virulence factors and variation". CLINICAL MICROBIOLOGY AND INFECTION 2010; 16(5): 411-418.
- 3. MIERNYK, KAREN M.; BULKOW, LISA R.; CASE, SAMANTHA L.; ET AL. "Population structure of invasive Streptococcus pneumoniae isolates among Alaskan children in the conjugate vaccine era, 2001 to 2013". DIAGNOSTIC MICROBIOLOGY AND INFECTIOUS DISEASE 2016; 86(2): 224-230.
- 4. BOGAERT, D; DE GROOT, R; HERMANS, PWM "Streptococcus pneumoniae colonization: the key to pnemococcal disease". Lancet Infectious Diseases 2004; 4: 144-154.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS

IVD	i <i>n vitro</i> diagnostic device	<del>*</del>	Keep dry
Ţ <u>i</u>	Consult instructions for use	1	Temperature limitation
$\subseteq$	Use by	***	Manufacturer
LOT	Batch code	$\sum_{n}$	Contains sufficient for <n> test</n>
DIL	Sample diluent	REF	Catalogue number
Control +	Positive Control	Control -	Negative Control

