# Vitassay qPCR

## Dengue

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- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- StepOne™ Real-Time PCR System (Applied Biosystems)
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTLite Real Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300/7500/7900HT Real-Time PCR
7500 Fast Dx Real-Time PCR System	ABI PRISM 7000
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 6 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 6 Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 3 Real-Time PCR
StepOne Plus™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
StepOne™ Real-Time PCR System	ViiA™ 7 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	<b>Bio-Rad</b>
<b>Bio-Rad</b>	CFX96 Touch™ Deep Well Real-Time PCR Detection
CFX96 Touch™ Real-Time PCR Detection System	iCycler iQ™ Real-Time PCR
Mini Opticon™ Real-Time PCR Detection System	iCycler iQ™5 Real-Time PCR
<b>Roche</b>	MyiQ™ Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™2 Real-Time PCR
LightCycler®96 Real-Time PCR System	<b>Eppendorf</b>
<b>Agilent Technologies</b>	Mastercycler™ep <i>realplex</i>
AriaMx Real-Time PCR System	<b>Stratagene / Agilent Technologies</b>
<b>Qiagen</b>	Mx3000P™ Real Time PCR System
Rotor-Gene®	Mx3005P™ Real Time PCR System
<b>Cepheid</b>	<b>Analytik Jena Biometra</b>
SmartCycler®	TOptical
	qTOWER 2.0
	<b>Abbott</b>
	Abbott m2000 RealTime System
	<b>BIONEER</b>
	Exicycler™ 96
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System*
	DTprime Real-time*
	<b>Qiagen</b>
	Rotor-Gene®
	<b>Cepheid</b>
	SmartCycler®

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## Intended use

Vitassay qPCR Dengue allows the qualitative detection of Dengue virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Dengue virus infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Dengue 4 x 8-well strip, low profile	7041002
Vitassay qPCR Dengue 4 x 8-well strip, high profile	7042002

## Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S002/ 7042S002	Dengue Virus strips high/low profile	-	4x8-well strip
7C002	Dengue Virus Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

## Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves
-

## Summary

Dengue is a mosquito-borne viral disease that has rapidly spread in recent years. Dengue virus is transmitted by female mosquitoes mainly of the species *Aedes aegypti* and, to a lesser extent, *Ae. albopictus*. One study, of the prevalence of dengue, estimates that 3.9 billion people, in 128 countries, are at risk of infection with dengue viruses.

Many people may experience no symptoms during a mild case of dengue fever. When symptoms do occur, they usually begin four to 10 days after you are bitten by an infected mosquito. Signs and symptoms of dengue fever most commonly include: high fever, headaches, muscle, bone and joint pain, pain behind your eyes. You might also experience widespread rash, nausea and vomiting and rarely, minor bleeding from your gums or nose. In some cases, symptoms worsen and can become life-threatening. Blood vessels often become damaged and leaky. And the number of clot-forming cells (platelets) in your bloodstream drops.

Dengue can be diagnosed by isolation of the virus, by serological tests (antigen and antibodies detection), or by molecular methods. Diagnosis of acute or recent dengue infection can be established by testing serum samples during the first 5 days of symptoms and/or early convalescent phase. Methods such as one-step real time RT-PCR are now widely used to detect dengue viral genes in acute-phase serum samples. This detection coincides with the viremia and the febrile phase of illness onset. It's important to consider that serological methods have cross reactivity with other flaviviruses including West Nile virus (WNV), St. Louis encephalitis virus (SLE), Japanese encephalitis virus (JEV) and yellow fever virus (YFV).

## Principle of the test

Vitassay qPCR Dengue is designed for the identification of Dengue virus in clinical specimens to aid in the assessment of infections caused by this virus.

Vitassay qPCR Dengue is based on the real-time amplification of the 3' non coding region encoded by the Dengue virus genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of Dengue virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Dengue is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an

internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

## **Precautions**

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and RNA extraction**

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (Qiagen).

Invisorb® Spin Universal Kit (Strattec).

## Positive control preparation

Reconstitute the lyophilized Dengue Virus Positive Control (red tube) in the 100 µL of PCR grade water (white tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

## Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative controls must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of resuspension buffer (green tube) into each well.
- Pipette 5 µL of RNA sample, negative and positive controls into each well
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Dengue virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the

Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer’s instructions.

For a valid diagnostic test run, the following control conditions must be met:

#### Positive control

The positive controls used in each run, must show an amplification curve for Dengue virus, which validates the reaction.

#### Negative control

The negative controls included in each run, must show the absence of signal for Dengue virus, which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Dengue Virus FAM	Internal Control HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	Virus Dengue Positive
-	+	-	+	Virus Dengue Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

+ **Positive:** Amplification signal

- **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

Overall, 34 serum clinical samples from symptomatic patients were tested by Real-time PCR using: Vitassay qPCR Dengue and RealStar® Dengue Virus RT-PCR Kit (Altona Diagnostics). Dengue virus was detected in 10 samples by Vitassay qPCR Dengue. The results were concordant between all tests. Vitassay qPCR Dengue shows a high sensitivity and specificity to detect Dengue virus.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Dengue virus template ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

### Analytical specificity

The analytical specificity for Dengue virus was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Dengue Virus		
Chikungunya virus strain S27 Petersfield	West Nile virus strain H160/99	West Nile virus Ug37
Zika virus strain MR 766	West Nile virus Heja	Yellow Fever virus strain 17D
St Louis Encephalitis virus strain 17D		

## Analytical reactivity

The reactivity of Vitassay qPCR Dengue was confirmed by the real-time amplification using Dengue 1 virus strain Hawaii, Dengue 2 virus strain New Guinea C, Dengue 3 virus strain H87 and Dengue 4 virus strain H241 as template.

## Compatibles real-time PCR equipment

Vitassay qPCR Dengue has been validated on the following equipments:

- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- StepOne™ Real-Time PCR System (Applied Biosystems)
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTLite Real Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of Dengue virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with serum samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.

- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Dengue virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers
<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 3 Real-Time PCR System
QuantStudio™ 5 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
<b>Bio-Rad</b>
CFX96 Touch™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
<b>Roche</b>
LightCycler® 480 Real-Time PCR System
LightCycler® 96 Real-Time PCR System
<b>Agilent Technologies</b>
AriaMx Real-Time PCR System
<b>Qiagen</b>
Rotor-Gene®
<b>Cepheid</b>
SmartCycler®

High profile Block Thermocyclers
<b>Applied Biosystems</b>
7300/7500/7900HT Real-Time PCR
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Real-Time PCR
QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Real-Time PCR
<b>Bio-Rad</b>
CFX96 Touch™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
MyiQ™ Real-Time PCR
MyiQ™2 Real-Time PCR
<b>Eppendorf</b>
Mastercycler™ep <i>realplex</i>
<b>Stratagene / Agilent Technologies</b>
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System
<b>Analytik Jena Biometra</b>
TOptical
qTOWER 2.0
<b>Abbott</b>
Abbott m2000 RealTime System
<b>BIONEER</b>
Exicycler™ 96
<b>DNA-Technology</b>
DTlite Real-Time PCR System*
DTprime Real-time*
<b>Qiagen</b>
Rotor-Gene®
<b>Cepheid</b>
SmartCycler®

\* See Attached III to configure exposure settings.

**Attached III: Optical measurement exposure setting**

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	150
	HEX	500
	ROX	500
	Cy5	100
<b>DTprime Real-time Detection Thermal Cycler (DNA- Technology)</b>	FAM	150
	HEX	3000
	ROX	2000
	Cy5	1500

## **Bibliography/Bibliografía**

1. Genetic and serologic properties of Dengue virus associated with an epidemic, Yap State, Micronesia, 2007. Lanciotti RS, Kosoy OL, Laven JJ, Velez JO, Lambert AJ, Johnson AJ, Stanfield SM, Duffy MR. Emerg Infect Dis. 2008 Aug;14(8):1232-9.
2. Quantitative real-time PCR detection of Dengue virus and evaluation with field-caught mosquitoes. Faye O, Faye O, Diallo D, Diallo M, Weidmann M, Sall AA. Virol J. 2013 Oct 22;10:311.
3. Detection of Dengue virus in urine. Gourinat AC, O'Connor O, Calvez E, Goarant C, Dupont-Rouzeyrol M. Emerg Infect Dis. 2015 Jan;21(1):84-6.
4. Centers for Disease Control and Prevention. Dengue Virus (<http://www.cdc.gov/Dengue/index.html>).
5. World Health Organization. Dengue virus and potential complications (<http://www.who.int/emergencies/Dengue-virus/en/>).

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.




Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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