
Vitassay qPCR

Zika

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- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- StepOne™ Real-Time PCR System (Applied Biosystems)
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300/7500/7900HT Real-Time PCR
7500 Fast Dx Real-Time PCR System	ABI PRISM 7000
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 6 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 6 Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 3 Real-Time PCR
StepOne Plus™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
StepOne™ Real-Time PCR System	ViiA™ 7 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	Bio-Rad
Bio-Rad	CFX96 Touch™ Deep Well Real-Time PCR Detection
CFX96 Touch™ Real-Time PCR Detection System	iCycler iQ™ Real-Time PCR
Mini Opticon™ Real-Time PCR Detection System	iCycler iQ™5 Real-Time PCR
Roche	MyiQ™ Real-Time PCR
LightCycler @480 Real-Time PCR System	MyiQ™2 Real-Time PCR
LightCycler @96 Real-Time PCR System	Eppendorf
Agilent Technologies	Mastercycler™ep <i>realplex</i>
AriaMx Real-Time PCR System	Stratagene / Agilent Technologies
DNA-Technology	Mx3000P™ Real Time PCR System
DTlite Real-Time PCR System	Mx3005P™ Real Time PCR System
DTprime Real-time Detection Thermal Cycler	Analytik Jena Biometra
Qiagen	TOptical
Rotor-Gene®	qTOWER 2.0
Cepheid	Abbott
SmartCycler®	Abbott m2000 RealTime System
	BIONEER
	Exicycler™ 96
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time
	Qiagen
	Rotor-Gene®
	Cepheid
	SmartCycler®

Intended use

Vitassay qPCR Zika allows the qualitative detection of Zika virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Zika virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Zika 4 x 8-well strip, low profile	7041003
Vitassay qPCR Zika 4 x 8-well strip, high profile	7042003

Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S003/ 7042S003	Zika Virus strips	-	4x8-well strip
7C003	Zika Virus Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Rehydration Buffer	green	1 vial x 1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Zika virus (ZIKV) is an emerging mosquito-borne pathogen (family *Flaviviridae*, genus *Flavivirus*) that was isolated in 1947 from a rhesus monkey in the Zika forest in Uganda. ZIKV is believed to be transmitted to humans by infected *Aedes spp.* mosquitoes. Several studies have demonstrated that ZIKV is endemic to Africa and Southeast Asia. In fact, in other regions, Zika virus infections have been reported in few cases until 2007 when an epidemic of ZIKV infection in humans occurred in Yap Island and Federated States of Micronesia, in the Pacific region.

In humans, ZIKV infection is characterized by mild fever (37.8°C–38.5°C); arthralgia (notably of small joints of hands and feet); myalgia, headache; retroorbital pain; conjunctivitis; and cutaneous maculopapular rash. While many cases of ZIKV infection are asymptomatic or have a mild medical history, not specific symptoms, and/or similar to that caused by the influenza virus. Reason why, ZIKV infection can be misdiagnosed during the acute (viremic) phase.

Confirmation of ZIKV infections, despite its low viremia, is based mostly on detection of virus RNA in serum by using reverse transcription PCR (RT-PCR). Besides the serological detection of IgM antibodies against ZIKV can be detected by ELISA, although few laboratories have this ability and Zika antibodies can present cross reactivity against other *flavivirus* (including Dengue). This can difficult the serologic diagnosis since a confirmation by neutralization assays is necessary. Consequently, the laboratory diagnosis is a challenge to be faced with the best possible approximation.

Principle of the test

Vitassay qPCR Zika is designed for the identification of Zika virus in clinical specimens to aid in the assessment of infections caused by this virus.

Vitassay qPCR Zika is based on the real-time amplification of specific conserved fragments of the *envelope* gene encoded by the Zika virus genome. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of Zika virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Zika is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of

the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pretreatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (Qiagen).

EZ1 Virus Mini Kit, EZ1 instrument (Qiagen).

High Pure Viral Nucleic Acid Kit (Roche).

Invisorb® Spin Universal Kit (Strattec).

Positive control preparation

Reconstitute the lyophilized Zika Virus Positive Control (red tube) in the 100 µL of PCR grade water (white tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative controls must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of resuspension buffer (green tube) into each well.
- Pipette 5 µL of RNA sample, negative and positive controls into each well
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Zika virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied

Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer’s instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for Zika virus, which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for Zika virus, which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Virus Zika FAM	Control Interno HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	Virus Zika Positive
-	+	-	+	Virus Zika Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

+ **Positive:** Amplification signal

- **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Overall, 61 clinical samples from symptomatic patients were tested by Real-time PCR using: Vitassay Zika qPCR, RealStar® Zika Virus RT-PCR Kit (Altona Diagnostics) and/or in-house RT-PCR assay (reported by Faye et al., 2013 and/or Lanciotti et al., 2008). Zika virus was detected in 6 samples by Vitassay Zika QPCR RT-PCR. The results were concordant between all test. Vitassay Zika qPCR shows a high sensitivity and specificity to detect Zika virus.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Zika virus template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Zika virus was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Zika Virus		
Chikungunya virus strain S27 Petersfield	Dengue 4 virus strain H241	West Nile virus Heja
Dengue 1 virus strain Hawaii	St Louis Encephalitis virus strain 17D	West Nile virus Ug37
Dengue 2 virus strain New Guinea C	West Nile virus strain H160/99	Yellow Fever virus strain 17D
Dengue 3 virus strain H87		

Analytical reactivity

The reactivity of Vitassay qPCR Zika was confirmed by the real-time amplification using Zika virus strain MR 766 as template.

Compatibles real-time PCR equipment

Vitassay qPCR Zika has been validated on the following equipments:

- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- StepOne™ Real-Time PCR System (Applied Biosystems)
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Zika virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with serum and urine samples. The use of other samples has not been established.

- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Zika virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 3 Real-Time PCR System
QuantStudio™ 5 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
Bio-Rad
CFX96 Touch™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Roche
LightCycler @480 Real-Time PCR System
LightCycler @96 Real-Time PCR System
Agilent Technologies
AriaMx Real-Time PCR System
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Qiagen
Rotor-Gene®
Cepheid
SmartCycler®

High profile Block Thermocyclers
Applied Biosystems
7300/7500/7900HT Real-Time PCR
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Real-Time PCR
QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Real-Time PCR
Bio-Rad
CFX96 Touch™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
MyiQ™ Real-Time PCR
MyiQ™2 Real-Time PCR
Eppendorf
Mastercycler™ ^{ep} <i>realplex</i>
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System
Analytik Jena Biometra
TOptical
qTOWER 2.0
Abbott
Abbott m2000 RealTime System
BIONEER
Exicycler™ 96
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time
Qiagen
Rotor-Gene®
Cepheid
SmartCycler®

Bibliography/Bibliografía

1. Genetic and serologic properties of Zika virus associated with an epidemic, Yap State, Micronesia, 2007. Lanciotti RS, Kosoy OL, Laven JJ, Velez JO, Lambert AJ, Johnson AJ, Stanfield SM, Duffy MR. *Emerg Infect Dis.* 2008 Aug;14(8):1232-9.
2. Quantitative real-time PCR detection of Zika virus and evaluation with field-caught mosquitoes. Faye O, Faye O, Diallo D, Diallo M, Weidmann M, Sall AA. *Viol J.* 2013 Oct 22;10:311.
3. Detection of Zika virus in urine. Gourinat AC, O'Connor O, Calvez E, Goarant C, Dupont-Rouzeyrol M. *Emerg Infect Dis.* 2015 Jan;21(1):84-6.
4. Centers for Disease Control and Prevention. Zika Virus (<http://www.cdc.gov/zika/index.html>).
5. World Health Organization. Zika virus and potential complications (<http://www.who.int/emergencies/zika-virus/en/>).

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.

Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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