Vitassay qPCR

Zika + Dengue + Chikungunya

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Intended use

Vitassay qPCR Zika+Dengue+Chikungunya allows the detection and differentiation of Zika virus, Dengue virus and Chikungunya virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Zika virus, Dengue virus and/or Chikungunya virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Zika+Dengue+Chikungunya 4x8 -well strip, low profile 7041005

Vitassay qPCR Zika+Dengue+Chikungunya 4x8-well strip, high profile 7042005

Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S005/ 7042S005	Zika+Dengue+Chikungunya Virus strips low/high profile	-	4x8-well strip
7C005	Zika+Dengue+Chikungunya Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1 mL
7003N	Negative control	yelow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Zika virus (ZIKV), Dengue (DENV) and Chikungunya virus (CHIKV) are transmitted by the sting of the same *Aedes* genus mosquito.

Besides, Zika virus was isolated in 1947 from a rhesus monkey in the Zika forest in Uganda. Typically, clinical manifestations of human ZIKV disease are similar to many others derived of arbovirus infections that occur without serious complications and may include a self-limiting febrile illness, arthralgia, myalgia, headache and maculopapular rash. Conjunctivitis, retroorbital eye pain, lymphadenopathy and diarrhea have also been reported.

Dengue is also an acute viral illness caused by RNA virus of the family Flaviviridae. Presenting features may range from asymptomatic fever to dreaded complications such as hemorrhagic fever and shock. Acute onset high fever, muscle and joint pain, myalgia, cutaneous rash, hemorrhagic episodes, and circulatory shock are the commonly seen symptoms.

Chikungunya virus was first isolated in Tanzania in 1953. Chikungunya virus causes severe, debilitating, often chronic arthralgia with high attack rates, resulting in severe morbidity and economic costs to affected communities.

These three arboviruses in the initial stages of the infection, cause similar clinical presentations, although the disease management is different, for that, an early detection and differentiation is crucial. Confirmation and identification of these infections is based mostly on detection of viral RNA in serum and urine by using reverse transcription PCR (RT-PCR). In fact, serology is not advised since this technique can cause cross reactivity between the antibodies of these arbovirus.

Principle of the test

Vitassay qPCR Zika+Dengue+Chikungunya test is based on the real-time amplification of specific conserved fragments of the *envelope* gene (Zika virus), of the 3' Non coding region (Dengue virus) and of the *NSP1* gene (Chikungunya virus). The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Zika+Dengue+Chikungunya test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In

addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Zika virus RNA target sequence is detected through the Cy5 channel, Dengue virus RNA target in FAM channel, Chikungunya virus RNA target in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (QIAGEN).

High Pure Viral Nucleic Acid Kit (Roche).

Invisorb® Spin Universal Kit (Strattec).

Positive control preparation

Reconstitute the lyophilized Zika+Dengue+Chikungunya Positive Control (red tube) in the 100 µL of PCR grade water (transparent tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative controls must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) into each well.
- Pipette 5 µL of RNA sample, negative and positive controls into each well
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the Cy5 (Zika virus), FAM (Dengue virus), ROX (Chikungunya) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached 2)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for Cy5 (Zika virus), FAM (Dengue virus) and ROX (Chikungunya), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Zika Virus	Dengue Virus	Chikungunya Virus	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Zika , Dengue and Chikungunya viruses Positive
-	-	-	+	-	+	Zika, Dengue and Chikungunya viruses Negative
+	-	-	+/-	-	+	Zika Virus Positive, Dengue and Chikungunya Viruses Negative
+	+	-	+/-	-	+	Zika and Dengue Viruses Positive, and Chikungunya Virus Negative
+	-	+	+/-	-	+	Zika and Chikungunya Viruses Positive, and Dengue Virus Negative
-	+	-	+/-	-	+	Dengue Virus Positive, Zika and Chikungunya Viruses Negative
-	+	+	+/-	-	+	Dengue and Chikungunya Viruses Positive, Zika Virus Negative
-	-	+	+/-	-	+	Chikungunya Virus Positive, Zika and Dengue Viruses Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 17 serum samples from Dengue symptomatic patients were tested by Real Time PCR using Vitassay qPCR Zika+Dengue+Chikungunya and RealStar® Dengue RT-PCR Kit (Altona Diagnostics). Dengue virus was detected in 5 samples by Vitassay qPCR Zika+Dengue+Chikungunya Real Time PCR.

A total of 4 clinical samples from INSTAND 2016 panel from Virus Genome Detection-Zikaviruses EQA Programme, were tested by Real Time PCR using Vitassay qPCR *Zika+Dengue+Chikungunya*. Zika virus was correctly detected in the 3 positive samples.

A total of 12 clinical samples from Dengue QCMD 2014 panel and 10 samples from Chikungunya QCMD 2015 panel were tested by Vitassay qPCR *Zika+Dengue+Chikungunya*. The results of Vitassay qPCR kit were concordant for all the samples.

The results show a high sensitivity and specificity to detect Zika, Dengue and Chikungunya viruses using Vitassay qPCR *Zika+Dengue+Chikungunya*.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Zika, Dengue and Chikungunya virus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Zika, Dengue and Chikungunya virus was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Zika	Dengue	Chikungunya
Chikungunya S27 Petersfield	Chikungunya S27 Petersfield	Zika MR 766
Dengue 1 Hawaii	Zika MR 766	Dengue 1 Hawaii
Dengue 2 New Guinea C	St Louis Encephalitis 17D	Dengue 2 New Guinea C
Dengue 3 H87	West Nile H160/99	Dengue 3 H87
Dengue 4 H241	West Nile Heja	Dengue 4 H241
St Louis Encephalitis 17D	West Nile Ug37	St Louis Encephalitis 17D
West Nile H160/99	Yellow Fever 17D	West Nile H160/99
West Nile Heja		West Nile Heja
West Nile Ug37		West Nile Ug37
Yellow Fever 17D		Yellow Fever 17D

Analytical reactivity

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR for Zika virus was confirmed by the real time amplification using Zika virus strain MR 766 as template.

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR for Dengue virus was confirmed by the real time amplification using Dengue 1 virus strain Hawaii, Dengue 2 virus strain New Guinea C, Dengue 3 virus strain H87 and Dengue 4 virus strain H241 as template.

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR for Chikungunya virus was confirmed by the real time amplification using Chikungunya virus S27 Petersfield as template.

Compatibles real-time PCR equipment

Vitassay qPCR Zika+Dengue+Chikungunya has been validated on the following equipments:

- LightCycler 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- CFX96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Zika virus, Dengue virus and Chikungunya virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with serum and urine samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Zika, Dengue and/or Chikungunya virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7500 Real-Time PCR
7500 Fast Dx Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 12K Flex 96-well Fast	QuantStudio™ 6 Flex 96-well
QuantStudio™ 6 Flex 96-well Fast	QuantStudio™ 7 Flex 96-well
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 5 Real-Time PCR
QuantStudio™ 5 Real-Time PCR System	ViiA™ 7 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	Bio-Rad
Bio-Rad	CFX96 Touch™ Deep Well Real-Time PCR Detection
CFX96 Touch™ Real-Time PCR Detection System	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™5 Real-Time PCR
LightCycler®480 Real-Time PCR System	Eppendorf
LightCycler®96 Real-Time PCR System	Mastercycler™ep <i>realplex</i>
Agilent Technologies	Stratagene / Agilent Technologies
AriaMx Real-Time PCR System	Mx3000P™ Real Time PCR System
DNA-Technology	Mx3005P™ Real Time PCR System
DTlite Real-Time PCR System	Analytik Jena Biometra
DTprime Real-time Detection Thermal Cycler	TOptical
Qiagen	qTOWER 2.0
Rotor-Gene®	Abbott
Cepheid	Abbott m2000 RealTime System
SmartCycler®	BIONEER
	Exicycler™ 96
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time
	Qiagen
	Rotor-Gene®
	Cepheid
	SmartCycler®

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Bibliography/Bibliografía

1. Genetic and serologic properties of Zika virus associated with an epidemic, Yap State, Micronesia, 2007 Lanciotti RS, Kosoy OL, Laven JJ, Velez JO, Lambert AJ, Johnson AJ, Stanfield SM, Duffy MR. *Emerg Infect Dis.* 2008 Aug;14(8):1232-9.
2. Quantitative real-time PCR detection of Zika virus and evaluation with field-caught mosquitoes. Faye O, Faye O, Diallo D, Diallo M, Weidmann M, Sall AA. *Virol J.* 2013 Oct 22;10:311.
3. Dengue virus: A global human threat: Review of literature. Hasan S, Jamdar SF, Alalawi M, Al Ageel Al Beajji SM. *J Int Soc Prev Community Dent.* 2016 Jan-Feb 6(1):1-6.
4. One-step real-time RT-PCR assays for serotyping dengue virus in clinical samples. Alm E, Lindegren G, Falk KI, Lagerqvist N. *BMC Infect Dis.* 2015 Nov2;15:493.
5. Pan-dengue virus detection by PCR for travelers returning from the tropics. Dumoulin A, Marti H, Panning M, Hatz C, Hirsch HH. *J Clin Microbiol* 2008 Sep;46(9):3104-6.
6. External quality assessment studies for laboratory performance of molecular and serological diagnosis of Chikungunya virus infection. Jacobsen S, Patel P, Schmidt-Chanasit J, Leparc-Goffart I, Teichmann A, Zeller H, Niedrig M. *J Clin Virol* 2016 Mar;76:55-65.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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