# Vitassay qPCR

## Campylobacter + Salmonella + Shigella

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## Intended use

Vitassay qPCR Campylobacter + Salmonella + Shigella allows the detection and differentiation of *Campylobacter*, *Salmonella* y *Shigella*/ enteroinvasive *Escherichia coli* (EIEC) by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Campylobacter*, *Salmonella* y/o *Shigella*/EIEC infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Campylobacter + Salmonella + Shigella 4x8 -well strip, low profile  
7041010

Vitassay qPCR Campylobacter + Salmonella + Shigella 4x8-well strip, high profile  
7042010

## Materials/reagents provided

Code	Reagent/Material	Colour	Quantity
7041S010/ 7042S010	Campylobacter + Salmonella + Shigella strips low/high profile	-	4x8-well strip
7C010	Campylobacter + Salmonella + Shigella Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

## Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)

- Filter tips
- Powder-free disposal gloves

## Summary

Gastroenteritis is a major cause of morbidity and mortality, worldwide, both in children 5 years old and in the general population.

*Campylobacter* are gram negative, slender, spirally curved, microaerophilic bacteria that live as commensal organisms in the gastrointestinal tract of many domestic and wild birds and mammals. *C. jejuni* and *C. coli* are by far the most important human pathogens in the genus and account for more than 95% of all clinical isolates worldwide.

The genus *Salmonella* contains two distinct species, designated *Salmonella enterica* (divided into six subspecies) and *Salmonella bongori*. *Salmonella* species are a leading bacterial cause of acute gastroenteritis.

Shigellosis produces inflammatory reactions and ulceration on the intestinal epithelium followed by bloody or mucoid diarrhea. It is caused by enteroinvasive *E. coli* (EIEC) as well as any species of the genus *Shigella*, namely, *S. dysenteriae*, *S. flexneri*, *S. boydii*, and *S. sonnei*. Shigellosis is endemic in many developing countries and also occurs in epidemics causing considerable morbidity and mortality.

## Principle of the test

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* test is based on the real-time amplification of a conserved region of the *16S rRNA* gene for the identification of *Campylobacter*, *invA* gene for *Salmonella* and *ipaH* gene for *Shigella*/EIEC. After DNA isolation, the presence of *Campylobacter*, *Salmonella* y/o *Shigella*/ enteroinvasive *Escherichia coli* (EIEC) is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Campylobacter* DNA target sequence is detected through the Cy5 channel, *Salmonella* DNA target in FAM channel and *Shigella*/EIEC DNA in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and DNA extraction**

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® RSC Blood DNA Kit, utilizando Maxwell® 16 instrument (Promega).

QIAamp DNA Mini kit (QIAGEN).

QIAamp DNA Stool Mini Kit (QIAGEN).

RIDA® Xtract (r-Biopharm).

Invisorb® Spin Universal Kit (Strattec).

### **Positive control preparation**

Reconstitute the lyophilized *Campylobacter* + *Salmonella* + *Shigella* Positive Control (red tube) in the 100 µL of PCR grade water (transparent tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) into each well.
- Pipette 5 µL of DNA sample, negative and positive controls into each well.
- Cover the wells with the caps provided. Spin down briefly(optional).
- Place the strips in the Real-time PCR instrument.

### Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the Cy5 (*Campylobacter*), FAM (*Salmonella*), ROX (*Shigella/EIEC*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached 2)

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer’s instructions.

For a valid diagnostic test run, the following control conditions must be met:

**Positive control**

The positive controls used in each run, must show an amplification curve for Cy5 (*Campylobacter*), FAM (*Salmonella*) and ROX (*Shigella*/EIEC), which validates the reaction.

**Negative control**

The negative controls included in each run, must show the absence of signal in Cy5, FAM and ROX which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

<b>Campylobacter (Cy5)</b>	<b>Salmonella (FAM)</b>	<b>Shigella/EIEC (ROX)</b>	<b>Internal Control</b>	<b>Negative Control</b>	<b>Positive Control</b>	<b>Interpretation</b>
+	+	+	+/-	-	+	<i>Campylobacter, Salmonella and Shigella/EIEC Positives</i>
-	-	-	+	-	+	<i>Campylobacter, Salmonella and Shigella/EIEC Negatives</i>
+	-	-	+/-	-	+	<i>Campylobacter Positive, Salmonella and Shigella/EIEC Negatives</i>
+	+	-	+/-	-	+	<i>Campylobacter and Salmonella Positives, Shigella/EIEC Negative</i>
+	-	+	+/-	-	+	<i>Campylobacter and Shigella/EIEC Positives, Salmonella Negative</i>
-	+	-	+/-	-	+	<i>Salmonella Positive, Campylobacter and Shigella/EIEC Negatives</i>
-	+	+	+/-	-	+	<i>Salmonella and Shigella/EIEC Positives, Campylobacter Negative</i>
-	-	+	+/-	-	+	<i>Shigella/EIEC Positive, Campylobacter and Salmonella Negatives</i>
+	+	+	+	+	+	<b>Invalid</b>
-	-	-	-	-	-	<b>Invalid</b>

(+) **Positive:** Amplification signal  
 (-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

A total of 400 faecal specimens from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Campylobacter + Salmonella + Shigella and RIDA®GENE Bacterial Stool Panel (R-biopharm) for *Salmonella* and *Campylobacter*.

Vitassay qPCR Campylobacter + Salmonella + Shigella detects *Campylobacter* in 58 positive samples. This test identified even 13 additional low positives that could be confirmed as positive by an additional commercial Real Time PCR Kit (Mericon-Campylobacter spp, QIAGEN).

Vitassay qPCR Campylobacter + Salmonella + Shigella detects *Salmonella* in 21 positive samples. This test identified even 1 additional low positive that could be confirmed as positive by culture for *Salmonella G*.

Vitassay qPCR Campylobacter + Salmonella + Shigella detects *Shigella*/EIEC in 2 positive sample. This result was confirmed by a commercial Real Time PCR Kit ("RIDA®GENE EHEC/EPEC" (r-Biopharm), which detects and differentiates EHEC, STEC (EHEC), EPEC, *Shigella dysenteriae* type 1 and *Shigella*/EIEC).

The specificity of Vitassay qPCR Campylobacter + Salmonella + Shigella was 100% for *Campylobacter*, *Salmonella* and *Shigella*/EIEC.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Campylobacter*, *Salmonella* and *Shigella*/EIEC templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  DNA copies per reaction.

## Analytical specificity

The analytical specificity for *Campylobacter*, *Salmonella* and *Shigella*/EIEC was tested within the panel of different microorganisms.

No cross-reactivity of *Campylobacter* was seen between any of the species:

<i>Shigella flexneri</i>	<i>Helicobacter heilmannii</i>	<i>Arcobacter butzleri</i>
<i>Shigella dysenteriae</i>	<i>Proteus vulgaris</i>	<i>Pseudomonas aeruginosa</i>
<i>Salmonella typhi</i>	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	<i>Enterococcus faecalis</i>
<i>Salmonella paratyphi A</i>	<i>Citrobacter freundii</i>	<i>Giardia intestinalis</i>
<i>Salmonella paratyphi B</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Yersinia enterocolitica</i> O:3
<i>Salmonella typhimurium</i>	<i>Serratia liquefaciens</i>	<i>Yersinia enterocolitica</i> O:9
<i>Salmonella bongori</i>	<i>Vibrio parahaemolyticus</i>	<i>Bacteroides fragilis</i>
<i>Salmonella enteritidis</i>	<i>Clostridium difficile</i>	<i>Cryptosporidium parvum</i>
<i>Salmonella enterica</i> subsp. <i>entérica</i>	<i>Clostridium perfringens</i>	<i>Entamoeba histolytica</i>
<i>Salmonella pullorum</i>	<i>Enterotoxigenic Escherichia coli</i>	Adenovirus serotipo 40
<i>Salmonella gallinarum</i>	<i>Enteropathogenic Escherichia coli</i>	Adenovirus serotipo 41
<i>Helicobacter pylori</i>	<i>Klebsiella oxytoca</i>	Rotavirus A
<i>Helicobacter hepaticus</i>	<i>Listeria monocytogenes</i>	Norovirus Genotipo I y II
<i>Helicobacter cinaedi</i>	<i>Candida albicans</i>	Astrovirus Genotipo I-VIII

No cross-reactivity of *Salmonella* was seen between any of the species:

<i>Helicobacter pylori</i>	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	<i>Arcobacter butzleri</i>
<i>Helicobacter hepaticus</i>	<i>Citrobacter freundii</i>	<i>Pseudomonas aeruginosa</i>
<i>Helicobacter cinaedi</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Enterococcus faecalis</i>
<i>Helicobacter heilmannii</i>	<i>Serratia liquefaciens</i>	<i>Bacteroides fragilis</i>
<i>Shigella flexneri</i>	<i>Vibrio parahaemolyticus</i>	<i>Yersinia enterocolitica</i> O:3 y O:9
<i>Shigella dysenteriae</i>	<i>Clostridium difficile</i>	<i>Cryptosporidium parvum</i>
<i>Campylobacter lari</i>	<i>Clostridium perfringens</i>	<i>Giardia intestinalis</i>
<i>Campylobacter fetus</i>	<i>Enterotoxigenic Escherichia coli</i>	<i>Entamoeba histolytica</i>
<i>Campylobacter coli</i>	<i>Enteropathogenic Escherichia coli</i>	Adenovirus serotipos 40/41
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	<i>Klebsiella oxytoca</i>	Rotavirus A
<i>Campylobacter upsaliensis</i>	<i>Listeria monocytogenes</i>	Norovirus Genotipos I y II
<i>Proteus vulgaris</i>	<i>Candida albicans</i>	Astrovirus Genotipos I-VIII

No cross-reactivity of *Shigella*/EIEC was seen between any of the species:

<i>Helicobacter pylori</i>	<i>Campylobacter coli</i>	<i>Candida albicans</i>
<i>Helicobacter hepaticus</i>	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	<i>Arcobacter butzleri</i>
<i>Helicobacter cinaedi</i>	<i>Campylobacter upsaliensis</i>	<i>Pseudomonas aeruginosa</i>
<i>Helicobacter heilmannii</i>	<i>Proteus vulgaris</i>	<i>Enterococcus faecalis</i>
<i>Salmonella typhi</i>	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	<i>Yersinia enterocolitica</i> O:3
<i>Salmonella paratyphi A</i>	<i>Citrobacter freundii</i>	<i>Yersinia enterocolitica</i> O:9
<i>Salmonella paratyphi B</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Bacteroides fragilis</i>
<i>Salmonella typhimurium</i>	<i>Serratia liquefaciens</i>	<i>Cryptosporidium parvum</i>
<i>Salmonella bongori</i>	<i>Vibrio parahaemolyticus</i>	<i>Giardia intestinalis</i>
<i>Salmonella enteritidis</i>	<i>Clostridium difficile</i>	<i>Entamoeba histolytica</i>
<i>Salmonella enterica</i> subsp. <i>enterica</i>	<i>Clostridium perfringens</i>	Adenovirus serotipos 40
<i>Salmonella pullorum</i>	<i>Enterotoxigenic Escherichia coli</i>	Adenovirus serotipos 41
<i>Salmonella gallinarum</i>	<i>Enteropathogenic Escherichia coli</i>	Rotavirus A
<i>Campylobacter lari</i>	<i>Klebsiella oxytoca</i>	Norovirus Genotipos I y II
<i>Campylobacter fetus</i>	<i>Listeria monocytogenes</i>	Astrovirus Genotipos I-VIII

### Analytical reactivity

Analytical reactivity of Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* for *Campylobacter* was evaluated against *Campylobacter jejuni* subsp. *jejuni*, *Campylobacter coli*, *Campylobacter lari*, *Campylobacter upsaliensis*, *Campylobacter fetus*, *Campylobacter concisus*, *Campylobacter hyointestinalis*, *Campylobacter gracilis*, *Campylobacter helveticus*, *Campylobacter curvus* and *Campylobacter rectus*, showing positive result.

Analytical reactivity of Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* for *Salmonella* was evaluated against *Salmonella enterica* subsp. *enterica*, *Salmonella paratyphi A*, *Salmonella paratyphi B*, *Salmonella typhimurium*, *Salmonella enteritidis*, *Salmonella typhi*, *Salmonella bongori*, *Salmonella pullorum*, *Salmonella gallinarum*, *Salmonella mbandaka* and *Salmonella braenderup*, showing positive result.

Analytical reactivity of Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* for *Shigella*/EIEC was evaluated against *Shigella dysenteriae*, *Shigella flexneri*, *Shigella boydii* and *Shigella sonnei*, showing positive result.

## Compatibles real-time PCR equipment

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of *Campylobacter*, *Salmonella* and/or *Shigella*/EIEC infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with human faecal samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different bacteria, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers
<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
<b>Bio-Rad</b>
CFX96 Touch™ Real-Time PCR Detection System
<b>Roche</b>
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
<b>Agilent Technologies</b>
AriaMx Real-Time PCR System
<b>DNA-Technology</b>
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
<b>Qiagen</b>
Rotor-Gene®
<b>Cepheid</b>
SmartCycler®

High profile Block Thermocyclers
<b>Applied Biosystems</b>
7500 Real-Time PCR
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Real-Time PCR
<b>Bio-Rad</b>
CFX96 Touch™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
<b>Eppendorf</b>
Mastercycler™ep <i>realplex</i>
<b>Stratagene / Agilent Technologies</b>
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System
<b>Analytik Jena Biometra</b>
TOptical
qTOWER 2.0
<b>Abbott</b>
Abbott m2000 RealTime System
<b>BIONEER</b>
Exicycler™ 96
<b>DNA-Technology</b>
DTlite Real-Time PCR System
DTprime Real-time
<b>Qiagen</b>
Rotor-Gene®
<b>Cepheid</b>
SmartCycler®

**Attached II: Detection channels of most common real time PCR equipment**

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

Real-time pcr thermocycler	Vitassay channel	Detection channel	Observations
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Passive reference option ROX is none
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene®Q Qiagen</b>	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

## Bibliography/Bibliografía

1. Campylobacter virulence and survival factors. Declan J. Bolton. Food Microbiology. 2015; 48: 99-108.
2. Automated 59 Nuclease PCR Assay for Identification of *Salmonella enterica*. J. Hoorfar *et al.* *Journal of Clinical Microbiology* 2000; 38: 3429–3435.
3. Detection of *Shigella* by a PCR assay targeting the *ipaH* gene suggests increased prevalence of shigellosis in Nha Trang, Vietnam. V. D. Thiem *et al.* *Journal of Clinical Microbiology* 2004; 42(5): 2031-2035.
4. Quantification of *Campylobacter* spp. in chicken carcass rinse by real-time PCR. N. Botteldoorn *et al.* *Journal of Applied Microbiology* 2008; 105(6):1909-1918.
5. Detection of *Shigella* by a PCR assay targeting the *ipaH* gene suggests increased prevalence of shigellosis in Nha Trang, Vietnam. V. D. Thiem *et al.* *Journal of Clinical Microbiology*
6. 2004; 42(5): 2031-2035.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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