Vitassay qPCR

Bordetella pertussis + B. parapertussis+ B. holmesii

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Intended use

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* allows the detection and differentiation of *Bordetella pertussis*, *Bordetella parapertussis* y/o *Bordetella holmesii* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Bordetella pertussis*, *Bordetella parapertussis* y/o *Bordetella holmesii* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* 4x8 -well strip, low profile 7041020

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* 4x8-well strip, high profile 7042020

Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S020/ 7042S020	<i>Bordetella pertussis</i> + <i>B. parapertussis</i> + <i>B. holmesii</i> strips low/high profile	-	4x8-well strip
7C020	<i>Bordetella</i> Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)

- Filter tips
- Powder-free disposal gloves

Summary

The genus *Bordetella* bacteria are gram negative bacilli which are isolated from the upper respiratory tract. *Bordetella pertussis* is the main causative agent of whooping cough, although *Bordetella parapertussis*, can also produce the disease, but with a milder symptomatology. Additionally, *Bordetella bronchiseptica* and *Bordetella holmesii* have been identified as causing respiratory tract diseases with clinical manifestations similar to that of the pertussis syndrome but mainly in immunocompromised individuals

The estimated incidence of this disease worldwide ranges from 30 to 50 million cases causing about 195,000 deaths a year. Most cases occur in developing countries where vaccination programs do not exist. In developed countries, despite the dramatic decrease in the incidence of this disease since the introduction of vaccination in the 1950s, a global resurgence of the disease has been observed in the last two decades.

Isolation of *B. pertussis* in culture is definitive for diagnosis, and although the low sensitivity (50%) remains the diagnostic reference method. PCR allows rapid diagnosis and enhances culture sensitivity.

Principle of the test

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* test is based on the real-time amplification of a conserved region of the insertion sequence *IS481* for the identification of *Bordetella pertussis/Bordetella holmesii*, *Bordetella holmesii* for the insertion sequence *hIS1001* and *Bordetella parapertussis* for the insertion sequence *pIS1001*. After DNA isolation, the presence of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Bordetella pertussis/Bordetella holmesii* DNA target sequence is detected through the FAM channel, *Bordetella holmesii* DNA target in ROX channel and *Bordetella parapertussis* DNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

EZ1 Virus Mini Kit (Qiagen)

RIDA® Xtract (r-Biopharm).

Invisorb® Spin Universal Kit (Strattec).

Positive control preparation

Reconstitute the lyophilized Bordetella Positive Control (red tube) in the 100 µL of PCR grade water (transparent tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) into each well.
- Pipette 5 µL of DNA sample, negative and positive controls into each well.
- Cover the wells with the caps provided. Spin down briefly(optional).
- Place the strips in the Real-time PCR instrument.

Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*Bordetella pertussis/Bordetella holmesii*), ROX (*Bordetella holmesii*), Cy5 (*Bordetella parapertussis*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached 2)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Bordetella pertussis/Bordetella holmesii*), ROX (*Bordetella holmesii*) and Cy5 (*Bordetella parapertussis*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

<i>Bordetella pertussis/Bordetella holmesii</i>	<i>Bordetella holmesii</i>	<i>Bordetella parapertussis</i>	Internal control	Negative control	Positive control	Interpretation
+	-	-	+/-	-	+	<i>Bordetella pertussis</i> Positive, <i>B. holmesii</i> and <i>B. parapertussis</i> Negatives
+	+	-	+/-	-	+	<i>Bordetella holmesii</i> Positive, <i>B. pertussis</i> and <i>B. parapertussis</i> Negatives
-	-	+	+/-	-	+	<i>Bordetella parapertussis</i> Positive, <i>B. holmesii</i> and <i>B. pertussis</i> Negatives
-	-	-	+	-	+	<i>B. pertussis</i> , <i>B. holmesii</i> and <i>B. parapertussis</i> Negatives
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 35 nasopharyngeal aspirates and pernasal swabs specimens from symptomatic patients were tested by Real Time PCR using Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* and RIDA@GENE *Bordetella* (R-biopharm).

Bordetella pertussis was detected in 14 samples by Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* but only 1 specimen could be confirmed as *Bordetella parapertussis* positive. 20 specimens resulted negative for *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii*.

The specificity of Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* was 100% for *Bordetella pertussis*, *Bordetella parapertussis* y *Bordetella holmesii*.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity for *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

<i>Legionella bozemanii</i>	<i>Haemophilus influenzae</i>	Influenza B/Phuket/3073/2013
<i>Legionella micdadei</i>	<i>Moraxella catarrhalis</i>	Respiratory syncytial virus (RSV)
<i>Legionella dumoffii</i>	Influenza A/California/7/2009(H1N1)	Human metapneumovirus A and B
<i>Legionella longbeachae</i>	Influenza A/Perth/16/2009(H3N2)	Human parainfluenza 1, 2, 3, 4
<i>Legionella pneumophila</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Human coronavirus 229E
<i>Mycoplasma pneumoniae</i>	Influenza A/Switzerland/9715293/2013	Human rhinovirus
<i>Streptococcus pneumoniae</i>	Influenza A/Turkey/Germany R2485+86/2014	Human Adenovirus 2 strain Adenoid 6
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza B/Brisbane/60/2008	Human Adenovirus 5
Methicillin-resistant <i>Staphylococcus aureus</i>	Influenza B/Florida/04/06	

Analytical reactivity

Analytical reactivity of Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* was evaluated against specific clinical strains of each species, showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR *Bordetella pertussis* + *B. parapertussis* + *B. holmesii* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Bordetella pertussis*, *Bordetella holmesii* and *Bordetella parapertussis* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with nasopharyngeal aspirates and pernasal swabs samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different bacteria, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
Bio-Rad
CFX96 Touch™ Real-Time PCR Detection System
Roche
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Agilent Technologies
AriaMx Real-Time PCR System
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Qiagen
Rotor-Gene®
Cepheid
SmartCycler®

High profile Block Thermocyclers
Applied Biosystems
7500 Real-Time PCR
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Real-Time PCR
Bio-Rad
CFX96 Touch™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
Eppendorf
Mastercycler™ep <i>realplex</i>
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System
Analytik Jena Biometra
TOptical
qTOWER 2.0
Abbott
Abbott m2000 RealTime System
BIONEER
Exicycler™ 96
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time
Qiagen
Rotor-Gene®
Cepheid
SmartCycler®

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

Real-time pcr thermocycler	Vitassay channel	Detection channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Passive reference option ROX is none
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Bibliography/Bibliografía

1. Novel Multitarget Real-Time PCR Assay for Rapid Detection of *Bordetella* Species in Clinical Specimens Kathleen M. Tatti, Kansas N. Sparks, Kathryn O. Boney, Maria Lucia Tondella *J Clin Microbiol.* 2011 Dec; 49(12): 4059–4066.
2. Detection and Differentiation of *Bordetella* spp. by Real-Time PCR. Christoph Koidl, Michael Bozic, Anja Burmeister, Markus Hess, Egon Marth, Harald H. Kessler. *J Clin Microbiol.* 2007 Feb; 45(2): 347–350.
3. Real-Time PCR Assay Targeting IS481 of *Bordetella pertussis* and Molecular Basis for Detecting *Bordetella holmesii*. Udo Reischl, Norbert Lehn, Gary N. Sanden, Mike J. Loeffelholz *J Clin Microbiol.* 2001 May; 39(5): 1963–1966).

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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