Vitassay qPCR

Flu A + Flu B

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- Applied Biosystems 7500 Fast Real-Time PCR System
- Applied Biosystems StepOne™ Real-Time PCR System
- Bio-Rad CFX96™ Real-Time PCR Detection System
- AriaMx Real-Time PCR System
- DNA-Technology DTPRime Real Time Detection Thermal Cyclers
- Rotor-Gene® Q
- SmartCycler®

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Intended use

Vitassay qPCR Flu A+ Flu B allows the detection and differentiation of Influenza A (Flu A) and Influenza B (Flu B) virus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Influenza virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Flu A+ Flu B 4x8-well strip, low profile	7041021
Vitassay qPCR Flu A+ Flu B 4x8-well strip, high profile	7042021

Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S021/ 7042S021	Flu A+ Flu B Virus strips low/high profile	-	4x8-well strip
7C021	Flu A+ Flu B Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Influenza A and B viruses cause the majority of viral lower respiratory tract infections; elderly and compromised individuals are especially at risk of developing severe illness and complications. Influenza viruses are a significant cause of morbidity and mortality worldwide.

Seasonal influenza epidemics impose a heavy burden on society, with 3–5 million cases and 250-500 000 deaths worldwide every year. The resulting economic impact is large and includes both direct and indirect costs.

After an incubation period of one to two days, the illness has an abrupt onset, with chills, fever, myalgia, headache, and anorexia. Pneumonia is the most important complication of influenza. There may be primary viral pneumonia, secondary bacteria pneumonia, or a combination of both.

Laboratory diagnosis of influenza has been accomplished using cell culture, rapid antigen testing, and, more recently, real-time PCR.

Principle of the test

Vitassay qPCR Flu A+ Flu B test is based on the real-time amplification of specific conserved fragments of the *M1* gene. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Flu A+ Flu B test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Influenza A RNA target sequence is detected through the FAM channel, Influenza B virus RNA target in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

RIDA® Xtract (r-Biopharm).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (ROCHE).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized Flu A+ Flu B Positive Control (red tube) in the 100 µL of PCR grade water (white tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) into each well.
- Pipette 5 µL of RNA sample, negative and positive controls into each well
- Cover the wells with the caps provided. Spin down briefly(optional).
- Place the strips in the Real-time PCR instrument.

Programm your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	63°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Influenza A), ROX (Influenza B) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3021P™ check that passive reference option ROX is none. (Attached 2)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (Influenza A) and ROX (Influenza B), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, and ROX which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Influenza A Virus	Influenza B Virus	Internal control	Negative control	Positive control	Interpretation
+	+	+/-	-	+	Influenza A and Influenza B viruses Positive
-	-	+	-	+	Influenza A and Influenza B viruses Negative
+	-	+/-	-	+	Influenza A Positive, Influenza B negative
-	+	+/-	-	+	Influenza B Positive Influenza A Negative
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 109 throat swabs from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Flu A+ Flu B and CLART® PneumoVir DNA array assay (Genomica).

Influenza A virus was detected in 75 samples by Vitassay qPCR Flu A+ Flu B Real Time PCR and in 77 samples by CLART® PneumoVir DNA array assay. For these 2 discordant samples, the low amount of template RNA is below the detection limit of the method used. For Influenza B virus, a total of 8 clinical samples were positives for both assays.

The specificity of Vitassay qPCR Flu A+ Flu B was 100% for both viruses (32 samples for Influenza A and 101 samples for influenza B).

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Influenza A and Influenza B virus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Influenza A and Influenza B virus was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Influenza A	Influenza B
<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>
<i>Legionella bozemanii</i>	<i>Legionella bozemanii</i>
<i>Legionella micdadei</i>	<i>Legionella micdadei</i>
<i>Legionella dumoffii</i>	<i>Legionella dumoffii</i>
<i>Legionella longbeachae</i>	<i>Legionella longbeachae</i>
<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>
<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
Methicillin-resistant <i>Staphylococcus aureus</i>	Methicillin-resistant <i>Staphylococcus aureus</i>
<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>
<i>Moraxella catarrhalis</i>	<i>Moraxella catarrhalis</i>
Influenza B/Brisbane/60/2008 virus	Influenza A/California/7/2009(H1N1) virus
Influenza B/Florida/04/06 virus	Influenza A/Perth/16/2009(H3N2) virus
Influenza B/Phuket/3073/2013	Influenza A/New Caledonia/20/99(H1N1) virus
Respiratory syncytial virus (RSV)	Influenza A/Switzerland/9715293/2013
Human metapneumovirus A and B	Influenza A/Turkey/Germany R2485+86/2014
Human coronavirus 229E	Human parainfluenza 1, 2, 3 and 4 viruses
Human rhinovirus	Respiratory syncytial virus (RSV)
Human Adenovirus 2 strain Adenoid 6	Human metapneumovirus A and B
Human Adenovirus 5	Human coronavirus 229E
Human parainfluenza 1, 2, 3 and 4 viruses	Human rhinovirus
	Human Adenovirus 2 strain Adenoid 6
	Human Adenovirus 5

Analytical reactivity

Analytical reactivity of Vitassay qPCR Flu A+ Flu B Virus Real Time PCR for Influenza A was evaluated against strains: A/California/7/2009(H1N1)pdm09-like virus, A/Perth/16/2009(H3N2)like virus, A/New Caledonia/20/99(H1N1), A/Switzerland/9715293/2013 A/Turkey/Germany R2485+86/2014, showing positive result.

Analytical reactivity of Vitassay qPCR Flu A+ Flu B Virus Real Time PCR for Influenza B was evaluated against strains B/Brisbane/60/2008-like virus (B/Victoria lineage), B/Florida/04/06, B/Phuket/3073/2013 (B/Yamagata lineage), showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Flu A+ Flu B has been validated on the following equipments:

Applied Biosystems 7500 Fast Real-Time PCR System
Applied Biosystems StepOne™ Real-Time PCR System
Bio-Rad CFX96™ Real-Time PCR Detection System
Agilent Technologies AriaMx Real-Time PCR System
DNA-Technology DTprime Real-time Detection Thermal Cycler
Rotor-Gene® Q (Qiagen)
SmartCycler® (Cepheid).

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Influenza virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swabs samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler in the list below, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300/7500/7900HT Real-Time PCR
7500 Fast Dx Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 12K Flex 96-well Fast	QuantStudio™ 6 Flex 96-well
QuantStudio™ 6 Flex 96-well Fast	QuantStudio™ 7 Flex 96-well
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 3 Real-Time PCR
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
StepOne Plus™ Real-Time PCR System	ViiA™ 7 Real-Time PCR
StepOne™ Real-Time PCR System	Bio-Rad
ViiA™ 7 Fast Real-Time PCR System	CFX96 Touch™ Deep Well Real-Time PCR Detection
Bio-Rad	iCycler iQ™ Real-Time PCR
CFX96 Touch™ Real-Time PCR Detection System	iCycler iQ™5 Real-Time PCR
Roche	Eppendorf
LightCycler @480 Real-Time PCR System	Mastercycler™ep <i>realplex</i>
LightCycler @96 Real-Time PCR System	Stratagene / Agilent Technologies
Agilent Technologies	Mx3000P™ Real Time PCR System
AriaMx Real-Time PCR System	Mx3021P™ Real Time PCR System
DNA-Technology	Analytik Jena Biometra
DTlite Real-Time PCR System	TOptical
DTprime Real-time Detection Thermal Cyclers	qTOWER 2.0
Qiagen	Abbott
Rotor-Gene®	Abbott m2000 RealTime System
Cepheid	BIONEER
SmartCycler®	Exicycler™ 96
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time
	Qiagen
	Rotor-Gene®
	Cepheid
	SmartCycler®

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/60	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3021P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Bibliography/Bibliografía

1. Methods for molecular surveillance of influenza. Wang R, Taubenberger JK. Expert review of anti-infective therapy. 2010;8(5):517-527.
2. WHO information for molecular diagnosis of influenza virus—update. World Health Organization. Available: http://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/. Accessed 2015 Dec 30.
3. Kuo RL, Yang SL, Liu YC, Chen LT, Mok CK, Kuo SM, Shih SR, Tsao KC. Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. J Virol Methods. 2014 Nov;208:41-6.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3021™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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