# Vitassay qPCR

## Metapneumovirus

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	<b>Eppendorf</b>
	Mastercycler™ ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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## Intended use

Vitassay qPCR Metapneumovirus allows the detection of Metapneumovirus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Metapneumovirus infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Metapneumovirus 4x8 -well strip, low profile 7041024

Vitassay qPCR Metapneumovirus 4x8-well strip, high profile 7042024

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S024/ 7042S024	Metapneumovirus strips low/high profile	-	4x8-well strip
7C024	Metapneumovirus Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

## Summary

First identified in 2001, Metapneumovirus is a newly discovered world-wide distributed pathogen, so limited epidemiological data is available. It is known that it causes both upper and lower respiratory tract infections, including bronchiolitis and pneumonia. These infections have been proven to be indistinguishable from the ones caused by RSV and influenza. Their symptoms having a lot in common and being: cough, fever, irritability, anorexia, wheezing and rhinitis. Metapneumovirus mainly affects children, but in adults appears with flu-like symptoms, and is a major cause of morbidity in elderly and immunocompromised patients. Metapneumovirus belongs to the *Paramyxoviridae* family, more concretely to the *Metapneumovirus* genus. It possesses two genotypes, A and B. Then, each genotype has two subtypes, these being: A1, A2 (which coincidentally is the most common among hosts), B1 and B2. Lastly, the A2 subtype is also divided into A2a and A2b sub-genotypes. The virus is characterized by a late winter-early spring epidemiologic pattern, even though summer infections cannot be excluded. As acute respiratory-tract infections are a major cause of disability in all age groups and isolation of the virus in cell culture requires prolonged incubation, identification through traditional techniques is difficult and slow. Real Time PCR has proven to increase the sensitivity of viral detection and is able to deal with traditional techniques problems. In order to do so, we target the *N* gene of the pathogen.

## Principle of the test

Vitassay qPCR Metapneumovirus test is based on the real-time amplification of specific conserved fragments of the *N* gene. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Metapneumovirus test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and RNA extraction**

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (QIAGEN).

High Pure Viral Nucleic Acid Kit (Roche).

Invisorb® Spin Universal kit (Stratec)

### **Positive control preparation**

Reconstitute the lyophilized Metapneumovirus Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

## Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) into each well.
- Pipette 5 µL of RNA sample, negative and positive controls into each well
- Cover the wells with the caps provided. Spin down briefly(optional).
- Place the strips in the Real-time PCR instrument

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Metapneumovirus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3022P™ check that passive reference option ROX is none. (Attached II)

## Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

### Positive control

The positive controls used in each run, must show an amplification curve for FAM (Metapneumovirus), which validates the reaction.

## Negative control

The negative controls included in each run, must show the absence of signal in FAM, which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Metapneumovirus	Internal control	Negative control	Positive control	Interpretation
+	+/-	-	+	Metapneumovirus Positive
-	+	-	+	Metapneumovirus Negative
+	+	+	+	Experiment fail
-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

A total of 86 respiratory samples from symptomatic patients were analyzed by two molecular diagnostic tests: Vitassay qPCR Metapneumovirus and CLART® PneumoVir DNA array assay (Genomica) + RealStar® hMPV RT-PCR (Altona Diagnostics).

Human Metapneumovirus was detected in all 45 samples using Vitassay qPCR Metapneumovirus.

These results show the high sensitivity and specificity to detect Human Metapneumovirus using Vitassay qPCR Metapneumovirus.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Bocavirus virus templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

### Analytical specificity

The analytical specificity for Metapneumovirus was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
<i>Bordetella pertussis</i>	<i>Staphylococcus aureus</i> resistente a la Meticilina	Virus Influenza B/Florida/04/06
<i>Legionella bozemanii</i>	<i>Haemophilus influenzae</i>	Virus Influenza B/Phuket/3073/2013
<i>Legionella micdadei</i>	<i>Moraxella catarrhalis</i>	Virus parainfluenza humano 1, 2, 3 y 4
<i>Legionella dumoffii</i>	Virus Influenza A/California/7/2009(H1N1)	Respiratory syncytial virus (RSV)
<i>Legionella longbeachae</i>	Virus Influenza A/Perth/16/2009(H3N2)	Metapneumovirus A y B humano
<i>Legionella pneumophila</i>	Virus Influenza A/New Caledonia/20/99(H1N1)	Coronavirus 229E humano
<i>Mycoplasma pneumoniae</i>	Virus Influenza A/Switzerland/9715293/2013	Rinovirus humano
<i>Streptococcus pneumoniae</i>	Virus Influenza A/Turkey/Germany R2485+86/2014	Adenovirus humano 2 cepa Adenoide 6
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Virus Influenza B/Brisbane/60/2008-like	Adenovirus humano 5

### Analytical reactivity

Analytical reactivity of Vitassay qPCR Metapneumovirus Virus Real Time PCR for Metapneumovirus was evaluated against Metapneumovirus A and B strains, showing positive result for both of them.

## Compatibles real-time PCR equipment

Vitassay qPCR Metapneumovirus has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- StepOne™ Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of Metapneumovirus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with asal, throat, tonsillar and pharyngotonsillar swabs, sputum and bronchoalveolare lavage samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by any of the virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	<b>Eppendorf</b>
	Mastercycler™ ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

**Attached II: Detection channels of most common real time PCR equipment**

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene®Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	150
	HEX	500
	ROX	500
	Cy5	100
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	150
	HEX	3000
	ROX	2000
	Cy5	1500

## **Bibliography/Bibliografía**

1. John V. Williams, et al. Human Metapneumovirus and Lower Respiratory Tract Disease in Otherwise Healthy Infants and Children. *The New England Journal of Medicine*, 2004; 350: 443-450.
2. Guy Boivin, et al. Virological Features and Clinical Manifestations Associated with Human Metapneumovirus: A New Paramyxovirus Responsible for Acute Respiratory-Tract Infections in All Age Groups. *The Journal of Infectious Diseases*, 2002; 186 (9): 1330-1334.
3. Wenhua Kong, et al. Circulation of Human Metapneumovirus Among Children With Influenza-Like Illness in Wuhan, China. *Journal of Medical Virology*, 2016; 88: 774-781.
4. Ann R. Falsey, et al. Human Metapneumovirus Infections in Young and Elderly Adults. *The Journal of Infectious Diseases*, 2003; 187 (5): 785-790.
5. Bernadette G, van den Hoogen, et al. Prevalence and Clinical Symptoms of Human Metapneumovirus Infection in Hospitalized Patients. *The Journal of Infectious Diseases*, 2003; 188 (10): 1571-1577.
6. Sai-Zhen Zeng, et al. Clinical Features of Human Metapneumovirus Genotypes in Children With Acute Lower Respiratory Tract Infection in Changsha, China. *Journal of Medical Virology*, 2015; 87: 1839-1845.
7. Amber K. Haynes, et al. Human Metapneumovirus Circulation in the United States, 2008 to 2014. *PEDIATRICS*, 2016; Volume 137: number 5

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3022™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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