
Vitassay qPCR

Pneumocystis jirovecii

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Termocicladores con bloque de bajo perfil	T. con bloque de alto perfil
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300/7500/7900HT Real-Time PCR
7500 Fast Dx Real-Time PCR System	ABI PRISM 7000/7700
QuantStudio™ 12K Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well Fast	QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 3 Real-Time PCR
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
StepOne Plus™ Real-Time PCR System	ViiA™ 7 Real-Time PCR
StepOne™ Real-Time PCR System	Bio-Rad
ViiA™ 7 Fast Real-Time PCR System	CFX96 Touch™ Deep Well Real-Time PCR Detection
Bio-Rad	iCycler iQ™ Real-Time PCR
CFX96 Touch™ Real-Time PCR Detection System	iCycler iQ™5 Real-Time PCR
Mini Opticon™ Real-Time PCR Detection System	MyiQ™ Real-Time PCR
Roche	MyiQ™2 Real-Time PCR
LightCycler®480 Real-Time PCR System	Eppendorf
LightCycler®96 Real-Time PCR System	Mastercycler™ep <i>realplex</i>
Agilent Technologies	Stratagene / Agilent Technologies
AriaMx Real-Time PCR System	Mx3000P™ Real Time PCR System
DNA-Technology	Mx3005P™ Real Time PCR System
DTlite Real-Time PCR System	Analytik Jena Biometra
DTprime Real-time Detection Thermal Cycler	TOptical
Qiagen	qTOWER 2.0
Rotor-Gene®	Abbott
Cepheid	Abbott m2000 RealTime System
SmartCycler®	BIONEER
	Exicycler™ 96
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time
	Qiagen
	Rotor-Gene®
	Cepheid
	SmartCycler®

Intended use

Vitassay qPCR *Pneumocystis jirovecii* allows the qualitative and quantitative detection of *Pneumocystis jirovecii* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *P. jirovecii* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Pneumocystis jirovecii 4 x 8-well strip, low profile 7041026
 Vitassay qPCR Pneumocystis jirovecii 4 x 8-well strip, high profile 7042026

Materials/reagents provided

Código	Reactivo/Material	Color	Cantidad
7041S026/ 7042S026	<i>Pneumocystis jirovecii</i> strips low/high profile	-	4x8-well strip
7C026	<i>Pneumocystis jirovecii</i> Positive Control	red	1 vial
7Q026	<i>Pneumocystis jirovecii</i> Quantitative standard	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x1 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, we recommend to separate in aliquots.
- Keep all reagents of in the dark.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Pneumocystis pneumonia (PCP) is a serious illness caused by the fungus *Pneumocystis jirovecii*. PCP is one of the most frequent and severe opportunistic infections in people with weakened immune systems, particularly people with HIV/AIDS. Although people with HIV/AIDS are less likely to get PCP today than in recent years, PCP is still a significant public health problem.

The following groups are at risk for *P. jirovecii* pneumonia: persons with HIV infection with low CD4 count, persons with primary immune deficiencies, persons receiving long-term immunosuppressive regimens, persons with hematologic and nonhematologic malignancies, persons with severe malnutrition.

The presentation of PCP in HIV-positive patients is well-known and consists of a triad of dyspnea, fever, and cough, whereas the presentation of PCP in HIV-negative patients is atypical and consists of a sudden outbreak, O₂ desaturation, and a rapid lethal outcome without therapy.

The laboratory diagnosis of Pneumocystis pneumonia (PCP), still relies on tinctorial and/or immunofluorescent staining of bronchoalveolar lavage (BAL) fluid samples. Nucleic acid amplification tests can overcome the difficulties of microscopic examination. Real-time quantitative PCR (qPCR) reduces the risk of false-positive results because of the closed-tube nature of the amplification process, and the resulting data are semi-quantitatives.

Principle of the test

Vitassay qPCR *Pneumocystis jirovecii* is based on the real-time amplification of specific conserved fragments of the *mt large-subunit rRNA* gene of *Pneumocystis jirovecii* genome. After DNA extraction, *P. jirovecii* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR *Pneumocystis jirovecii* is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence is detected through the FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used). In addition, quantification of *P. jirovecii* DNA can be achieved by generating a standard curve using *P. jirovecii* quantitative standard provide with the kit.

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the laboratories.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

EZ1 Mini Kit, EZ1 instrument (Qiagen)

Invisorb® Spin Universal Kit (Strattec).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega)

Positive control and Quantitative Standard preparation

Reconstitute the lyophilized *Pneumocystis jirovecii* Positive Control and the Quantitative Standard (red tubes) in the 100 µL of PCR grade water (white tube) supplied. To ensure a complete resuspension, vortex the tube thoroughly.

For *Pneumocystis jirovecii* Positive Control, after first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Standard Curve Preparation

To perform a quantitative assay, a standard curve must be prepared by serial dilution of *Pneumocystis jirovecii* Quantitative Standard (which contains approximately 2×10^7 copies/ μl^*). Prepare the standard curve dilution series as follows:

- Pipette 90 μl of PCR grade Water in 6 microcentrifuge tubes.
- Add 10 μl of *Pneumocystis jirovecii* Quantitative Standard to the first tube to get a standard with about 2×10^6 copies/ml. Mix by vortex and centrifuge briefly.
- Add 10 μl of standard with $\sim 2 \times 10^6$ copies/ μl to the second tube to get a standard with approximately 2×10^5 copies/ μl . Mix by vortex and centrifuge briefly.
- Repeat the previous step sequentially to complete the dilution series for standards with about 2×10^4 , 2×10^3 , 2×10^2 and 2×10^1 copies/ μl .

* For more details see “Certificate of Analysis”, where lot number and DNA copy number of *Pneumocystis jirovecii* Quantitative Standard are detailed.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. Peel off protective aluminium seal from the strips.
- Pipette 15 μL of resuspension buffer (green tube) into each well.
- Pipette 5 μL of DNA sample, negative, positive controls and standard curve dilutions into each well
- Cover the wells with the optical caps provided. Spin down briefly.
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler.

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*Pneumocystis jirovecii*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none.

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Pneumocystis jirovecii*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal for FAM (*Pneumocystis jirovecii*), which validates the reaction.

The experiment seems to be fail if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

Qualitative analysis

The result interpretation is summarized in the following table:

<i>Pneumocystis jirovecii</i> FAM	Control Interno HEX	Negative Control	Positive Control	Interpretation
+	+/-	-	+	<i>Pneumocystis jirovecii</i> Positive
-	+	-	+	<i>Pneumocystis jirovecii</i> Negative
+	+	+	+	Invalid
-	-	-	-	Invalid

+ Positive: Amplification signal
-Negative: No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quantitative analysis

Using standard curve dilution series results a standard curve for quantitative analysis can be generated:

$$Ct = m \log (Q) + b \quad \text{Where: } \begin{array}{l} Ct = \text{Threshold Cycle} \\ m = \text{Slope} \\ Q = \text{Concentration} \\ b = \text{Intercept} \end{array}$$

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{(Ct-b)/m}$$

The concentration of the samples is displayed in copies/μl and refers to the concentration in the Eluted DNA not to the original clinical sample. To determine the target concentration of the original sample, take into account the dilutions of the extraction procedure and PCR Set-up.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The clinical performance of Vitassay qPCR *P. jirovecii* was evaluated using 14 bronchoalveolar lavage from patient's with suspicious of *Pneumocystis* pneumonia. We obtained a positive result for 11 samples and a negative result for 3 of them. The results were compared with other real-time PCR commercial kit (FTD *Pneumocystis jirovecii*, fast-track DIAGNOSTICS). The concordance between both tests was 100%.

The results show a high sensitivity and specificity to detect *Pneumocystis jirovecii* using VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Pneumocystis jirovecii* template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity for *Pneumocystis jirovecii* was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

<i>Staphylococcus aureus</i>	Influenza A/California/7/2009(H1N1) virus
<i>Bordetella pertussis</i>	Influenza A/Perth/16/2009(H3N2) virus
<i>Fluoribacter bozemaniae</i>	Influenza A/New Caledonia/20/99(H1N1) virus
<i>Legionella micdadei</i>	Influenza A/Switzerland/9715293/2013 virus
<i>Legionella dumoffii</i>	Influenza A/Turkey/Germany R2485+86/2014 virus
<i>Legionella longbeachae</i>	Influenza B/Brisbane/60/2008-like virus
<i>Legionella pneumophila</i>	Influenza B/Florida/04/06 virus
<i>Mycoplasma pneumoniae</i>	Influenza B/Phuket/3073/2013 virus
Methicillin-resistant <i>Staphylococcus aureus</i>	Human Adenovirus 5
<i>Haemophilus influenzae</i>	Human Adenovirus 2 strain Adenoid 6
<i>Moraxella catarrhalis</i>	Human coronavirus 229E
<i>Streptococcus pneumoniae</i>	Human parainfluenza 1, 2, 3 and 4 viruses
Respiratory syncytial virus (RSV)	Human metapneumovirus A and B
Human rhinovirus	

Analytical reactivity

The reactivity of Vitassay qPCR *Pneumocystis jirovecii* was confirmed by the real-time amplification using *Pneumocystis jirovecii* strain as template.

Compatibles real-time PCR equipment

Vitassay qPCR *Pneumocystis jirovecii* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- SmartCycler® (Cepheid)

For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Pneumocystis jirovecii* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- Although this assay can be used with other types of samples it has been validated with bronchoalveolar lavage (BAL) specimens.
- The quality of the test depends on the quality of the sample; proper DNA from respiratory clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Pneumocystis jirovecii*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached 1: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you don't find your thermocycler in the list, please contact with your supplier.

Low profile Block Thermocyclers
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 3 Real-Time PCR System
QuantStudio™ 5 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
Bio-Rad
CFX96 Touch™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Roche
LightCycler @480 Real-Time PCR System
LightCycler @96 Real-Time PCR System
Agilent Technologies
AriaMx Real-Time PCR System
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Qiagen
Rotor-Gene ®
Cepheid
SmartCycler®

High profile Block Thermocyclers
Applied Biosystems
7300/7500/7900HT Real-Time PCR
ABI PRISM 7000/7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Real-Time PCR
QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Real-Time PCR
Bio-Rad
CFX96 Touch™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™5 Real-Time PCR
MyiQ™ Real-Time PCR
MyiQ™2 Real-Time PCR
Eppendorf
Mastercycler™ <i>ep realplex</i>
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System
Analytik Jena Biometra
TOptical
qTOWER 2.0
Abbott
Abbott m2000 RealTime System
BIONEER
Exicycler™ 96
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time
Qiagen
Rotor-Gene ®
Cepheid
SmartCycler®

Bibliography/Bibliografía

1. *Pneumocystis jirovecii* pneumonia in patients with and without human immunodeficiency virus infection. Su YS, Lu JJ, Perng CL, Chang FY. *J Microbiol Immunol Infect*. 2008 Dec. 41(6):478-82.
2. *Pneumocystis* pneumonia in patients with HIV infection: clinical manifestations, laboratory findings, and radiological features. Fujii T, Nakamura T, Iwamoto A. *J Infect Chemother*. 2007 Feb. 13(1):1-7.
3. Detection of *Pneumocystis jirovecii* by two staining methods and two quantitative PCR assays. Rohner P, Jacomo V, Studer R, Schrenzel J, Graf JD. *Infection*. 2009 Jun;37(3):261-5.
4. Cryptococcal vs *Pneumocystis (carinii) jiroveci* pneumonia (PCP): Clinical & Microbiology differential diagnostic considerations. Cunha BA, Schoch PE, Berbari N. *Infect Dis Pract*. 2006. 30:514-517.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.



Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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