Vitassay qPCR

Tick Borne Disease

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Intended use

Vitassay qPCR Tick Borne Disease allows the detection and differentiation of viral RNA or genomic DNA specific for Tick Borne Encephalitis Virus (TBEV), *Rickettsia* spp., *Babesia microti*, *Babesia divergens*, *Ehrlichia chafeensis*, *Ehrlichia muris*, *Borrelia burgdorferi* sensu lato (*s.l.*), *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* in blood, serum, tissue samples and microbiological culture from ticks, biopsy skin, cerebrospinal fluid and synovial fluid from patients with signs and symptoms of Tick Borne diseases. This test is intended for use in the Tick-Borne Disease diagnosis alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR Tick Borne Disease 12x8-well strip, low profile	7041062
Vitassay qPCR Tick Borne Disease 12x8-well strip, high profile	7042062

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S062A/ 7042S062A	<i>Borrelia</i> , <i>Anaplasma</i> & <i>Coxiella</i> strips low/high profile	-	4 x 8-well strip
7041S062B/ 7042S062B	<i>Rickettsia</i> , <i>Babesia</i> & <i>Ehrlichia</i> strips low/high profile	-	4 x 8-well strip
7041S062C/ 7042S062C	TBEV strips low/high profile	-	4 x 8-well strip
7C062	Tick Borne Disease Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.

- Keep all reagents in the darkness.

Additional equipment and material required

- RNA/DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Tick populations are increasing, and their geographic ranges are expanding, as the suitable habitats for these arthropod vectors and the pathogens that they carry. Ticks can transmit bacterial, parasitic, and viral pathogens and often harbor more than one agent simultaneously. They also can transmit more than one pathogen when taking a blood meal. Tick Borne Disease affecting humans include Lyme disease, relapsing fever, babesiosis, anaplasmosis, ehrlichiosis, Q fever, Tick Borne encephalitis and spotted fever among others.

***Rickettsia* spp.**

Members of the genus *Rickettsia* (family *Rickettsiaceae*; order *Rickettsiales*) may be classified into spotted fever group (SFG) *rickettsiae*, typhus group *rickettsiae*, the *Rickettsia bellii* group, and the *Rickettsia canadensis* group. Spotted fever is caused by bacteria of the genus *Rickettsia* and is widely distributed by different geographical areas, being able to transmit depending on it by different ticks such as: *Amblyomma variegatum* (Sub-Saharan Africa), *Hyalomma anatolicum excavatum* and *Ixodes ricinus* (Europe), *Demacentor andersoni* and *Demacentor nitens* (North and Central America), *Carios capensis* (Japan, Asia) and *Amblyomma cajennense* (South America). Some of the spotted fever symptoms are fever, headache, malaise, myalgias, nausea, rash, vomiting and abdominal pain.

Anaplasma

Anaplasma phagocytophilum is a zoonotic tick-borne pathogen transmitted by infected *Ixodes* ticks. White-tailed deer and the white-footed mouse are its primary reservoirs. *A. phagocytophilum* is distributed over large portions of Europe and parts of central Asia, where it is vectored by *Ixodes ricinus* and *Ixodes persulcatus* ticks. In North America, is now found on the East Coast and in midwestern regions of the United States, transmitted by *Ixodes scapularis*, and on the Pacific Coast, transmitted by *Ixodes pacificus*. Some of the anaplasmosis symptoms are high fever, chills, malaise,

generalized myalgia and rash, severe headache, and gastrointestinal and respiratory symptoms.

Borrelia

The *Borrelia* genus is divided into the Lyme disease group (*B. burgdorferi* sensu lato complex) and relapsing fever *Borrelia* species (including *B. hermsii*, *B. turicatae*, *B. parkeri*, and *B. miyamotoi*). The spirochaete *B. burgdorferi* is a zoonotic tick-borne pathogen transmitted by the bite of an infected *Ixodes* tick. It is highly invasive and infects more humans in Europe, Asia, Russia and North America than any other tick-borne bacteria. Common signs and symptoms of Lyme disease are fever, chills, headache, fatigue, muscle and joint aches and one or more erythema migrans rashes. If patients get worse, the symptoms will be neurological and cardiac such as cranial nerve palsies or chest pain.

B. miyamotoi and *B. hermsii* are genetically more closely related to the relapsing fever group. While spirochete *B. hermsii* is transmitted by its tick vector, *Ornithodoros hermsi*, spirochete *B. miyamotoi* is transmitted by infected *Ixodes* ticks. Relapsing fever common symptoms are similar to anaplasmosis such as fever, chills, myalgia, fatigue, arthralgia, lymphadenopathy and possible erythema migrans rash.

Ehrlichia

Ehrlichioses are known as important emerging tick-borne diseases in humans, as well as in domestic animals, and are caused by infection with *Ehrlichia* spp. The most common *Ehrlichia* spp. to infect humans are *E. chaffeensis* and *E. muris*. These pathogens are transmitted by infected *Amblyomma americanum* ticks, and its major reservoir is the white-tailed deer. Ehrlichiosis common symptoms are similar to anaplasmosis such as fever, headache, chills, malaise, myalgia and nausea. Other symptoms are cough, arthralgia or rash. Severe complications may include renal failure, pneumonia, acute respiratory distress syndrome, neurologic disorder and intravascular coagulation.

Coxiella

Coxiella burnetii is the agent of Q fever, or “query fever,” a zoonosis first described in Australia in 1937. The main reservoirs of *C. burnetii* are cattle, sheep, and goats and transmission to humans is most frequently due to inhalation of aerosolized bacteria and urine, feces or milk of infected animals, and by tick bites. The tick species that are the most frequent vectors of *C. burnetii* belong to the genera *Ixodes*, *Rhipicephalus*, *Amblyomma*, and *Dermacentor*. The Q fever presentation is extremely variable, and infection may lead to asymptomatic, acute disease or chronic infection. The mild symptoms are flu-like such as fever, chills, headache or cough, and the severe cases of patients may develop pneumonia or hepatitis.

Babesia

Babesiosis is a parasitic disease caused by protozoa of the genus *Babesia*. In humans, the most prevalent are infections caused by *B. microti* and less frequently *B. divergens*. The pathogen vector is the ticks of the genus *Ixodes*, and the infections occur by tick bites. Disease is manifested by general malaise, fatigue and flu-like symptoms including fever, chills, sweating, joint and muscle pain. Because of the excessive erythrocyte lysis, babesiosis can lead to hemolytic anemia, intravascular coagulopathy, or thrombocytopenia, among others.

Virus Tick Borne Encephalitis

Tick-borne encephalitis virus (TBEV), which is a member of the genus *Flavivirus*, within the family *Flaviviridae*, is the causative agent of tick-borne encephalitis (TBE), a potentially fatal neurological infection affecting humans in Europe and Asia. The majority of TBEV infections occur through a tick bite, although a small number of infections occur through consuming infected unpasteurized milk. During the initial short febrile period can include fatigue, headache and pain in the neck, shoulders, and lower back, accompanied by high fever and vomiting. The patients can progress to a second phase, characterized by acute central nervous system symptoms with a high fever, which can cause meningitis, encephalitis, myelitis or radiculitis.

Test principle

Vitassay qPCR Tick Borne Disease test is based on the real-time amplification of specific conserved fragments of the 3'UTR region for TBEV, the 23S rRNA gene for *Rickettsia spp.*, the CCT-eta rRNA for *Babesia microti*, the hsp70 gene for *Babesia divergens*, the GroEI gene for *Ehrlichia chafeensis* and *Ehrlichia muris*, the 23S rRNA gene for *Borrelia burgdorferi* s.l., *Borrelia miyamotoi* and *B. hermsii*, the msp2 gene for *Anaplasma phagocitophilum* and the IS1111 gene for *Coxiella burnetii*. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in by polymerase chain reaction. After nucleic acids' extraction, the presence of the targets (RNA/DNA) is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Tick Borne Disease test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. Each kit includes three kind of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of *Borrelia burgdorferi* (s.l.), *Borrelia miyamotoi* and/or *B. hermsii*, *Anaplasma phagocitophylum* and *Coxiella burnetii* (*Borrelia*, *Anaplasma* & *Coxiella* 8-well strips). After the amplification reaction, *Borrelia burgdorferi* s.l. and *Borrelia miyamotoi* DNA targets are detected in FAM channel, *Anaplasma phagocitophylum* DNA target is detected in ROX channel, *Coxiella burnetii* DNA target is detected in Cy5 channel and the internal control (IC) is detected in HEX, VIC or JOE channels (depending on the equipment used). The second strip contains the multiplex reaction mix for the detection of *Rickettsia* spp., *Babesia microti* and/or *Babesia divergens* and *Ehrlichia chafeensis* and/or *Ehrlichia muris* (*Rickettsia*, *Babesia* & *Ehrlichia* 8-well strips). After the amplification reaction, *Rickettsia* spp. DNA target is detected in HEX, VIC, or JOE channels (depending on the equipment used), *Babesia microti* and *Babesia divergens* are detected in ROX channel and *Ehrlichia chafeensis* and *Ehrlichia muris* are detected in FAM channel. The third strip contains the monoplex reaction mix for the detection of TBEV (TBEV 8-well strips), whose RNA is amplified and detected in FAM channel.

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Make sure to use a well for determining TBEV, another well for determining *Rickettsia* spp., *Babesia microti*/*Babesia divergens* and *Ehrlichia*

chafeensis/Ehrlichia muris and another well for determining *Borrelia burgdorferi* s.l., and *Borrelia miyamotoi*, *Anaplasma phagocitophilum* and/or *Coxiella burnetii*. Pay attention to not mix them throughout the process.

- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA/DNA extraction

For the isolation, identification, and characterization of pathogens, the proper collection and transport of clinical specimens is critical. Specimens (ticks, blood, serum, tissue samples, biopsy skin, cerebrospinal fluid and synovial fluid) should be properly collected in a clean area and be processed as soon as possible to avoid loss of viability of the etiological agents for suitable microbiological culture and/or to prevent nucleic acid degradation, and to guarantee the test quality.

For longer samples' storage, the samples should be frozen at -20°C. In this case, the sample should be totally thawed and brought to room temperature before testing. Homogenize sample as thoroughly as possible prior to preparation. Freezing and thawing cycles before isolating nucleic acids should be avoided. It is recommend using fresh samples or immediately frozen.

Given that the *B. burgdorferi* spirochetes number in infected tissues or body fluids of patients is very low, for yielding reliable and consistent PCR results, appropriate procedures for sample collection and transport and DNA preparation from clinical samples are critical.

Blood samples might be collected using vacutainer system (tubes stabilized with EDTA) and be stored at 4°C for up to one week. After blood centrifugation, blood components (serum and plasma) can be obtained and should be stored at 4°C for up to one week or at -20°C indefinitely. Nucleic acids isolation from whole blood, serum or plasma could be performed using 100 µL and/or following the manufacturer's recommendations.

Tissue samples should be collected in a clean container and stored immediately at -20°C or -80°C until use.

Ticks could be collected directly from patients using a mechanical method (forceps). For bacterial growth, tissues from each tick could be pooled and cultivated in primary cell line.

For biopsy skin, the cutaneous lesion peripheral border should be identified and under sterile conditions 4-mm-diameter punch biopsy specimen could be obtained from the lesion peripheral aspect. The specimens should be stored immediately at -20°C or -80°C. For microorganism isolation, each skin biopsy specimen could be placed in a polystyrene tube containing 6 ml of isolation/broth medium without antibiotics and incubated at the corresponding temperature during the corresponding time. If the tissue is kept in BSK medium for over 24h, some spirochetes will have migrated from the skin biopsy to the culture medium. In this case, DNA should be prepared from both the skin biopsy and the medium.

Cerebrospinal fluid specimen should be processed in a microbiology laboratory within 1 hour after collection. If processing within 1 hour is not feasible, cerebrospinal fluid specimen should be inoculated into Trans-Isolate (T-I) medium or similar for transport to the laboratory. Refrigeration is not recommended. These samples should be maintained at room temperature prior testing. If samples can be processed within 1 hour, centrifuge the cerebrospinal fluid specimen for 15 min at 1000 x g. Take the sediment and seed it in primary media. If not feasible, incubate the T-I medium over-night.

Synovial fluids samples from patients with Lyme arthritis could be collected without additives and stored overnight initially at 4°C and then transferred to -20°C or -80°C until use.

After specimen microbiology cultivation (ticks, biopsy skin, cerebrospinal fluid and synovial fluid), nucleic acid could be extracted using 1 mL of the culture after being pelleted by centrifugation (2 min 8000 x g). To perform an extraction of bacterial DNA from gram positive Bacteria the lysozyme addition is needed.

To perform more suitable sample collection, transport, and storage, we suggest following the manufacturer's recommendations appearing in the instructions for use of the extraction kit used.

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Invisorb® Spin Universal Kit (Invitex Molecular)

MagDEA Dx SV kit, utilizando el magLEAD® 6gC instrument (Precision System Science Co.)

QIAamp Viral RNA Mini Kit (QIAGEN)

Positive control preparation

Reconstitute the lyophilized Tick Borne Disease Positive Control (red tube) in 300 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA/DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) in the first strip through the FAM (*Borrelia burgdorferi* s.l. and *Borrelia miyamotoi*), HEX, JOE or VIC (Internal Control (IC)), ROX (*Anaplasma phagocitophilum*), and Cy5 channels (*Coxiella burnetii*); in the second strip through the FAM (*Ehrlichia chafeensis* and *Ehrlichia muris*), HEX, JOE or VIC (*Rickettsia* spp.), and ROX (*Babesia microti* and *Babesia divergens*) and in the third reaction mix through the FAM (TBEV). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (See Attached II).

Analysis and interpretation of results

The results' analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve in the pathogen panel channels.

Negative control

The negative control included in each run must show the signal absence in the pathogen panel channels.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive well. The assay should be repeated.

Results' interpretation for *Borrelia*, *Anaplasma* & *Coxiella* 8-well strips:

- A sample is considered positive for *Borrelia burgdorferi* s.l./*Borrelia miyamotoi* if there is an amplification signal in FAM channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal*.
- A sample is considered positive for *Anaplasma phagocitophilum* if there is an amplification signal in ROX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal.
- A sample is considered positive for *Coxiella burnetii* if there is an amplification signal in Cy5 channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal.

*Sometimes, the detection of internal control is not necessary because a high copy number of targets can cause preferential amplification of target-specific nucleic acids.

- A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the internal control amplification.
- In case of internal control signal absence in sample wells we recommend repeating the assay by diluting the sample 1:10 or to repeat the extraction to check for possible inhibition problems.

Results' interpretation for *Rickettsia*, *Babesia* & *Ehrlichia* 8-well strips:

- A sample is considered positive for *Rickettsia* spp. if there is an amplification signal in HEX, VIC or JOE channel, and the Ct value obtained is less than 40.
- A sample is considered positive for *Babesia microti/Babesia divergens* if there is an amplification signal in ROX channel, and the Ct value obtained is less than 40.
- A sample is considered positive for *Ehrlichia chafeensis/Ehrlichia muris* if there is an amplification signal in FAM channel, and the Ct value obtained is less than 40.
- A sample is considered negative, if the sample shows no amplification signal in the detection system, but the internal control of *Borrelia*, *Anaplasma* & *Coxiella* assay is positive.

Results' interpretation for TBEV 8-well strips:

- A sample is considered positive for TBEV if there is an amplification signal in FAM channel, and the Ct value obtained is less than 40.
- A sample is considered negative, if the sample shows no amplification signal in the detection system, but the internal control of *Borrelia*, *Anaplasma* & *Coxiella* assay is positive.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to review all the parameters, and the sigmoid shape of the curve. It is also recommended to repeat the assay, preferably in duplicate.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in the first strip. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The Vitassay qPCR Tick Borne Disease kit's clinical performance for *Borrelia*, *Anaplasma* & *Coxiella* 8-well strips was tested using 95 DNA samples extracted from microbiological culture from ticks, biopsy skin, cerebrospinal fluid, and synovial fluid. In this study, a total of 17 well characterized *Borrelia* strains comprising 9 different *Borrelia burgdorferi* sensu lato genospecies were included. Additionally, two relapsing fever control strains *B. hermsii* and *B. miyamotoi* were included as well as the potentially

cross-reactive spirochaetes *Leishmania* spp. and *Treponema* spp. were tested. The Vitassay qPCR Tick Borne Disease successfully detected all tested *Borrelia* sensu lato genospecies and the relapsing fever group strains *B. hermsii* and *B. miyamotoi*. No cross reactivity was observed with DNA from *Leptospira* and *Treponema* species for Vitassay assay.

Vitassay qPCR Tick Borne Disease was evaluated by INSTAND and QCMD panels from 2017 to 2021. A total of 186 samples from QCMD and INSTAND programs were tested. Vitassay qPCR Tick Borne Disease found 40 samples positive for *Borrelia* spp., 14 samples positive for TBEV and 10 samples positive for *Coxiella*. The results were compared with the EQA program final reports. All samples were correctly detected. Sensitivity and specificity values for each target are shown in the following table:

Target	OA	TP	TN	FP	FN	SE	SP	PPV	NPV
Borrelia	1 (0.96-1)	40	56	0	0	1 (0.91-1)	1 (0.93-1)	1 (0.91-1)	1 (0.93-1)
Coxiella	1 (0.96-1)	10	86	0	0	1 (0.69-1)	1 (0.95-1)	1 (0.69-1)	1 (0.95-1)
TBEV	1 (0.98-1)	14	120	0	0	1 (0.76-1)	1 (0.97-1)	1 (0.76-1)	1 (0.97-1)

OA= Overall agreement, TP = True Positive, TN = True Negative, FP = False Positive, FN = False Negative, SE = Sensibility, SP = Specificity, PPV = Positive Predictive Value, NPV = Negative Predictive Value.

In conclusion, the results show a high sensitivity and specificity to detect those pathogens using Vitassay qPCR Tick Borne Disease.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of ≥ 10 RNA/DNA copies per reaction.

Analytical specificity

The analytical specificity for Tick Borne Disease assay was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species except the targeted pathogens of each assay:

Cross-reactivity testing			
<i>Anaplasma marginale</i>	-	<i>Borrelia garinii</i>	-/+
<i>Bartonella henselae</i> strain Houston-1	-	<i>Borrelia japonica</i>	-/+
<i>Borrelia hermsii</i>	-/+	<i>Borrelia miyamotoi</i>	-/+
<i>Borrelia lusitanae</i>	-/+	<i>Borrelia spielmanii</i>	-/+
<i>Borrelia valaisiana</i>	-/+	<i>Coxiella burnetii</i> strain Nine Mile Q	-/+
<i>Borrelia azfeli</i> strain P-Ko/1984	-/+	<i>Theileria annulata</i>	-
<i>Borrelia bavariensis</i>	-/+	<i>Rickettsia conorii</i> . strain Moroccan	-/+
<i>Borrelia bisetti</i>	-/+	<i>Leptospira</i>	-
<i>Borrelia burgdorferi</i> sensu stricto strain IRS	-/+	<i>Treponema phagedenis</i>	-
<i>Borrelia burgdorferi</i> sensu stricto strain B31	-/+	<i>Trypanosoma cruzi</i>	-
<i>Japanese encephalitis</i>	-	<i>Plasmodium falciparum</i> (strain 3D7)	-
West Nile Virus (Heja)	-	Tick Borne Encephalitis Virus (TBEV) strain Neudorfl	-/+
West Nile Virus (NY99)	-	Zika virus strain FB-GWUH-2016	-
West Nile Virus (Ug37)	-	Zika Virus (African strain)	-
Yellow Fever Virus (17D strain)	-	Zika Virus (Asiatica strain PF13/251013-18)	-
Yellow Fever Virus French Neurotropic	-	Zika Virus strain 11474/16 (French Polynesia)	-
Chikungunya virus strain S27 Petersfield	-	Zika Virus strain 11468/16(French Polynesia)	-
Chikungunya virus Martinique isolate	-	Dengue 1 virus strain Hawaii	-
Chikungunya virus strain F24	-	Dengue 2 virus strain New Guinea C	-
Chikungunya virus WHO IS (R91064)	-	Dengue 3 virus strain H87	-
St Louis Encephalitis virus	-	Dengue 4 virus strain H241	-

Analytical reactivity

The Vitassay qPCR Tick Borne Disease kit's analytical reactivity was evaluated against the following strains showing positive results:

TBEV strain Neurdofl	<i>Borrelia bisettiae</i>
<i>Borrelia garinii</i> (strains PHei, PWudII, PRef and PLa)	<i>Borrelia burgdorferi</i> sensu stricto (strains B31, IRS and PBre)
Synthetic sequence of <i>Babesia microti</i>	<i>Rickettsia coronii</i> strain Moroccan
Synthetic sequence of <i>Babesia divergens</i>	<i>Borrelia garinii</i> strains OspA Typ3
Synthetic sequence of <i>Ehrlichia chafeensis</i>	<i>B. japonica</i>
Synthetic sequence of <i>Ehrlichia muris</i>	<i>B. lusitanae</i> (strain Poti B2)
<i>Borrelia azfeli</i> (P-Ko/1984 and PVPM strains)	<i>B. spielmanii</i> (strain PSigII)
<i>Borrelia bavariensis</i> (strain PBi)	<i>Borrelia valaisiana</i> (strain VS116)
<i>Borrelia bissettii</i> (strain PGeb)	<i>B. hermsii</i>
<i>B. miyamotoi</i>	<i>Anaplasma phagocitophylum</i>
<i>Coxiella burnetii</i> strain Nine Mile Q	

Compatibles real-time PCR equipment

Vitassay qPCR Tick Borne Disease has been validated on the following equipment:

- 7500 Fast Real-Time PCR System (Applied Biosystems)¹
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Cobas LightCycler480 System (Roche)
- Step One Real Time PCR System (Applied Biosystems)

I: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All obtained results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with RNA/DNA extracted from blood, serum, tissue samples and microbiological culture from ticks, biopsy skin, cerebrospinal fluid, and synovial fluid. The use of other samples has not been established.
- The test quality depends on the sample's quality; proper RNA/DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Tick Borne Disease, either samples containing high concentrations of target RNA/DNA or contamination due to PCR products from previous reactions.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Excycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Excycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™ 5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cycler	DTprime Real-time Detection Thermal Cycler ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®48 0II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliografía

- Eldin C, Mélenotte C, Mediannikov O. et al. (2017). From Q fever to *Coxiella burnetii* Infection: A Paradigm Change. *Clinical Microbiology Reviews*. 30(1):115-190.
- Gayle A and Ringdahl E. (2001). Tick-Borne diseases. *American Family Physician*. 64(3):461-466.
- Guzman N, Yarrarapu SNS and Beidas SO. (2021). *Anaplasma Phagocytophilum*. In: StatPearls Publishing, Treasure Island (FL). PMID: 30020713.
- Johnson DK, Schiffman, EK, Davis JP, et al. (2015). Human Infection with *Ehrlichia muris*-like Pathogen, United States, 2007-2013. *Emerging Infectious Diseases*. 21(10), 1794-1799.
- Kawahara M, Ito T, Suto C, et al. (1999). Comparison of *Ehrlichia muris* Strains Isolated from Wild Mice and Ticks and Serologic Survey of Humans and Animals with *E. muris* as Antigen. *Journal of Clinical Microbiology*. 37(4):1123-1129.
- Kazar J. (2005). *Coxiella burnetii* Infection. *Annals of the New York Academy of Science*. 1063(1):105-114.
- Madison-Antenucci S, Kramer LD, Gebhardt LL and Kauffman E. (2020). Emerging Tick-Borne Diseases. *Clinical microbiology Reviews*. 33(2):e00083-18.
- Mansfield KL, Johnson N, Phipps, LP, et al. (2009). Tick-borne encephalitis virus – a review of an emerging zoonosis. *Journal of General Virology*. 90(8):1781-1794.
- Merhej V and Raoult D. (2011). Rickettsial evolution in the light of comparative genomics. *Biological Reviews*. 86(2):379-405.
- Parola P, Paddock CD, Socolovschi C, et al. (2013). Update on Tick-Borne Rickettsioses around the world: a Geographic Approach. *Clinical Microbiology Reviews*. 26(4):657-702.
- Raoult D, Marrie TJ and Mege JL. (2005). Natural history and pathophysiology of Q fever. *The Lancet Infectious Diseases*. 5(4): 219-226.
- Rożej-Bielicka W, Stypułkowska-Misiurewicz H and Gołąb E. (2015). Human Babesiosis. *Przegl Epidemiol*. 69(3):489-494.
- Schwan, TG, Raffel, SJ, Schrumf ME and Porcella SF. (2007). Diversity and Distribution of *Borrelia hermsii*. *Emerging Infectious Diseases*. 13(3):436-442.
- Wood H and Artsob H. (2012). Spotted Fever Group Rickettsiae: A Brief Review and a Canadian Perspective. *Zoonoses Public Health*. 59(2):65-79.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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