

VA Vitassay

identifying pathogens worldwide



VA Vitassay

OUR COMPANY

ABOUT US

Vitassay Healthcare


More than ten years ago, we set out with a simple mission: to help improve public health around the world.

Since 2015, we have partnered with healthcare professionals by providing reliable diagnostic solutions that make a real difference in people's lives, while striving to reduce our environmental footprint every step of the way.

Every Vitassay product is the result of a shared commitment to innovation, quality, and user experience. Developed by R&D experts and manufactured under strict quality standards and **ISO 13485** certification, our solutions are designed to deliver confidence when it matters most.

Because **in diagnostics, every result counts**. That's why our easy-to-use tests provide fast, robust, and accurate answers, helping healthcare professionals make informed decisions and improve patient outcomes worldwide.

OUR JOURNEY



PRODUCTS FOR HELPING GLOBAL HEALTH

RAPID TEST

- Faecal samples: Gastrointestinal infections and tumour markers
- Urine samples • Respiratory samples

REAL-TIME PCR

- Respiratory • Gastrointestinal tract • Vector - Borne and Tropical diseases
- Sexual diseases • Antimicrobial resistance • Immunocompromised

TURBIDIMETRY

- Quantitative detection of inflammatory and tumour markers



OUR PURPOSE

To provide the finest solutions to improve health worldwide.



OUR TARGET

To be a company of reference and global reach, specialised in products for rapid and accurate diagnosis.



OUR STRENGTHS

To offer personalised attention and join efforts to identify needs and offer support in every step.

VITASSAY RAPID TEST



Faecal

Bacteria Single Card X25

7355020	H. pylori
7355026	E. coli
7355028	Campylobacter
7355029	Salmonella
7355030	Salmonella typhi
7355022	Clostridium difficile antigen GDH

Bacteria Double Card X25

7455021	H. pylori + Transferrin
7455023	Clostridium difficile Toxin A + B
7455027	EHEC VT1 + VT2
7455031	Salmonella typhi + paratyphi
7455032	Yersinia enterocolitica O:3 + O:9

Parasites Single Card X25

7355033	Crypto
7355034	Giardia
7355036	Entamoeba

Parasites Double Card X25

7455035	Crypto + Giardia
---------	------------------

Virus Single Card X25

7355009	Adenovirus
7355008	Rotavirus
7355011	Astrovirus
7355013	Norovirus
7355018	Enterovirus

Virus Double Card X25

7455010	Rotavirus + Adenovirus
7455014	Norovirus GI + GII
7455015	Rotavirus + Norovirus

Urine

Single Card X25

7355047	Strep. pneumoniae
7355044	Legionella

Double Card X25


7455049	Strep. pneumoniae + Legionella
---------	--------------------------------

FAECAL TEST PROCESS

01 Sample collection



02 Close and shake



03 Break the cap



04 Dispense 4 drops



05 Wait 10 minutes



TEST RESULTS

Positive



Negative



Invalid



Tumour MarkersSingle Card  X25

7355001 FOB

7355004 Calprotectin

7355046 Lactoferrin

Double Card  X25

7455003 FOB + Transferrin

7455005 Calprotectin 50 + 200

7455006 Calprotectin + Lactoferrin

RespiratorySingle Card  X25

7355038 Strep A

7355039 RSV

7355040 Adenovirus Resp.

7355052 SARS-CoV-2

Double Card  X25

7455041 RSV + Adenovirus Resp.

7455042 Influenza A + B

ALL IN ONE**Faecal**Bacteria Multitest  X10

7715024 Clostridium difficile GDH + Toxin A + B (3 in 1)

Parasites Multitest  X10

7715037 Crypto + Giardia + Entamoeba (3 in 1)

Virus Multitest  X10

7715012 Rotavirus + Adenovirus + Astrovirus (3 in 1)



7715016 Rotavirus + Adenovirus + Norovirus (3 in 1)

7715017 Rota + Adeno + Astro + Noro (4 in 1)

7715050 Rota + Adeno + Astro + Noro + Enterovirus (5 in 1)

Respiratory X10

7715043 Influenza A + B + RSV + Adenovirus (4 in 1)

7715053 SARS-CoV-2 + Influenza A + B + RSV
+ Adenovirus (5 in 1)**FAECAL TEST PROCESS****01 Sample collection****02 Close and shake****03 Insert the vial****04 Screw the lid****05 Wait 10 minutes****TEST RESULTS****+ Positive****- Negative****✗ Invalid**

VITASSAY REAL-TIME PCR

Vitassay real-time PCR assays are your best partners for the molecular diagnosis.

- Lyophilised master mix with all the ingredients for a successful PCR/RT-PCR.
- Storage and transport at room temperature.
- Long shelf-life: > 12 months.
- Compatibility with the majority of real-time PCR detection systems.



Complete your assays with:

Vitassay Microorganism-NA isolation kit.

Nucleic acids isolation kit intended to be used in automated systems based on magnetic-particle technology.

Ref. 7071001

TEST PROCESS

01 Resuspend



02 Add the sample




03 Fit the cover



04 Place in the thermal cycler



05 Interpret the results



Respiratory Bacteria

	Vitassay qPCR	Target
7041020	Bordetella pertussis + B. parapertussis + B. holmesii	<i>Bordetella pertussis</i> , <i>B. parapertussis</i> , <i>B. holmesii</i>
7041048	MTBC/NTM	<i>M. tuberculosis</i> , <i>M. tuberculosis</i> Complex, Non-tuberculous mycobacteria
7041049	Primary Atypical Pneumonia	<i>C. pneumoniae</i> , <i>M. pneumonia</i> , <i>L. pneumophila</i>
7041065	HSM	<i>H. influenzae</i> , <i>S. pneumoniae</i> , <i>M. catarrhalis</i>
7041071	MTB complex	<i>M. tuberculosis</i> Complex

Respiratory Virus

7041021	Flu A + Flu B	Influenza A, Influenza B
7041022	RSV A + RSV B	RSV A, RSV B
7041025	Parainfluenza 1/3+2/4	Parainfluenza virus types 1, 2, 3, 4
7041027	Flu A + Flu B + RSV	Influenza A, Influenza B, Respiratory Syncytial virus A/B
7041028	AMB	Adenovirus, Metapneumovirus, Bocavirus
7041031	Type I / Flu	Influenza A (H1N1)pdm09, H3N2
7041032	Type II / Flu	Influenza A (H1N1)pdm09, H3N2, H5N1, H7N9
7041040	HCoV	Coronavirus 229E, NL63, OC43, HKU1
7041041	Rhinovirus & Enterovirus	Rhinovirus, Enterovirus
7091046	SARS-CoV-2	SARS-CoV-2 (<i>ORF1ab</i> and <i>N</i> genes)
7091052	Flu + RSV + SARS CoV-2 N1/N2	Influenza A/B, RSV A/B, SARS-CoV-2 (<i>N1/N2</i> genes)
7091054	SARS-CoV-2 + Flu A + Flu B	SARS-CoV-2 (<i>N1/N2</i>), Influenza A, Influenza B
7091060	SARS-CoV-2 (<i>ORF1ab</i> , <i>E</i> & <i>N</i>)	<i>ORF1ab</i> , <i>E</i> , <i>N</i>

Respiratory Fungi

7041026	Pneumocystis jirovecii	<i>Pneumocystis jirovecii</i> (quantitative)
7041076	Aspergillosis	<i>A. fumigatus</i> , <i>A. flavus</i> and <i>A. terreus</i>

Vector - Borne & Tropical Diseases

7041001	Chikungunya	Chikungunya virus
7041002	Dengue	Dengue virus
7041003	Zika	Zika virus
7041004	West Nile	West Nile virus
7041005	Zika+Dengue+Chikungunya	Zika virus, Dengue virus, Chikungunya virus
7041006	Dengue 1+2+3+4	Dengue virus type 1, 2,3,4
7041050	Malaria 5	<i>P. falciparum</i> , <i>P. vivax</i> , <i>P. ovale</i> , <i>P. malariae</i> , <i>P. knowlesi</i>
7041038	Japanese Encephalitis	Japanese Encephalitis virus
7041062	Tick Borne Disease	<i>Borrelia</i> , <i>Coxiella</i> , <i>Anaplasma</i> + <i>Babesia</i> , <i>Rickettsia</i> , <i>Ehrlichia</i> +TBEV

Gastrointestinal Tract Bacteria

	Vitassay qPCR	Target
7041008	<i>Clostridium difficile</i>	<i>Clostridium difficile</i>
7041009	<i>Clostridium difficile</i> toxins A/B	<i>Clostridium difficile</i> toxin A/B
7041042	<i>Clostridium difficile</i> toxins A+B	<i>Clostridium difficile</i> toxin A, toxin B
7041010	<i>Campylobacter</i> + <i>Salmonella</i> + <i>Shigella</i> /EIEC	<i>Campylobacter</i> , <i>Salmonella</i> , <i>Shigella</i> /EIEC
7041011	<i>Campylobacter</i> + <i>Salmonella</i> + <i>Yersinia enterocolitica</i>	<i>Campylobacter</i> , <i>Salmonella</i> , <i>Yersinia enterocolitica</i>
7041059	EHEC, EPEC & EIEC	Enterohemorrhagic <i>E. coli</i> (EHEC), Shiga toxin-producing <i>E. coli</i> (STEC), Enteropathogenic <i>E. coli</i> (EPEC), Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i>
7041064	<i>E. coli</i> pathotypes	<i>E. coli</i> (EHEC), Shiga toxin-producing <i>E. coli</i> (STEC), Enteropathogenic <i>E. coli</i> (EPEC), Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i> , Enterotoxigenic <i>E. coli</i> (ETEC)

Gastrointestinal Tract Virus

7041014	Astrovirus	Astrovirus
7041017	Rotavirus	Rotavirus
7041018	Sapovirus	Sapovirus
7041070	Noro (I&II)	Norovirus GI, Norovirus GII

Gastrointestinal Tract Parasites

7041012	<i>Cryptosporidium</i> + <i>Entamoeba histolytica</i> + <i>Giardia lamblia</i>	<i>Cryptosporidium</i> spp., <i>Entamoeba histolytica</i> , <i>Giardia lamblia</i>
7041029	<i>Dientamoeba fragilis</i>	<i>Dientamoeba fragilis</i>
7041030	<i>Blastocystis hominis</i> + <i>Dientamoeba fragilis</i>	<i>Blastocystis hominis</i> , <i>Dientamoeba fragilis</i>

Sexual Diseases

7041035	Sexually Transmitted Infections	<i>Neisseria gonorrhoeae</i> , <i>Chlamydia trachomatis</i> , <i>Mycoplasma genitalium</i> , <i>Trichomonas vaginalis</i> , <i>Ureaplasma urealyticum</i> , <i>Ureaplasma parvum</i> , <i>Mycoplasma hominis</i>
7041036	Genital ulcer	HSV-1, HSV-2, <i>Treponema pallidum</i>
7041037	HPV (16+18)	HPV 16, HPV 18
7041051	Vaginosis	<i>C. albicans</i> , <i>G. vaginalis</i> , <i>T. vaginalis</i>
7041063	High-Risk HPVs	HPV 16, HPV 18, HPV 35/58/66, HPV 33/45/51, HPV 52/59/68, HPV 31/39/56
7041080	MonkeyPox	Mpox virus

Antimicrobial Resistance

	Vitassay qPCR	Target
7041034	H. pylori ClariRes	<i>H. pylori</i> Clarithromycin resistant, <i>H. pylori</i> wild type
7041043	Vancomycin Resistance	<i>VanA</i> , <i>VanB</i>
7041057	Carbapenemase-generating Enterobacteriaceae	<i>NDM</i> , <i>VIM</i> , <i>OXA-48</i> , <i>OXA-48-like</i> , <i>KPC</i> , <i>IMP</i>
7041066	MRSA	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA), methicillin-sensitive <i>S. aureus</i> (MSSA), methicillin-resistant coagulase-negative <i>Staphylococci</i> (MRCoNS)
7041068	NG CiproRes	<i>N. gonorrhoeae</i> , <i>gyrA</i> wild type, <i>gyrA</i> mutant (ciprofloxacin resistant)
7041069	β -lactamases + Colistin resistance	<i>blaCTX</i> , <i>blaTEM</i> , <i>blaSHV</i> , <i>mcr-1</i>
7041077	Macrolide Res- 23S rRNA Mutations	<i>Mycoplasma genitalium</i> 23S rRNA Wild type, <i>M. genitalium</i> 23S rRNA mutant (macrolide resistant)

Immunocompromised

7041045	Viral Meningitis	HSV-1, HSV-2, Varicella Zoster Virus
7041047	BK & JC	BK Virus, JC Virus
7041067	Human Herpes Virus	HHV-6, HHV-7, HHV-8
7041072	BK (q)	BK Virus (quantitative)
7041075	Entero + Parecho + Parotitis	Enterovirus, Parechovirus, Mumps virus
7041073	Gram-negative Nosocomial Agents I	<i>E. coli</i> , <i>A. baumannii</i> , <i>Enterobacter</i> spp.
7041074	Gram-negative Nosocomial Agents II	<i>K. pneumoniae</i> , <i>P. aeruginosa</i> , <i>P. mirabilis</i>

DISCOVER OUR GROWING PORTFOLIO

Latest developments for Research Use Only

RUO

7041078R	Toxoplasmosis RUO	<i>Toxoplasma gondii</i>
7041081R	Cytomegalovirus RUO	Cytomegalovirus (quantitative)
7041082R	Andes virus RUO	Hantavirus Andes strain
7041083R	Measles & Rubella RUO	Measles virus, Measles vaccine strains, Rubella virus

VITASSAY TURBIDIMETRY

Turbidimetry is the perfect quantitative technique for analysing multiple samples in the routine diagnosis or patients' follow-up.

Biomarker screening and early detection is key for a favourable patient prognosis, and a quick treatment onset. Vitassay turbidimetry: the smart solution for global health.

Turbidimetry Kits

- 7115001** FOB Turbidimetric Assay.
7115001-100V FOB Turbidimetry vial (100 pc).
- 7115002** Calprotectin Turbidimetric Assay.
7115001-100V Calprotectin Turbidimetry vial (100 pc).
- 7115004** Transferrin Turbidimetric Assay.
7115001-100V Transferrin Turbidimetry vial (100pc).
- 7115005** Elastase Turbidimetric Assay.
7115001-100V Elastase Turbidimetry Kit vial (100pc).



TEST PROCESS

01 Open the sample vial




02 Sample collection




03 Homogenise & Centrifuge



04 Run the Analyser



05 Interpret the results



ADVANTAGES

Fully automated system with intuitive interpretation.




At Vitassay we are at the customer's service. Since 2015, our products have reached every corner of the world thanks to our distributors.

VA RAPID TEST

- Room temperature
- Fast result
- Easy interpretation
- Unique format: 5-in-1
- Long shelf-life

VA REAL-TIME PCR

- Room temperature
- Lyophilised & Ready-to-use
- User-friendly
- Highest accuracy
- Long shelf-life



VA TURBIDIMETRY

- Quantitative
- High stability
- Multiple samples



Vitassay Healthcare, S.L.U

Parque Tecnológico Walqa
Ctra. N-330 Km. 566
22197 Huesca (Spain)
Ph. (+34) 974 001 193

info@vitassay.com
www.vitassay.com

ISO 13485 Certified

VITASSAY/GEN/Sep26
F09-83 EN Ed01